
Dispergierwerkzeuge / Dispersion tools / Outils dispersants

S10D - 7G-KS-65
S10D - 7G-KS-110

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OPERATING INSTRUCTIONS

MODE D'EMPLOI

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Safety instructions

- No liquid may get into the drive unit of the disperser.
- Adhesion of the disperser equipment or the rotating accessories can get into the medium you are working on.
- Follow the health and safety regulations accident prevention regulations.
- Do not use any damaged dispersion tools e.g. with hairline cracks or splits.
- Note the operating instructions of the disperser unit.
- Only operate in a vertical working position.

Product information

- The packaging unit endorses 12 or 25 pieces.
- The dispersion tool consists of shaft tube and rotor shaft.
- Dispersion tools made of plastics are not sterile and not free of DNA or RNA.
- The tools are ideal for applications, with which no cross contamination is permissible.
- The tools are ideally appropriate for a single-use. However they can be used, depending upon the application, several times too.
- In this case please consider the cleaning references.

Application instruction

- Attention!** Never run the dispersion tool dry. The bearing will be destroyed without cooling by the medium. Always ensure that the dispersing element does not run dry in some cases should spouts develop.
- The fill level of the medium must be between the MIN and MAX marks.
- Medium may leak out of the overflow hole around the MAX mark.
- The dispersion tool is not suitable for continuous operation.
- Dispersion causes the medium to heat up.
- Store the dispersing elements in a dry place, in their packaging and at room temperature.
- The solid bodies in samples must not be frozen or too hard.
- The optimal rotating frequency, usually 20,000-30,000 rpm, must be determined by attempts.
- Likewise the optimal dispersion duration must be determined by attempts. Usually a few seconds are sufficient, maximum duration is one minute.
- Longer application times bring no improvements, increase however the sample temperature substantially.
- Dispose of the tools after use according to the regulations. Observe the laboratory regulations.

Designation	Application	Max. permissible number of revolution (1/min)	Optimal speed range (1/min)	Working volume (ml)
S10D-7GKS-55	Liquid/ liquid or soft fabric	30.000	20.000-30.000	0,25-25
S10D-7GKS-110	Liquid/ liquid or soft fabric	30.000	20.000-30.000	0,25-25

Materials and chemical resistance

- The shaft tubes are made of polycarbonate (PC), the rotor shafts are made of polysulfone (PSU).
- The plastic materials used are approved by FDA, (Food and Drug Administration Approval).
- The dispersion tool has a good stability against weak acids, chlorides, hypochlorides and many other chemicals.

Cleaning, disinfection and sterilization

- The plastic dispersion tools are not sterile, but you can sterilize them.

Method of sterilisation	Not allowed procedures
UV-irradiation	UV-irradiation
Hot-air sterilisation 160-190 °C	Hot-air sterilisation 160-190 °C
Autoclave up to 125 °C with 2 bar	Autoclave up to 125 °C with 2 bar
Ethylen oxide	Ethylen oxide
Formaldehyde (4)	Formaldehyde (4)
Glycerin (2%)	Glycerin (2%)
Gammabestrahlung bis zu 5 Mrad	Gammabestrahlung bis zu 5 Mrad
Alkohol (70%)	Alkohol (70%)
Hypochloride (6%)	Hypochloride (6%)
Formaldehyd (40%)	Formaldehyd (40%)
Radialwash	Radialwash
Count-off	Count-off
Water/ Ethanol/ SDS	Water/ Ethanol/ SDS

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