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Warranty

In accordance with IKA® warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs. The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

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

Safety instructions

- Read the operating instructions in full before starting up and follow the safety instructions.
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- The dispersion tool may heat up during operation. In particular in the processing of hot medium, the tool must be cooled before it can be removed from the drive.
- When in operation the dispersion tool can become hot. If the dispersion tool is not inserted into the drive flange correctly or the medium is hot, it can become extremely hot.
- Don't touch rotating parts during operation.
- No liquid may get into the drive unit of the dispenser.
- The work has to be interrupted immediately if you notice unusual noise and/or emission of liquid from the upper rinsing drilling. The seal must be controlled, and replaced if necessary.
- While working with the dispenser, the user must wear his personal protective equipment, in particular safety goggles.
- Follow the health- and safety regulations accident prevention regulations.
- Note the operating instructions of the dispenser unit.
- Note the operating instructions of the rotating accessories can get into the medium you are working on.
- Abuse of the dispenser equipment or the rotating accessories can get into the medium the bearing will be destroyed.
- Attention!** The dispersion tool must not be used for continuous operation.
- Dispensing causes the medium to heat up.
- Use the use of hot medium to heat up and relating frequency must be determined by attempts.
- The function of the dispenser tool depends on the condition of the sharp edges of rotor and stator. Abrasive media can round off these edges fast, whereby the dispense effect diminishes.

Product information

The shaft bearing of the dispersion tool consists of a two ball bearings and a single-thrust slide ring seal (SGC). Content of package: shaft and dispersion head. The plastic material used are approved by FDA (Food and Drug Administration Approval).

Application instructions

Application: Wet crushing, suspension, emulsifying (Batch-operation).

Spare parts list

Item	Designation	Item	Designation
1	Cylindrical pin	12	O-ring 9,25 x 1,78
4	Deepgroove ball bearing	14	Deepgroove ball bearing
6	Seeger circlip ring	15	Seeger circlip ring
8	Stator	2001	Rotor shaft, kompl.
9	Rotor	2002	Friction ring seal, kompl.
10	O-Ring 4,7 x 1,42	2003	Shaft, tube, kompl.

für S 25 KV 2802: 8* – Protection ring, 9* – Dissäher

Maintenance and cleaning

Allowed procedures
 By germ-killing disinfectants alcohol, phenol formalin...
 Chemical procedures
 Disinfectant remainders are to be removed with germ-free water.
Hot air
 Germ-killing by hot air with 160 to 190 °C (approx. 30 min). Life span of the bearings shortens.
 For pre-cleaning the dispersion tool operate it in a dissonant for dissolving the residues of the substance. Due to the large flow rate, rotor and stator are cleaned to a large extent. The dispersion tool must be divided and cleaned immediately after working.
 For cleaning purposes disassemble the dispersion tool as described below (shaft key included in the tool set). The rotor is held fast by means of the generator key and screwed off from above by turning the shaft using the shaft key.
 • Unscrew the stator (8) from the shaft (2) left-hand thread. Insert transverse rib of the shaft key into the front-end slot of the shaft tube. By using the generator key the stator may now be screwed off in a clockwise direction. In doing so, the dispersion tool and the stator should be held up in vertical position, thus preventing the spring of the slide ring sealing from falling out.
 • Tap the shaft (2001) and the ball bearing (4) downward in the direction of the stator thread out of the shaft tube.
 • In this position the parts of the slide ring seal can be removed from the shaft (2001).
 • The slide rings and O-rings must be replaced in the case of wear. The sliding surfaces of the file and counter ring must be clean.
 The assembly of the dispersion tool takes place in reverse order. Make certain that the O-ring (10) and the thrust collar (3004) are properly fitted during assembly.
Attention! The rotor must not be tightened too strongly, because otherwise the thread will be damaged.

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