
## IKA® RO 5/10/15

## IKA® RT 5/10/15


Betriebsanleitung	DE	3
Ursprungssprache		
Operating instructions	EN	9
Mode d'emploi	FR	15
Руководство пользователя	RU	21
使用说明	ZH	27


Indicaciones de seguridad	ES	33	Bezpečnostní pokyny	CS	42
Veiligheidsinstructies	NL	34	Biztonsági utasítások	HU	43
Norme di sicurezza	IT	35	Varnostna opozorila	SL	44
Säkerhetsanvisningar	SV	36	Bezpečnostné pokyny	SK	45
Sikkerhedshenvisninger	DA	37	Ohutusjuhised	ET	46
Sikkerhetsanvisninger	NO	38	Drošības norādes	LV	47
Turvallisuusohjeet	FI	39	Nurodymai dėl saugumo	LT	48
Instruções de segurança	PT	40	Инструкции за безопасност	BG	49
Wskazówki bezpieczeństwa	PL	41	Indicații de siguranță	RO	50
			Υποδείξεις ασφάλειας	EL	51

IKA® RO 5/10/15

IKA® RT 5/10/15

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	Page		
Control panel and display	02	Setting the safe temperature limit	11
Contents	09	Stirring function	11
Declaration of conformity	09	Eco mode	12
Warranty	09	Operating modes	12
Explication of warning symbols	09	Direction of rotation reversal	12
Safety instructions	10	Function heating	12
Unpacking	11	Maintenance and cleaning	12
Correct use	11	Error codes (RT 5/10/15)	13
Commissioning	11	Accessories	13
Safe temperature limit	11	Technical data	14

## Declaration of conformity

We declare under our sole responsibility that this product is in compliance with the regulations 2014/35/EU, 2014/30/EU and 2011/65/EU and conforms to the standards or normative documents: EN 61010-1, EN 61010-2-010, EN 61010-2-051, EN 61326-1, EN 60529 and EN ISO 12100.

## Warranty

In accordance with **IKA**® warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

## Explication of warning symbols

General hazard.

This symbol identifies information **that is of absolute importance to ensure health and safety**. Failure to observe this information may be detrimental to health or may result in injuries.

This symbol indicates information **which is important for ensuring that the appliance functions without any technical problems**.

Failure to observe this information could damage the appliance.

This symbol indicates information **which is important for proper use of the appliance and / or ensuring that the appliance functions correctly**.

Failure to observe this information can lead to inaccurate results.

ATTENTION - Risk of damage due to magnetism.

DANGER - Reference to the endangerment by a hot surface

## Safety instructions

### For your protection

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Socket must be earthed (protective ground contact).

#### Attention - Magnetism!

Effects of the magnetic field have to be taken into account (e.g. data storage media, cardiac pacemakers ...).

#### RT

#### Risk of burns!

Exercise caution when touching the housing parts and the heating plate.

The heating plate can reach temperatures in excess of 125 °C. Pay attention to the residual heat after switching off.

#### RT

Ensure that the mains power supply cable does not touch the heating base plate.

- Wear your personal protective equipment in accordance with the hazard category of the media to be processed. Otherwise there is a risk from:
  - splashing and evaporation of liquids
  - ejection of parts
  - release of toxic or combustible gases.
- Set up the appliance in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- The feet of the appliance must be clean and undamaged.
- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
- Gradually increase the speed.
- Reduce the speed if
  - the medium splashes out of the vessel because the speed is too high
  - the appliance is not running smoothly
  - the container moves on the base plate.

#### RT

**Caution!** Only process and heat up any media that has a flash point higher than 180 °C.

The safe temperature limit must always be set to at least 25 °C lower than the fire point of the media used.

Beware of hazards due to:

- flammable materials
- combustibles media with a low boiling temperature
- glass breakage
- incorrect container size
- overfilling of media
- unsafe condition of container.
- The appliance may heat up when in use.
- The base plate can heat up due to the action of the drive magnets at high motor speeds, even if the heater is not operational.

- Process pathogenic materials only in closed vessels under a suitable extractor hood. Please contact IKA® if you have any questions.

**Do not** operate the appliance in explosive atmospheres, with hazardous substances or under water.

Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.

- Please observe the operating instructions for any accessories used.
- Safe operation is only guaranteed with the accessories described in the „Accessories“ chapter.
- Accessories must be securely attached to the device and cannot come off by themselves. The centre of gravity of the assembly must lie within the surface on which it is set up.
- Always disconnect the plug before fitting accessories.
- The appliance can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- The socket for the mains cord must be easily accessible.
- It may be possible for wear debris from rotating accessory parts to reach the material being processed.
- When using PTFE-coated magnetic bars, the following has to be noted:

*Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300 °C - 400 °C. Only elementary fluorine, chlorotrifluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.*

(Source: Römpps Chemie-Lexikon and „Ullmann“ Volume 19)

### For protection of the equipment

- The appliance may only be opened by experts.
- The voltage stated on the type plate must correspond to the mains voltage.
- Do not cover the device, even partially e.g. with metallic plates or film. This results in overheating.
- Protect the appliance and accessories from bumps and impacts.
- Ensure that the base plate is kept clean.
- Observe the minimum distances between the devices, between device and wall and minimum distances (min. 800 mm) above the assembly, see Fig. 2.


Fig. 2

## Unpacking

### • Unpacking

- Please unpack the device carefully
- In the case of any damage a detailed report must be sent immediately (post, rail or forwarder)

### • Delivery scope

#### RO

- Magnetic stirrer  
RO 5/10 or 15
- power supply unit
- operating instruction

#### RT

- Heating magnetic-stirrer RT 5/10 or 15
- mains cable
- operating instruction

## Correct use

### • Use

- For mixing and / or heating liquids

### • Range of use (inside only)

- Laboratories    - Schools
- Pharmacies    - Universities

The device is suitable for use in all areas except:

- Residential areas
- Areas that are connected directly to a low-voltage supply network that also supplies residential areas.

The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly contrary to the manufacturer's specifications or if the appliance or the printed circuit board are modified parties.

## Commissioning

### RO/RT

Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.

The device is switched on and off via the switch (A) located on the right-hand side of the device. On the front panel a small arrow labelled Power above the display (B) indicates the location of the switch.

After the device is switched on a display self-test is performed. The display cycles through the following information:

- All segments light up
- The software version and the selected mode are displayed
- The selected "Eco Mode" is displayed  
(E when ECO mode is active)
- If the drive has not been started the display reads „OFF“. If the drive has been started the selected speed is displayed.

## Safe temperature limit

### RT

The max. attainable heating plate temperature is restricted to 130 °C by a permanently set safety temperature limit. Once this limit has been attained, the device switches off the heating.

The safe temperature limit must always be set at least 25 °C lower than the flash point of the media to be processed!

## Setting the safe temperature limit

### RT

The safety temperature is set as follows:

Switch the device on using the On/Off switch on the right side of the device.

As soon as the display shows SAFE, press and hold down the Temp membrane key and use

the +/- membrane keys (**H** or **I**) to set the desired safety temperature.

The set safety temperature is displayed.

Once the Temp membrane key is released, the safety temperature is saved.

## Stirring function

### RO/RT

The stirring function is started and stopped by pressing the „Start/Stop“ button (**C**).

Pressing the +/- buttons (**E** or **D**) increases or reduces the speed.

## ECO mode

### RO/RT

The heat generated in the device can be reduced by lowering the power of the input coils.

Pressing the "ECO Mode" (**F**) key to reduce the current power setting. ECO mode is indicated by an "E" preceding the speed. When the power is set to ECO mode, the speed is limited to 600 rpm.

To make the stirring power more independent from the speed range, press the "ECO Mode" membrane key (**F**) and then select one of the following setting options.

- P50 (Stirring power 50% independent of the speed range)
- P75 (Stirring power 75% independent of the speed range)
- P100 (Stirring power 100% independent of the speed range)
- (Standard stirring performance with the stirring power regulated by the speed)

## Operating modes

### RO/RT

To toggle between modes, press the Start/Stop key (**C**) when switching on.

#### Operating Mode "A"

All the parameter settings are retained when the device is switched off or disconnected from the mains.

The "Start/Stop" (**C**) key starts the device with the set parameters.

#### Operating Mode "B"

All parameter settings, including "Start/Stop", are saved when the device is switched off.

#### Operating Mode "D"

In the D mode the device operates as per the A mode but with the following exception:

- Switch on and then press the "Temp" membrane key to confirm the safety temperature.
- The value flashes until it is confirmed.

*Factory setting: mode A*

## Direction of rotation reversal

### RO/RT

The (**G**) key activates the "Reverse Rotation" function.

When the „Reverse Rotation“ function is switched on, a decimal point (**K**) appears in the LCD display.

The running time (30, 60, 90 or 120sec) is set by pressing the (**G**) key followed by the +/- keys.

After the preset running time, the drive is switched off and restarted at minimum speed in the opposite direction after a pause specific to the speed.

In this way the magnetic rod is captured and any liquid still rotating is gradually slowed down. Then the drive accelerates to the target speed.

Note: If the "Reverse Rotation" key (**G**) is pressed again, the drive will keep rotating in the same direction.

When the device is started up again, the drive will rotate in the direction last set.

## Function Heating

### RT


The heating function is started and stopped by pressing the „Temp“ button (**J**).

Pressing the +/- buttons (**H** or **I**) increases or reduces the temperature.

## Maintenance and cleaning

The equipment is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

### Cleaning

For cleaning disconnect the main plug..

Use only cleansing agents which have been approved by **IKA**® to clean devices.

Use to remove:

- |                        |   |
|------------------------|---|
| Dyes                   | Isopropyl alcohol                             |
| Construction materials | water containing ten side / isopropyl alcohol |
| Cosmetics              | water containing ten side / isopropyl alcohol |
| Foodstuffs             | water containing tenside                      |
| Fuels                  | water containing tenside                      |
- Do not allow moisture to get into the appliance when cleaning.
  - Wear protective gloves when cleaning the devices.
  - Please consult with **IKA**® before using any cleaning or decontamination methods, other than those recommended here.

### Ordering spare parts

When ordering spare parts, please give:

- Device type
- Manufacturing number, see type plate
- Software version
- Item number and designation of the spare part, see [www.ika.com](http://www.ika.com).

### Repair

**Please send instrument in for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.**

For this you should request the “**Decontamination Clearance Certificate**” from IKA®, or use the download printout of it from the IKA® website [www.ika.com](http://www.ika.com).

Return the instrument in its original packaging. Storage packaging is not sufficient. Also, please use suitable shipping package materials.

## Error Codes (RT 5/10/15)

Error-code	Cause	Effect	Solution
Er 3	Temperature inside device too high	Heating off	- Switch off device and allow to cool down
Er 6	Break in safety circuit	Heating off	- Plug in contact plug - Plug in contact thermometer or PT 1000 temperature sensor - Replace faulty connection cable, plug, or contact thermometer
Er21	Safety relay doesn't open during test	Heating off	- Switch off device <i>- Warning! Only to be carried out by authorised service personnel: Check PCB</i>
Er22	Safety temperature is too low during test	Heating off	- Switch off device <i>- Warning! Only to be carried out by authorised service personnel: Check PCB and safety temperature sensor</i>
Er25	Heating and switching element monitoring	Heating off	- Check the external temperature controller - Switch off device <i>- Warning! Only to be carried out by authorised service personnel: Carry out an internal test on the device to check the plug-in connector for the heating element, the safety temperature sensor or the PCB</i>
Er26	Difference between temperature of safety sensor and temperature of control sensor <b>control temperature</b> > (safety temperature + 40 K)	Heating off	- Switch off device <i>- Warning! Only to be carried out by authorised service personnel: Carry out an internal test on the device to check the plug-in connector for the temperature sensor</i>
Er31	Triac error detection doesn't work	Heating off	- Switch off device <i>- Warning! Only to be carried out by authorised service personnel: Check PCB</i>
Er46	Difference between temperature of safety sensor and temperature of control sensor <b>control temperature</b> > (control temperature + 40 K)	Heating off	- Switch off device <i>- Warning! Only to be carried out by authorised service personnel: Carry out an internal test on the device to check the plug-in connector for the temperature sensor</i>

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- Please contact the service department;
- Send the device for repair, including a short description of the fault.

## Accessories

- Stirring bar:            ∅ 8 mm; length 30 mm
- RSE                        Stirring bar remover

## Technical data

### RO

#### Power supply unit

Input

**V** 100 - 240

**A** 1,67

Output

**Vdc** 24; 40 W.LPS; (limited power source)

**Hz** 50/60

Protection class

2 (double insulated) 

Operating voltage

**Vdc** 24

	RO 5	RO 10	RO 15
<b>A</b>	0,5	1,0	1,5
<b>W</b>	12	24	36
<b>W</b>	2		
<b>W</b>	17		
<b>rpm</b>	0 - 1200 in 10 rpm-steps adjustable		
<b>rpm</b>	0 - 600 in 10 rpm-steps adjustable		
<b>mm</b>	120 x 450	180 x 450	270 x 450
<b>mm</b>	120 x 570 x 60	190 x 570 x 60	280 x 570 x 60
<b>kg</b>	3,0	4,0	7,0

Power consumption (max.)

Power consumption, standby operation

Power output

Speed ranges

P50, P70, P100, ----  
ECO-Mode

Set-up surface

Dimension:

Weight:

### RT

Unit

Operating voltage

**Vac** 230 ± 10%

**Vac** 115 ± 10%

Frequency

**Hz** 50 / 60

Speed ranges

P50, P70, P100, ----  
ECO-Mode

**rpm** 0 - 1000 in 10 rpm-steps adjustable

**rpm** 0 - 600 in 10 rpm-steps adjustable

Power consumption (max.)

Set-up surface

Dimension:

Weight:

Adjustment and display resolution

Surface temperature

	RT 5	RT 10	RT 15
<b>W</b>	185	395	600
<b>mm</b>	110 x 495	180 x 495	270 x 495
<b>mm</b>	120 x 610 x 60	190 x 610 x 60	280 x 610 x 60
<b>kg</b>	4,0	6,5	9,4
<b>K</b>	1		
<b>°C</b>	120*		
<b>°C</b>	50 - 150		

**W** 185

**mm** 110 x 495

**mm** 120 x 610 x 60

**kg** 4,0

**K** 1

**°C** 120\*

Safety temperature limit (adjustable)

**°C** 50 - 150

Fuse:

230 V

100 V / 115 V

	RO/RT 5	RO/RT 10	RO/RT 15
<b>A</b>	T 2 (5x20)	T 4 (5x20)	T 6,3 (5x20)
<b>A</b>	T 4 (5x20)	T 6,3 (5x20)	T 10 (5x20)

**A** T 2 (5x20)

**A** T 4 (5x20)

### RO/RT

Setting resolution

**rpm** ±5

Deviation between the stirrer points

**%** 0

Permissible duration of operation

**%** 100

Permissible ambient temperature

**°C** +5 to +40

Permissible relative humidity

**%** 80

EN 60529 protection class

IP 40

Operation at a terrestrial altitude

**m** max. 2000

Stirred quantity max. (water)

400 ml per stirrer point in the 600 ml beaker glass

in other vessels

Stirrer points

Distance from stirrer point to stirrer point

	RO/RT 5	RO/RT 10	RO/RT 15
<b>ltr</b>	2	4	6
<b>mm</b>	5	10	15
<b>mm</b>	90	90 x 90	90 x 90

**ltr** 2

**mm** 5

**mm** 90

 **\*Caution! The specified maximum hotplate temperature refers to instruments with a full load. Without a load, the (setup) plate temperature can reach 130 °C (+10 °C / -5 °C). Load to determine the specified values: 250 ml Erlenmeyer flask filled with 200 ml water.**

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DOMINIQUE DUTSCHER SAS
**IKA®-Werke GmbH & Co.KG**