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Contents

	Page
Explication of warning symbols	5
Safety instructions	5
Unpacking	6
Correct use	6
Set-up and Commissioning	6
Maintenance and cleaning	7
Accessories	7
Technical data	7
Warranty	7

Explication of warning symbols

General hazard

This symbol identifies information **that is of vital importance for safeguarding your health and safety**. Disregarding this information can lead to health impairment and injuries.

This symbol identifies information and notes that are of importance for correct handling of the unit.

Safety instructions

• **Read the operating instructions in full before starting up and follow the safety instructions.**

• Keep the operating instructions in a place where they can be accessed by everyone.

- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- The stand pillar must be inserted into the baseplate so that the guide groove is at the back. Other alignments are not permissible, as otherwise the stand can topple over when a unit is mounted on it.
- Set up the stand in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.

• Any small unevennesses that cause the stand to rock must be compensated for with the height-adjustable feet.

- The T 65 D ULTRA-TURRAX® dispersion unit must be lifted on to the guide head of the telescopic stand with the aid of an assistant, and secured with the two screws and toothed lock washers (supplied), in accordance with the T 65 D ULTRA-TURRAX® operating instructions. The unit must be mounted in front, vertically above the baseplate; a different alignment is not permissible, as otherwise the assembly could topple over.

- The overhead stirrer IKA® RW 47 D must be lifted on to the guide head of the telescopic stand with the aid of an assistant, and secured with the four screws (supplied), in accordance with the IKA® RW 47 D operating instructions. The unit must be mounted in front, vertically above the baseplate; a different alignment is not permissible, as otherwise the assembly could topple over.

- This telescopic stand must be fitted with a safety disconnection device, which stops the unit if the stand is moved upwards when the unit is operating.
- Switch the unit off before moving the telescopic stand upwards.
- When lowering the telescopic stand, take care that no parts of the body get between the safety switches (danger of crushing).
- Only switch the unit on when the tool (agitator or dispersion unit) is in the vessel.
- The handwheel bolt on the stand can come loose due to vibration, which would lead to the unit moving down. So, check the tightness of the handwheel bolt regularly when using the unit, and retighten it if necessary.
- Raise or lower the stand (after slackening the handwheel bolt) by applying the manual force as close as possible to the stand pillar. The further the force is applied from the pillar, the more force you need to apply, plus there is an increased risk of the equipment toppling over.
- Secure the vessel against spinning round, by using a clamp, for example.
- Beachten Sie die Betriebsanleitung des Dispergiergerätes T 65 D ULTRA-TURRAX® oder des Rührgerätes IKA® RW 47 D!
- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
- Protect the appliance and accessories from bumps and impacts.

Unpacking

• Unpacking

- Please unpack the device carefully
- In the case of any damage a detailed report must be set immediately (post, rail or forwarder)

• Delivery scope

Telescopic stand

- Telescopic stand corresponding to the type ordered T 653 oder R 474
- Operating instructions

For technical reasons associated with shipping, the telescopic stand is supplied stripped down, as pre-assembled modules..

Correct use

Use


The T 653 telescopic stand is used solely for mounting **T 65 D ULTRA-TURRAX®** dispersion units.

The R 474 telescopic stand is used solely for mounting **IKA® RW 47 D** agitator units.

The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly contrary to the manufacturer's specifications or if the appliance is modified by third parties.

Set-up and commissioning

Setting up the stand

The stand pillar is inserted into the baseplate so that the guide groove is at the back. The black handwheel bolt is then on the right. The stand pillar is then secured from below with a hexagon nut, toothed lock washer and a large shaped washer.

All four rubber feet are height-adjustable, so the stand can be set up so it does not rock, even on uneven floors.

In addition, a hole is provided in the baseplate of the R 474 telescopic stand, for screwing on a 16 mm diameter support rod. A clamp can be fixed on this rod.

Mounting the unit

The **IKA® RW 47 D** agitator unit or the **T 65 D ULTRA-TURRAX®** dispersion unit are placed on the arm of the stand and bolted on from below with the bolts supplied. The help of an assistant is needed for screwing on the T 65 D ULTRA-TURRAX® dispersion unit.

The cable clips included in the delivery package are required for fixing the power cable of the unit.

Safety disconnection device

The height-adjustable telescopic stand must be fitted with a safety disconnection device. We recommend installation of the SI 400 safety limit switch together with the SI 474 clamp support. You will find information on their installation in the operating instructions of the units.

Commissioning

The stand can be raised and lowered by loosening the handwheel bolt. The integral gas shock absorber helps to keep the manual force needed for this as light as possible. The manual force must be applied close to the stand pillar, so that the stand can be raised and lowered with only a slight use of force. If the force for raising or lowering is applied at a distance from the stand pillar, a turning moment results, leading to increased friction on the pillar guide. As a result, the stand can only be moved with great difficulty, or cannot be moved at all. Besides this, the danger of the stand toppling over is increased.

When the unit is at the correct height, the handwheel bolt is tightened by hand. Check the tightness of the handwheel bolt regularly when using the unit, and retighten it if necessary.

Maintenance and Cleaning

The stand pillar must always be lightly oiled, so that raising and lowering is not made more difficult due to increased friction.

Apart from this, the telescopic stand is maintenance-free.

Spare parts order

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item and designation of the spare part, see www.ika.com, spare parts diagram and spare parts list

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Only clean **IKA®** appliances using these **IKA®** approved cleaning agents:

Dirt	Cleaning agent
Dyes	Isopropanol
Building materials	Water containing detergent, Isopropanol
Cosmetics	Water containing detergent, Isopropanol
Food	Water containing detergent,
Fuels	Water containing detergent,
Other materials	Please consult IKA®

- Wear protective gloves when cleaning the devices.
- Please consult **IKA®** before using any cleaning or decontamination methods, other than those recommended here..

Accessories

• SI 400 Safety limit switch

• SI 474 Support holder

Warranty

In accordance with **IKA®** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual..

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DOMINIQUE DUTSCHER SAS
GmbH & Co.KG