IKA®

20000005394

M20_072015

Betriebsanleitung Ursprungssprache	DE	3
Operating instructions	EN	6
Mode d'emploi	FR	9
Instrukcja obsługi	PL	12

IKA®-Werke, Germany
Reg. No. 004343

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CE-DECLARATION OF CONFORMITY

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 60 010-1; EN 50 081; EN 50 082; EN 55 014 and EN 60 555.

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In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.

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Warranty

You have purchased an original **IKA®** laboratory machine which meets the highest engineering and quality standards.

In accordance with **IKA®** guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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Contents

	Page
Safety instructions	6
Correct use	6
Unpacking	6
Commissioning	6
Motor Circuit Breaker	7
Description	7
Grinding	7
Grinding while cooling	7
Replacing the grinding tool	8
Maintenance and cleaning	8
Accessories	8
Technical data	8

grinders operate without creating dust or lost. Beaters, blades and grinding are replaceable. The useful volume of the grinding chamber is 250ml for the M20 grinder.

Hard and brittle ground material is crushed by means of beaters or blades rotating at high speed.

Oil substances or those that are given to smearing and which are not capable of being ground at room temperature have to be cooled for grinding (e.g. embrittled with liquid nitrogen).

Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The M20 delivery quantity includes: one M20 grinder, one M21 stainless steel beater, one set of tools and operating instructions.

Safety instructions

The M20 device may not be used in areas subject to explosion hazards.

Protective glasses must be worn when performing experimental cuts. In addition, suitable protective gloves must be worn when adding coolants.

Experiments involving comminution must be performed behind a protective wall (Plexiglas®-wall, exhaust hood or the like) in order to prevent injury caused by the sudden occurrence of excess pressure (possible bursting of the cover).

Please comply with all of the relevant safety regulations.

Correct use

IKA universal grinder M20 is suitable for carrying out low-loss grinding of hard and brittle substances as of grain size 6 - 7mm. The

Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged.

Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.

NOTE: The grinder is designed for short term operation.

Maximum operating time - 7 minutes.

After that, the grinder must cool off for approximately 10 minutes. This condition does not apply if the grinder is connected to a cooling unit.

Motor Circuit Breaker

The grinder is protected against overload by means of an overcurrent protective device. In the event of overload the overcurrent protective device interrupts the current supply. Restarting is carried out by pressing in the black button at the rear of the grinder. The motor is switched on via push-button at the front.

Description

The device consists of three basic units (drive unit, middle part and cover) which are connected to each other by means of quick-release locks. When assembling the three units make certain that the hose coupling for cooling the middle part is located at the rear of the device. When putting on the cover make certain that the safety pin in the middle part is pushed down by the pressure pin on the cover. Otherwise, a tight and reliable connection is not ensured and the grinder will not start.

The device is started by pressing the push-button at the front.

Grinding

Fineness is mainly dependent on how long the grinding operation is performed.

If materials are to be ground that are smaller than a certain grain size then it should be sieved and the residue in the sieve is put back into the grinder. For sharp-edged material which is thus very abrasive, it is recommended that the entire quantity is left in the grinder in order to additionally exploit the properties of the material.

Grinding while cooling

Plastic, soft or smeary materials have to be cooled if they are to be ground.

There are two hose couplings located at the rear of the grinder in order to connect the cooling brine. Cooling with tap water is usually sufficient for substances which become hot as a result of the grinding operation and are thus given to smearing or those that produce oil.

The grinding chamber may only be filled halfway especially when drugs containing oils and similar materials are ground.

If more refrigeration is required then the ground matter should be cooled prior to being put into the grinding chamber.

The following procedure shall apply if liquid nitrogen is used in order to embrittle the material to be ground:

Safety information:

Suitable protective gloves must be worn when filling and transferring liquid nitrogen!

Protective glasses must be worn!

- The sample to be ground is put into a suitable vessel (e.g. evaporating dish) and covered with enough liquid nitrogen until it is completely covered.

When using liquid nitrogen make certain that the sample is not capable of reacting with oxygen.

- Cool grinding chamber using liquid nitrogen.
- Put the deep-cooled sample into the grinder and grind approximately 15 seconds.
- Then open the grinder and add liquid nitrogen for only as long as it takes for the sample to become deep-cooled again (embrittled).

Safety information:

Make certain that the liquid nitrogen has evaporated before beginning the grinding operation.

NOTE: Danger of overpressure!

- Repeat the grinding operation for approximately 15 seconds.
- Repeat the sequence of cooling and grinding until desired fineness is reached.

ess has been achieved.

- If the required fineness has not been achieved after 5 minutes, then a further reduction in grain size is unlikely.

The respective grinding chambers are sealed by means of labyrinth sealings. Though they are dust-proof, they are not leak-proof. This should be taken into account for substances where liquids may escape.

Replacing the grinding tools

The grinding tool is fixed by means of the retaining tube and the cap nut is unscrewed by using the SW10 spanner. The grinding tool, which rests on a double-edge track, may then be pulled out.

Maintenance and cleaning

The IKA M20 is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices.

Electrical equipment may not be put into the cleansing agent for cleaning.

The operator shall consult the manufacturer before using a method of cleaning or decontamination other than the one which is recommended by the manufacturer in order to make certain that the method which the operator intends to employ will not damage the equipment.

Accessories

M 20.1	spar grinding chamber with cover
M 21	spare beater, stainless steel
M 22	beater, stainless steel
M 23	star-type blade, stainless steel

Technical data

Design voltage:	VAC	230±10%
or	VAC	115±10%
Design frequency:	Hz	50/60
Permissible ambient temperature:	°C	+5 bis +40
Permissible relative humidity:	%	80
Degree of protection		IP 21
as in DIN 60 529:		20 000
Rated speed:	1/min	universal motor
Motor:		
Maximum motor input:	W	550
Grinding chamber useful volume:	ml	250
Permissible operating time (without cooling):	min	7
Interference suppression:		as in VDE 0871
Dimensions:	mm	170 x 170 x 350
Weight:	kg	6,6

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DOMINIQUE DUTSCHER SAS
IKA®-Werke GmbH & Co.KG