



illustra **MicroSpin** G-50 Columns

Product Booklet

Table of Contents

| | | |
|---|-------------------------------------|----|
| 1 | Introduction | 3 |
| 2 | Components | 4 |
| 3 | Description | 5 |
| 4 | Protocol | 9 |
| 5 | Appendices | 12 |
| 6 | Related products | 16 |
| 7 | Quick Reference Protocol Card | 21 |

1 Introduction

Product codes

27533001 (50 purifications)

27533002 (250 purifications)

About

For rapid buffer exchange or desalting, dye terminator or primer removal and removal of labeled nucleotides from labeling reactions.

Important

Read these instructions carefully before using the products.

Intended use

The products are intended for research use only, and shall not be used in any clinical or *in vitro* procedures for diagnostic purposes.

It is the responsibility of the user to verify the use of the illustra™ MicroSpin™ G-50 Columns for a specific application, as the performance characteristics of this product have not been verified for any specific organism.

Safety

All chemicals should be considered as potentially hazardous. For use and handling of the products in a safe way, refer to the Safety Data Sheets.

Storage



All kit components should be stored at room temperature (20°C-25°C). Do not freeze.

Expiry

For expiry date please refer to outer packaging label.

2 Components

Kit contents

| Identification | Pack size | 50 purifications | 250 purifications |
|---|------------------------------------|------------------|-------------------|
| Product code | 27533001 | 27533002 | |
|  | illustra MicroSpin G-50 columns | 50 | 250 |
|  | Collection tubes | 50 | 250 |

Refer to the Certificate of Analysis for a complete list of kit components.

Materials to be supplied by user

Disposables:

1.5 mL DNase-free microcentrifuge tubes.

Equipment needed

Microcentrifuge that accomodates 1.5 mL microcentrifuge tubes
Vortex mixer.


3 Description

Background


illustra MicroSpin G-50 Columns contain Sephadex™ G-50 DNA grade resin. They allow DNA purification by the process of gel filtration. Molecules larger than the largest pores in the Sephadex are excluded from the gel and elute first. Intermediate size molecules penetrate the matrix to varying extents, depending on their size and the resin used. Penetration of the matrix retards progress through the column; very small molecules elute last. The volume required to elute these small molecules is dependent on the volume available both inside and outside the pores i.e. the bed volume. Gel filtration resins do not exhibit a fixed exclusion limit when used in a spin-column format. Exclusion limits of gel filtration resins are only meaningful in continuous flow processes where the molecules being purified have sufficient time to reach equilibrium between the time spent in the gel filtration medium and the time spent in the eluent stream. In spin-column chromatography, the observed exclusion properties that allow the product to pass through the gel while the smaller impurities are retained depends on experimental factors, such as: the resin used, sample volume, product size, and the g forces used in the purification process.

The basic principle

Illustration



Step procedure

| Step | Comments | Component |
|--------------------|--|---|
| Column Preparation | The resin is re-suspended and excess storage buffer removed by centrifugation. | illustra MicroSpin G-50 column  |
| Sample Application | The sample is applied to the column. | |
| Elution | Purified sample is eluted by centrifugation. | |

Product specifications

| Sample Type: | Automated sequencing reactions |
|---------------------------------|---|
| Principle | Gel filtration |
| Column matrix | Sephadex G-50 DNA grade F |
| Column storage buffer | TE buffer (10 mM Tris/HCl, 1 mM EDTA) containing 0.05% Kathon™ CG/ICP Biocide as preservative. |
| Input sample volume | 12 to 50 µL |
| Percent sample recovery | Variable—depends on input sample |
| Maximum column loading capacity | 10 µg |
| Length of labeled DNA recovered | > 20 bases (N.B. there is no maximum length of probe that can be purified) |
| Nuclease Testing | Column components are tested in nickase, single and double-stranded exonuclease and RNase assays. |

| Sample Type: | Automated sequencing reactions |
|-------------------------------|---|
| Major subsequent applications | Dependent on input sample, but includes blotting and sequencing applications. |

When to use an illustra MicroSpin G-50 column

The illustra MicroSpin G-50 column is designed for the rapid purification of DNA for use in a wide range of applications, including desalting, buffer exchange, removal of dye terminators from cycle sequencing reactions and removal of labeled nucleotides from DNA labeling reactions. Good product yield and purity is obtained with sample volumes from 12 to 50 µL. It is suitable for any DNA greater than 20 bases in length and will not remove or denature enzyme. For guidelines to consider for use of illustra MicroSpin G-50 Columns, please see [Guidelines for use of an illustra MicroSpin G-50 column, on page 14](#).

Cytiva provides a wide range of nucleic acid purification products, some of which may be better suited to your application. These products and the application for which they have been optimized are summarized in the table shown below. illustra AutoSeq G-50 Columns and illustra ProbeQuant™ G-50 Micro Columns are provided pre-equilibrated in the optimal buffer for the application for which they are designed.

Table 1. Optimized Cytiva products for specific applications

| Application | Product | Product code | Pack size |
|---|---|--------------|-------------------|
| PCR reaction and enzymatic DNA reaction purification | illustra GFX™ PCR DNA and Gel Band Purification Kit | 28903470 | 100 purifications |

| Application | Product | Product code | Pack size |
|---|---|--|--|
| 50 bp–10 kbp size range Extraction of DNA from agarose gels | | 28903471 | 250 purifications |
| Dye terminator removal from automated sequencing reactions | illustra AutoSeq G-50 | 27534001 27534002 27534003 | 50 purifications 250 purifications 1000 purifications |
| Unincorporated labeled nucleotide removal from a DNA labeling reaction (> 20 mers) | illustra ProbeQuant G-50 Micro Columns | 28903408 | 50 purifications |
| Purification of oligonucleotides following synthesis, buffer exchange and de- salting. Gravity format, 500 mL loading volume | illustra NAP™-5 Columns | 17085301 | 20 purifications |
| Spin column format, 150 µL loading volume | illustra MicroSpin G-25 Columns | 27532501 | 50 purifications |

4 Protocol

Note: Columns are NOT transferable between Cytiva kits, e.g., the composition of the MicroSpin G-50 Columns is not the same as the composition of the ProbeQuant G-50 Micro Columns.

For Materials & Equipment to be supplied by user see *Materials to be supplied by user, on page 4* and *Equipment needed, on page 5*.

Protocol for purification of a range of sample types

Column Preparation

| Step | Action |
|------|--------|
|------|--------|

- | | |
|---|---|
| 1 | Re-suspend the resin in the column by vortexing. |
| 2 | Loosen the cap one-quarter turn and twist off the bottom closure. |
| 3 | Place the column in the supplied Collection tube for support. |
| 4 | For removal of labeled nucleotides from DNA labeling reactions, spin for 1 minute at $735 \times g$. For removal of dye terminators following cycle sequencing reactions, spin for 1 minute at $2\,000 \times g$. |

Note:

See *RPM calculation from RCF, on page 12* for RPM calculation from RCF.

Step Action



NOTICE

Use columns immediately after preparation to avoid drying out of the resin. If the column resin appears dry, displaced or cracked after the first spin, this is usually indicative of over-centrifugation (too fast or too long). Re-hydrate the column with 250 μ L of TE buffer, vortex and re-centrifuge, checking the settings. Spin speed can be reduced by 20% if necessary. Do not use the pulse button on the microcentrifuge as this may over-ride the speed setting.

- 5 Proceed immediately to the next part of the protocol.
-

Sample Application

Step Action

- 1 Place the column into a fresh DNase-free 1.5 mL microcentrifuge tube (user supplied)
- 2 Slowly apply 12–50 μ L sample to the topcenter of the resin, being careful not to disturb the resin bed.

| Step | Action |
|------|--------|
|------|--------|

Note:

The resin will have come away from the column slightly to form a pillar. It is essential that the sample being purified is applied slowly and is not allowed to run down the sides of the resin bed. Avoid touching the resin bed with the pipette tip.

- | | |
|---|---|
| 3 | Proceed to the next part of the protocol. |
|---|---|

Elution

| Step | Action |
|------|--------|
|------|--------|

- | | |
|---|---|
| 1 | For removal of labeled nucleotides from DNA labeling reactions, spin for 2 minutes at $735 \times g$. For removal of dye terminators following cycle sequencing reactions spin for 2 minutes at $2\,000 \times g$. The purified sample is collected in the bottom of the 1.5 mL microcentrifuge tube. |
| 2 | Cap the microcentrifuge tube. |
| 3 | Store the purified probe at -20°C . |

5 Appendices

RPM calculation from RCF

The appropriate centrifugation speed for a specific rotor can be calculated from the following formula:

$$\text{RPM} = 1000 \times \sqrt{(\text{RCF}/1.12r)}$$

Where RCF = relative centrifugal force; r = radius in mm measured from the center of the spindle to the bottom of the rotor bucket; and RPM = revolutions per minute.

E.g. if an RCF of $735 \times g$ is required using a rotor with a radius of 73 mm, the corresponding RPM would be 3000.

Table below shows appropriate RPM for various microcentrifuges.

Table 2. Appropriate RPM for an RCF of $735 \times g$ and $2000 \times g$

| Microcentrifuge | Appropriate RPM for an RCF of $735 \times g$ | Appropriate RPM for an RCF of $2000 \times g$ |
|----------------------------|--|---|
| Heraeus Biofuge 15 | 2800 | 4600 |
| Beckman GS15R | 2100 | 3600 |
| Hettich Mikro 24–48 | 2630 | 4300 |
| Hettich Mikro EBA12 | 2700 | 4400 |
| Eppendorf Centrifuge 5415C | 3000 | 4900 |
| Eppendorf Centrifuge 5417C | 2700 | 4400 |

Guidelines for use of an illustra MicroSpin G-50 column

illustra MicroSpin G-50 columns can be used for a wide variety of DNA purification applications. The DNA to be purified must be at least 20 bases in length. When using these columns, consider the following guidelines:

20× rule The best results will be obtained when the product being purified is at least 20 times larger than the largest impurity. If the difference in size is less than 20-fold, either purity or yield may be compromised.

Purity versus yield In general, purity is inversely proportional to yield. Larger sample volumes will provide higher yield but lower purity, and vice-versa.

Non-specific binding The non-specific binding exhibited by the illustra MicroSpin G-50 column is relatively insignificant, allowing purification of samples in the nanogram range. There will be a uniform proportional loss of sample which is due to the nature of spin column chromatography.

Retention For a given sample volume, product retention is relative to molecular size. As the size of the product increases, its relative retention decreases.

Loading volumes Load 12–25 µL onto a column for dye terminator removal. We recommend use of illustra AutoSeq G-50 columns for this application as they have been optimized for sequencing reaction clean-up, and for salt-sensitive analyzers that utilize capillary loading.

Load 50 µL for removal of unincorporated labeled nucleotides from DNA labeling reactions. We recommend use of illustra ProbeQuant G-50 Micro-Columns for this application, especially when handling small (ng range) quantities of DNA. Recovery of DNA from illustra MicroSpin G-50 columns is at least 10% less than that from illustra ProbeQuant G-50 micro-columns.

For larger sample volumes, either use more than one column or reduce the sample volume by drying or precipitation. For smaller sample volumes, dilute the sample to improve product recovery.

Enzyme Removal For purification of DNA fragments 50 bp–10 kbp in length, following an enzymatic reaction, we recommend using the illustra GFX PCR DNA and Gel Band Purification Kit, as the enzyme will be removed during the spin column process. If using an illustra MicroSpin G-50 column, you must Phenol Chloroform extract prior to loading onto the column to ensure enzyme removal.

Troubleshooting guide

This guide may be helpful in the first instance. However, if problems persist or for further information, please contact Cytiva technical services. Telephone numbers are on the back page.

Alternatively log onto [cytiva.com/illustra](https://www.cytiva.com/illustra)

Table 3. Problem: Poor sample purity

| Possible cause | Suggestions |
|--------------------|---|
| Poor sample purity | <ul style="list-style-type: none">• Ensure the sample volume was within acceptable range prior to loading (see Guidelines for use of an illustra MicroSpin G-50 column, on page 14).• Ensure sample is CAREFULLY pipetted into center of resin. Do not disturb the column. Do not allow the sample to run into the sides of the resin bed.• Use the column immediately after completing Column Preparation step. Do not allow the resin to become dried out or cracked. |

6 Related products

A full range of Molecular Biology reagents can be found in the Cytiva catalog and on the web site

cytiva.com/illustra

A full range of Detection Products and available pack sizes can be found in the Cytiva catalog and on the web site cytiva.com/newhyperfilm

| Application | Product | Product code | Pack size |
|-----------------|------------------------|--------------|-----------|
| Blotting | Hybond™-N+(82 mm) | RPN82B | 50 discs |
| | Hybond-N+ (15 × 20 cm) | RPN1520B | 10 sheets |
| | Hybond-NX (82 mm) | RPN82T | 50 discs |

| Application | Product | Product code | Pack size |
|-----------------------------|--------------------------------------|--------------|--------------|
| | Hybond-NX (15 × 20 cm) | RPN1520T | 10 sheets |
| | Hybond-N (82 mm) | RPN82N | 50 discs |
| | Hybond-N (15 × 20 cm) | RPN1520N | 10 sheets |
| | Hybond-XL (82 mm) | RPN82S | 50 discs |
| | Hybond-XL (15 × 20 cm) | RPN1520S | 10 sheets |
| | Hybond blotting paper (20 × 20 cm) | RPN1520S | 100 sheets |
| Radioactive labeling | Rediprime II DNA Labeling System | RPN1633 | 30 reactions |
| | Rediprime DNA Labeling Beads (-dCTP) | 27924001 | 1 kit |
| | Megaprime DNA Labeling System, dNTP | RPN1604 | 30 reactions |
| | Megaprime DNA Labeling System, dCTP | RPN1606 | 30 reactions |
| | Nick Translation Kit, dNTP | N5500 | 20 reactions |
| | Nick Translation Kit, dCTP | N5000 | 20 reactions |
| | 5'-End Labeling Kit | RPN1509 | 20 reactions |

| Application | Product | Product code | Pack size |
|--|--|--------------|---------------------|
| | Rapid-Hyb Buffer | RPN1635 | 125 mL |
| Detection | Hyberfilm™ MP (18 × 24 cm) | 28906843 | 50 sheets |
| | Hyberfilm MP Enveloped (18 × 24 cm) | 28906850 | 50 sheets |
| | Hypercassette™ | RPN11642 | 1 |
| Purification of DNA probes and oligonucleotides | illustra MicroSpin G-25 Columns | 27532501 | 50 purifications |
| | illustraProbeQuant G-50 Micro Columns | 28903408 | 50 purifications |
| | illustra NICK Columns | 17085502 | 50 purifications |
| | illustra NAP-5 Columns | 17085302 | 50 purifications |
| Purification of DNA from PCR, agarose gel bands and enzymes | illustra GFX PCR DNA & Gel Band Purification Kit | 28903470 | 100 purifications |
| | illustra GFX 96 PCR Purification Kit | 28903445 | 10 × 96 well plates |
| Preparation of PCR products for automated sequencing | ExoSAP-IT™ | US78200 | 100 reactions |

| Application | Product | Product code | Pack size |
|---|---|--------------|--|
| Dye terminator removal from automated sequencing reactions | illustra AutoSeq G-50 Columns | 27534001 | 50 purifications |
| Kits containing ready-to-use mix for PCR amplification | illustra Hot Start Master Mix | 25150001 | 100 reactions |
| | illustra PuReTaq Ready-To-Go™ PCR Beads | 27955701 | 96 reactions in 0.2 mL tubes/plate |
| | illustra PuReTaq Ready-To-Go PCR Beads | 27955702 | 5 × 96 reactions in 0.2 mL tubes/plate |
| | FidelityTaq™ PCR Master Mix Plus (2 ×) | E71182 | 100 reactions |
| | FidelityTaq Master Mix Plus | E71183 | 100 reactions |
| Premixed nucleotides for PCR amplification | illustra DNA Polymerization Mix dNTP Set (A,C,G,T) 20 mM each | 28406557 | 10 µmol |
| | illustra DNA Polymerization Mix dNTP Set (A,C,G,T) 20 mM each | 28406558 | 40 µmol (4 × 10 µmol) |

| Application | Product | Product code | Pack size |
|-------------|---|--------------|-------------|
| | illustra PCR Nucleotide Mix dNTP Set (A,C,G,T) 25 mM each | 28406560 | 500 μ L |
| | illustra PCR Nucleotide Mix dNTP Set (A,C,G,T) 2 mM each | 28406562 | 1 mL |

7 Quick Reference Protocol Card

Cue card

Quick Reference Protocol Card

illustra MicroSpin G-50 Columns

27533001 (50 purifications)

27533002 (250 purifications)

A. Protocol for purification of a range of sample types

⊕ :Add ⌚ :Spin

1. Column preparation

- Re-suspend the resin in the column by vortexing
- Loosen the cap one-quarter turn and twist off the bottom closure
- Place the column in the supplied Collection tube
- ⌚ 1 minute 735 × g for removal of labeled nucleotides from DNA labeling reactions OR
- ⌚ 1 minute 2 000 × g for removal of dye terminators from cycle sequencing reactions



2. Sample application

- Place the column into a fresh DNase-free 1.5 ml microcentrifuge tube (user supplied)
- ⊕ 12–50 µl of sample to the top-center of the resin with care



3. Elution

- ⌚ 2 minutes 735 × g for removal of labeled nucleotides from DNA labeling reactions OR
- ⌚ Spin 2 minutes 2 000 × g for removal of dye terminators from cycle sequencing reactions
- Retain eluate
- Store the purified sample at -20°C





cytiva.com

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

illustra, GFX, Hybond, Sephadex, Hyperfilm, MicroSpin, ProbeQuant, NAP and Ready-To-Go are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

FideliTaq and ExoSAP-IT are trademarks of Affymetrix, Inc. Kathon is a trademark of Dow Chemical Co.

All other third-party trademarks are the property of their respective owners.

© 2020–2021 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

27533001PL AB V:6 03/2021