

## **Vivaspin protein concentrator spin columns**



Vivaspin™ ultrafiltration spin columns are designed for fast, non-denaturing concentration and desalting of proteins and other biological samples.

Vivaspin™ ultrafiltration spin columns are designed for fast, non-denaturing concentration and desalting of

proteins and other biological samples.

- **Ease-of-use:** One-step sample concentration in a single tube, for minimal sample handling and reduced sample loss, as well as pH compatibility over a broad range – pH 1 to 9.
- **Sample protection:** Patented Vivaspin dead-stop technology ensures samples cannot be concentrated to dryness and enables direct concentrate recovery.
- **Ease-of-use:** One-step sample concentration in a single tube, for minimal sample handling and reduced sample loss.
- **High flow rates:** And minimal membrane blockage through the vertical polyethersulfone (PES) membrane design.
- **Minimal protein adsorption:** Thanks to the Vivaspin™ polycarbonate tube design.
- **Wide selection:** Vivaspin sample concentrators are suitable for sample volumes ranging from 100 µL to 20 mL, with a range of molecular weight cut-off (MWCO) values from Mr 3000 to 100 000.
- Download the Vivaspin protein concentrator data file to learn more and to determine the right MWCO (molecular weight cut-off) and volume for your enrichment study.

## Spin column protein concentration – how does it work?

Membrane spin columns such as Vivaspin combine the principle of size exclusion with centrifugation to achieve rapid concentration and high recovery of target biomolecules in a single tube. An upper compartment containing sample is separated from the lower compartment by a semipermeable membrane with an appropriate MWCO value chosen to exclude molecules of interest. Centrifugation is applied to force solvent through the membrane, leaving a more concentrated sample in the upper chamber. The Vivaspin column design permits protocols for one-step biomolecule concentration, rapid sample recovery (reverse spin protocol), desalting, and buffer exchange.

## Product Specifications

---

Parameter	Vivaspin 500 MWCO 10 000
Material (Membrane)	Polyethersulfone (PES)

<b>Parameter</b>	<b>Vivaspin 500 MWCO 10 000</b>
Material	Polycarbonate (PC)
Filtrate Vessel Material	Polycarbonate (PC)
Storage conditions	4 - 30°C
Dead-Stop Volume	5 µl
Hold-up Volume of Membrane	<5 µl
Wetted Part	No
Includes	Instructions supplied with the concentrators
Molecular Weight [MW] Cut-Off	10000
Volume	100-500 µl