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Contents

1 Safety	16
1.1 Safety Icons.....	16
1.2 Risks that the operator can be subjected to.....	16
1.3 Qualified users	16
2 General information	17
3 Technical specifications	17
4 Contents list and accessories	18
4.1 Standard contents	18
4.2 Accessories	18
5 Introduction	19
5.1 Definitions.....	19
5.1.1 Sterilization	19
5.1.2 Desinfection	19
5.1.3 Culture media.....	19
5.1.4 Purge.....	19
5.1.5 Saturated steam	19
5.1.6 Atmospheric purge.....	19
6 Equipment description	20
6.1 Equipment	20
6.2 Features	20
6.3 Control panel.....	21
6.4 Elements position	21
7 Installation	22
8 Operation	23
9 Error messages	26
10 Spare parts	26
11 Lid emergency opening	26
12 Recommendation to achieve a perfect sterilization	27

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1 Safety

The equipment comes complete with safety features.

This manual indicates areas of possible risk.

1.1 Safety Icons

Identifies risky situations and safety measures that should be adhered to. The icons make reference within the paragraph, marked as a lined grey.

Risk danger

Danger risk

Beware and follow the instructions as indicated.

Electrical risk

Risk of electrical shock if accessing zones are shown with this icon.

Beware and follow the instructions as indicated.

Risk of burns in high temperature areas

The temperature in the zone indicated with this icon can exceed 60°C. Use thermal protective gloves to complete described functions.

Beware and follow the instructions as indicated.

Important information

- Important information in obtaining the best results or optimum performance of the equipment.
- Important information to extend the life of the equipment and maintain its optimum performance.

1.2 Risks that the operator can be subjected to

- Possibility of touching surfaces that are over 60°C.
- Possible escape of steam vapours.
- Electrical risk.


1.3 Qualified users

This equipment should only be used by personnel that are qualified users.

This equipment should only be used by personnel that have understood this manual or have previous adequate instruction on the use of steam sterilizers.

2 General information

2.1 Handle the equipment with care. Unpack and check that all items coincide with the delivery note. If you see any discrepancy or damage of any kind, notify your distributor as soon as possible.

2.2 Read this manual before using the equipment.

2.3 These instructions are an important part of the equipment and should be made available to all users.

2.4 If in any doubt or clarity on the use of this instrument, contact your local distributor or the service department of J.P. Selecta, s.a.u.

2.5 **ATTENTION, EQUIPMENT NOT CLEAN AND DISINFECTED SHALL NOT BE REPAIRED.**

2.6 Any unauthorised modification, lack of maintenance or deliberate damage to the machine, transgresses the directive 89/655/CEE where the manufacturer will not be responsible for damages that could occur due to undue care and attention.

2.7 Do not use the equipment for fluids that can produce vapours that can form inflammable or explosive mixtures.

DO NOT OPEN THE LID OF THE AUTOCLAVE IF THE MACHINE IS NOT CONNECTED TO THE MAINS SUPPLY AND THE MAIN SWITCH IS ILLUMINATED. IF NOT CONNECTED, DO NOT FORCE THE DOOR HANDLE.

3 Technical specifications

Voltage supply 115-230V 50/60Hz according to the characteristics plate indications.


Part number:		4001758	4001759
Capacity:		50L	80L
Chamber dimensions:	diameter:	30cm	40cm
	height:	70cm	60cm
External dimensions:	height:	117cm	107cm
	width:	48cm	58cm
	depth:	62cm	72cm
Range:		115-134°C	
Power:		3000W	3600W
Weight:		75Kg	100Kg

4 Contents list and accessories

4.1 Standard contents

The standard equipment consists of the following items:


Model Presoclave III	4001758	4001759
Shelf	0005601	0005600
Instructions manual	80350	80350

5600-5601

4.2 Accessories

- Wired baskets made in AISI 304 stainless steel.

Model:	4001758	4001759	
Part number:	1000495	1000496	
Dimensions:	diameter:	28cm	38cm
	height:	20cm	28cm
Number of baskets:	3 units	2 units	

- Basket 4 compartments made in AISI 304 stainless steel for 1000495 code: 1001217
- Basket 6 compartments made in AISI 304 stainless steel for 1000496 code: 1001218

5 Introduction

The Presoclave III autoclaves are versatile and suitable for a wide range of applications, sanitation, industrial processes and quality control. These autoclaves can sterilise solids, with or without being wrapped, liquids, culture media.

5.1 Definitions:

5.1.1 Sterilization

Sterilization is the destruction or elimination of all forms of microorganisms, including spores present in inanimate objects.

5.1.2 Disinfection

Disinfection is the process of the destruction of microorganism agents. The temperature used is less than that of sterilization.

5.1.3 Culture media

A growth medium or culture medium is any substance in which microorganisms or cells can grow and multiply. Microorganisms can therefore be isolated and used in susceptibility testing. Generally it is produced as dry powder or granular form, but can also be produced as hydrates and pre-prepared. Should not be used unless they have been sterilized.

5.1.4 Purge

The purge cycle is when air is eliminated from inside the vessel with the objective of producing saturated steam vapour.

5.1.5 Saturated steam

Water vapour at a temperature corresponding to the ebullition point of the liquid origin.

5.1.6 Atmospheric purge

The term atmospheric purge is when air is gravimetrically expelled by vapour for a determined period.

6 Equipment description

The Presoclave III autoclave belongs to a series of saturated vapour autoclaves that are distinguished for the following features:

1. Arm / joist locking system.
2. Temperature and time electronic control circuit.
3. Safety locking system that only allows lid opening with pressure 0 and temperature under 95°C.
4. Maximum pressure limiter that disconnects the heating element in case of over-pressure.
5. Limiter thermostat to protect the heating element.

6.1 Equipment


- Safety system:

Safety over pressure valve.
Safety pressure switch.
Safety thermostat.
Lid closed detector.
Protective thermal lid cover.

- Chamber and lid made in stainless steel.
- External case made in stainless steel.
- Manual emptying valve.
- Waste water tank connection.

6.2 Features

- Sterilization temperature from 115°C to 134°C.
- Temperature resolution 1°C.
- Sterilization time from 3 minutes to 99 minutes.
- Time resolution 1 minute.
- Suitable to work with solids and liquids atmospheric purge from 115°C to 134°C.

6.3 Control panel

1. Pressure gauge
2. Open lid indicator lamp
3. Open valve indicator lamp
4. End cycle indicator lamp
5. Overtemp alarm indicator lamp
6. State indicator lamp
7. Temperature and error display
8. Time and zero pressure display
9. Temperature set key
10. Time set key
11. Increase value key
12. Decrease value key
13. Stop key
14. Starting up key
15. Depressurization button


6.4 Elements position

16. Emptying manual valve control
17. Open lid handle

18. Main switch

19. Safety valve

20. Power supply cord

21. Thermostat reset button

22. Condenser outlet and chamber emptying

7 Installation

7.1 Locate the autoclave close to an appropriate power source.

7.2 The autoclave needs to be placed on a flat stable horizontal surface, leaving a minimum space of 10 cm around the instrument.

7.3 Adjust the front feet to level the autoclave and set the autoclave in position.

7.4 Fix a plastic or rubber hose to the pump output connection (22) with a receptacle to collect any condensate produced by the atmospheric purge.

**¡WARNING! DON'T TRY TO OPEN THE LID
UNLESS READING
«8 OPERATION» SECTION**

8 Operation

8.1 OPENING THE LID

- The autoclave must be connected to mains supply and the main switch turned on.
- Pressure gauge: 0 bar.
- Temperature in the internal chamber must be under 95°C.
- Press the stop button (13) to unlock the locking system of the lid opening lever.


ONLY WITH THIS CONDITIONS LID IS READY TO BE OPEN

8.2 START

- Start on using the main switch (18)
- Open the lid as we can see on 8.1 section, using the handle (17)
- Place the shelf at the bottom of chamber.
- Be sure the valve (16) is closed.
- Fill of decalcified water until shelf level.
- Put the material to sterilize and close the lid.
- Start the cycle as shows the section 8.3

BEFORE STARTING THE CICLE CHECK ALWAYS THERE IS ENOUGH WATER IN THE CHAMBER. IF NOT, THE HEATING ELEMENT COULD BE DAMAGED.

8.3 TEMPERATURE CYCLE SELECTION

- To select the working temperature (from 115^a to 134^oC), press «°C» key (9) and, simultaneously, the increase value key (11) or decrease value key (12). Temperature is shown on display.

8.4 TIME CYCLE SELECTION

- To select the cycle time (from 3 min to 99 min), press the set time key (10) and, simultaneously, the increase value key (11) or decrease value key (12). Time is shown on display.

8.5 START WORKING CYCLE

- Once time finished previous operations on sections 8.2, 8.3 and 8.4, press the start on key (14) and autoclave will work with the following indications:

Lamp (6): Lighting

Display (7): Real temperature is shown.

Display (8): 6 blinking horizontal lines.


8.6 CYCLE WORKING STEPS

- 1 Purge: The autoclave expulses the chamber internal air to get saturated steam. The Presoclave III autoclaves makes a process called atmospheric purge. The temperature increases up to reach 103^oC. The purge electrovalve holds opened during 3 minutes once time 103^oC have been reached.

Lamp (6): Lighting

Display (7): Show temperature below 103^oC.

Display (8): 6 blinking horizontal lines.


- 2 Temperature increasing until set temperature: The temperature continues increasing from 103^oC to the set temperature, with all valves closed.

Lamp (6): Lighting

Display (7): Real temperature is shown between 103^oC and set temperature.

Display (8): 6 blinking horizontal lines.

- 3 Sterilization: 20 seconds after set temperature is reached, the autoclave starts the sterilization during all the time that we have set.

Lamp (6): Permanent lighted on.

Display (7): Real temperature is shown, the same set temperature.

Display (8): Remain minutes to finish the cycle.

DO NOT OPEN THE VALVE 16 UNTIL THE CYCLE IS COMPLETELY FINISHED TO AVOID POSSIBLE BURNS WITH HEAT WATER OR STEAM.

- End cycle: When sterilization is finished, cold process starts. The lamp (6) is blinking and the lamp (4) lights on. The display (7) shows the internal real temperature.
In case of solids sterilization, fast steam output can be made using button (15), the steam output is made with output (22).

¡WARNING!

- **THE DIRECT CONTACT WITH PRESSURE STEAM CAN PRODUCE SERIOUS BURNS.**
- **IF YOU USE A HOSE TO MAKE THE STEAM OUTPUT, DON NOT HOLD IT WITH HANDS TO AVOID BURNS.**

YOU CAN NOT OPEN THE LID UNTIL THE PRESSURE IS NOT 0 BAR AND THE INTERNAL TEMPERATURE UNDER 95°C. SEE SECTION 8.1

When the lid is ready to be opened, the display (7) shows the word «End» and the lamp (6) is off.

If liquids or culture media are sterilized must leave cooler until the pressure will be 0 bar. During the cooling, is the same case of manual depressurize, the display shows the internal temperature.

When the pressure is 0, the temperature is under 95°C, the display (7) shows the word «End». Therefore, you can open the lid with its handle (17).

9 Error messages

Error 1 Pressure sensor failure or overpressure.

Error 2 PT100 probe failure or overtemp.

Error 3 External alarm. Pressure switch or safety thermostat triggered.
To reset the safety thermostat, press its button number 21.

The thermostat is adjusted in the factory, it must not be operated except to reset it.

Error 4 Open door interlocking during the cycle.

Error 5 Opened door during a cycle.

Error 6 Pressure over +0.05bar with the autoclave stopped.

Error 7 E2PROM failure.

Error 8 Out of limits temperature during the sterilization.

Error 9 Mains supply failure during cycle.

10 Spare parts

07292	RS-232 communication wire
07315	26 pins communication wire
10044	R4 closing
13006	Contactora
15410	Fuse 16 A
16055	Solenoid electromagnet
16057	Electrovalve NA 3 mm 24 V
16212	Pressure switch
20096	Bipolar switch
29358	Frontal circuit
29459	Power circuit
36027	Solid relay 25A
39001	Heating element 2800W/230V
43031	Thermostat
43071	Pt100 probe

11 Lid emergency opening

In case of power supply failure, it is possible to open the lid manually.

Attention: Before opening the cover by using this system makes sure that the gauge pressure is 0 bar.

Remove the top drawer from the display by pulling hard (it is not screwed down, but snapped). A yellow cap will appear, gently pull it until you hear a click. The closure will be opened.

12 Recommendation to achieve a perfect sterilization

The material to be sterilized needs to be clean, free from any encrusted or residual material, we recommend to wash it with a good detergent and distilled water, then rinse it with distilled water.

Do not overload the baskets, drums or trays, always leave a space of 1 or 2 cm between items to allow a free passage of vapour and to ease drying.

12.1 Sterilization of instruments in bulk:

Place the instruments on a tray with waterproof paper leaving a space between each instrument.

Try to keep similar type material instruments together on a tray.

12.2 Tubing:

The tubing needs to be clean and drained. Place the tubing so that the openings on both ends are open. Make sure that the tubes are not touching the autoclave wall, blocked or twisted and that nothing can block the inside.

12.3 Vessels:

Never put hermetically sealed vessels in the autoclave.
Place the contents open end down to prevent any residual water being left inside.

12.4 Liquids:

Place the receptacle with the liquid inside on a drip tray, making sure that the receptacle can withstand the temperature.

Do not fill the receptacle more than 2/3rds of its maximum capacity, DO NOT CLOSE THE LID IF IT HAS A HERMETICAL SEALED, only use a lid if the pressure can be released.

Try to avoid using narrow neck vessels.

At the end of the cycle, let the autoclave cool down slowly (until the pressure is 0 kg/cm²).

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