**AUTOCLAVES PARA ESTERILIZACIÓN A VAPOR
AUTOCLAVES FOR STEAM STERILIZATION**

**PRESOCLAVE III 50L
PRESOCLAVE III 80L**

**4001758
4001759**

INDICE

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1 Safety

The equipment comes complete with safety features.

This manual indicates areas of possible risk.

1.1 Safety Icons

Identifies risky situations and safety measures that should be adhered to. The icons make reference within the paragraph, marked as a lined grey.

Risk danger

Danger risk

Beware and follow the instructions as indicated.

Electrical risk

Risk of electrical shock if accessing zones are shown with this icon.

Beware and follow the instructions as indicated.

Risk of burns in high temperature areas

The temperature in the zone indicated with this icon can exceed 60°C. Use thermal protective gloves to complete described functions.

Beware and follow the instructions as indicated.

Important information

- Important information in obtaining the best results or optimum performance of the equipment.
- Important information to extend the life of the equipment and maintain its optimum performance.

1.2 Risks that the operator can be subjected to

- Possibility of touching surfaces that are over 60°C.
- Possible escape of steam vapours.
- Electrical risk.

1.3 Qualified users

This equipment should only be used by personnel who are qualified users.

This equipment should only be used by personnel that have understood this manual or have previous adequate instruction on the use of steam sterilizers.

2 General information

2.1 Handle the equipment with care. Unpack and check that all items coincide with the delivery note. If you see any discrepancy or damage of any kind, notify your distributor as soon as possible.

2.2 Read this manual before using the equipment.

2.3 These instructions are an important part of the equipment and should be made available to all users.

2.4 If in any doubt or clarity on the use of this instrument, contact your local distributor or the service department of J.P. Selecta, s.a.u.

2.5 **ATTENTION, EQUIPMENT NOT CLEAN AND DISINFECTED SHALL NOT BE REPAIRED.**

2.6 Any unauthorised modification, lack of maintenance or deliberate damage to the machine, transgresses the directive 89/655/CEE where the manufacturer will not be responsible for damages that could occur due to undue care and attention.

2.7 Do not use the equipment for fluids that can produce vapours which can form inflammable or explosive mixtures.

DO NOT OPEN THE LID OF THE AUTOCLAVE IF THE MACHINE IS NOT CONNECTED TO THE MAINS SUPPLY AND THE MAIN SWITCH IS ILLUMINATED. IF NOT CONNECTED, DO NOT FORCE THE DOOR HANDLE.

3 Technical specifications

Voltage supply 115-230V 50/60Hz according to the characteristics plate indications.

Part number:		4001758	4001759
Capacity:		50L	80L
Chamber dimensions:	diameter:	30cm	40cm
	height:	63cm	60cm
External dimensions:	height:	117cm	119cm
	width:	48cm	58cm
	depth:	62cm	72cm
Range:		115-134°C	
Power:		3700W	3700W
Weight:		75Kg	100Kg

4 Contents list and accessories

4.1 Standard contents

The standard equipment consists of the following items:

Model Presoclave III	4001758	4001759
Cover resistance shelf	0009958	0009959
Instructions manual	80350	80350

4.2 Accessories

- Wired baskets made in AISI 304 stainless steel.


Model:	4001758	4001759
Part number:	1000495	1000496
Dimensions:		
diameter:	28cm	38cm
height:	20cm	28cm

Number of baskets:	3 units	2 units
--------------------	---------	---------

- BLIND baskets in AISI 304 stainless steel with folding handle.

Model:	4001758	4001759
Part number:	1000490	1000491
Dimensions:		
diameter:	27cm	36cm
height:	19cm	28cm

Number of baskets:	3 units	2 units
--------------------	---------	---------

- Basket 4 compartments made in AISI 304 stainless steel for 1000495 code: 1001217
- Basket 6 compartments made in AISI 304 stainless steel for 1000496 code: 1001218

5 Introduction

Presoclave III autoclaves are versatile and suitable for a wide range of applications, sanitation, industrial processes and quality control. These autoclaves can sterilise solids, with or without being wrapped, liquids, culture media.

5.1 Definitions:

5.1.1 Sterilization

Sterilization is the destruction or elimination of all forms of microorganisms, including spores present in inanimate objects.

5.1.2 Disinfection

Disinfection is the process of destruction of microorganism agents. The temperature used is less than that of sterilization.

5.1.3 Culture media

A growth medium or culture medium is any substance in which microorganisms or cells can grow and multiply. Microorganisms can therefore be isolated and used in susceptibility testing. Generally it is produced as dry powder or granular form, but can also be produced as hydrates and pre-prepared. Should not be used unless they have been sterilized.

5.1.4 Purge

The purge cycle is when air is eliminated from inside the vessel with the objective of producing saturated steam vapour.

5.1.5 Saturated steam

Water vapour at a temperature corresponding to the ebullition point of the liquid origin.

5.1.6 Atmospheric purge

The term atmospheric purge is when air is gravimetrically expelled by vapour for a determined period.

6 Equipment description

Presoclave III autoclave belongs to a series of saturated vapour autoclaves that are distinguished for the following features:

1. Arm / joist locking system.
2. Temperature and time electronic control circuit.
3. Safety locking system that only allows lid opening with pressure 0 bar and a temperature under 95°C.
4. Maximum pressure limiter that disconnects the heating element in case of over-pressure, set at 2.5 bar.
5. Limiter thermostat to protect the heating element, set at 160°C.


6.1 Equipment

- Safety system:
 - Safety over pressure valve.
 - Safety pressure switch.
 - Safety thermostat.
 - Lid closed detector.
 - Protective thermal lid cover.
- Chamber and lid made in stainless steel.
- External case made in stainless steel.
- Manual emptying valve.
- Waste water tank connection.

6.2 Features

- Sterilization temperature from 115°C to 134°C.
- Temperature resolution 1°C.
- Sterilization time from 3 minutes to 99 minutes.
- Time resolution 1 minute.
- Suitable to work with solids and liquids atmospheric purge from 115°C to 134°C.

6.3 Control panel

1. Pressure gauge
2. Open lid indicator lamp
3. End cycle indicator lamp
4. Overtemp alarm indicator lamp
5. State indicator lamp
6. Temperature and error display
7. Time and zero pressure display
8. Temperature set key
9. Time set key
10. Increase value key
11. Decrease value key
12. Stop key
13. Starting up key
14. Depressurization button


6.4 Elements position

15. Emptying manual valve control
16. Open lid handle

17. Main switch

- 18. Safety valve
- 19. Power supply cord
- 20. Thermostat reset button
- 21. Condenser outlet and chamber emptying

7 Installation


1. Locate the autoclave close to an appropriate power source.
2. The autoclave needs to be placed on a flat stable horizontal surface, leaving a minimum space of 10 cm around the instrument.
3. Adjust the front feet to level the autoclave and set the autoclave in position.
4. Fix a plastic or rubber hose to the pump output connection (21) with a receptical to collect any condensate produced by the atmospheric purge.

You must ensure that the purge condensate tube is not clogged or submerged in water after the cycle, since the chamber will empty and you will not be able to open the lid.

Note that during the purge phase, the equipment draws steam through this tube and this steam is in form of water inside the container. Make sure the tube is not in contact with the bottom of the container to ensure that at the end of the cycle it does not absorb water from it.

7.1 Unblocking the door

In case of open lid blocked by a safety system failure make as follows:

- 1 Stop the machine.
- 2 Disconnect the machine from the mains.
- 3 Place an outlet and a tube in the aeration valve (11).
- 4 Carefully open the valve to avoid steam leakages. Be careful to avoid possible burns.
- 5 Wait for the total steam emptying.
- 6 Extract the upper cover where display is located in order to put it aside (it is not necessary to disconnect any wire). It does not use screws, pull tightly.
- 7 Remove the plastic cap and pull slightly until the lever is released.

**¡WARNING! DON'T TRY TO OPEN THE LID
UNLESS READING
«8 OPERATION» SECTION**

8 Operation

8.1 Opening the lid

- The autoclave must be connected to mains supply and the main switch turned on.
- Pressure gauge: 0 bar.
- Temperature in the internal chamber must be under 95°C.
- Press the stop button (12) to unlock the locking system of the lid opening lever.

ONLY WITH THIS CONDITIONS, LID IS READY TO BE OPEN





8.2 Start

- Start on using the main switch (17)
- Open the lid as we can see on 8.1 section, using the handle (16)
- Place the shelf at the bottom of chamber.
- Be sure the valve (15) is closed.
- Fill of decalcified water until shelf level.
- Put the material to sterilize and close the lid.
- Start the cycle as shown on section 8.3
- When closing the lid, make sure that the open door indicator (2) is off, otherwise the device will not start.

**BEFORE STARTING THE CICLE,
ALWAYS CHECK IF THERE IS ENOUGH WA-
TER IN THE CHAMBER. IF NOT, THE HEATING
ELEMENT COULD BE DAMAGED.**

8.3 Temperature cycle selection

- To select the working temperature (from 115^a to 134°C), press «°C» key (8) and, simultaneously, the increase value key (10) or decrease value key (11). Temperature is shown on display.

8.4 Time cycle selection

- To select the cycle time (from 3 min to 99 min), press the set time key (9) and, simultaneously, the increase value key (10) or decrease value key (11). Time is shown on display.

8.5 Start working cycle

- Once previous operations finished indicated on sections 8.2, 8.3 and 8.4, press the start ON key (13) and autoclave will work with the following indications:

Lamp (5): Lighting

Display (6): Real temperature is shown.

Display (7): 6 blinking horizontal lines.

8.6 Cycle working steps

- 1 Purge: The autoclave expulses the chamber internal air to get saturated steam. **Presoclave III** autoclaves make a process called atmospheric purge. The temperature increases up to reach 103°C. The purgue electrovalve holds opened during 3 minutes once time 103°C have been reached.

Lamp (5): Lighting

Display (6): Show temperature below 103°C.

Display (7): 6 blinking horizontal lines.

- 2 Temperature increasing until set temperature: The temperature continues increasing from 103°C to the set temperature, with all valves closed.

Lamp (5): Lighting

Display (6): Real temperature is shown between 103°C and set temperature.

Display (7): 6 blinking horizontal lines.

- 3 Sterilization: 20 seconds after set temperature is reached, the autoclave starts the sterilization during all the time that we have set.

Lamp (5): Permanent lighted on.

Display (6): Real temperature is shown, the same set temperature.

Display (7): Remain minutes to finish the cicle.

DO NOT OPEN THE VALVE 16 UNTIL THE CYCLE IS COMPLETELY FINISHED TO AVOID POSSIBLE BURNS WITH HEAT WATER OR STEAM.

- End cycle: When sterilization is finished, cold process starts. The lamp (5) is blinking and the lamp (3) lights on. The display (6) shows the internal real temperature.
In case of solids sterilization, fast steam output can be made using button (14), the steam output is made with output (21).


¡WARNING!

- **THE DIRECT CONTACT WITH PRESSURE STEAM CAN PRODUCE SERIOUS BURNS.**
- **IF YOU USE A HOSE TO MAKE THE STEAM OUTPUT, DON NOT HOLD IT WITH HANDS TO AVOID BURNS.**

YOU CAN NOT OPEN THE LID UNTIL THE PRESSURE IS NOT 0 BAR AND THE INTERNAL TEMPERATURE UNDER 95°C. SEE SECTION 8.1

When the lid is ready to be opened, the display (6) shows the word «End» and the lamp (5) is off.

If liquids or culture media are sterilized, you must leave it cool until the pressure will be 0 bar. During the cooling, it is the same case of manual depressurize, the display shows the internal temperature.

When the pressure is 0, the temperature is under 95°C, the display (6) shows the word «End». Therefore, you can open the lid with its handle (16).

9 Maintenance

9.1 Temperature calibration

If a correction needs to be done, do as follows:


- Press the button (8) «°C» and (9) «clock» simultaneously for more than 5 seconds.

- While continuing to press these two keys, press keys (10) and (11) to modify this parameter. For fast forward press and hold the arrow (10 or 11) for a few seconds.

Note: Correction range from -9.9 to + 9.9°C. The modified value is added or restored to the measured temperature by correcting it.

9.2 Safety thermostat

8. Safety thermostat regulation screw (grey).
9. Safety thermostat reset button (white).

9.2.1 Rearming the safety thermostat

Safety thermostat is an element to prevent the temperature from exceeding a value. If by a fortuitous event the thermostat is triggered, that is, it opens the power circuit, the heating element will stop working. Usually there is an indicator light, or an error message, that indicates that the thermostat is in this position.

If the thermostat is triggered, check that the equipment looks good and it is not overheating at any point.

To change its status, simply press the white button (9), and you will hear a "click" that indicates that the thermostat has been reset. If the equipment has not cooled down sufficiently, the thermostat cannot be reset, and you must wait for the temperature to drop before it can be reset.

9.2.2 Adjusting the safety thermostat

In some equipment such as ovens or baths and with certain processes, it may be useful to adjust the thermostat trigger value, but for the rest of equipments this operation should NEVER be carried out.

From factory, the thermostat is adjusted to its maximum temperature.

9.3 Gasket

It is necessary to check that the gasket is in good condition. With the door open, remove the gasket and check that it is not deformed, has cracks or is broken. To check it, remove the gasket from the cover, dry it with a cloth and dry the groove of the cover where the gasket is located. Reassemble and do a test cycle. To ensure a good seal with a use of 1 cycle a day, it is recommended to replace the gasket every two years.

9.4 Cleaning the chamber

When working with steam, it is possible that traces of materials are deposited on the walls of the chamber. It is recommended from time to time to clean the inside of the chamber with a clean cloth. Also remove the heater cover and clean the bottom of the chamber.

10 Error messages

Error 1 Pressure sensor failure or overpressure.

Error 2 PT100 probe failure or overtemp.

Error 3 External alarm. Pressure switch or safety thermostat triggered.
To reset the safety thermostat, press its button number 20.

The thermostat is adjusted in the factory, it must not be operated except to reset it.

Error 4 Open door interlocking during the cycle.


Error 5 Open door during a cycle.

Error 6 Pressure over +0.05bar with the autoclave stopped.

Error 7 E2PROM failure.

Error 8 Out of limits temperature during sterilization.

Error 9 Mains supply failure during cycle.

11 Spare parts

07292	RS-232 communication wire
07315	26 pins communication wire
10044	R4 closing
13006	Contactora
15410	Fuse 16 A
16055	Solenoid electromagnet
16057	Electrovalve NA 3 mm 24 V
16212	Pressure switch
20096	Bipolar switch
29358	Frontal circuit
29459	Power circuit
36027	Solid relay 25A
39001	Heating element 2800W/230V
43031	Thermostat
43071	Pt100 probe

12 Lid emergency opening

In case of power supply failure, it is possible to open the lid manually.

Attention: Before opening the cover by using this system, make sure that the gauge pressure is 0 bar.

Remove the top drawer from the display by pulling hard (it is not screwed down, but snapped). A yellow cap will appear, gently pull it until you hear a click. The closure will be opened.

13 Recommendation to achieve a perfect sterilization

The material to be sterilized needs to be clean, free from any encrusted or residual material. We recommend to wash it with a good detergent and distilled water, then rinse it with distilled water.

Do not overload the baskets, drums or trays, always leave a space of 1 or 2 cm between items to allow a free passage of vapour and to ease drying.

13.1 Sterilization of instruments in bulk

Place the instruments on a tray with waterproof paper leaving a space between each instrument.

Try to keep similar material instruments type together on a tray.

13.2 Tubing

The tubing needs to be clean and drained. Place the tubing so that the openings on both ends are open. Make sure that the tubes are not touching the autoclave wall, blocked or twisted and that nothing can block the inside.

13.3 Vessels

Never put hermetically sealed vessels in the autoclave. Place the contents open end down to prevent any residual water being left inside.

13.4 Liquids

Place the receptacle with the liquid inside on a drip tray, making sure that the receptacle can withstand the temperature.

Do not fill the receptacle more than 2/3rds of its maximum capacity, DO NOT CLOSE THE LID IF IT HAS A HERMETICAL SEALED, only use a lid if the pressure can be released.

Try to avoid using narrow neck vessels.

At the end of the cycle, let the autoclave cool down slowly (until the pressure is 0 kg/cm²).

Notice to customers:

The product is made up of various components and various materials that must be recycled or, failing that, deposited in the corresponding debris removal sites when the product's life has been completed or when otherwise it is necessary to dispose of it. To do this, the end user who acquires the product must know the current regulations of each municipality and / or locality based on the waste electrical and electronic equipment. The user who acquires this product must be aware of and responsible for the potential effects of the components on the environment and human health as a result of the presence of hazardous substances. Never place the product in a conventional container of citizen scope if a previous dismantling and knowledge of the components that incorporates. If you do not know the procedure to follow, consult with the city council for more information.