Making Testing Easier

Eclipse Farm^{3G}

ZE/EF25

ZE/EF50

Single tube tests for detection of
inhibitory substances in milk

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www.zeulab.com

SCOPE

ECLIPSE Farm^{3G} is a qualitative test designed for the detection of antibiotics and inhibitors in milk. This test can be used in raw (with or without azidiol), skimmed, heated or powder milk from cow, sheep, goat or buffalo. The analysis of raw cow milk can be performed in individual or commingled milk. For the other species, the analysis should be performed only with commingled milk.

PRINCIPLE

ECLIPSE Farm^{3G} is based on the inhibition of microbial growth. The tubes provided contain a specific culture medium spread with *Geobacillus stearothermophilus* spores and a pH indicator. When the tubes are incubated in absence of antibiotics, spores germinate and cells grow (negative). If samples contain inhibitors at higher concentrations than the detection limit, microorganisms will not grow (positive). The system has been optimized by combining the kit with COMET8 equipment, which integrates the processes of incubation, time control and interpretation of results.

ADDITIONAL MATERIAL (NOT PROVIDED)

- COMET8 (ref. ZE/COM8) or FX incubator (ref. ZE/FX) or FS incubator (ref. ZE/FS).
- Negative control sample (without antibiotics).
- Positive control. Freeze dried Penicillin G (ref. ZE/PG5), Oxytetracycline (ref. ZE/OXITETRA5), Sulfathiazole (ref. ZE/SULFA5).

NOTES

- ✓ To achieve the best assay performance, it is recommended to analyse samples by using COMET8, which incubates, stops the assay automatically and interprets the results. If it is not possible, use an FX or FS incubator.
- ✓ Before starting, install and switch on COMET8 equipment following the manual provided.
- ✓ A new pipette should be used for each sample. To use the provided pipettes: squeeze the bulb, put the stem in the sample and release the bulb to draw up liquid. Excess liquid is siphoned automatically into overflow tanks. Squeeze the bulb again to release 100 µl.
- ✓ It is recommended to analyse a negative control sample of the same species and in the same conditions than the tested samples.
- ✓ Colours change over time, therefore results should be read at the time of the end of the assay.
- ✓ This test is extremely sensitive to antibiotics and other antibacterial substances such as detergents and disinfectants. Any contamination with those substances should be prevented.

- ✓ Milk contains natural inhibitors which do not interfere with the result. However, these substances appear in high concentrations in colostrum, at the end of lactation and in mastitis cases, and they may alter the test results.
- ✓ The range of colours obtained using COMET 8 is different to that seen with other incubation systems. Colours should not be compared with colour card shown on the last page.
- ✓ Comet8 was formerly called e-Reader8.
- ✓ If a negative control sample is not applied when using an incubator, the assay should be stopped at 2h 55min of incubation. Please, be aware that sensitivity of the test can be affected.

TEST PROCEDURE (Flowchart Procedure on the last page)

Please pay attention to "Notes" in page 2.

Sample application

1. Cut the tubes needed with scissors. Take care not to damage the foil on the remaining tubes to prevent the agar from drying out and stored immediately at 4-12°C.
2. Pierce or remove the foil covering each tube and add 100 µl of milk sample using the pipettes provided in the kit. (See "Notes" for pipettes use).
 - To detect antibiotics in **cow and buffalo milk**: Seal carefully the test tubes with the adhesive film provided.
 - To detect antibiotics in **goat or sheep milk, and samples containing azidiol***: Seal the test tubes with tape. Test tubes must be incubated at room temperature for 1 h. Afterwards, wash the wells 2-3 times with distilled water. To do this, empty the tubes, fill with distilled water and empty again, removing excess water by gently tapping the tube against absorbent paper. Seal with the adhesive film provided.

*Azidiol can be purchased from PanReac AppliChem (Product Code 176131).

Incubation and reading

There are two options: A) with COMET8 and B) with incubator and visual reading.

A) With COMET8

3. Place the sample tubes in COMET8. Follow COMET8 instructions to begin the assay. Choose the COMET8 method depending on your analysis:

Eclipse Single	Cow	Each test is started individually. No negative control is required.
Eclipse NC	Cow, buffalo	The samples are analysed at the same time with a negative control sample at position 1.
Eclipse NC SGA	Sheep, goat and samples with azidiot	The samples are analysed at the same time with a negative control sample of the same characteristics at position 1.

4. COMET8 will stop automatically the assay and results will be shown on the screen.

B) With incubator and visual reading:

3. Incubate the sealed test tubes in a pre-warmed incubator at $63 \pm 1^\circ\text{C}$.
4. Check the tubes at the time indicated on the product certificate. This time is given as a guide. The incubation should be stopped when the negative control has turned yellow as shown in the colour card (around $2\text{h}40 \pm 15\text{min}$). If the negative control sample is not yellow, incubate the tube again for further 10-15 min to get a yellow colour.

Results: (see colour chart)

- A yellow colour indicates the absence of antibiotics (negative).
- A blue-purple colour indicates the presence of antibiotics (positive).
- A yellow-green colour indicates presence of antibiotics in a concentration close to the detection limit (doubtful). Analysis should be repeated.

STABILITY AND STORAGE

The kit components should be stored at $4\text{-}12^\circ\text{C}$ and protected from light. The kit is provided with a minimum shelf life of 9 months. See the kit expiring date on the package.

SAFETY

Good laboratory practices are recommended when using this kit. A Material Safety Sheet (MSDS) is available from your local distributor or ZEULAB upon request.

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ECLIPSE is an *in vitro* diagnostic kit for antibiotics screening in milk. Suspected non-compliant results obtained with screening methods for official control purposes are liable to physicochemical confirmation (Commission Decision 2002/657/EC). ZEULAB does not assume any legal responsibility.

To the best of our knowledge, the information contained herein is accurate and complete. However, nothing herein shall be construed to imply any warranty or guarantee.

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