# IC BOVINO IC CAPRINO

ZE/ICB      ZE/ICB100  
ZE/ICC      ZE/ICC100

Immunochemical test for milk  
species identification

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**SCOPE**

IC BOVINO kit is a qualitative test to detect the presence of cow's milk in goat or sheep's milk. The test is an immunoassay based on the detection of bovine immunoglobulins (IgG) in milk and cheese samples.

IC CAPRINO kit is a qualitative test to detect the presence of goat's milk in sheep's milk. The test is an immunoassay based on the detection of caprine immunoglobulins (IgG) in milk and cheese samples.

For quantitative testing of up to 2% mix, IC BOVINO QUANTI (ZE/ICBQ, ZE/ICBQ100) and IC CAPRINO QUANTI (ZE/ICCQ, ZE/ICCQ100) are available. Ask your local distributor or ZEULAB for more information.

KIT COMPONENTS	ZE/ICB & ZE/ICC	ZE/ICB100 & ZE/ICC100
Strip tests	25	100
Dilution solution	1	2
Disposable pipettes	25	100
Test tubes	25	-
Divisible strip of 8 wells	-	13
Tube/well rack	yes	yes

Whole, skimmed, fresh and pasteurised milk; whey samples and cheese can be used with these kits. Nevertheless, be aware that heat treatments over 72°C/15-20 sec can begin to denaturize milk immunoglobulins. The higher the heat treatment, the higher the immunoglobulins denaturalization will be. Therefore, if possible, test samples prior to being processed.

**MILK ASSAY PROCEDURE 1+3**

1. Take a test tube or a well from the divisible strip. Place it into the rack.
2. Apply 1 drop of milk and 3 drops of dilution solution. Stir gently.
3. Place the strip test into the test tube/well following the arrows direction.
4. Read the results at 3, 5 or 10 minutes depending on the Limit of Detection desired, (see table 1+ 3).

Table 1+ 3  
LOD visual reading

Assay time	3 min	5 min	10 min
LoD - IC Bovino	1%	0,5%	0,25%
LoD - IC Caprino	-	1%	0,5%

**MILK ASSAY PROCEDURE 2 + 3**

1. Take a test tube or a well from the divisible strip. Place it into the rack.
2. Apply 2 drops of milk and 3 drops of dilution solution. Stir gently.
3. Place the strip test into the test tube/well following the arrows direction.
4. Read the results at 3, 5 or 10 minutes depending on the limit of detection desired (see table 2 + 3).

Table 2 + 3  
LOD visual reading

Assay time	3 min	5 min	10 min
LoD - IC Bovino	0,5%	0,25%	0,1%
LoD - IC Caprino	1%	0,5%	0,25%

**ASSAY PROCEDURE FOR CHEESE**

1. Weigh 5g of an homogeneous cheese sample and add 10 mL of distilled water. Mix with a pestle and mortar to obtain a fine homogenised sample.
2. Centrifuge at 3000 x g for 10 minutes. Alternatively, the sample can be filtered using a filter paper.
3. Take a test tube or a well from the divisible strip. Place it into the rack.
4. If a micropipette is available, it is recommended to mix 30 µL of the cheese extract with 130 µL of the dilution buffer. If not, apply 1 drop of cheese extract and 4 drops of dilution solution. Stir gently.
5. Place the strip test into the test tube/well following the arrows direction.
6. Read the result at 10 minutes.

**INTERPRETATION OF THE RESULTS****Visual reading**

**NEGATIVE:** A BLUE line is observed in the white central zone of the strip test. This line indicates that the test has been run properly.

**POSITIVE:** Two lines, BLUE and RED, are observed in the white central zone of the strip test.

**Objective reading - IRIS (Assay procedure 1+3, reading at 10 min)**

Place the strip in the specific tray for IC in the correct position (see figure 1 on the last page) and carry out the reading following the procedure described in the IRIS user guide. Results will be displayed on the screen and on the app as:

**NEG** (negative), **POS** (positive) or **POS++** (high positive).

For more information about levels of mixtures, please check the product certificate.

**Note:** There is a specific reading method for each batch. The method is transferred

to IRIS via the smartphone app using a QR code that is available on the product certificate, or using IRIS Lite. When different batches are used, pay special attention to the method uploaded on IRIS.

## NOTES

- Use a new disposable pipette for each sample. The pipettes included in the kit provide the correct volume of sample.
- Milk samples can be stored under refrigeration for 1 or 2 days before the assay. For longer period of time, samples should be frozen at -20°C.
- The kit is designed to detect traces of milk mixture. As example, 0.125% of mixture implies that a tank of 10 000 L is contaminated with 12.5 L of other milk specie.
- When pure cow milk is tested with IC Bovino the result is positive.
- When pure goat milk is tested with IC Caprino the result is positive.
- The intensity of the line might change over time. Do not read results after 12 minutes since the beginning of the assay.
- To keep records of the result use IRIS or take pictures of the strips when the assay is finished.
- If the blue line does not appear on the white central zone, the test has not been run properly or the reagents have been damaged. The result must be considered INVALID (INV) and repeated. If it happens again, please contact ZEULAB.
- During cheese making, the lost of soluble proteins is variable, depending on the variety of cheese and the production system of each cheese maker. For information about test sensitivity in cheese, please contact ZEULAB or your distributor.

## STABILITY

Store the test under dry conditions and temperatures between 15 and 25°C. The strips are very sensitive to humidity and must be kept inside the container provided with the desiccant bag. The kit has a shelf life of approximately 12 months. Check the expiration date on the package.

## SAFETY

This kit should be used following good laboratory practices. A Safety Data Sheet (SDS) is available on request.

**IC KIT** is an in vitro diagnostic test. In analysis implicating legal processes, the results should be re-evaluated with an official reference method. ZEULAB, S.L. does not assume any legal responsibility.

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