IC OVINO

ZE/ICO

Immunochromatographic test for detection of sheep's milk in cow's, goat's and buffalo's milk

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info@zeulab.com • www.zeulab.com

SCOPE

IC OVINE kit is a qualitative test to detect the presence of sheep's milk in cow's, goat's and buffalo's milk. The test is based on the detection of ovine immunoglobulins (IgG) in milk samples.

KIT COMPONENTS

	ZE/ICO
Strip tests	25
Dilution solution	1
Disposable pipettes	25
Test tubes	25
Tube rack	1

SAMPLE INFORMATION

Whole, skimmed, fresh or pasteurised milk samples can be analysed with this kit. Nevertheless, be aware that heat treatments over 72°C/15-20 seconds can begin to denaturise milk immunoglobulins. The higher the heat treatment, the higher the immunoglobulin denaturalization will be. Therefore, if possible, test samples prior to being processed.

Samples can be stored under refrigeration for 1 or 2 days before the assay. For longer period of time, samples should be frozen at -20°C.

No milk sample preparation is needed.

DETECTION LIMITS

Limit of Detection in milk: 1.0% of sheep's milk in milk of cow, goat and buffalo.


TEST PROCEDURE

1. Take a test tube and place it into the rack provided.
2. Apply 1 drop of milk and 3 drops of dilution solution. Stir gently.
3. Place the strip test into the test tube following the arrows direction.
4. Wait for 10 minutes and read the results.
5. Reading of results:

Visual reading

NEGATIVE: One red line (C) is observed in the white central zone of the strip test. This control line indicates that the test has been run properly.

POSITIVE: Two red lines (C and T) are observed in the white central zone of the strip test.

NOTES

- Use a new disposable pipette for each sample. The pipettes included in the kit provide the correct volume of sample.
- The test line (T) might not appear when testing samples of pure sheep milk or mixtures of more than 70% of sheep's milk due to "prozone or Hook effect".
- Results should be read within the next 2 minutes after the end of the analysis.
- The intensity of the line might change over time. Cut the absorbant pad after reading the result.
- If the control red line does not appear on the white central zone, the test has not been run properly or the reagents have been damaged. The result must be considered INVALID. The test must be repeated and contact ZEULAB if it happens again.

STABILITY

Store the test under dry conditions and temperatures between 15 and 25°C. The strips are very sensitive to humidity and must be kept inside the container provided with a desiccant bag. The kit has a shelf life of approximately 12 months. Check the expiration date on the package.

SAFETY

This kit should be used following good laboratory practices. A Material Safety Data Sheet is available on request.

DOMINIQUE DUTSCHER SAS

IC KIT is an *in vitro* diagnostic test. In analysis implicating legal processes, the results should be re-evaluated with an official reference method. ZEULAB, S.L. does not assume any legal responsibility.

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