PROTEON GLUTEN EXPRESS

ZE/PR/GL10

Immunochromatographic test for
detection of gluten in food and surfaces

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SCOPE AND TEST PRINCIPLE

The celiac disease is a disorder caused by an inappropriate immune response to prolamins from gluten found in some wheat, rye, barley and, to a lesser extent, in oats.

PROTEON GLUTEN EXPRESS is a qualitative and immunochromatographic rapid test for detection of gluten traces in food and working surfaces. The assay is based on the reaction of the prolamins with a specific antibody coated with red particles. The appearance of two lines, red and blue, on the strip test is a positive result.

Limit of detection in food: 3 ppm of gluten.

Limit of detection in surfaces: 1 µg of gluten.

KIT COMPONENTS

Component	ZE/PR/GL10	How to use the components
STRIP TESTS	10	
EXTRACTION BUFFER	1	Take it with one of the 10 ml syringe.
ANALYSIS BUFFER	1	Take it with the 1 ml syringe.
FILTER BAGS	10	To manipulate the sample during the analysis.
PIPETTES	10	To take the sample.
TUBES	20	To make the dilutions.
SYRINGES (1 ml and 10 ml)	3	To take Extraction and Analysis Buffer
SWABS	10	To analyse surfaces.
PRODUCT CERTIFICATE	Yes	

SAFETY

A SAFETY DATA SHEET is available from your local distributor or ZEULAB. Good laboratory practices would be advised.

STABILITY AND STORAGE

Store the kit under dry conditions and temperatures between 15 and 25°C.

The kit has a shelf life of 12 months. Check the expiration date on the package.

GENERAL INFORMATION BEFORE YOU BEGIN

- The strips are very sensitive to moisture, keep them inside the supplied container and fully closed.
- The test should be performed between 15 and 25 °C.
- Samples must be representative (at least 50 g) and homogeneous.

- The provided syringe should be used to add the analysis buffer to each sample. Extreme care is needed not to touch the bag content to avoid contamination.
- To use the provided pipettes, squeeze the bulb, put the stem in the sample and release the bulb to draw up liquid. The central capillary will fill up and the excess liquid will be siphoned automatically into overflow tanks. Squeeze the bulb again to release 100 µl from the central capillary.
- The kit is designed for detecting gluten traces. In the analysis of food, very high concentrations of gluten can lead to content misinterpretation. This is due to "prozone or hook effect". For this reason, gluten cereal samples can not be directly analysed. Dilutions higher than 1/1000 are recommended for the analysis.
- It is recommended to make an initial evaluation to test matrix, specificity and sensibility.
- For manipulating the strip, hold it from the extreme opposite end of the arrows area.
- Please, contact ZEULAB for further information.

NOTES for food analysis and rinse water

- Additional material
 - Balance
 - Beater/Grinder
- Use one filter bag, pipette, tube and swab for each sample.
- Make sure that the sample does not reach the area of the test strip arrows.
- To use the filter bag, add the sample and the buffer inside the inner bag. Close it with the help of the wire, folding it upon itself.
- For **samples containing polyphenols, tannins or antioxidants** special procedure must be followed:
0.4 g of Polyvinylpyrrolidone (Sigma PVP-10) and 1 g of gelatine (fish gelatine) to the extraction buffer and continue the procedure in step 2. These samples must be incubated at 50°C for 40 minutes. The minimum level of detection for these samples is 10 ppm.

NOTES for surfaces analysis


- The area should be visually cleaned.
- Swab a representative area of at least 100cm².
- The result obtained from the test is the presence or absence of the allergen on the analysed surface, and it cannot be extrapolated to ppm values.
- The test limit of detection on surfaces is subjected to the use of the swabs provided by ZEULAB.

ANALYSIS OF FOOD

Please, read **INFORMATION PRIOR TO BEGIN THE ASSAY**.

The following procedure can be applied to solid and liquid samples. For liquid samples such as rinsing water or juices the analysis can start in step 4.

1- Weigh 1 g (1 ml) of the sample inside the provided bag and add 10 ml of EXTRACTION BUFFER.

2- Rub the bag with your hands (1-2 min approximately) to obtain a homogeneous mixture.

3- The filtered sample (placed between the filter and the outer bag) is ready to be analysed.

4- For **detecting 20 ppm of gluten** (two consecutive dilutions are needed):

- Dilution 1: take 100 μ l of the filtered sample and add 10 ml of ANALYSIS BUFFER.

- Dilution 2: take 2 ml of the Dilution 1 and add 2 ml of ANALYSIS BUFFER.

5- Pour the diluted sample into the test tube cap.

6- Introduce the strip test into the tube cap following the arrows' direction.

7- Wait for 10 min and read the test result following indications on "Interpretation of the results". Wait for 10 min and read the test result following indications on "**Interpretation of the results**".


Dilutions for levels of detection of 3 and 10 ppm

Level of detection	3 ppm	10 ppm
Dilution	1/30	1/100
Sample Volume	100 μ l	100 μ l
ANALYSIS BUFFER Volume	3 ml	10 ml

ANALYSIS OF SURFACES

(Please, read NOTES before running the assay)

1. Add 0.5 ml of ANALYSIS BUFFER into a tube, use the 1 ml syringe.
2. Introduce a swab into the tube and soak it in the ANALYSIS BUFFER. The swab must be held and manipulated by its opposite end in order to avoid cross contamination.
3. Rub the swab along in all possible directions on the analysed surface. Rotate from time to time the swab's point in contact with the surface.
4. Place the swab in the initial tube with the ANALYSIS BUFFER and stir gently into the solution.
5. Dispose of the swab.
6. Introduce the strip test into the tube following the arrows' direction.
7. Wait for 10 min and read the test result following indications on "Interpretation of the results".


INTERPRETATION OF RESULTS

NEGATIVE: Only a blue line is observed in the white central zone of the strip test. This is the control line and it indicates that the test works properly.

POSITIVE: Two lines, red and blue, are observed in the white central zone.

When no blue line is shown, the result must be considered **INVALID**. This may occur when the assay has not been performed properly, when the matrix interferes or when the reagents are damaged. Repeat the assay with a new strip.

PROTEON is an *in vitro* diagnostic kit. In analysis implicating legal processes, the results should be reevaluated with an official reference method. ZEULAB does not assume any legal responsibility.

To the best of our knowledge, the information contained herein is accurate and complete. However, nothing herein shall be construed to imply any warranty or guarantee.

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0.5 ml of AB
en 0,5 ml de AB*