PROTEON DUO MILK EXPRESS

ZE/PR/LBC10
ZE/PR/LBC10SW

Immunochromatographic test
for milk detection

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SCOPE AND TEST PRINCIPLE

PROTEON DUO MILK EXPRESS is a qualitative immunochromatographic test, in strip format, to detect milk proteins in food and working surfaces. The test is based on the detection of casein and β -lactoglobulin (β -LG), which react with specific antibodies coated to red particles in the strip.

The double detection of these main proteins is particularly useful because only one test is required to identify the presence of milk, independently of its origin. In addition, the detection of both groups of proteins is an advantage when testing highly processed foods, as caseins are more resistant to high temperatures.

- Limit of detection in food:
Casein: 4.2 ppm (5.0 ppm of milk proteins or 14.0 ppm powder milk MoniQA)*
 β -LG: 0.1 ppm (1.0 ppm of milk proteins or 2.8 ppm powder milk MoniQA)*
- Limit of detection in surfaces:
0.7 μ g of milk proteins (2.0 μ g powder milk MoniQA)*

Consult the technical sheet or contact ZEULAB for more information.

*References: Bobe et al., 2007; Ng-Kwai-Hang, 1986; MoniQA reference material Certificate

SAFETY

A SAFETY DATA SHEET is available from your local distributor or ZEULAB. This kit should be used following good laboratory practices.

STABILITY AND STORAGE

Store the test under dry conditions and temperatures between 4 and 12°C. Check the expiry date on the package.

PROTEON is an *in vitro* diagnostic kit. In analysis implicating legal processes, the results should be re-evaluated with an official reference method. ZEULAB does not assume any legal responsibility. To the best of our knowledge, the information contained herein is accurate and complete. However, nothing herein shall be construed to imply any warranty or guarantee.

KIT COMPONENTS

	ZE/PR/LBC10	ZE/PR/LBC10SW	How to use the components
ANALYSIS BUFFER	1		
ANALYSIS BUFFER TUBE (AB) (0.5 ml)		10	
FILTER BAG	10		To filter the sample during the analysis
SYRINGE (10 ml)	1		To take the analysis buffer for food and rinse water analysis
PIPETTES	10		To take the sample
TUBES	10		To put the sample in and incubate with the strip
RACK	1		Support for tubes
SYRINGE (1 ml)	1		To take the analysis buffer for surface analysis
SWABS	10	10	To swab the working surfaces
TEMPLATES		10	
STRIP TESTS	10	10	

- With reference ZE/PR/LBC10; food, rinse water and surfaces can be analysed.
- With reference ZE/PR/LBC10SW only surfaces can be analysed.

NOTES


- The strips are very sensitive to moisture, keep them inside the supplied container and fully closed.
- The test should be performed between 15 and 25 °C.
- 15 minutes before performing the analysis, remove the container with the strips from the refrigerator. Temper it without opening it to avoid condensation.
- Use a filter bag, a pipette, tube and swab for each sample.
- For manipulating the strip, hold it from the opposite end of the arrows area.
- The test has been validated in a large number of representative matrices. It is recommended to verify that the test works with your specific

matrix before using it. To do this, analyze at least one positive sample of your matrix and one negative to ensure that the test works correctly. It is also recommended to re-verify it if there are changes in the ingredients or in the production processes of the matrix. Please, contact ZEULAB for further information.

- **WARNING:** The kit detects allergen traces. Due to the high sensitivity of the Proteon Duo Milk Express kit, high concentrations of milk may lead to an invalid or negative result due to the “prozone or Hook effect”. Dilutions higher than 1/1000 are recommended for milk analysis.

NOTES for food analysis and rinse water (Ref. ZE/PR/LBC10)

- Additional material:
 - Beater/grinder
 - Balance
- For using the filter bag, add sample and buffer inside the inner bag. Close it with the help of the wire, folding it upon itself.
- Make sure that the sample does not reach the area of the test strip arrows.

NOTES for surfaces analysis (Ref. ZE/PR/LBC10 and ZE/PR/LBC10SW)

- The area should be visually clean.
- Swab a representative area of at least 100 cm².
- To delimit the area, you can use the plastic template provided with the kit in ref. ZE/PR/LBC10SW.
- The result obtained from the test is the presence or absence of the allergen on the analysed surface, and it cannot be extrapolated to ppm values.
- The test limit of detection on surfaces is subjected to the use of the swabs provided by ZEULAB.

PROCEDURE

(Please, read NOTES before running the assay)

► Sample Preparation**A. Analysis of food and rinse water (Ref. ZE/PR/LBC10)**

1. Sampling
 - For solid samples, take a representative sample (50 g) and mash or crush to homogenise. Weigh 1 g inside of the filter bag and add 10 ml of ANALYSIS BUFFER using the 10 ml syringe.
 - For liquid samples and rinse water, place 1 ml of sample (or weigh 1 g) into de filter bag and add 9 ml of ANALYSIS BUFFER using the 10 ml syringe.
2. Rub the bag with your hands to obtain a homogeneous mixture (1-2 min approximately).
3. The filtered sample is ready to be analysed (go to section Assay).

B. Analysis of surfaces

1. Prepare the ANALYSIS BUFFER tube:
 - Reference ZE/PR/LBC10. Add 0.5 ml of ANALYSIS BUFFER into a tube using the 1 ml syringe.
 - Reference ZE/PR/LBC10SW. Open a tube of ANALYSIS BUFFER.
2. Introduce a swab into the tube and soak it in the ANALYSIS BUFFER. The swab must be held and manipulated by its opposite end in order to avoid cross contamination.
3. Rub the swab along in all possible directions on the analysed surface. Rotate, from time to time, the swab's point in contact with the surface.
4. Place the swab in the initial tube with the ANALYSIS BUFFER and stir gently into the solution.
5. Dispose of the swab.

► Assay

1. Take the necessary amount of the prepared sample for the assay:
 - Reference ZE/PR/LBC10:
 - Surface analysis. The sample is already prepared for the assay.

- Other matrices. Add 6 drops of filtered sample using the disposable pipette in the provided tube. Use a new disposable pipette and tube for each sample.
 - Reference ZE/PR/LBC10SW. The sample is already prepared for the assay.
2. Introduce the strip test into the tube following the arrows' direction.
 3. Wait for 10 minutes and read the test result following indications on "Interpretation of Results".

INTERPRETATION OF RESULTS

(See the "Results" figure on the last page).

NEGATIVE: Only a blue line is observed in the white central zone of the strip test. This is the control line and it indicates that the test works properly.

POSITIVE: For a positive test, in addition to the blue line, one or two red lines appear in the white central zone.

INVALID: When no blue line is shown. This may occur when the assay has not been performed properly, when the matrix interferes or when the reagents are damaged. Repeat the assay with a new strip and contact ZEULAB for further information.

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1mL / 9 ml AB
0.5 ml of AB
ml de AB