## PROTEON FISH EXPRESS

ZE/PR/P10

ZE/PR/P10SW

Immuno-chromatographic test  
for fish detection

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[www.zeulab.com](http://www.zeulab.com)

**SCOPE AND TEST PRINCIPLE**

PROTEON FISH EXPRESS is a qualitative immunochromatographic test, in strip format, to detect fish proteins using a combination of antibodies that detects a common antigen in all vertebrate fish. The cartilaginous fishes (shark and sturgeons) can be also detected but with less sensitivity. Considering fish proteins is one of the main responsible of food allergic reactions, this target can be use as an indicator of presence of fish proteins in food, rinse water and working surfaces. The test is based on the detection of an antigen present in fish muscle, which reacts with a specific antibody coated to red particles in the strip. The appearance of two red lines, on the strip test is a positive result.

Limit of detection in food matrix: 2 ppm fish proteins.

Limit of detection in surfaces: 0.25 µg of fish proteins.

Consult the technical note or contact ZEULAB for more information.

**KIT COMPONENTS**

	<b>ZE/PR/P10</b>	<b>ZE/PR/P10SW</b>	<b>How to use the components</b>
ANALYSIS BUFFER	1		
ANALYSIS BUFFER TUBE (AB) (Tubes 0.5 ml)		10	
FILTER BAG	10		To filter the sample during the analysis
SYRINGE (10 ml)	1		To take the analysis buffer for food analysis and rinse water analysis
PIPETTES	10		To take the sample
TUBES	10		To put the sample in and incubate with the strip
RACK	1		Support for tubes
SYRINGE(1 ml)	1		To take the analysis buffer for surface analysis
SWABS	10	10	To swab the working surfaces
TEMPLATES		10	
STRIP TESTS	10	10	
PRODUCT CERTIFICATE	YES	YES	

- With reference ZE/PR/P10; food, rinse water and surfaces can be analysed.
- With reference ZE/PR/P10SW only surfaces can be analysed.

## SAFETY

A SAFETY DATA SHEET is available from your local distributor or ZEULAB. This kit should be used following good laboratory practices.

## STABILITY AND STORAGE


Store the test under dry conditions and between 15 and 25°C. The kit has a shelf life of 12 months. Check the expiry date on the package.

## NOTES

- The strips are very sensitive to moisture, keep them inside the supplied container and fully closed.
- The test should be performed between 15 and 25 °C.
- The kit is designed for detecting allergen traces. In the analysis of food, very high concentrations of the allergen can lead to content misinterpretation. This is due to "prozone or Hook effect".
- It is recommended to make an initial evaluation to test matrix, specificity and sensibility. Please, contact ZEULAB for further information.

### NOTES for food analysis and rinse water (Ref. ZE/PR/P10.)

- Additional material:
  - Beater/grinder
  - Balance
- Use a filter bag, a pipette, tube and swab for each sample.
- For using the filter bag, add sample and buffer inside the inner bag. Close it with the help of the wire, folding it upon itself.
- Make sure that the sample does not reach the area of the test strip arrows. For manipulating the strip, hold it from the opposite end of the arrows area.

NOTES for surfaces analysis (Ref. ZE/PR/P10 and ZE/PR/P10SW)

- In the analysis of surfaces, the area should be visually clean.
- It is recommended to swab a representative area of at least 100 cm<sup>2</sup>.
- To delimit the area, you can use the plastic template provided with the kit in ref. ZE/PR/P10SW.
- The result obtained from the test is the presence or absence of the allergen on the analysed surface, and it cannot be extrapolated to ppm values.
- The test limit of detection on surfaces is subjected to the use of the swabs provided by ZEULAB.
- Please, contact ZEULAB for further information.

**PROCEDURE**

ANALYSIS OF FOOD AND RINSE WATER (Ref. ZE/PR/P10)

*(Please, read NOTES before running the assay)*

A. Sample preparation

A.1 Solid samples

1. Take a representative sample (50 g) and mash or crush to homogenise. Weigh 1 g inside of the filter bag and add 10 ml of ANALYSIS BUFFER using the 10 ml syringe.
2. Rub the bag with your hands to obtain a homogeneous mixture (1-2 min approx.)
3. The filtered sample is ready to be analysed (paragraph B "Assay").

A.2 Liquid samples

1. Take 1 ml of sample (or weigh 1 g) and 9 ml of ANALYSIS BUFFER using the 10 ml syringe. Add them inside of the filter bag.
2. Rub the bag with your hands (1-2 min approx).

B. Assay

1. Add 6 drops of filtered sample using the disposable pipette in the provided tube. Use a disposable pipette and tube for each sample.
2. Introduce the strip test into the tube following the arrows' direction.

3. Wait for 10 min and read the test results following indications on "Interpretation of Results".

#### ANALYSIS OF SURFACES (Ref. ZE/PR/P10)

##### A. Sample preparation

1. Add 0.5 ml of ANALYSIS BUFFER into a tube using the 1 ml syringe.
2. Introduce a swab into the tube and soak it in the ANALYSIS BUFFER. The swab must be held and manipulated by its opposite end in order to avoid cross contamination.
3. Rub the swab along in all possible directions on the analysed surface. Rotate, from time to time, the swab's point in contact with the surface.
4. Place the swab in the initial tube with the ANALYSIS BUFFER and stir gently into the solution.
5. Dispose of the swab.

##### B. Assay

6. Introduce the strip test into the tube following the arrows' direction.
7. Wait for 10 minutes and read the test result following indications on "Interpretation of Results".

#### ANALYSIS OF SURFACES (Ref. ZE/PR/P10SW)

##### A. Sample preparation

1. Open a tube of ANALYSIS BUFFER.
2. Introduce a swab into the tube and soak it in the ANALYSIS BUFFER. The swab must be held and manipulated by its opposite end in order to avoid cross contamination.
3. Rub the swab along in all possible directions on the analysed surface. Rotate, from time to time, the swab's point in contact with the surface.
4. Place the swab in the initial tube with the ANALYSIS BUFFER and stir gently into the solution.
5. Dispose of the swab.

## B. Assay


- Introduce the strip test into the tube following the arrows' direction.
- Wait for 10 minutes and read the test result following indications on "Interpretation of Results".

**INTERPRETATION OF RESULTS**

**NEGATIVE:** Only a red line is observed in the strip test. This is the control line and it indicates that the test works properly.

**POSITIVE:** Two red lines are observed in the white central zone.

**INVALID:** When no line or either only bottom red line is shown. This may occur when the assay has not been performed properly, when the matrix interferes or when the reagents are damaged. Repeat the assay with a new strip and contact ZEULAB for further information.

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**PROTEON** is an *in vitro* diagnostic kit. In analysis implicating legal processes, the results should be re-evaluated with an official reference method. ZEULAB does not assume any legal responsibility. To the best of our knowledge, the information contained herein is accurate and complete. However, nothing herein shall be construed to imply any warranty or guarantee.

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1g/10 ml AB
1 mL/9 ml AB
0,5 ml of AB
de AB*