COMET8

Comet8 was formerly named e-Reader. *Comet8 se llamaba anteriormente e-Reader.*

ZEULAB, S.L.

C/ Bari, 25 dpdo. • 50197 Zaragoza (SPAIN)

Tel.: +34 976 731 533

info@zeulab.com

www.zeulab.com

DESCRIPTION AND SCOPE

COMET8 is an autonomous system that coupled to Eclipse Farm and Explorer 2.0 tests allows to simplify and standardize the screening of a broad range of antibiotics in milk, meat, egg, kidney and blood.

Eclipse and Explorer are bacteria inhibition growth tests, where the presence and absence of antibiotics is detected by changes in the media colour.

COMET8 incubates the tests and measures the colour in real time to monitor the bacteria growth kinetic. An internal software integrates colour and time parameters to determine the endpoint of the assay. COMET8 also interprets the results as qualitative and numerical values. It can be used with the touchscreen or with the COMET8 app and a smartphone or tablet. The results are displayed on the screen and in the app and stored for traceability.

CONTENTS AND REFERENCES

	ZE/COM8
COMET8	1
Power cable	1
USB cable	1

Reference	Model	N° of methods
ZE/COM8/XXX	Basic	1 method ^a
ZE/COM8/CTE	Complete	All methods

^a Any method can be chosen.

SPECIFICATIONS

Outer dimensions (cm): 17 x 7.5 x 23 cm

Weight: 0.97 kg

Working temperature: 65°C ± 0.5°C / 63°C ± 0.5°C

Screen type: Touch and colour screen

Screen size: 9.6 cm x 5.5 cm

Electrical requirement: 12V / 15V \equiv 1.2A / 1.6A

Power supply: Input 100-240 V ~ 50/60 Hz , 0.6 A ; Output 15V \equiv 1.6A , 24.0W



Light source: White diode LED 5700k

Measuring frequency: 2 min

Preheating time: 30 min

Computer interface: USB connection

Operating PC system: Windows 7, 8, 10

SAFETY PRECAUTIONS

- Pay attention to the following symbols and warnings:

Metallic surfaces can be hot during use. Please, do not touch them during operation.

- Only use the device under the specified conditions.
- Do not attempt to disassemble or alter the device
- ZEULAB is not responsible of an inadequate usage of the device.

GENERAL CONSIDERATIONS

- COMET8 must be connected to electrical power.
- COMET8 works with the touch screen or the COMET8 app.
- COMET8 needs approximately 30 min to reach and maintain the incubation temperature.
- In case of a power cut, COMET8 will continue running for a short time while the assay temperature does not drop 5°C below the set point temperature.
- Incubation temperature can be verified by introducing a calibrated temperature probe in the probe holder located in the heating block. Contact your distributor or Zeulab for more information.
- COMET8 stores more than 20.000 analyses. When the memory is full, new results will replace the previous ones.

AVAILABLE METHODS

Information on the available methods is shown in the table below.

METHOD	ABBREVIATION (XXX)	DESCRIPTION	CUT OFF
ECLIPSE SINGLE	ECS	For analysis of antibiotics in cow milk. Individual assays are performed. Each assay can be started individually. No negative control is needed.	60
ECLIPSE NC ECLIPSE NEGATIVE CONTROL	ENC	For analysis of antibiotics in cow milk. The assay is performed using a negative control (NC). Samples and NC are tested in one assay at the same time.	56
ECLIPSE NC SGA ECLIPSE NEGATIVE CONTROL SGA	SGA	For analysis of antibiotics in sheep milk, goat milk, and samples containing azidol. The assay is performed using a NC. Samples and NC are tested in one assay at the same time.	66
EXPLORER	EXP	For analysis of antibiotics in meat and egg. The assay is performed using a NC. Samples and NC are tested in one assay at the same time.	56
EXPLORER KIDNEY	EXK	For analysis of antibiotics in kidney. The assay is performed using a NC. Samples and NC are tested in one assay at the same time.	76
EXPLORER BLOOD	EXB	For analysis of antibiotics in blood. The assay is performed using a NC. Samples and NC are tested in one assay at the same time.	56

COMET8 SETTINGS

In "complete" devices, "Settings" is accessed from the main menu. On "basic" devices, this option is available in the lower right corner.

From the "Settings" option, date and time can be set and also access the "About" and "WiFi" options.

In "About" language can be set and the acoustic warning activated or deactivated. In addition, the configuration of your device, firmware version and serial number are displayed. "Demo" option should be in "off" position.

In "WiFi" the available WiFi networks can be scanned, the appropriate one selected and the connection checked. In addition, the firmware version can be updated by tapping on "Update". This also can be set using Comet8 app.

Changes to the configuration can only be made when no test is in progress.

COMET8 CONNECTIVITY

COMET8 communicates via Bluetooth with your mobile phone or tablet.

COMET8 sends the data to the cloud www.test4all.es via Wi-Fi and provides the results in real time wherever you are and without the need to be near the device. If COMET8 does not have Wi-Fi communication, the results will be uploaded to the cloud when the device and mobile phone are connected.

If the area where the assay is performed does not have internet coverage, COMET8 can be used with the touch screen or the app in offline mode. When COMET8 has internet access using Wi-Fi or mobile internet coverage, results will be synchronized in the cloud and emails will be sent.

TEST4ALL ACCOUNT

The Test4all account allows the management of devices, users and results. It is required for the use of the device with the Comet8 app and for downloading results from the device.

There are two categories of users: "Admin" and "User":

- "Admin" allows registering and modifying devices and users and accessing all the registered results and devices of the company.
- "User" can only use the assigned devices and see the results to which authorization has been given.

In the "Results" section, the results of each of the registered devices can be checked and downloaded all the data in CSV format for opening in excel.

Test4all account creation

Account creation can be accessed from www.test4all.es by clicking on "Sign up" or from the app by selecting "Create account".

Several fields must be completed and the device key, that appears on your Comet8 certificate, must be entered. Click on "Confirm account" to receive an email and complete the process.

COMET8 APP

The COMET8 application (app) is used for operation from the mobile phone or tablet. It can be downloaded from the Google Play Store (Android) and Apple App Store (iOS). Enter "Comet8 Zeulab" for your search.

Connecting to COMET8

Open the app, log in with your Test4all account credentials, tap on "Connect" and select the COMET8 device from the list. Wait for the connection.

App configuration

To configure the different parameters of the app, it must be connected to a COMET8.

Language: allows you to select the language

User: allows you to log out the user. The next time the app is opened, it will ask for the user and password.

WiFi: allows you to configure the WiFi of the Comet8 device.

About: allows you to view the firmware version loaded in the device and update it. Here you can also see the version of the app.

PROCEDURE

Assays can be performed by using the touch screen or the Comet8 app.

Before starting an assay

- Switch on COMET8 and wait until the instrument reaches incubation temperature.
- Review the test protocol and follow the kit instructions.

Start an assay

Using the touch screen

- Select the method of analysis if using COMET8 "complete".
- Select sample position 1 to 8 if using Eclipse Single. If using a method with negative control, press any of the position numbers.
- Enter the "Batch Code" indicated on the product certificate of the kit. The Batch Code is the only information required to start the assay. It is recommended to add more information using the "Assay data" option to ensure traceability, although it is not mandatory.
- Place the tubes with the samples in the selected positions by pressing them to the bottom. For methods with negative control, remember to insert the negative control in position 1. Seal the tube with the adhesive film.
- Touch "Start" to start the test.

Using the Comet8 app

- Open the app, log in with your Test4all account credentials, tap on "Connect" and select the COMET8 device from the list. Wait for it to connect.
- Tap on "Launch test". If you are using a COMET8 "complete" version, select the method of analysis.

ENGLISH

- Select the position (1 to 8) that the sample will occupy. If it is a Negative Control (NC) method, all samples will be selected at once.
- Tap "Batch code" to enter the code that appears on the Eclipse or Explorer product certificate. This is the only information required to start the assay. It is recommended to add the kit lot number and sample identification under "Test details" to ensure traceability, although it is not mandatory.
- Place the tubes with the samples in the selected positions by pressing them to the bottom. For methods with negative control, remember to insert the negative control in position 1. Seal the tube with the adhesive film.
- Touch "Start test" to start the test.

Canceling an assay

Tests can be cancelled at any time during the current assay by clicking on the test position and confirming the cancellation. If a negative control method is being used, the entire assay will be cancelled.

RESULTS

Once the test is finished (2h - 3h15, depending on the test used) the results will be displayed on the Comet8 screen and on the mobile to which the equipment is connected in the "status" and/or "results" section.

The results are assigned a numerical value and a colour: green for negative and red for positive. The cut-off point has been calculated for each method (see table page 3). Results equal to or higher than the cut-off point are identified as positive.

If Comet8 has Wi-Fi configured and active, the results are uploaded to the Test4all.es cloud where they can be consulted at any time, and emails are automatically sent to the configured accounts.

If Comet8 does not have Wi-Fi configured or active, the results can be uploaded to the cloud by connecting the device to the app. Emails will then also be sent to the configured accounts.

When Comet8 is connected to the app, the last 100 results are displayed in the "results" section. By clicking on the reference number of the test you can see the complete information and share the result to another email from the "Share" button.

The results can be exported to a .csv file for opening with Excel from the results section of www.test4all.es. It is recommended to save the results periodically.

MAINTENANCE AND CLEANING

70% ethanol may be used to clean the device.

It is recommended that the optics and temperature of the device be checked annually. Contact your distributor or Zeulab for more information.

MESSAGES AND TROUBLESHOOTING

The following table provides a brief description of the most common error messages and troubleshooting.

MESSAGE	DESCRIPTION	TROUBLESHOOTING
Waiting for assay temperature	Temperature is below 65 / 63°C	Wait until COMET8 reaches the assay temperature.
Batch code not valid	Batch Code has not been entered or is incorrect.	The Batch Code can be found in the product certificate of every test.
NTB	<ul style="list-style-type: none"> No tube in the assay position chosen Tube is not placed correctly 	Place the test tube in the heating block or push the tube down to the bottom. Start a new assay.
NPS	<p>There has been a power cut for a few minutes.</p> <p>The power supply has not recovered in time and the test temperature has dropped 5°C below the setpoint temperature.</p>	Check the power supply of the equipment

ENVIRONMENTAL OPERATING CONDITIONS

Location	Indoor use. No direct sunlight, avoid dusty atmospheres.
Altitude	< 2000m
Relative humidity	≤ 80%, 30°C; ≤ 50%, 31 ~ 40°C
Temperature	5 ~ 40°C

REFERENCES

- Mata L., Sanz D., Razquin P. (2015). *Performance of Eclipse Farm Test Coupled to e-Reader for Antibiotic Residues Detection in Raw Milk*, Food Analytical Methods, pp 1-9
- Reybroec W., Ooghe S., ILVO (2015). *Validation of the Eclipse Farm 3G + e-Reader*.
- Mata L., Sanz D., Razquin P. (2014) *Validation of the Explorer® 2.0 test coupled to e-Reader® for the screening of antimicrobials in muscle from different animal species*, Food Additives & Contaminants: Part A, 31:9, 1496-1505.
- Razquin, P. et al. Validation of the Eclipse Farm 4G & COMET for Detection of Antibiotics in Raw Bovine Milk: AOAC Performance Tested MethodSM 022101. J. AOAC Int. 2021
- Explorer 2.0 kit insert, G-COM-EX.08.
- Eclipse Farm kit insert, G-COM-EF.01.
- Eclipse Farm 4G kit insert, G-COM-EF.12.

INPUT PANEL

State	ON-OFF
Assay ref.	Number assigned by COMET8. It is essential to track the result.
Method/Met	Method used by COMET8 to analyse the sample
Sample ID	Sample identification. Optional. With the exception of Eclipse single. Position 1 is reserved for negative control sample
Time	Elapsed time (minutes).
Value	Arbitrary units assigned by COMET8 to each finished assay.
Results	Qualitative results: positive (POS), negative (NEG), without tube (NTB), cancelled (CAN), no power supply (NPS).
<	To go back to previous screen.
Batch code	Code that provides to COMET8 the necessary information to run the assay. It can be found in the kit product certificate of each test tube. Compulsory.
Start	To begin the assay
Assay details	To enter extra Information about the assay. Optional.
Batch No.	To enter the batch number. It can be found at the label box of the kit. Optional.
User	To enter user name or code. Optional.
Date	Date is assigned by the COMET8.
Save	To save data.
< (keyboard)	To delete data.
OK (keyboard)	To accept and continue.
See Data / Results	To see the information about finished assays. Data analysis are listed by the assay reference number in descending order and data can be exported in a PDF and Excel file.
▲ ▼ ⏮ ⏭	To move up or down in data list.
🔊	To turn audio warning on.
🔇	To turn audio warning off.
+/- </>	To modify device settings parameters.

COMET8 is an in vitro diagnostic system for antibiotic screening. In analysis implicating legal processes, the results should be validated with an official reference method. ZEULAB, S.L. does not assume any legal responsibility.

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

2 *Select position tapping the numbers and enter the "Batch code" indicated on the kit product certificate.

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3 Place the tubes with the samples, seals de tubes and Tap "Start". (If the method requires negative control sample, place it in position 1).

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4 Numerical and qualitative results (+ or -) will be displayed on the screen at the end of the assay.

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*Users of COMET8 complete should select the method previously.

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