A : GRIGLIA - GRID - GRILLE - CUADRÍCULA - RASTER - ROOSTER
 B : QUADRATO - SQUARE - CARRÉ - CUADRADO - QUADRAT - VIERKANT
 C : SETTORE - SECTOR - SECTEUR - SECTORES - SEKTOR - SEKTORN

ITALIANO

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$$[\text{Cellule/ml}] = (67 / 5) \times 10^2 \times 10^4 = 13,4 \times 10^6$$

ENGLISH

CALCULATION METHOD FOR CELLS / μL IN URINARY SEDIMENT

FAST READ 102® is a disposable plastic device composed of 10 counting chambers. With each device you can analyze 10 samples.

Each room contains a GRID composed of 10 SQUARES, each of which is divided into 16 smaller squares (called SECTORS).

One of the advantages of employing FAST-READ is the ease in determining the cells per μL in the specimen.

Chamber volume reading: 7 μL


Grid size: 2 mm x 5 mm

Volume grid: 1 μL

Volume of the square: 0.1 μL

Each grid is composed of: 10 squares (1mm x 1mm)

Each square consists of: 16 sectors

- 1) Perform the count on a centrifuged fresh urine sample, after having decanted.
- 2) Gently resuspend the sediment.
- 3) Using a pipette, introduce the sample into the well and examine under the microscope the area of the grid
- 4) Count the number of cellular elements within N squares

$$\text{Cell.} / \mu\text{L} = \left(\frac{\sum \text{cells counted in square N}}{N} \right) \times \text{concentration factor} \times 10$$

Concentration factor = Volume of sediment / Volume centrifuged urine

10: conversion from 0,1 μl to 1 μl

For uncentrifuged urine, don't multiply the number of cells counted for the concentration factor.

METHOD OF CELL COUNTS FOR DILUTED SAMPLES (CELLS / ML).

After filling the counting chamber with the sample, proceed to the counting of cells distributed in N squares.

Considering that the grid consists of 10 squares, each square has a dimension of 1 x 1 mm, a depth of 0.1 mm and a volume of 0.1 μl, the formula for determining the concentration of cells (cells / ml) is:

$$\text{Cells} / \text{ml} = \left(\frac{\sum \text{cells counted in square N}}{N} \right) \times \text{dilution factor} \times 10^4$$

10⁴ = conversion 0,1 μl in 1 ml

Attention to the cells on the edges, you should only count those on either side, to avoid the risk of over or under.

In the example shown we perform the cell count of a sample diluted 100 times:

N = 5 (number of squares considered for counting)

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