



# LP ITALIANA SPA

Via C. Reale, 15/4 - 20157 Milano - ITALY  
Tel. +39 02 3933061  
Fax +39 02 39313484  
www.lpitaliana.com

Capitale Sociale € 309.600,00  
R.E.A. MI 882798  
Reg. Imp. MI C.F. e P.I. 01794050151  
info@lpitaliana.com

## TECHNICAL DATA SHEET

### WATER SAMPLING BOTTLES

Dedicated to optimise the sampling of water or food for microbiology analyses, LP bottles are designed and manufactured to assure functionality, safety and user-friendliness.

LP Sampling Bottles come in four sizes: 125 ml, 250 ml, 500 ml 1.000 ml.

Key features, common for all models:

- Manufactured with virgin raw material, in conformity of the environmental norms: PETG – PP - HDPE
- Strong and virtually unbreakable.
- Lightweight, well balanced and easy to handle. Its square cross section and rounded edges assure the best transport and storage optimization, while minimizing their packaging volume. All types are available in single package, too.
- Two types of opening : "narrow neck" - "wide neck"
- The volume of the inserted liquid can be easily controlled by means of the clear moulded graduation (see table); accuracy  $\pm 2\%$ .
- Closure: all bottles are closed by one leakproof screw cap with a special inert inner gasket and with a safety seal.
- Caps: all made in HDPE
- All closures are designed with vertical ribs for easy use and openness, even when wearing gloves
- Traceability : The product label on each box includes the lot number and expiration date. In addition, each bottle bears a pre-pasted label that allows the recording of the main sample identification data :


. - Batch number,

. - Expiry date

. - A unique number, in clear and in form of bar code.

All these allow full traceability of each individual sample.

- Sterilization: based on ISO 11137 - SAL 10-6 performed by means of ionizing radiation.
- Shelf life: see product tables.
- The whole manufacturing process, packaging included, is checked according applicable norms and LP's proprietary quality criteria
- Thiosulfate content, where not otherwise specified: 20 mg / l

Firma: 	Distribuito a: DG – DVE – VEE – VEI - UA	
Documento: RSSG 099 Rev. U	Pg. 1 of 6	Data Emissione: 23.07.2015

# LP ITALIANA SPA


## WIDE NECK BOTTLES FEATURES


produced in: PETG – PP – amber PP

Models with 50 mm nominal opening, facilitate the filling, minimize the risk of contamination during sampling (in accordance with the recommendations of the ISO 19458:2006 - sampling techniques):

Common characteristic:

Nominal Capacity (ml)	<b>125</b>	<b>250</b>	<b>500</b>	<b>1.000</b>
Brim Capacity (ml)	205	330	590	1.150
Graduation (ml)	25 – 125 (150)	50 - 250	100 - 500	200 - 1.000
Bottle Square (mm)	55 x 55	60 x 60	70 x 70	90 x 90
Bottle Height w/o cap (mm)	88	110	143	170
Bottle Height with cap (mm)	93	115	145	175
Cap Diameter (mm)	50			
Shelf life (months)	With Na Thiosulfate	24		
	w/o Na Thiosulfate	60		
Case quantity	350	216	120	72



Firma: 	Distribuito a : DG – DVE – VEE – VEI - UA
Documento: RSSG 099 Rev. U	Pg. 2 of 6
	Data Emissione: 23.07.2015

# LP ITALIANA SPA

**Square bottles in virgin PETG** Glass-like clarity for complete inspection during filling and pouring.

Code	Capacity (ml)	Single pack	Thiosulfate
<b>292155<sup>1</sup></b>	<b>250</b>	no	yes
<b>292158</b>	<b>250</b>	no	yes
<b>292258</b>	<b>250</b>	no	no
<b>292358</b>	<b>250</b>	yes	yes
<b>292458</b>	<b>250</b>	yes	no
<b>295155<sup>1</sup></b>	<b>500</b>	no	yes
<b>295158</b>	<b>500</b>	no	yes
<b>295198<sup>2</sup></b>	<b>500</b>	no	yes
<b>295258</b>	<b>500</b>	no	no
<b>295298<sup>2</sup></b>	<b>500</b>	no	no
<b>295358</b>	<b>500</b>	yes	yes
<b>295458</b>	<b>500</b>	yes	no
<b>291158</b>	<b>1.000</b>	no	yes
<b>291258</b>	<b>1.000</b>	no	no
<b>291358</b>	<b>1.000</b>	yes	yes
<b>291458</b>	<b>1.000</b>	yes	no



Transparency: \_\_\_ glass like  
Cap colour: \_\_\_ natural white

Norms:

- .- CE marking: \_\_\_ N/A
- .- food contact: \_ raw material complies with EC regulations 1935/2004 and 10/2011

Fields of application:

- .- industry
- .- microbiology
- .- environmental sampling

<sup>1</sup>Thiosulphate content : 100 mg / l ; <sup>2</sup>Resistant up to 100°C

**Square bottles in natural, virgin PP**

Code	Capacity (ml)	Single pack	Thiosulfate
<b>299148</b>	<b>125</b>	no	yes
<b>299248</b>	<b>125</b>	no	no
<b>299348</b>	<b>125</b>	yes	yes
<b>299448</b>	<b>125</b>	yes	no
<b>292148</b>	<b>250</b>	no	yes
<b>292248</b>	<b>250</b>	no	no
<b>292348</b>	<b>250</b>	yes	yes
<b>292448</b>	<b>250</b>	yes	no
<b>295148</b>	<b>500</b>	no	yes
<b>295248</b>	<b>500</b>	no	no
<b>295348</b>	<b>500</b>	yes	yes
<b>295448</b>	<b>500</b>	yes	no
<b>291148</b>	<b>1.000</b>	no	yes
<b>291248</b>	<b>1.000</b>	no	no
<b>291348</b>	<b>1.000</b>	yes	yes
<b>291448</b>	<b>1.000</b>	yes	no



Transparency: \_\_\_ translucent  
Cap colour: \_\_\_ natural white

Norms:

- .- CE marking: \_\_\_ N/A
- .- food contact: \_ raw material complies with CE regulation 1935/2004

Fields of application:

- .- industry
- .- microbiology
- .- environmental sampling

**Square bottles in amber virgin PP** To be used in those in cases, for example with drinking water, where it is appropriate to protect from light the potentially photosensitive elements eventually contained in the sample.

Code	Capacity (ml)	Single pack	Thiosulfate
<b>292648</b>	<b>250</b>	no	yes
<b>292748</b>	<b>250</b>	no	no
<b>292848</b>	<b>250</b>	yes	yes
<b>292948</b>	<b>250</b>	yes	no
<b>295648</b>	<b>500</b>	no	yes
<b>295748</b>	<b>500</b>	no	no
<b>295848</b>	<b>500</b>	yes	yes
<b>295948</b>	<b>500</b>	yes	no
<b>291648</b>	<b>1.000</b>	no	yes
<b>291748</b>	<b>1.000</b>	no	no
<b>291848</b>	<b>1.000</b>	yes	yes
<b>291948</b>	<b>1.000</b>	yes	no



Transparency: \_\_\_ no  
Cap colour: \_\_\_ natural white

Norms:

- .- CE marking: \_\_\_ N/A
- .- food contact: \_ N/A

Fields of application:


- .- industry
- .- microbiology
- .- environmental sampling


# LP ITALIANA SPA

## NARROW NECK BOTTLES - in natural, virgin HDPE.

The nominal diameter of the neck of 31 mm, facilitate both pouring precision and opening.  
 Body colour: matt, semi-translucent.  
 Suitable for water sampling, they are mainly used in food industry.

Nominal Capacity (ml)	<b>250</b>	<b>500</b>	<b>1.000</b>
Brim Capacity (ml)	290	550	1.090
Graduation - every 50ml - (ml)	50 - 250	100 - 500	100 - 1.000
Body Transversal section (mm)	52 x 52	70 x 70	80 x 90
Bottle Height w/o cap (mm)	135	151	195
Bottle Height with cap (mm)	138	154	198
Cap Diameter (mm)	34		
Shelf life (months)	With Na Thiosulfate	24	
	w/o Na Thiosulfate	60	
Case quantity	231	120	72



Firma: 	Distribuito a : DG – DVE – VEE – VEI - UA
Documento: RSSG 099 Rev. U	Pg. 4 of 6
	Data Emissione: 23.07.2015

# LP ITALIANA SPA

## NARROW NECK BOTTLES technical data

Code	Capacity (ml)	Single pack	Thiosulfate	Cap color
292168	250	No	yes	red
292268	250	No	no	blue
295168	500	No	yes	red
295268	500	No	no	blue
295768	500	No	no	blue
291168	1.000	No	yes	red
291268	1.000	No	no	blue
291768	1.000	No	no	blue
292368	250	Yes	yes	red
292468	250	Yes	no	blue
295368	500	Yes	yes	red
295468	500	Yes	no	blue
291368	1.000	Yes	yes	red
291468	1.000	Yes	no	blue



Transparency: \_\_\_\_\_  
translucent

Norms:

.- CE marking: \_\_ N/A

.- food contact: \_\_ raw material complies with CE regulation 1935/2004

Fields of application:

- .- industry
- .- microbiology
- .- environmental sampling


## SODIUM THIOSULFATE BUFFERING

- Chlorides: with chlorinated water samples, buffering of free or residual chlorine is suggested, to inhibit its bactericide action during transportation and storing. All LP bottles are available pre-dosed with Sodium Thiosulfate (see "Suggestions") in conformity with:
  - ISO 19458:2006 (same as French Standard NFT 90-40): \_\_20 mg/l or
  - International Standard ISO 5667-3: \_\_\_\_\_80 mg/l
 (Bottles filled to nominal capacity)

In case of highly chlorinated water, by customer request, bottles with higher quantity of Na Thiosulfate are available.

Note: in these cases the standard buffer amount is 100 mg/l, which can be increased up to 120 mg/l (typically used for swimming pool water sampling) with no influence on sample's quality (see "Suggestions"). Sodium Thiosulfate has no influence on the sample, that is why it is possible to use bottles containing Sodium Thiosulfate with non chlorinated samples.

- Sterility and individual wrapping: LP Water Sampling Bottles, with and without Na Thiosulfate, are all irradiated to a SAL (sterility assurance level) of  $10^{-6}$ . The sterilization is guaranteed until seal is broken and the cap is opened (see table for shelf life). For those applications where it is important to avoid any possible contamination induced by the bottle itself (i.e. sampling by immersion), sterility is necessary for the external surface of the bottle as well as the internal one. For this requirement, all LP Water Sampling Bottles can be supplied individually wrapped (flow pack).
- Guarantee: screw caps have a tamper evident safety seal. If the cap seal is unbroken closure and inner sterility are guaranteed.

Firma: 	Distribuito a: : DG – DVE – VEE – VEI - UA
Documento: RSSG 099 Rev. U	Pg. 5 of 6
	Data Emissione: 23.07.2015

# LP ITALIANA SPA

## SUGGESTIONS FOR THE BEST USE

- Microbiological analysis of water samples must always be performed within the shortest possible time after sampling. If the sample has aged too long, many factors can influence its bacterial contents, all of them related to the quality of the water: presence of toxic or nutritive substances for bacterial flora, saltiness, pH and so on.

As a general guideline, LP ITALIANA suggests to transport and store samples at a temperature of (+4 to +10)°C and to analyse them within 24 hours.

- 2) Buffering ratio

As a bactericide, sporicide, fungicide and viroicide, a sodium salt (hypochlorite NaClO and/or chlorite NaClO<sub>2</sub>) is usually dissolved in water or, more frequently, one mixture of the two. It is not possible to know neither the composition of the mixture, nor the amount of each dissolved salt, therefore it is not possible to supply any indication of how much free chlorine is to be buffered by a predetermined amount of sodium thiosulfate (Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>).

Furthermore, depending on inactivation dynamics, it is hard to predict which amount of Sodium Thiosulfate is required to inactivate even a known amount of free residual chlorine.

We suggest considering the following indications:

Buffering ratio between thiosulfate and hypochlorite → 1 Mole: 1 Mole

Buffering ratio between thiosulfate and Chlorite → 4 Moles: 1 Mole


In order to know the actual ratio in weight, such relationships must be related to the respective molecular weights, but quantitative indications are not set here because they could be not pertinent or misleading to each single case. We limit ourselves only to give evidence that in one limit case (all chlorite) the amount of Thiosulfate required is approximately quadruple than the one required for the opposite case ("all hypochlorite").

We suggest considering that 18 mg Sodium Thiosulfate will inactivate at least 2 mg/l and up to 5 mg/l of free chlorine residual, which is sufficient for the majority of samples.

## DISPOSAL

All components of LP bottles (bottles, caps and gaskets) are manufactured using virgin raw materials that are completely recyclable and environment friendly; they can be disposed:

- Emptied and cleaned are fully recyclable; or
- Emptied, burned with energy recovery; or
- By incineration, at an appropriate temperature and with an adequate air insufflations. They give residuals of H<sub>2</sub>O and CO<sub>2</sub>, only; or
- Delivered to a controlled landfill, in accordance with local laws on waste disposal.

Firma: 	Distribuito a: : DG – DVE – VEE – VEI - UA	
Documento: RSSG 099 Rev. U	Pg. 6 of 6	Data Emissione: 23.07.2015