



Disposable Hemocytometer

System Neubauer Improved

DHC-N01

NanoEnTek

## NanoEnTek Inc.

Website : www.nanoentek.com  
E-mail : sales@nanoentek.com



**NanoEnTek Inc.**  
851-14, Seohaeh-ro, Paltan-myeon, Hwaseong-si, Gyeonggi-do, 18531, Korea  
Tel: +82-2-6220-7940 / Fax: +82-2-6220-7999

EC REP

**MT Promedt Consulting GmbH**  
Altenhofstrasse 80, 66386 St. Ingbert, Germany

US Branch  
**NanoEnTek USA Inc.**  
5627 Stoneridge Drive Suite 304, Pleasanton, CA 94588, USA  
Tel : +1-925-225-0108, +1-888-988-0108(Toll free) / Fax : +1-925-225-0109

All the materials in this manual are protected by Korean and international copyright laws. This manual may not be reproduced, translated, published or distributed by any means without the permission of the copyright owner.

The information in this manual is described as correctly as possible, but may be changed without prior consent or notification.

Documentation : NESMU-DHC-001E (V.0.5)

Copyright 2013, by NanoEnTek Inc.  
All rights reserved. Published in Korea.

## Unpacking

When you receive the C-Chip (DHC-N01) for the first time, you will find the following components in your package.

Disposable hemocytometer  
Instruction manual

## Safety Precautions

For analyzing hazardous or potential infectious materials:

- Take necessary precautions
- Handle with care
- Dispose in an appropriate way

Long exposure to solvents will cause the slide to warp. Xylene and toluene based mounting media should be avoided. Glycerol, gelatin, and other aqueous-based media are recommended.

## Safety Symbols

The safety symbols on the C-Chip (DHC-N01) are intended to inform you of potential danger or a particular caution. Before use, please read and consult the guide for the symbols and their meanings.

Batch code (Lot Number) Use by

Do not reuse Manufactured by

Consult instructions for use

**NOTE :** The C-Chip (DHC-N01) is for **single use only. Do not reuse.** It should be used immediately after unsealing. The warranty on the C-Chip included in the conditions of supply is valid for 24 months from the date of manufacturing. The **expiration date** is indicated on the front side of outer box.

## Introduction

The C-Chip (DHC-N01) is a disposable plastic hemocytometer used for manual cell counting. It consists of surface-patterned two enclosed chambers with two ports for sample injection (Fig. 2).

The DHC-N01 grid pattern is exactly same as the Neubauer Improved. It consists of 9 large squares, each measuring 1 x 1 mm, and the depth of the chamber is 0.1 mm. Each square has a total volume of 0.1 mm<sup>3</sup> or 10<sup>-4</sup> cm<sup>3</sup> (Fig. 1).

The central square is divided into 25 small squares with triple lines and four corner squares are divided into 16 small squares.





Figure 1 Grid pattern

Figure 2 DHC-N01



## Counting with C-Chip

## A. General Methods

1. Mix the samples well.
2. Load 10 μL of sample into the sample injection area in Fig. 2, so that it fills the chamber by capillary action.
3. Count the cells under the microscope.

$$\text{Cells per mL} = \text{average count per square} \times \text{dilution factor} \times \text{volume factor}$$

## B. Mammalian cell counting

1. Treat the cell samples with trypsin-EDTA.
2. Carefully remove the supernatant with a pipette tip without disturbing the pellet.
3. Add an appropriate volume of growth media or PBS to dilute to a final concentration of 5 x 10<sup>3</sup> cells/mL to 5 x 10<sup>6</sup> cells/mL.
4. Thoroughly re-suspend the cell pellet with a pipette
5. Check visually if there are any cell clumps or agglomerates.
6. Load 10 μL of sample into the sample injection area in Fig. 2.
7. Count the cells in 5 large squares.

$$\text{Cells per mL} = \frac{\text{cells in 5 large squares}}{5} \times \text{dilution factor} \times 10^4 (\text{volume factor})$$

## C. Erythrocyte counting (1:200 dilution)

1. Dilute blood using accepted laboratory methods.
2. Load 10 μL of diluted sample into the sample injection area in Fig. 2.
3. Count the erythrocytes in the 5 small squares (four small corner squares and one small middle square) of the large center square.

$$\text{RBCs per mL} = \frac{\text{cells in 5 corner squares}}{4} \times 20 (\text{dilution factor}) \times 10^4 (\text{volume factor})$$

## D. Leukocyte counting (1:20 dilution)

1. Dilute blood using accepted laboratory methods.
2. Load 10 μL of diluted sample into the sample injection area in Fig. 2.
3. Count the leukocytes in the 4 large corner squares.

$$\text{WBCs per mL} = \frac{\text{cells in 4 corner squares}}{4} \times 20 (\text{dilution factor}) \times 10^4 (\text{volume factor})$$

## Trouble shooting

In case of poor visibility results,

Carefully load samples into the C-Chip in order to prevent the introduction of **air bubbles**. Observe after removing the dust from samples. Adjust the focus of the microscope.

# C-Chip

## INSTRUCTIONS

Grid pattern of Bürker  
**DHC-B01**

Grid pattern of Bürker-Türk  
**DHC-B02**

Grid pattern of Fuchs-Rosenthal  
**DHC-F01**

Grid pattern of Malassez  
**DHC-M01**



NanoEntek


### Grid pattern of Bürker **DHC-B01**

The C-Chip (DHC-B01) grid pattern is exactly same as the Bürker. It consists of 9 large squares, each measuring 1 x 1 mm, and the depth of the chamber is 0.1 mm. Each square has a total volume of 0.1 mm<sup>3</sup> or 10<sup>-4</sup> cm<sup>3</sup>.

The large squares are subdivided into 16 group squares with 0.2 mm sides. The group squares correspond in size to the Neubauer grid, but have no further subdivisions.

$$\text{Loading volume} = 10\mu\text{L}$$

$$\text{Cells per mL} = \text{average count per large square} \times \text{dilution factor} \times 10^4 \text{ (Volume factor)}$$



Grid pattern of DHC-B01

### Grid pattern of Bürker-Türk **DHC-B02**


The C-Chip (DHC-B02) grid pattern is exactly same as the Bürker-Türk. It consists of 9 large squares, each measuring 1 x 1 mm and the depth of the chamber is 0.1 mm.

Each square has a total volume of 0.1 mm<sup>3</sup> or 10<sup>-4</sup> cm<sup>3</sup>.

The large squares are subdivided into 16 group squares with 0.2 mm sides. In the central large square, each group is subdivided into 16 mini squares with 0.05 mm sides (= 0.0025 mm<sup>2</sup>).

$$\text{Loading volume} = 10\mu\text{L}$$

$$\text{Cells per mL} = \text{average count per large square} \times \text{dilution factor} \times 10^4 \text{ (Volume factor)}$$



Grid pattern of DHC-B02


### Grid pattern of Fuchs-Rosenthal **DHC-F01**

The C-Chip (DHC-F01) grid pattern is exactly the same as the Fuchs-Rosenthal. It consists of 16 large squares orientated by triple lines, each measuring 1 x 1 mm, giving a total area 4 x 4 mm.

The depth of each chamber is 0.2 mm, giving a on large square with triple line has a volume of 0.2 μL, total volume for counting area of 3.2 μL (3.2 mm<sup>3</sup>).

$$\text{Loading volume} = 20\mu\text{L}$$

$$\text{Cells per mL} = \text{average count per large square} \times \text{dilution factor} \times 2 \times 10^5 \text{ (Volume factor)}$$



Grid pattern of DHC-F01


### Grid pattern of Malassez **DHC-M01**

The C-Chip (DHC-M01) grid pattern is exactly same as the Malassez. It consists of 25 large squares, each measuring 0.2 x 0.25 mm and the depth of the chamber is 0.2 mm. Each square has a total volume of 0.01 mm<sup>3</sup> or 10<sup>-5</sup> cm<sup>3</sup>.

The large squares are subdivided into 20 small squares with an area of 0.0025 mm<sup>2</sup>.

$$\text{Loading volume} = 20\mu\text{L}$$

$$\text{Cells per mL} = \text{average count per large square} \times \text{dilution factor} \times 10^5 \text{ (Volume factor)}$$



Grid pattern of DHC-M01