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MAINTENANCE  
&  
USERS MANUAL

MANTELLI  
RISCALDANTI  
MM-MA

MM-MA  
HEATING  
MANTLES

FALC INSTRUMENTS s.r.l. Treviglio (BG) Italy

## INDICE

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## INDEX

1.0 SECURITY.....	page 1
1.1 GENERAL.....	page 1
1.2 OPERATORS QUALIFICATION.....	page 1
1.3 RESIDUAL RISKS.....	page 1
2.0 PRELIMINARY OPERATIONS.....	page 1
2.1 CHECK FOR DAMAGE DURING TRANSPORT.....	page 1
2.2 IN THE EVENT OF DAMAGE.....	page 1
2.3 CLEANING THE INSTRUMENT.....	page 1
3.0 INSTRUMENTS CHARACTERISTICS.....	page 2
4.0 COMMANDS.....	page 3
5.0 INSTALLATION AND USE PROCEDURES.....	page 4
5.1 PHYSICAL CHARACTERISTICS.....	page 4
5.2 FUNCTIONING.....	page 4
6.0 INTENDED AND UNINTENDED USES.....	page 5
7.0 MAINTENANCE.....	page 6
7.1 INSTRUMENTS MAINTENANCE.....	page 6
7.1.1 ORDINARY MAINTENANCE.....	page 6
7.1.2 PROGRAMMED MAINTENANCE.....	page 6
7.2 CLEANING.....	page 6
7.2.1 CLEANING PRODUCTS.....	page 6
7.3 TRANSPORT.....	page 7
7.4 STORING.....	page 7
7.4.1 STORING CHARACTERISTICS.....	page 7
8.0 ELECTRICAL WIRING DIAGRAM.....	page 8
9.0 DEMOLITION.....	page 11
9.1 OPERATORS QUALIFICATIONS.....	page 11
9.2 DEACTIVATING AND DEMOLITION.....	page 11
9.3 DEACTIVATING AND DEMOLITION PROCEDURES.....	page 11
9.4 RESIDUAL RISKS AFTER DEMOLITION.....	page 11

This manual should be conserved for the whole technical life of the instrument so that it can be easily consulted in case of need.

Conserve the manual in a dry and reachable place for a rapid consultation. Do not stain the manual to avoid creating stains that could prevent the lecture of this text.

## 1.0 SECURITY

### 1.1 GENERAL

Criteria and notions have been used in the construction of this instrument that guarantee that the instrument complies with the basic safety requirements under Low Voltage Directive 2006/95/EC.

The accurate risk analysis carried out by the manufacturer means that the majority of risks connected with the use of the instrument have been eliminated.

Full documentation of the adopted safety measures can be found in the technical file held by the manufacturer.

The manufacturer recommends that the instructions, procedures and recommendations contained in this manual, together with the security regulations in the working area, be carried out scrupulously even for the use of the intended uses of protections.

*FALC INSTRUMENTS cannot be held liable for any eventual damage to property, personal injury or injury to pets arising from the use of the instrument by unqualified operators and/or from an improper use of the instrument.*

### 1.2 OPERATOR QUALIFICATIONS

This instrument must only be used by trained, qualified and authorised personnel fully familiar with the information contained in this manual.

Residual risks during use can be contact of the user's limbs with the heating element or with the heated substances. We advise to wear gloves or use appropriate tools to move the hot containers.

#### **ATTENTION**

The non use and/or tampering with the protections can cause risks to the instrument user or other persons exposed; therefore *FALC INSTRUMENTS cannot be held liable for any eventual damage to property personal injury or injury to pets arising from these procedures*

## 2.0 PRELIMINARY OPERATIONS

### 2.1 CHECK FOR DAMAGE DURING TRANSPORT

Carefully remove the packaging and visually check the instrument for damage. Any dents in the visual parts indicate that the instrument has been knocked during transport, which could compromise its normal operation. Check specially for the good conditions of the following parts:

Check specially for the good conditions of the following parts:

- Instruments case.
- Controls and pilot lights.
- Check that all accessible screws and nuts are tight.

### 2.2 IN THE EVENT OF DAMAGE

Damage due to transport must be attributed to the carrier and immediately notified to the manufacturer or his representative.

### 2.3 CLEANING THE INSTRUMENT

Remove any dust or external dirt accumulated during transport

### 3.0 INSTRUMENTS CHARACTERISTICS

Model		MM 100	MM 250	MM 500	MM 1000	MM 2000
Tank capacity	ml	100	250	500	1000	2000
Temperature	°C	+350				
Precision	°C	± 5				
Heating Power	W					
External Dimensions ØxH	mm					
Weight	Kg					
Item		625.0500.01	625.0500.02	625.0500.05	625.0500.10	625.0500.20

**ATTENTION** → The maximum indicated temperatures refer to the heating element.

Power V/Hz 230/50-60	Security class 0	Protection class IP 30
-------------------------	---------------------	---------------------------

The MM heaters are the result of a long experience. Steel external casing painted with anti acid epoxy paint. The semi rigid nests which adapt themselves to the flask and like so reduce the breakage of the same, the thermal insulation is made of ceramic fibre. The temperatures are controlled by means of a numbered scale thermo regulator. All heaters are doted with an ON/OFF green light switch and spy light which indicates the functioning of the same.

Model		MA 3	MA 4	MA 6	MA 10	MA 20	MA 50
Flask capacity	ml	3000	4000	6000	10000	20000	50000
Temperature	°C	+350					
Precision	°C	± 5					
Heating power	W	600	900	900	1500	2100	3000
External dimensions ØxH	mm	310 x 290	380 x 310	380 x 310	460 x 330	540 x 510	700 x 650
Weight	Kg	6,0	8,0	8,0	11,0	16,5	20,0
Item		625.0520.03	625.0520.04	625.0520.06	625.0520.10	625.0520.20	625.0520.50

**ATTENTION** → The maximum indicated temperatures refer to the heating element.

Power V/Hz 230/50-60	Security class 0	Protection class IP 44
-------------------------	---------------------	---------------------------

The MA heating mantles for large capacities are constructed in order to have a higher security in respect with the MM and are doted with a liquid recuperation system, in case of flask breakage. The external casing and the nest, coated with fibre glass cloth, are completely in STAINLESS STEEL. The resistances are placed underneath the nests and are insulated with rock wool in order to utilize the maximum heating effect and minimally reducing heat dispersion on the external casing. The temperature is controlled by a numbered scale thermo regulator. All heaters are supplied with a green ON/OFF switch and spy light which indicates functioning.

## 4.0 COMMANDS

Heating Element	Space dedicated to the housing of the flasks.
Luminous Spy Light	Indicates the functioning of the resistance
Thermo regulator	Serves to regulate temperature.
General IG Switch	Serves to switch on and off the instrument.
Drain tap	Serves to drain the liquid deposited in the heating element in case of flask breakage. ( Drain tap is available only on MA models).

## 5.0 INSTALLATION AND USE PROCEDURES

### 5.1 PHYSICAL CHARACTERISTICS

Position the instrument in a suitable for use (norm D.L. 626/94-c risk valuation).

Other than the overall dimensions, the instrument also requires the following conditions:

- In the immediate vicinity of the instrument you must provide power sources. The instrument has only one electrical power source placed on the back of the same. Check that the electrical distribution line is the same as the power required by the instrument and that it is earth grounded (for tension data please control the label placed on the back of the instrument).
- Provided a sufficient flat space for the normal use and maintenance of the instrument, and any peripheral equipment.

The instrument must be placed in a covered area, not exposed to the weather. Good lighting conditions for safe normal operation and maintenance of the instrument. The room where the instrument is used must be lighted sufficiently to allow for operation without risks due to areas in shadow. The accepted temperatures for normal functioning is of a minimum +5°C to a maximum of +40°C with maximum relative humidity of 80% for temperatures up to +31°C with linear reduction 50 % at +40°C.

#### ATTENTION

FALC INSTRUMENTS declines any responsibilities for inconveniences, damages or malfunctioning that should arise from the non respect of the above indications.

### 5.2 FUNCTIONING

#### 1) ELECTRICAL CONNECTIONS

The instrument has a single external electrical supply point placed on the back of the same.

#### 2) TRIAL START UP

Before regular and continuous use of the instrument, check that it works perfectly by conducting at least one full trial cycle .

#### 3) SWITCHING ON THE INSTRUMENT

- Plug in the cable's plug in the socket placed on the back of the instrument.
- Plug in the supply line in at the mains.
- Turn to position "1" the general IG switch.

#### 4) SET PHASE

- Make sure that the drain tap is closed (only for MA models).
- Turn clockwise the thermo-regulator knob and stop when in correspondence with the required position.
- Turning the knob anticlockwise will decrease the temperature.

MA POSITION	°C
1	55
2	90
3	160
4	200
5	320

MM POSITION	°C	MM POSITION	°C
1	50	6	240
2	90	7	275
3	125	8	315
4	160	9	350
5	200		

**IMPORTANT →** The above temperatures refer to the heating element's temperatures.

**IMPORTANT →** The above temperatures are only significant and are not to be taken absolutely as a value.

## 5) NORMAL STOP PROCEDURE

- Turn thermo regulator to position OFF.
- The luminous spy light, if on, turns off.
- Turn to position “0” the general IG switch.
- Remove flask from instrument.  
**ATTENTION** → movement of hot objects.
- Unplug the instrument from the mains.

## 6) LIQUID RECUPERATION PROCEDURE

*Operation practicable only for MA models in case of breakage of the heating element in the flask..*

- Set to OFF the thermo regulator knob.
- The luminous spy light if on, switches off.
- Turn to position “0” the general IG switch.
- Remove the large pieces of broken glass from the heating element.  
**ATTENTION** → movement of hot objects.
- Connect the one end part of a tube to the drain pipe and place the other end part in a container capable of containing all the liquid that was inside the flask before its breakage.
- Open drain tap.
- When the liquid is completely drained remove all pieces of glass.  
**ATTENTION** → movement of hot objects.
- Close drain tap and remove the tube.

### **ATTENTION**

FALC INSTRUMENTS does not assume any responsibilities for damage to property, personal injury or injury to pets arising from the non respect of the prescribed precautions in the DPI.

(See chapter “SECURITY” on page 1).

## 6.0 INTENDED AND UNINTENDED USES

### ☺ INTENDED USES

- 👍 The instruments are projected and manufactured for heating and maintaining the temperature of solutions.
- 👍 The instruments are projected and manufactured for treating substances normally used in laboratories.
- 👍 The instruments are projected and manufactured to be used in laboratories for professional, industrial and teaching purposes.
- 👍 The instruments are to be used in the biochemical, microbiological, haematological and diagnostics field.

### ☹ UNINTENDED USES

- 👎 All uses are unintended unless specifically indicated in *INTENDED USES*
- 👎 Using the instrument in environment conditions other than those indicated
- 👎 The use of the instrument in explosive atmospheres
- 👎 Electro-medical use
- 👎 Washing the instrument with water jets
- 👎 The use of acids and solvents in direct contact with the instrument

### **ATTENTION**

*FALC INSTRUMENTS does not assume any responsibilities for eventual damage to property, personal injury or injury to pets arising from not using the instrument as specified in the above mentioned uses.*

## 7.0 MAINTENANCE

### 7.1 INSTRUMENT MAINTENANCE

Carry out the maintenance operations with the instrument switched off and unplugged. Make sure that the instrument is cold before starting any type of operation.

Maintenance can be carried out in safety by expert personnel trained in preparation and maintenance of the instrument and authorised for the use, after having carefully read and fully understood the recommendations and instructions provided in this section.

#### 7.1.1 ORDINARY MAINTENANCE

By ordinary maintenance we mean all operations that can be carried out by the user. These consist in cleaning operations, periodic and preventive inspections which allow the use of the instrument in complete safety.

#### 7.1.2 PROGRAMMED MAINTENANCE

By programmed maintenance we mean all operations that can be carried out only by authorised personnel. These consist in periodic and preventive inspections of the instrument.

Nevertheless if the authorised and trained personnel should not be able to resolve a determined problem, contact the technical assistance.

### 7.2 CLEANING

#### ATTENTION

Carry out cleaning operations with the instrument switched off and unplugged from the mains.

At the end of each working day, thoroughly clean the instrument, specially the following parts:

- Plate and superior surface of the instrument.
- Frontal part of commands.
- Surface surrounding the instrument.
- The floor around the working area.

#### 7.2.1 CLEANING PRODUCTS

For cleaning operations it is sufficient to have a normal detergent and industrial water. In case of stubborn dirt use only isopropyl alcohol.

#### ATTENTION

Avoid the use of solvents that could damage the synthetic material or the paint work. In particular avoid the use of benzene, nitroperchlorite and trichloroethylene.

### 7.3 TRANSPORT

If the instrument is transported by opened motor vehicles, the packaged instrument and its fittings must be covered completely by a plastic tarpaulin or equivalent to protect against the weather.

### 7.4 STORING

The indications contained in this section must be observed during temporary storage.

- ◆ Installing the instrument non immediately after being supplied.
- ◆ Deactivating the instrument, its storing whilst waiting for relocation.

#### 7.4.1 STORING CHARACTERISTICS


- Up to 2000m above sea level.
- Accepted temperatures: from +5°C to +40°C.
- Maximum relative humidity 80% at temperatures up to 31°C with linear reduction to 50% at +40°C.
- Adequate weather proofing.
- Sufficient room to manoeuvre for safe and easy lifting and transport.
- Horizontal surface capable of bearing more than the weight of the instrument.
- NEVER PLACE OR REST ANY OBJECT ON THE INSTRUMENT.


## 8.0 ELECTRICAL WIRING DIAGRAMM

Model: MM


NR	COMPONENTS	NR	COMPONENTS
A	FUSES	D	CORAL RESISTANCES
B	RIVET FASTENING SWITCH	E	THERMO REGULATOR
C	TWO POLE SWITCH	F	LUMINOUS SPY LIGHT

NR	COMPONENTS	NR	COMPONENTS
A	FUSES	D	CORAL RESISTANCE
B	PANEL FUSE HOLDER	E	THERMO REGULATOR
C	TWO POLE SWITCH	F	LUMINOUS SPY LIGHT

NR	COMPONENTS	NR	COMPONENTS
A	FUSES	D	CORAL RESISTANCE
B	PANEL FUSE HOLDER	E	THERMO REGULATOR
C	TWO POLE SWITCH	F	LUMINOUS SPY LIGHT

## 9.0 DEMOLITION

### 9.1 OPERATORS QUALIFICATIONS

A qualified operator who has read and fully understood the contents of chapter “*SECURITY*” on page 1 and the contents of chapter “*DEMOLITION*” on page 11.

### 9.2 DEACTIVATING AND DEMOLITION

Once the instrument has reached the end of its technical operation life, the instrument has to be deactivated. Putting the same out of function and putting the instrument in a condition not to be used any more for the scopes which it had initially been projected and constructed, it should, however be possible for the raw materials of the same to be re-used.

It is necessary for a deactivating and demolition to be carried out in security of the instrument.

#### ATTENTION

*FALC INSTRUMENTS* cannot be held liable for personal injury or injury to pets arising from the re-use of the single parts of the instrument or installation situations different to the original ones.

### 9.3 DEACTIVATING AND DEMOLITION PROCEDURES

The deactivating and demolition procedures of the instrument have to be carried out only by fully qualified and trained personnel.

-  Switch-off the instrument.
-  Disconnect the power sources.
-  Dismantle the following parts:
  - Electrical parts.
  - Plastic Parts.
  - Resin Support Feet.
-  In case of movement of the instrument refer to chapter *INSTALLATION AND OPERATING PROCEDURE* page 4, *TRANSPORT* and *STORING* page 7.

### 9.4 RESIDUAL RISKS AFTER DEMOLITION

If the indications of section *DEACTIVATING AND DEMOLITION* page 11, are accurately carried out, there are no residual risks after demolition.

#### ATTENTION

The instrument has been constructed with non biodegradable material. Dispose of the instrument in an authorised deposit for the disposal.

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FALC INSTRUMENTS s.r.l.

#### SUBJECT: GUARANTY

For a period of 12 months from the date of delivery to the carrier. With the word "guaranty" we engage ourselves to substitute any and every faulty instrument, product components or part, due to faults or defect of our manufacturing process, not due to transport or misuse, as soon as we are informed and the item sent back to us for verification. The item should be sent back to us at your charge with a report explaining and describing the fault. When the problem has been verified and confirmed, the item shall be substituted or changed repaired and sent back to the customer free of charge. In a case the item has been tampered with, misused, or used not following the instructions, or in case wrong repair trials have been done on it, there shall be no guaranty whatsoever.

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The following unit complies with the safety requirements outlined by the EC DIRECTIVE concerning the guidelines for ELECTROMAGNETIC COMPATIBILITY (2004/108/EC) and LOW VOLTAGE REGULATIONS (2006/95/EC).

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The units are manufactured in compliance with the regulations

*CEI EN 61000-6-2*

*CEI EN 61000-6-4*

*CEI EN 61010-1*

*CEI EN 60355-1*

*CEI EN 61326-1*

*Il responsabile legale / legal responsible*

**Scotti Francesco**