**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**  
**操作指南**

## AM 4 Multi-positions Heating Magnetic Stirrer

F20500590

### RC 2 Heating Plates

F20700620, F20710620

#### General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息

Before using the unit, please read the following instruction manual carefully.

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Attention, surface chaude!

Prudencia, superficie caliente!

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.

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The product can be used with flammable liquids.

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**This unit must be used for indoor laboratory applications only.**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a way not specified by the manufacturer or with non-specified accessories, product's safety may be compromised.

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**This unit has been designed and manufactured in compliance with the following standards:**

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**IEC/EN 61010-1  
IEC/EN 61010-2-051  
IEC/EN 61010-2-010**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improve their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power cable provided with the instrument.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The vessel must be made of a suitable material to withstand the foreseen temperature.

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Magnetic field doesn't interfere with the function of cardiac pacemakers or data media if they are farther than 20 cm to the instrument.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 625°C and a flash point higher than 600°C.

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It is responsibility of the user to appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The AM 4 multi-position heating magnetic stirrer is suitable for heating and stirring liquids inside a suitable vessel placed on the aluminium plate. Stirring is possible by means of a magnetically driven stirring bar inserted into the vessel.

The RC 2, 2-position heating plate, is used to heat liquids inside a suitable vessel placed on the aluminium plate.

The instruments are used in laboratory for general purpose and for all those applications that require precise temperature control of the heating plate and, for AM 4 also the mixing of liquids.

The Aluminum alloy heating plate ensures:

- Optimal heat distribution and a high specific power thanks to the circular configuration
- Temperature homogeneity
- High resistance to thermal stress and thermal shock

For the AM 4, the powerful brushless motor is able to stir from 100 up to 1500 rpm. SpeedServo technology maintains the speed constant as viscosity changes.

**Note:** it's important to choose the most suitable magnetic stirring bar related to the type and the quantity of the liquid to be stirred as well as to the type of the beaker. Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. The stirring bar which satisfies most applications is code A00000356 (8 x 40 mm).

## 2. Installation

- Unpacking
  - Check the integrity of the unit after unpacking
- The box includes
  - AM 4 multi-position heating magnetic stirrer or RC 2, 2-position heating plate
  - Instruction manual
  - N°4 Magnetic stir bar Ø 8 x 40 mm (A00000356) (only for AM 4)
- Installation
  - Place the unit on a non-flammable surface
  - Ensure that the technical specification of the unit and of the power supply line are the same
  - Ensure that the electrical network is grounded
  - Ensure that the main switch is on position "0" (OFF)
  - Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power cable

**Note:** keep the power cable far away from the hot plate.

- Place the container with the liquid and the stir bar (only for AM 4) on the heating plate

## 3. Operation

- Power-up
  - Turn the unit on using the main switch
  - The right led between the two knobs lights up
- Stirring (only for AM 4)
  - Rotate the right knob to set the stirring speed, the led over the speed knob lights up
  - The real speed increases until setpoint achievement
  - A microprocessor ensures constant speed even if the viscosity changes (SpeedServo)
  - Turn the knob to the 0 positions to stop stirring
- Heating
  - Rotate the left knob to set the heating plate temperature
  - If the set temperature is higher than the actual one, the led over the temperature knob lights up
  - When the heating plate temperature exceeds 50 °C, the left led between the two knobs lights up until the temperature drops below 50 °C

**Note:** This warning is active even if the main switch is off but not if the instrument is not powered.

- Turn the knob to the 0 positions to stop heating

## 4. Error and warning messages

When the control panel shows one of the following error messages, the instrument functions automatically stop.

| Display indicator                          | Malfunction:  |
|--|---|
| Flashing Heating led and Hot plate led ON  | Heating plate overtemperature (T > 330°C)             |
| Flashing Heating led and Hot plate led OFF | Heating plate safety probe overtemperature (T >370°C) |
| Flashing Stirring led (only for AM 4)      | The stirring system doesn't run correctly             |

To remove the error message, restart the instrument from the main switch. For the overtemperature of the safety probe is necessary to disconnect the instrument from the power supply.

If the alarm persists on the control panel, please contact Velp Scientifica's technical assistance service.

## 5. Maintenance

No routine or extraordinary maintenance is required except periodic cleaning.

### Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

### Repair

Repairs must be carried out by authorized Velp personnel only.

The transport of the instrument by freight forwarders, couriers or others must be carried out using the original shockproof packaging. Follow any instructions on the original packaging (e.g. palletizing).

It is the responsibility of the user to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

## 6. Technical data

|                         |                                     | AM 4                                | RC 2                             |                |
|-------------------------|-------------------------------------|-------------------------------------|----------------------------------|----------------|
| <b>General features</b> | Model                               | F20500590                           | F20700620                        | F20710620      |
|                         | Voltage                             | 230 V<br>50/60 Hz                   | 230 V<br>50/60 Hz                | 115 V<br>60 Hz |
|                         | Power input                         | 2550W                               | 1200W                            |                |
|                         | Current consumption                 | 11 A                                | 5,2 A                            | 10.4 A         |
|                         | Dimensions (WxHxD)                  | 715x90x270 mm<br>(28.1x3.5x10.6 in) | 340x90x270<br>(13.4x3.5x10.6 in) |                |
|                         | Weight                              | 5.8 kg (12.8 lb)                    | 2.6 kg (5.8 lb)                  |                |
|                         | Construction material               | Aluminum                            |                                  |                |
|                         | Working in continuous               | Admitted                            |                                  |                |
|                         | Maximum load on the plate           | 20 kg (44 lb)                       |                                  |                |
|                         | Noisiness                           | << 80 dBa                           |                                  |                |
|                         | Environmental temp. admitted        | +5...+40 °C                         |                                  |                |
|                         | Storage temp. admitted              | -10...+60 °C                        |                                  |                |
|                         | Max humidity                        | 80%                                 |                                  |                |
|                         | Overvoltage category                | II                                  |                                  |                |
|                         | Pollution degree CEI EN61010-1      | 2                                   |                                  |                |
| Max altitude            | 4000 m                              |                                     |                                  |                |
| <b>Heating</b>          | Heating plate power output          | 600 W                               |                                  |                |
|                         | Heating plate dimensions            | Ø 135 mm                            |                                  |                |
|                         | Temperature range                   | 50 ÷ 310 °C                         |                                  |                |
|                         | Heat control type                   | Analogic                            |                                  |                |
|                         | Heating plate construction material | Aluminum                            |                                  |                |
|                         | Safety circuit                      | Separated with dedicated probe      |                                  |                |
|                         | Hot plate alarm                     | Over 50°C                           |                                  |                |
| Overtemperature alarm   | Over 330°C                          |                                     |                                  |                |
| <b>Stir</b>             | Stirring capacity                   | 15 l H <sub>2</sub> O x 4           |                                  |                |
|                         | Speed range                         | 100 ÷ 1500 rpm                      |                                  |                |
|                         | Stir control type                   | Analogic                            |                                  |                |
|                         | Motor type                          | BLDC                                |                                  |                |
| Motor rating output     | 10 W                                |                                     |                                  |                |

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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|                |   |
|----------------|---|
| EN 61010-1     | Safety requirements for electrical equipment for measurement, control, and laboratory use<br>Part 1: General requirements |
| EN 61010-2-010 | Particular requirements for laboratory equipment for the heating of material  |
| EN 61010-2-051 | Particular requirements for laboratory equipment for mixing and stirring  |
| EN 61326-1     | Electrical equipment for measurement, control and laboratory use - EMC requirements -<br>Part 1: General requirements     |

and satisfies the essential requirements of the following directives:

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|                        |  |
|------------------------|--|
| 2006/42/EC             | Machinery directive                                    |
| 2014/30/EU             | Electromagnetic compatibility directive                |
| 2015/863/EU (RoHS III) | Restriction of the use of certain hazardous substances |
| 2012/19/EU (WEEE)      | Waste of electric and electronic equipment             |

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung /  
加修改.

## 11. Declaration of conformity <sup>UK</sup><sub>CA</sub>

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

- S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016
- S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

- EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use  
Part 1: General requirements
- EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material
- EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring
- EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -  
Part 1: General requirements

and satisfies the essential requirements of regulations:

- S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008
- S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic  
Equipment Regulations 2012
- S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO 14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortexers and Shakers  
Dispersers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Open Circulating Baths  
Pumps

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