Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
操作指南

AREC.X Heating Magnetic Stirrer

F20500554, F20510554

AREC Heating Magnetic Stirrer

F20500551, F20510551

REC Heating Plate

F20700553, F20710553

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de instrucciones.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch.
使用前请仔细阅读以下使用说明书。

Caution, hot surface!
Attenzione, superficie calda!
Attention, surface chaude!
Prudencia, superficie caliente!
Vorsicht, heiße Oberfläche!
注意：加热面高温

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2012/19/UE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2012/19/UE.
No tirar el equipo en los residuos urbanos, como exige la Directiva 2012/19/UE.
Dieses Gerät unterliegt der Richtlinie 2012/19/UE und darf nicht mit dem normalen Hausmüll entsorgt werden.
按照EEC指令2002/96/CE，不要将该设备作为城市垃圾处理。

The product can be used with flammable liquids.
Il prodotto può essere utilizzato con liquidi infiammabili.
Le produit peut être utilisé avec des liquides inflammables.
El producto puede utilizarse con líquidos inflamables.
Das Produkt kann mit brennbaren Flüssigkeiten verwendet werden.
本产品可与易燃液体混合使用。

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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**IEC/EN 61010-1
IEC/EN 61010-2-051
IEC/EN 61010-2-010**

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power supply cable provided with the instrument.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The vessel must be made of a suitable material to withstand the foreseen temperature.

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Magnetic field doesn't interfere with the function of cardiac pacemakers or data media if they are farther than 20 cm to the instrument.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use explosive and dangerous materials for which the equipment is not designed for. The stirrer must not be used in explosive atmospheres, in bain-marie and to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 935 °C and a flash point higher than 910 °C.

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The unit is fitted with two fuses (2xT5 A L 250 V (for 230V), 2xT8 A 250 V (for 115 V)), found in the socket on the back. To replace one or more disconnect the mains cable and, using a screwdriver, lift up the small cover on the fuse box.

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(2xT5 A L 250 V (for 230V), 2xT8 A 250 V (for 115 V)),

It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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DOMINIQUE DUTSCHER SAS

1. Introduction

The AREC.X / AREC / REC with ceramic plate is used to heat and mix (not available on REC) liquids inside a suitable container placed on the ceramic plate. The stirring (not available on REC) is possible thanks to a magnetic drive stir bar placed into the vessel.

The instrument is used in the laboratory for general use and for all those applications that require a precise regulation of the stirring speed (not available on REC) and of the heating plate temperature.

The ceramic heating plate allows to reach higher temperature than a traditional aluminium heating plate; furthermore, the ceramic is an inert and very hard material, so it offers great resistance to almost any type of chemical or mechanical aggression keeping the surface characteristics unaltered over time.

The heating plate works up to 550 °C.

The motor (not available on REC) is able to stir up to 1500 rpm.
The electronic board allows to maintain constant speed as viscosity changes.

The AREC.X can work in combination with the Pt100 external probe or the VTF external thermoregulator for precise and accurate control of the sample temperature.

Note: it's important to choose the most suitable magnetic stirring bar related to the type and quantity of the liquid to be stirred as well as to the type of the beaker. Size and shape of the magnetic stir bar determine the stirring efficiency at any given speed. The stirring bar which satisfies most applications is code A00000356 (8 x 40 mm).

2. Installation

- Unpacking
 - Check the integrity of the unit after unpacking

- The box includes
 - AREC.X / AREC / REC
 - Instruction manual
 - Magnetic stir bar (Not for REC)
 - Power supply cable

- Installation
 - Place the unit on a non-flammable surface
 - Make sure that the technical specification of the unit and the power supply line are the same
 - Make sure the electrical network is grounded
 - Make sure that the main switch is on position "0" (OFF)
 - Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power supply cable
 - **Note:** keep the power supply cable far away from the hot plate.
 - Place the container with the liquid and the stir bar (not for REC) on the on the stirring plate

3. Operation

Power-up

- Turn the unit on using the main switch
- The display shows the software version, the set starting mode (A or B) and the last temperature setpoint

If the set starting mode is A (see chap. 4), the heating and the stirring (not available on REC) functions are not active and the display shows OFF. Otherwise, the unit starts working at the last set temperature and speed (not available on REC) setpoints.

Stirring (not available on REC)

- Rotate the right knob to set the stirring speed
- The dedicated led is on when the stirring function is active

Heating

- Rotate the left knob to set the heating plate temperature
- The display shows the current heating plate temperature for 3 seconds and the setpoint for 1 second

Note: when the display shows the temperature setpoint, the dedicated led "Set" is on.

Note: if the heating function is not active and the heating plate temperature is more than 50 °C, the display shows the blinking message "Hot". This warning is not active if the instrument is not powered.

External probe connection (only for AREC.X)

- Turn off the instrument using the main switch
- Screw the threaded support rod into its seat on the back of the instrument (optional)
- Fasten the clamp onto the support rod (optional)
- Place the external temperature probe into the clamp (optional); then place it into the sample contained inside the flask
- Plug the probe into the dedicated socket on the back of the instrument
- Turn the unit on using the main switch

When the external probe is connected, the dedicated led "Probe" is on.

The heating function works as described above but the settable and showed temperature will be that read by the external probe.

Note: make sure that the set "heating plate temperature limit" (see chap. 4) is suitable for the liquid in use. The temperature of the heating plate will not exceed the set value. By decreasing the temperature limit, it's possible to minimize potential overshoot and the time for reaching the working temperature will be longer.

External thermoregulator connection (only for AREC.X)

- Turn off the instrument using the main switch
- Screw the threaded support rod into its seat on the back of the instrument
- Fasten the VTF thermoregulator onto the support rod
- Place the temperature probe into the sample contained inside the flask
- Plug the VTF connection cable into the dedicated socket on the back of the instrument
- Turn the unit on using the main switch

When the external thermoregulator is connected, the dedicated led "Probe" is on and the display shows "--".

Rotate the temperature control knob to the maximum and set the working temperature on the VTF.

Note: make sure that the set "heating plate temperature limit" (see chap. 4) is suitable for the liquid in use. The temperature of the heating plate will not exceed the set value. By decreasing the temperature limit, it's possible to minimize potential overshoot and the time for reaching the working temperature will be longer.

4. Setting mode

- At the power-up, when the display is showing the software version, turn the heating control knob from the minimum to the maximum twice to access the setting mode
- Rotate the heating control knob to browse through the available menus, listed below
- For AREC.X and AREC
 - Rotate the speed control knob from the minimum to the maximum and then to the minimum again to access the selected menu
 - The display shows the saved value for 3 seconds; then it shows the new value according to the position of the speed control knob
 - The new value is saved if the speed control knob is not touched for 4 seconds
The display will show "End"; then it shows the setting mode menus
- For REC
 - Do not touch the knob for 2 seconds to access the selected menu
 - The display shows the saved value for 3 seconds; then it shows the new value according to the position of the heating control knob
 - The new value is saved if the knob is not touched for 4 seconds;
The display will show "End"; then it shows the setting mode menus
- To exit the setting mode, restart the instrument

Text shown		Default value	Range	Menu's name and description	
CAL	P10	0.0 °C	-10 ÷ +10 °C 0,1 °C steps	<u>Probe calibration</u> It allows to align the external temperature probe reading to a reference thermometer.	Only for AREC.X
FS	tCJ	550	50 ÷ 550 °C 10 °C steps	<u>Heating plate temperature limit</u> It limits the maximum value of temperature setpoint for the heating plate.	
OPt	tyPE	b	A – b	<u>Starting mode</u> <ul style="list-style-type: none"> ➢ A: when the instrument is powered, the heating and the stirring (not available on REC) functions are not active and the display shows OFF ➢ b: when the instrument is powered, it starts working at the last temperature and speed (not available on REC) setpoints 	
OPt	P10	FAS	FIn - FAS	<u>Thermoregulation mode with external probe</u> <ul style="list-style-type: none"> ➢ FIn: good and optimized temperature control, minimized overshoot and oscillations, slow rise in temperature ➢ FAS: fast rise in temperature, increased overshoot 	Only for AREC.X
ALL	P10	On	On - OFF	<u>External probe safety alarm</u> It allows to enable / disable AL 6 and AL 7 error messages (see chap. 5).	Only for AREC.X

5. Error and warning messages

When the display shows an error message, the instrument's functions are stopped automatically.

AL1	Heating plate overtemperature	
AL2	Excessive heating time	
AL3	The stirring system doesn't run correctly	Not available on REC
AL4	External probe overtemperature or probe disconnected when the instrument is ON	Only for AREC.X if the external probe is connected
AL6	Slow temperature increase read by the external probe	Only for AREC.X if the external probe is connected and if ALL P10 is On
AL7	Fast temperature decrease read by the external probe	

To remove the error message, restart the instrument.

If the alarm persists on the display, please contact VELP Scientifica's technical service department.

6. Maintenance

No routine or extraordinary maintenance is necessary except periodically cleaning.

Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

Repair

Repairs must be carried out by authorized Velp personnel only.

The instrument must be transported using the original shockproof packaging when freight forwarders, couriers or other are taking care of the unit transport. Follow the indications on the original packaging (e.g. place the instrument on a pallet).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

7. Technical data

		AREC.X		AREC		REC	
General features	Model	F20500554	F20510554	F20500551	F20510551	F20700553	F20710553
	Voltage	230 V 50/60 Hz	115 V 60 Hz	230 V 50/60 Hz	115 V 60 Hz	230 V 50/60 Hz	115 V 60 Hz
	Power input	800 W		800 W		780 W	
	Current consumption	3.5 A	7 A	3.5 A	7 A	3.4 A	6.8 A
	Dimensions (WxHxD)	203 x 94 x 344 mm (8 x 3.7 x 13.5 in)					
	Weight	3.5 kg (7.7 lb)				2.7 kg (5.9 lb)	
	Construction material	Technopolymer					
	Working in continuous	Admitted					
	Maximum load on the plate	25 kg					
	Noisiness	<< 80 dBa					
	Environmental temperature admitted	+5...+40 °C					
	Storage temperature admitted	-10...+60 °C					
	Max humidity	80%					
	Overvoltage category	II					
	Pollution degree CEI EN61010-1	2					
	Max altitude	4000 m					
Heating	Heating plate power output	770 W					
	Heating plate dimensions	180 x 180 mm (7 x 7 in)					
	Temperature range	25 ÷ 550 °C (5 °C steps)					
	Temperature resolution	1 °C					
	Heating plate construction material	Ceramic					
External probe (only for AREC.X)	Type	Pt100 Class A – Ø 3mm		-		-	
	Temperature range	0 ÷ 250 °C (1 °C steps)		-		-	
	Temperature resolution	1 °C		-		-	
	Accuracy	± 1 °C *		-		-	
	Electrical data	5 VDC – 1 W (max)		-		-	
VTF (only for AREC.X)	Temperature range	- 10 ÷ 300 °C (1 °C steps)		-		-	
	Temperature resolution	0.2 °C		-		-	
	Accuracy	± 0.5 °C *		-		-	
Stir (not available on REC)	Stirring capacity	15 L H ₂ O				-	
	Speed range	0 ÷ 1500 rpm				-	
	Motor rating output	10 W				-	

* in the following conditions: 800 ml water in 1 liter glass beaker (diameter 105 mm), stirring bar 8 x 40 mm, 600 rpm, 50 °C.

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
11. Declaration of conformity ^{UK}_{CA}

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-010	Particular requirements for laboratory equipment for the heating of material
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements

and satisfies the essential requirements of the following directives:

2006/42/EC	Machinery directive
2014/30/EU	Electromagnetic compatibility directive
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (RAEE)	Waste of electric and electronic equipment

plus modifications.

It may therefore be presumed that the product is in conformity with the technical provisions of the following UK Statutory Instruments (and their amendments):

- The Electrical Equipment (Safety) Regulations 2016 (S.I. 2016/1101)
- The Electromagnetic Compatibility Regulations 2016 (S.I. 2016/1091)

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

DOMINIQUE DUTSCHER SAS

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO 14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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VELP Scientifica Srl