LP-70

UV protective glasses

INFORMATION FOR USERS SAFETY FACE-SHIELD

GENERAL INFORMATION

These Glasses have been designed for your safety and comfort.

It complies with the European Directives RG 425/2016, DI 2014/30/CE and DI 2014/35/CE et satisfies all the essential requirements detailed within standards NF EN ISO 16321-1 and EN 170.

Our UV Filter Glasses meet OSHA requirements and ANSI-Z87.1-1979

SPECIFIC INFORMATION

Our UV light filter glasses are made from polycarbonate with additives for effective protection against ultraviolet radiation. They offer the best protection against 'UV light' between 254 and 400 Nm. UV light causes damage to the eyes, particularly when taking photosensitising medication.

Use: for use with sources that emit predominantly UV light and where glare is not a significant factor. Use with low-pressure mercury vapour lamps, such as those used to stimulate fluorescence or 'black lights'.

MARKING

OSHA : Occupational Safety and Health Administration

ANSI-Z87.1-1979: practice for occupational and educational eye and face protection

UVPS : manufacturing

ADJUSTMENT

The GLASSES have a unique size

IMPORTANT: all materials used in the manufacture of the glasses have been specially selected so as not to cause an allergic reaction. The wearer must report any suspicion of allergy, even though such reaction is unlikely.

MAINTENANCE

The glasses should be washed in tepid, soapy water, then wiped dry with a soft and clean cloth. Do not use abrasive or chemical aggressive cleaning preparations as solvents and detergents can alter the mechanical properties of the face-shield.

DESINFECTATION : Ultraviolet Germicidal Lamp or Ultra-Sonic cleaning machine

STORAGE AND TRANSPORT


It is recommended that storage and transport of glasses should be in the manufacturer's original packing. No undue strain should be exerted on the packaging.

OBSOLESCENCE

The glasses must be replaced in these circumstances:

- . Glasses showing scratches, cracks, distortion or dull spots in centre zone
- . Glasses when broken or perforated

Decontamination, collection and elimination of waste for France only

The buyer ensures and finances the decontamination, the collection and the disposal of waste electrical and electronic equipment (WEEE) under the conditions provided in the Articles 21 and 22 of the Decree No. 2005-829 dated of 20 July 2005.

Improper disposal may be harmful to the environment and human health.

UVItec Ltd
Unit 3.05, St John's
Innovation Centre Cowley
Road
CAMBRIDGE CB4 OWS,
United Kingdom
Tel : 44 (0) 1223 421 270
Fax : 44 (0) 1223 421 271

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