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CertoClav MultiControl

12/18 L

Laborautoklav

(User manual English: page 23)

Inhalt

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
| Step | Time[d:h:m:s] | T[deg C] |
|------|--|----------|
| S01A | 00:00:00:00 | 89.5 |
| S02H | 00:00:04:42 | 120.1 |
| S03H | 00:00:04:44 | 120.1 |
| S04H | 00:00:04:46 | 120.2 |
| S05H | 00:00:04:48 | 120.2 |
| S06H | 00:00:04:50 | 120.2 |
| S07H | 00:00:04:52 | 120.3 |
| S08H | 00:00:04:54 | 120.4 |
| S09E | 00:00:04:56 | 120.5 |
| 40 | 14 Tage 121 Grad finished successfully | |

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(German version: page 1)

Content

| | |
|--|-----------|
| Introduction | 25 |
| Declaration of conformity MultiControl2 | 26 |
| Safety remarks | 27 |
| Technical Data | 29 |
| <i>Overview - Autoclave</i> | <i>28</i> |
| <i>Directives.....</i> | <i>28</i> |
| <i>Overview - Controller.....</i> | <i>29</i> |
| <i>Stickers.....</i> | <i>29</i> |
| Scope of Delivery..... | 31 |
| Setup | 32 |
| <i>Autoclave</i> | <i>31</i> |
| <i>Controller</i> | <i>31</i> |
| Operation | 32 |
| <i>Before first use</i> | <i>32</i> |
| <i>Loading the autoclave</i> | <i>32</i> |
| <i>Close and Locking</i> | <i>32</i> |
| <i>Starting a program</i> | <i>33</i> |
| <i>Venting.....</i> | <i>33</i> |
| <i>Sterilization / Autoclaving.....</i> | <i>33</i> |
| <i>End of cycle - Cooling down.....</i> | <i>33</i> |
| <i>Opening.....</i> | <i>34</i> |
| <i>Autoclave - cleaning and maintenance.....</i> | <i>34</i> |
| <i>Spare parts.....</i> | <i>35</i> |
| Controller MC 2 | 36 |
| <i>Menu Hierarchy.....</i> | <i>36</i> |
| <i>Screens.....</i> | <i>37</i> |
| <i>Cleaning and maintenance</i> | <i>40</i> |
| Error messages | 40 |
| Usage of protocol printer | 41 |
| Electric schematic | 42 |
| BRIEFOPERATIONGUIDE | 43 |

Introduction

Thank you for choosing CertoClav MultiControl Alu and for the trust you placed in us by purchasing this autoclave.

The MultiControl Alu autoclave is equipped with microprocessor-controlled technology and a high precision temperature sensor, and it speeds up work in the laboratory. With a free choice of time and temperature, the MultiControl Alu offers a versatile range of applications in sterilization, quality control and material testing.

The basic version of the PC Software allows visualization and storage of your process data. It is available as a free download at www.certoclav.com. The Professional version of the software (Art. 8200210) additionally enables you to create multi-level programs with varying time and temperature settings.

For documentation purposes, the data can be printed on your CertoClav protocol printer (Art. 1300187) via the printer interface.

The device was manufactured according to modern manufacturing methods and 100% tested by our quality control department. We grant a 2-year warranty, except for parts which are subject to wear and tear, e.g. seals, from the date of purchase. If you are running long term cycles, the seals may wear out faster than usually and may need to be replaced after 6 months already. Deposits and discoloration are a result of usage. No guarantee is given in the event of improper use, e.g. if the load contains corrosive material.

The CertoClav MultiControl Alu is a technical device and requires a certain amount of specific knowledge to operate and achieve reliable results. Please read the instruction manual carefully to learn about the usage and all options for operation.

Please mind the safety remarks.

If you have any questions, please do not hesitate to contact us.

By e-Mail: support@certoclav.com or telephone: +43 732 674 278-0

Your CertoClav Team

Declaration of conformity MultiControl

This is to certify, that the **CERTOCLAV** Product Range MultiControl complies with the following regulations of the European Community:

2014/35/EU Low Voltage Directive

Associated Standards:

EN 61010-2-040:2015 Safety requirements for electrical equipment for measurements, control and laboratory Part 2-040: Particular requirements for sterilizers and washer-disinfectors used to treat medical materials.

EN 61010-1:2010 Safety requirements for electrical equipment for measurements, control and laboratory use – Part 1: General requirements

The Certoclav pressure equipment is constructed according Low Voltage Directive. 2014/30/EU Electromagnetic compatibility

Associated Standards:

EN 61326-1:2013 Electrical Equipment for measurements, control and laboratory use - EMC Requirements-

Part 1: General requirements (provided replacement: EDIN EN 61326-1:2018)

EN 61000-3-2:2014 Electromagnetic Compatibility (EMC) — Part 3- 2: Limits —

Limits for harmonic current emissions (equipment input current equal to or less than 16 A per phase)

EN 61000-3-3:2013 Electromagnetic compatibility — Part 3- 3: Limits — Limits for voltage fluctuation and flicker on the public low-voltage system (equipment having an input current equal to or less 16 A per phase) and not subject to conditional connection.

2011/65/EU Restriction of the use of certain hazardous substances (RoHS) Directive

Associated Standards:

EN 50581:2012 Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances

Manufacturing quality:

CertoClav manufactures according to modern manufacturing methods using high-quality materials. Decades of experience and proven employees guarantee high product quality and maximum customer service. All devices are 100% tested for function and quality. CertoClav products are manufactured exclusively in Austria.

Manufacturer:

CERTOCLAV Sterilizer GmbH
A-4060 Leonding / Austria

Michael Dirix CEO

Leonding, January 2019

CertoClav Sterilizer GmbH
Peintner Strasse 10
A-4060 Leonding / Austria
www.certoclav.com support@certoclav.com

CERTOCLAV

Safety instructions

-  Only use in accordance with this instruction manual. Keep the manual in a place accessible to all users of the device.
-  Use the autoclave chamber only with the MultiControl 2 Controller. Do not plug in directly to a power supply. The chamber will heat in an uncontrolled way.
-  Never manipulate on safety devices. Operate only controller and chamber of the same series. Do not mix lid and chamber of different types of autoclaves.
-  Surfaces are hot while in operation (metal and plastic parts). Wear protective gloves.
-  Must not be installed in a space that is subject to the danger of explosions.
-  Do not install the controller close to high power electronics.
-  Do not install closer than 30cm to a CertoClav autoclave in order to protect the controller from heat.
-  Use only appropriate material for wall mounting
-  Do not open the case.
-  Install the controller only in a full regulating circle where the heat probe PT100 can gauge the activity of the heater attached.
-  Use only with rated voltage.
-  Only for use in laboratories. Operate only by trained personnel. Do not use in domestic settings.
-  Never operate without water.
-  Target water amount: 1 – 1,5 liters (min. 0,7L, max. 2,5L)
-  The CertoClav MultiControl is made of aluminum alloy. Aggressive media and saline solutions will cause corrosion.
-  Do not autoclave tightly sealed containers (especially glass bottles) since there is no pressure compensation.
-  Only use water to generate steam. We recommend distilled water. DO NOT use demineralized water from an ion exchanger, this will cause pitting. If only demineralized water is available, mix with tap water (approx. 1/3)!
-  Do not release steam in the direction of a person – danger of burning. We recommend connecting a temperature resistant hose (inner diameter 10mm) to release hot steam into a suitable container or sink. We offer this hose and container as accessory: Art. no. 8582016 – hose 1.0m
-  Do not open the autoclave before the acoustical signal shows the end of the cycle and the autoclave is cooled down to 80°C.

- ⚠ If the pressure is released quickly when autoclaving liquids, there is a possibility that the liquid will boil. This should be avoided we recommend cooling the autoclave to 70°C before opening the valve.
- ⚠ When running long time cycles danger of cooking dry and overheating. The autoclave could be damaged. Check regularly water level.
- ⚠ Do not transport when filled. Do not open with force.
- ⚠ When not in use, always remove the water from the chamber (to avoid deposits and pitting!)
- ⚠ This autoclave must only be maintained by CertoClav or authorized retailers. Always use original spare parts.

Autoclave – Safety devices

Secure Locking mechanism and steam release valve within the lid

Are providing safety against opening under pressure and will avoid a pressure build up, if the unit is not closed the correct way. By turning the red knob to “open”, steam will be released. As long as the unit is not securely locked, the steam release valve remains open and so prevents any build-up of pressure.

Pressure Control valve:

When pressure/temperature is exceeded by about 10% over the highest working pressure (2,7 bar), steam is released through a safety mechanism. The temperature sinks.

Safety valve:

At 3,6 to 4 bar a safety valve opens in the event the pressure control valve fails. This safety valve closes automatically after pressure dropped. In any case than the pressure control valve has to be changed, as it is defect.

Overpressure safety recess

Should all safety devices blocked at an excess pressure at around 4,7 bar the lid gasket would be squeezed through an overpressure safety recess at the rim of the chamber with a loud noise. If this occurs do not operate the unit anymore and return it to the dealer or the manufacturer

- ⚠ If this occurs do not operate the unit anymore and return it to the dealer or the manufacturer for examination and check. It is dangerous to work with such a unit.

Technical Data

Overview - Autoclave

| | |
|---|------------------------------------|
| Ambient temperature | -10°...+40°C |
| Volume of the autoclave chamber | 12 L or 18 L |
| Maximum charge | 9 kg with max. 3 L water inside |
| Sterilization cycle | as programmed |
| Venting | Flow process with fractionation |
| Max. Working pressure | 2,7 Bar |
| Test pressure | 4,1 Bar |
| Chamber designed according to | Low-voltage-directive 2014/35/EU |
| Voltage | 230V ~10% 50 Hz |
| Power | Controller: 5,5 W, Heating: 1900 W |
| Interference suppression classification | N |

Directives

- **Low-voltage-directive 2014/35/EU**

CertoClav MultiControl constructed according **Low-Voltage-Directive 2014/35/EU**.

Tested in accordance with:

EN61010-1:2001 (IEC 1010-1:2001) Safety requirements for electrical equipment for measurement, control and laboratory use – Part 1: General requirements

EN61010-2-040:2005 (IEC 1010-2-040:2005) Particular requirements for sterilizers and washer-disinfectors used to treat medical materials

- **Electromagnetic Compatibility Directive 2014/30/EU**

Associated standards:

EN 61326-1:2013 Electrical Equipment for measurements, control and laboratory use - EMC-Requirements Part 1: General requirements

EN 61000-3-2:2006 + A1:2009 + A2:2009 Electromagnetic Compatibility (EMC) — Part 3- 2: Limits - Limits for harmonic current emissions (Equipment input current equal to or less than <= 16A per phase)

EN 61000-3-3:2013 Electromagnetic Compatibility — Part 3- 3: Limits — Limits for voltage fluctuation and flicker on the public low-voltage system (equipment having an input current to or less <= 16 A per phase), and not subject to conditional connection.

EN 61010-1:2010

- **RoHS 2011/65/EU: Directive on the restriction of the use of certain hazardous substances**

Microbiological test certificate and air removal test by Österreichische Bundesstaatliche Bakteriologisch-Serologische Untersuchungsanstalt, Vienna (Federal Austrian Bacteriological and Serological Research Institute), Institute for Applied Hygiene, Graz (Institute of Applied Hygienic, Graz) and ANAWA, Munich.

CERTOCLAV MultiControl is **not** suitable for sterilization of tight sealed bottles and loads with aggressive media. For porous loads the **CERTOCLAV** MultiControl can't be recommended. Smaller amounts of textiles with low density can be sterilized. Liquids can only sterilized in small amounts. Large amounts will not reach the required temperature. Please contact CertoClav for special firmware versions if you want to sterilize liquids.

Overview - Controller

| | |
|------------------------------|----------------|
| Protection | IP 20 |
| Enclosure material | ABS |
| Power connection | 230V 50/60Hz |
| Heater connection | 230V 50/60Hz |
| Maximum output | 2000W |
| Control algorithm | PID |
| Temperature sensor interface | PT100, class A |

Stickers

| | |
|---|---|
|  | <p>Directive 2012/19/EU (WEEE): The supplier will take the supplied commodity back on the expense of customer after use completion and ensure best available treatment, recovery and recycling according to the laws duly</p> |
|  | <p>ATTENTION: Consult manual before use</p> |
|  | <p>ATTENTION: Hot surface! Do not touch!</p> |
|  | <p>CE marked according to directives NSR 2014/35/EC, EMC 2014/30EC and EN 61010-1:2010.</p> |

Scope of Delivery

After opening the packaging, please inspect the product for eventual transport damage. In case of damaged items contact your dealer and do not operate the unit for safety reasons.

The following parts are delivered. Please check them carefully. Should one of these parts be missing, please contact your dealer or call us directly.

The following parts are delivered with the autoclave, please check:

- Autoclave Chamber MultiControl Alu and lid
- Controller unit MC2
- Manometer with gasket and nut
- Steam release cock with venting nozzle, gasket and nut
- Instrument tray with water level indicator
- Cleaning needle for venting nozzle (at steam release cock)
- The basic version of the PC software - available on www.certoclav.com at or contact us support@certoclav.com .
- Mains cable
- Power connection cable from controller to autoclave
- Manual and warranty certificate

If something is missing please contact CertoClav Sterilizer GmbH.

Additional we offer the following accessories:

- Art. no. 8583010 Wire basket diam. 23, stainless steel
- Art. no. 8582010 Exhaust air filter (0,01 micron) complete with condensation bucket 5L and hose (1,0m)
- Art. no. 8582011 Exhaust air filter cartouche
- Art. no. 8582020 Condensation bucket (5L) and hose (1,0m)
- Art. no. 8582016 Hose (1,0m) for mounting to steam release cock (diam. 10mm)
- Art. no. 1250499 Service Set 1 with all gaskets
- Art. no. 8500482 Water distiller. We recommend using distilled water
- Art. no. 1300187 protocol printer with connection cable
- Art. no. 8500496 thermal paper for protocol printer
- Art. no. 8500210 PC Software professional version.


For further information about our product range visit: www.certoclav.com

Setup

Autoclave

- Set up is easy. Either put the autoclave on a table or the floor, so that you get a convenient work height.
- Unpack the manometer and the steam release cock. Mount them by putting them with gaskets in the holes and fix with nut from inside of the lid. These appliances have been packed separately to avoid damage in transportation.
- We recommend mounting a hose on the steam release valve, inner diameter 10 mm. (Art. no. 8582012)
- Plug in the connection cable to the socket on the back side. Connect the other end to the MultiControl 2 controller.

Controller

- 1) USB port for PC
- 2) Connection for CertoClav protocol printer
- 3) Port for a 4-wire PT100 temperature sensor
- 4) Port for connection to a heater or autoclave (220V, 10A max)
- 5) Mains connection (mains cable) with switch and fuse holder

1. Connect the controller (port 4) with a heater/autoclave using connection cable.
2. Connect the controller (port 5) with mains cable to a power socket.
3. Screw in the cable of the PT100 probe (port 3).
4. Optionally connect a printer (Port 2)
5. The controller can be mounted on the wall using the holes on the side of the enclosure.

Operation

Before first use

- READ INSTRUCTIONS BEFORE FIRST USE!
- The first run should be without load; only with water.

Loading the autoclave

- Fill in 1.5 liters of distilled water.
- It is possible to insert more water, yet this will increase heating up and cooling down time.
- In order to avoid calcification of your appliance, we recommend that you use distilled water and decalcify the appliance regularly.
- Put in instrument tray and place loading. It is important that steam can circulate freely, so you must have a distance between the loading and the chamber wall.

 Attention: Do not autoclave tightly sealed containers (especially glass bottles) since there is no pressure compensation.

Close and Locking

- Check if the gasket is inserted in the lid.
- Check if the orange turning knob is on position "OPEN".
- We recommend mounting a hose to the steam release cock to avoid danger of hot steam. (Art. 8582016)
- Place the lid on the chamber and turn it clockwise until the handles are in line with each other.
- Turn the orange knob to position "Locked"
- Plug in the power cable. Plug in the heater cable into the pot and controller. Plug in the temperature sensor cable and screw it into the controller.

 ATTENTION: on Longtime cycles avoid cooking dry and overheating (could damage the autoclave). Check the water level. The steam release cock has to be closed– remove the seal for this.

Starting a program

Select a sterilization program in the control unit and push start.

For detailed description about using the controller please read chapter "Controller MC2".

Venting

If venting is activated, then the sterilization chamber is vented at 115°C for 10 minutes. Venting is necessary for even temperature distribution in the chamber and optimal temperature transmission.

Sterilization / Autoclaving

For long time cycles close the cock. Remove the seal to enable this. This does NOT lead to a loss of warranty. The seal is only for process quality control in sterilization processes where the cock must not be closed.

End of cycle – Cooling down

As the last step the controller cools down the autoclave to 80°C or 70°C for liquid program. Thereafter an acoustic signal sounds and the controller displays "program finished". The autoclave can now be opened.

Opening

- Turn the orange knob to position “OPEN” to unlock. Steam release valve will be opened. Rest steam escapes.
- Turn lid anticlockwise and remove.
- To get better drying of the loading remove the basket from the chamber and let it cool down to room temperature before taking goods out.

 Should the cycle interrupted caused by switch off or electrical failure or other reason – the sterilization is not completed correct and has to be repeated. For opening please follow the safety remarks for cooling down.

Autoclave - cleaning and maintenance

For CLEANING the appliance wipe inside and outside with a damp cloth. Never immerse the lid and chamber of the unit in water.

Before start of cleaning or disinfection plug off mains cable and make sure that appliance has cooled down. For DISINFECTION of metal surface disinfection liquids with following substances are suitable: Ethanol, Propyl alcohol, Isopropyl alcohol, Glutaraldehyde, quaternary ammonium compounds. Please follow the application notes of the manufacturer for the cleaning and disinfection detergents.


Discoloration and stains occur when not using distilled water. These stains may block the pressure control valve. The steam release valve nozzle should be cleaned regularly with a needle. Plastic surfaces shall be cleaned with a dry cloth.

To remove discoloration, add some citric acid to the aqua and boil for about 10 minutes, do not close the autoclave. Empty the chamber immediately. We recommend working only with distilled water.

Repairs on valves, safety features and electric must only be carried out by an appointed CertoClav partner or by the manufacturer. Use original CertoClav spare parts only.

Every two years the unit should be checked by an appointed retailer or CertoClav.

Replacement of the lid gasket

1. Put out the old gasket and clean the groove.
2. Oil the groove and the gasket with cooking oil
3. Press gasket into groove ribbed side pointing downward (letters show to the middle of the lid) Do not stretch the gasket when inserting, as it will get too long. Rather try to squeeze it a little bit, it will seal better.

WE RECOMMEND THE **CERTOCLAV SERVICE SET ART. No.: 1250499**. It contains all spares which are under wear and tear and must be replaced.

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
Art. no. 1250499 Service Set 1: includes all gaskets and valves which are necessary for maintain the lid. Gaskets should be replaced yearly. **Long time cycles cause a higher wear and tear**; replacements could be necessary already after 6 months.


Controller MC 2

1. Down / minus button
2. OK / Enter button
3. Up / plus button

- Use the up/down buttons (1, 3) to select a point in the menu or to increase/decrease a value.
- Use the ok button (2) to enter a menu point or to confirm an edited value.

Menu Hierarchy

Screens

Main menu

```
Main Menu
>Programs
  Settings
```

Programs: Access to programs


Settings Access to settings

Programs

```
Programs
>Standard
  Prions
  Liquids
```

Shows all available programs

| | |
|-----------------|---|
| Standard | 121°C / 15 min. with 10 min venting (115°C), cooling down to 90°C |
| Prions | 134°C / 20 min. with 10 min venting (115°C), cooling down to 80°C |
| Liquids | 121°C / 32 min., cooling down to 80°C |
| Quick | 134°C / 4 min. with 10 min venting (115°C), cooling down to 80°C |
| Soft | 115°C / 32 min. with 10 min venting (115°C), cooling down to 80°C |
| Boil | 105°C / 32 min., cooling down to 80°C |
| Back | Back to Main Menu |

Process example. Temperature curve of the standard program

Program – recurring menu points in program menu:

```
Standard
>Start
  Info
  New Program
```

- **Start** Starts the selected program
- **Info** Displays program steps
- **New program** Overwrites existing program to create a new, single step program
- **Back** Back to Programs Menu

Settings

```
Settings
>Cycles : 0
  Version : 1.9.4
  Serial : MC2.0305
```

- **Version** Version of firmware
- **Serial** Serial no. of the controller
- **Contrast** Display contrast. Selectable from 0-5
- **Factory reset** Resets controller to factory setting, including programs
- **Language** Language selector: German / English
- **Back** Back to Main Menu

Install a new program

To create a new program, you have to overwrite one of the pre-installed programs! The new specific programs are shown in the menu with temperature and time (eg. 125deg16d18h08m). It is not possible to give the new program a specific name.

- Use the \wedge/V – buttons to go to ‚Programs‘ and confirm with OK.
- Choose with the \wedge/V – buttons the program you want to overwrite and confirm with OK.
- Choose with the \wedge/V – buttons ‚New program‘ and confirm with OK.

```
Program
  Start
  Info
  > New Program
```

- Vent: choose with the OK-button if you want venting YES/NO.
- Then, go on with the \wedge/\vee – buttons to temperature. Confirm with ‚OK‘. Then you can program the required temperature with the \wedge/\vee – buttons. To fix it, you have to confirm with OK.
- Now, go on with the \wedge/\vee – buttons to program the time (days/hours/minutes). With the OK-button you enter the submenu. With the \wedge/\vee – buttons you can select the respective value (days/hours/minutes) and confirm it with ‚OK‘.

```
New Program
  Vent.   :   No
  Temp.   :  121 degC
> Days   :   00
```

```
New Program
  Hours   :   00
  Minutes :   00
> -----
```

```
New Program
  -----
  Save
> Cancel
```

- Finally, you choose ‚Save‘ with the \wedge/\vee – buttons and confirm with ‚OK‘.
- You will automatically return to the main menu.

Status display while in operation:

```
Start
>Time    :   00:00:01
  Temp.   :   28.3 deg
  STOP
```

- **1 Achieve xxx.x deg** Number and description of current step
- **Time** Time elapsed since start (hh:mm:ss) or (dd:hh)
- **Temperature** Current temperature
- **Stop** Option to manually stop the running program

Cleaning and maintenance

If necessary, the controller can be wiped with a moist, not wet, cloth. Maintenance must only be done by the manufacturer or an authorized body. No regular maintenance necessary.

Error messages

Error messages are shown on the display as follows:

| Message | Possible reason | Troubleshooting |
|--|--|--|
| "Program stopped" | The program has been stopped manually by user. | Let the autoclave cool down before you start again. |
| "Power loss" | Program has been stopped because of power loss. | Ensure that the fuse holder has been assembled tightly and that the power cable is mounted correctly. |
| "Delta-T error" | The allowed temperature range of T-1/+2°C has been exceeded during holding time. | Check if the temperature cable has been plugged and screwed correctly into the controller. Check if the autoclave chamber is not out of water. |
| "Sensor error" | Error in temperature sensor. | Most often this means the sensor is not plugged in into controller. Other possible reason: cable break. |
| "Heater error" | Maximum allowed temperature exceeded (143°C) OR target temperature is not reached. | Please check if temperature sensor is working properly or autoclave chamber is not out of water. |
| "Profiler error" | | Please navigate to Settings -> Factory reset and press the OK button. This is necessary to solve software issues. |
| Display background light is turned off permanently | Fuses defective | Please check if the power switch is turned on. Please replace both 5mm 10A slow acting fuses placed at the shelf of the black power connector. |
| Other errors | | Navigate to Settings -> Factory reset and press the OK button. This is necessary to solve software issues. |

| Description | Possible reason | Troubleshooting |
|--|--|---|
| Display contrast is bad | | Navigate to Settings -> Contrast and press the OK button until the contrast is well. |
| Autoclave runs out of water during sterilization | Lid gasket defective Lid gasket insert wrong Autoclave not locked correct Steam release cock is open ... | Check if the gasket has been placed correctly into the lid. The text "Place This Side Facing Down" on the gasket must be readable for the user. Check if the red knob on the handle is in locking position during sterilization. For cycles with sterilization time longer than 60 minutes, close the steam release cock. |
| Water is dropping out of the gasket of the lid | On heat-up, the gasket of the lid will be compressed. If the gasket is wet also, this can cause a few drops coming out of the gasket during heat-up. | If the dropping doesn't stop after heat-up, please remove the gasket and clean it with distilled water. If the problem is not resolved, replace the lid gasket. |

If the errors cannot be resolved contact your dealer or the manufacturer
E-Mail: support@certoclav.com / Phone: +43 (0) 732 / 674 278

Usage of protocol printer

After connecting a CertoClav protocol printer (Art. No. 1300187) the following protocol is printed:

Header

10 Type of autoclave
15 Program name
30 Cycle number

Program steps


SxxA Heat up
SxxH Hold
SxxC Drying
SxxE End

Footer

40 Control message.
Sterilization successful?

Electric schematic

BRIEFOPERATIONGUIDE

MULTICONTROL Alu

Loading

1. Insert the CertoClav bottom plate into the sterilization chamber.
2. Fill the chamber with 1-1.5 liters of distilled water.
For CertoClav Essential, add 1/3 tap water
3. Place the items to be sterilized on the bottom plate.

Closing and locking

1. Check if the gasket is inserted in the lid.
2. Turn the orange knob to "OPEN".
3. Place the lid on the chamber and turn it clockwise until the handles are in line with each other.
4. Turn the orange knob to position "CLOSE" for locking the lid.
5. Connect a tube to the steam release cock, usual it is opened and sealed.
6. Plug in the power cable. Plug in the heater cable into the pot and controller. Plug in the temperature sensor cable and screw it into the controller.

Starting a program

Select a sterilization program in the control unit and push start.

Venting

If venting is activated, then the sterilization chamber is vented at 115°C for 10 minutes.

Sterilization

In case of sterilization times longer than 60 minutes, close the black steam release cock after the sterilization temperature has been reached. Turn knob clockwise.

Opening

Wait until the temperature of the Chamber is below 80°C. Turn red knob to OPEN and open the lid of autoclave. If liquids have been sterilized, wait until the chamber temperature drops below 70°C°. Otherwise there is the risk of burns from delayed boiling.

Sterilization timetable

If you want to choose a new program in the control unit, consider the following sterilization holding time for different volume of liquid media.

| Volume of loading liquid | Sterilization time until 121°C |
|---------------------------------|---------------------------------------|
| 3000 ml | 30 Minutes |
| 2000 ml | 25 Minutes |
| 1000 ml | 20 Minutes |
| 500 ml | 15 Minutes |

DOMINIQUE DUTSCHER SAS