MODE D'EMPLOI

Référence 432033

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Code	Status	Possible cause	Troubleshooting
E-1	The program was terminated manually by the user.	User pressed the "STOP" button during sterilization.	-
E-2	The communication between external touch screen and mainboard is interrupted.	The serial cable (COM 4) between touch screen and autoclave has loosened.	Check the cable connection between autoclave and touch PC. Otherwise restart the device and the touch PC. If the problem persists, navigate to Settings->Device-> Factory Reset to reset the touch PC to the factory settings.
E+2	The sterilization pressure was exceeded by more than 40kPa.	The vacuum pump is defective.	Check the vacuum pump and then perform a vacuum test.
E31	The temperature of the chamber is higher than 150 °C.	The steam generator sensor is defective.	Check, whether the temperature sensor (PT1000) in the chamber is properly inserted. If necessary, replace the temperature sensor.
E32	The temperature of the chamber heating is higher than 280 °C.	The steam generator sensor is defective.	Check, whether the temperature sensor (PT1000) of the chamber heater is properly inserted. If necessary, replace the temperature sensor.
E51	The temperature in the chamber is less than 0 °C.	The steam generator sensor is defective or the temperature in the chamber is too low.	Check, whether the chamber temperature is above 0 °C. Check that the temperature sensor (PT1000) in the chamber is properly inserted. If necessary, replace the temperature sensor.

Code	Status	Possible cause	Troubleshooting
E52	The temperature of the chamber heating is less than 0 °C.	The steam generator sensor is defective or the temperature in the chamber is too low.	Check, whether the chamber temperature is above 0 °C. Check that the temperature sensor (PT1000) in the chamber is properly inserted. If necessary, replace the temperature sensor.
E63	The temperature of the steam generator is less than 0 °C. The temperature of the steam generator is higher than 230 °C.	The steam generator temperature sensor is defective. Steam generator temperature is higher than 230 °C. Steam temperature control is unstable, causing steam generator temperature to become hotter than 230 °C.	Check, whether the steam generator temperature sensor is properly installed. Check steam generator temperature sensor (PT1000). If necessary, replace the temperature sensor.
E61	The temperature control is unstable. The control temperature inside the chamber was exceeded by 6 °C.	The chamber temperature sensor is defective.	Check, whether the temperature sensor (PT1000) is properly installed. If necessary, replace the temperature sensor.
E62	The temperature of the chamber heating is higher than 133 °C. The temperature control is unstable.	The chamber temperature sensor is damaged. The mainboard is damaged.	Check, whether the chamber temperature sensor is properly inserted. Check the temperature sensor (PT1000). Check the mainboard.
E41	During the preheating period, a chamber temperature of 100 °C was not reached after 8 minutes.	Chamber heating is damaged.	Check the chamber heating. Check the resistance of the heating element. Check the heating rod wiring.
E42	During the preheating period, a chamber temperature of 110 °C was not reached after 8 minutes.	Chamber heating is damaged.	Check the chamber heating. Check the resistance of the heating element. Check the heating rod wiring.
E5	During the venting phase, the pressure in the chamber is still greater than 0.3 bar after 10 minutes.	Solenoid valve is blocked.	Check the water drain valve. Clean or replace if necessary.

Code	Status	Possible cause	Troubleshooting
E52	The temperature of the chamber heating is less than 0 °C.	The steam generator sensor is defective or the temperature in the chamber is too low.	Check, whether the chamber temperature is above 0 °C. Check that the temperature sensor (PT1000) in the chamber is properly inserted. If necessary, replace the temperature sensor.
E63	The temperature of the steam generator is less than 0 °C. The temperature of the steam generator is higher than 230 °C.	The steam generator temperature sensor is defective. Steam generator temperature is higher than 230 °C. Steam temperature control is unstable, causing steam generator temperature to become hotter than 230 °C.	Check, whether the steam generator temperature sensor is properly installed. Check steam generator temperature sensor (PT1000). If necessary, replace the temperature sensor.
E61	The temperature control is unstable. The control temperature inside the chamber was exceeded by 6 °C.	The chamber temperature sensor is defective.	Check, whether the temperature sensor (PT1000) is properly installed. If necessary, replace the temperature sensor.
E62	The temperature of the chamber heating is higher than 155 °C. The temperature control is unstable.	The chamber temperature sensor is damaged. The mainboard is damaged.	Check, whether the chamber temperature sensor is properly inserted. Check the temperature sensor (PT1000). Check the mainboard.
E41	During the preheating period, a chamber temperature of 100 °C was not reached after 8 minutes.	Chamber heating is damaged.	Check the chamber heating. Check the resistance of the heating element. Check the heating rod wiring.
E42	During the preheating period, a chamber temperature of 110 °C was not reached after 8 minutes.	Chamber heating is damaged.	Check the chamber heating. Check the resistance of the heating element. Check the heating rod wiring.
E5	During the venting phase, the pressure in the chamber is still greater than 0.3 bar after 10 minutes.	Solenoid valve is blocked.	Check the water drain valve. Clean or replace if necessary.

Code	Status	Possible cause	Troubleshooting
E6	The door was opened during operation or the door sensor is damaged.	Door sensor is damaged.	Check the door sensor and micro-switch.
E7	The measured absolute atmospheric pressure is less than 70kPa.	Incorrect air pressure. The measured absolute atmospheric air pressure is too low for this location and must be corrected.	Pressure sensor calibration. See chapter 7.5.
E8	In the pre-vacuum phase, the temperature in the chamber rises by less than 3 °C within 5 minutes.	No water in the tank or air escapes from the chamber.	Check the water level in the fish water tank. Check the chamber heater and the mainboard.
E9	During the holding phase during sterilization, the measured pressure deviates by 0.3 bar from the intended pressure.	No water in the tank. Steam generator heating defective or water pump defective.	Check the steam generator heating, water pump and water level.
E10	The electronic interlock has the wrong end position.	Electronic interlock defective. Mainboard defective.	Check the electronic interlock and the mainboard.
E11	The electronic interlock has the wrong end position.	Electronic interlock defective. Mainboard defective.	Check the electronic interlock and the mainboard.
E12	The vacuum in the sterilization chamber does not reach less than -70kPa twice in a row with a program that has at least 3 vacuum phases.	Vacuum pump defective, fan defective.	Check the vacuum pump and fan.
E99	The communication between front display and mainboard is interrupted.	Ribbon cable not connected. CPU IC has loosened.	Check the flat ribbon cable and the connection plugs between front board and main board. Check whether the microcontroller IC on the mainboard is still firmly in the socket.

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