

VACUETTE® Luer Adapter



Intended use

The **VACUETTE®** Luer Adapter is a sterile, single-use device to be used with Safety/Blood Collection sets or Luer needles in routine blood collection. This device is to be used by trained healthcare professionals only in accordance with these instructions.

Product Description

The **VACUETTE®** Luer Adapter is manufactured from stainless steel and polypropylene and is fitted with a safety valve at one end. The perforated label not only serves to simplify identification, but also acts as a seal of integrity.


General Instructions for Use

- Grasp both ends of the luer case and carefully twist.
- Remove the short end of the cover.
- The perforation of the seal will break as the cap is twisted.

NOTE: If the paper label is broken, discard and use another luer adapter

- Thread the needle end of the luer adapter into a holder.
 - Make sure the Luer Adapter is securely seated.
- NOTE: It is recommended to use the **VACUETTE®** Luer Adapter with **VACUETTE®** Holders.*

- Remove the shield from the Luer Adapter and attach a Safety/Blood Collection Set or Luer Needle.
- Perform Venipuncture. (see Luer needle and/or Blood Collection Set Instruction for Use for more details)
- After venipuncture, the Holder + Luer + puncture device should be disposed of together as a unit in a sharps container according to your facility's protocol.



Storage of VACUETTE® Luer Adapter prior use

Recommended storage temperature: 4–36°C (40–97° F).









*NOTE: Avoid exposure to direct sunlight. Exceeding the maximum recommended storage temperature may lead to impairment of the **VACUETTE®** Luer Adapter quality.*

Precautions

- The handling of **VACUETTE®** Luer Adapter must be in accordance to the policies and procedures of your facility.
- HBV, HCV, HIV and other infectious diseases can be transmitted through contact with biological samples. In the event of exposure to biological samples, ensure that appropriate medical attention is obtained.
- Discard all **VACUETTE®** Luer Adapters in biohazard containers approved for their disposal.
- Do not use **VACUETTE®** Luer Adapters after their expiry date.
- Wear gloves during venipuncture and when handling blood collection tubes to minimise exposure hazard.

Label Information

References

	Item number		Temperature limit	ISO 6710 *Single – use containers for venous blood specimen collection * Clinical Laboratory Standards Institute (CLSI) H1-A4 Evacuated Tubes and Additives for Blood Specimen Collection- Fourth Edition; Approved Standard H3-A6 Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fourth Edition H2-A4 Methods for the Erythrocyte Sedimentation Rate (ESR) Test-Third Edition; Approved Standard H21-A3 Collection, Transport, and Processing of Blood Specimens for Coagulation Testing and General Performance of Coagulation Assays; Approved Guideline-Third Edition
	LOT number: Batch number		Do Not Reuse	
	Expiry Date. Use by the end of the month indicated		Consult Instructions For Use	
	Sterilization Using Ethylene Oxide		Manufacturer	