# Bedienungsanleitung Deutsch Operating instructions English Mode d'emploi Français

BOLA-Star-Exsikkatoren V 1871-07 / V 1873-07 / V 1875-07 / V 1879-07  
V 1899-07 / V 1910-07

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
## 2.0 Operating instructions

### BOLA-Star-Desiccators

Please read the information listed below thoroughly. We urgently ask you to observe the information about security, use and maintenance. BOLA-Desiccators are made of high-quality and transparent plastic material. Their task is to „dry“ humidity-sensitive products as well as optical and microelectronic parts securely, safe and in a dustfree atmosphere. This happens by adding SILICA GEL which considerably sinks the air humidity into the desiccator by its absorbing effect.

## 2.1 Supplied with

- 1 Star-Desiccator, transparent
- 4 Shelves made of PMMA
- 8 Rails for shelves made of aluminium
- 1 Hygrometer (fitted into the door)
- 1 Desiccant tray
- 1 Bag filled with 500 ml Silica Gel
- 1 Detailed operating instructions
- 1 Solid and ecologically friendly packing

## 2.2 Security instruction

Please do not take a damaged desiccator into operation. To avoid possible injuries, this desiccator should be checked for broken plates before starting to work.

## 2.3 Taking into operation

This desiccator is ready for operation as soon as the rails and shelves are installed and the tray filled with the desiccant is entered.

### 2.3.1 Desiccators with magnetic lock

It is recommended to hold the desiccator when opening (extremely high magnetic force which assures a safe closing).



Due to the strong magnetic lock, it is harder to open the door the longer it was closed before. This is normal and can be compensated by pulling stronger, preferably on the upper side of the grip.

Please consider that the magnet always has to be free of contamination, especially from metal parts.



### 2.3.2 Installation of rails and shelves

Before installing the rails and shelves, the door should be opened. The rails are hooked into the lateral mounted suspension arrangements. Please note that the flat side of the rails should sit close to the suspension arrangements. After a removal of the protecting foil, the shelves can be inserted into the rails.

### 2.3.3 Before entering the tray into the desiccator, it should be filled with the SILICA-GEL. Please note that it has got the shape of pearls and easily rolls out of the tray during filling it in.

## 2.4 Mode of function

### 2.4.1 Thermometer/Hygrometer

Limits: Temperature from  $-10^{\circ}\text{C}$  to  $+60^{\circ}\text{C}$   
Humidity from 10% to 99% R. H.

Max. error:  $\pm 1^{\circ}\text{C}$  within  $-10^{\circ}\text{C}$  to  $+50^{\circ}\text{C}$   
 $\pm 5\%$  R. H. within 25% to 75% R. H. and  $0^{\circ}\text{C}$  to  $50^{\circ}\text{C}$

Battery: 1,5 Volt AAA (for Super Star)  
LR44 or AG 14

Use: Open the battery compartment by pushing the cover on the back of the instrument downwards as indicated by the arrow. Remove the battery safety strip and the protection foil on the display. The instrument is now ready for use.

The temperature can be changed from  $^{\circ}\text{C}$  to  $^{\circ}\text{F}$  by the sliding switch at the back of the instrument (only for Super Star)

#### *Displaying maximum/minimum values:*

- Pressing the MAX/MIN button, the highest values since the unit was last re-set appear.
- Pressing the same button again gets you the lowest values since the unit was last re-set.
- To return to the current values, press the MAX/MIN button again.

#### *Resetting the MAX/MIN memory to the present values:*

Press the RESET/CLEAR button for 3 seconds

If the battery is nearly used up, a low battery flag is shown at the bottom-left corner of the lower display. Replace the battery to ensure the unit works normally.

*Waste disposal:*

The hygrometer as well as the packing are made of recyclable materials. Please put used batteries to a special waste site

#### 2.4.2 Drying with Silica Gel

The regulation of the air humidity into the desiccator is made by adding different quantities of SILICA-GEL (one pack is included in delivery, others are available as accessories). The saturation of the SILICA-GEL can be seen by a change of the colour (The gel can be reactivated into an external oven at approx. 100° to 180°C). The actual degree of atmospheric humidity is visible on the hygrometer fitted in the door. To prevent the air humidity getting into the desiccator when putting in or taking out parts, the door should be opened only for a very short period of time.

#### 2.4.3 Swapping with gas (only Super Star Desiccator V 1875-07)

The Super Star Desiccators allow to swap the stored products with gas. The corresponding bores to lead the gas in are already made. The tubing and the connecting pieces are supplied with the desiccators, the installation is easy.

The ports are protected by means of non-return valves, i. e. without a connection to the suitable counterpart there is not connection to the outside. This connection only exists after putting on the movable couplings.

**Function:** Cut the tubing with scissors or a knife to the requested length and press it on the movable couplings (tip: put some drops of water on the coupling before). Please **always** mount both couplings – otherwise the mixture of air and gas could burst the door open due to overpressure. We recommend - according to the weight of the gas in comparison with the air – to lead in the gas from above (for gas which is lighter than air) or from below (for gas which is heavier than air). Only this way an effective swapping can be effected. Please also consider that there must not be any tensile or compressive strength on the tubing since this could lead to a damage of the panels.

To disconnect the tubing, please dismount the couplings.

The non-return valves mounted on the desiccator and on the tubing keep the connections tight. Please consider that the sealing surfaces have to be stored carefully.

## 2.5 Service / Maintenance

2.5.1 Only use soft clothes and no detergents which contain solvents or which scrub.

2.5.2 Only use the desiccator in clean and dustfree conditions. Never expose it to direct insolation or other sources of heat.

- 2.5.3 After a long-time use the desiccant (SILICA-GEL) is exhausted – visible by a discoloration of the pearls – and can no longer absorb any humidity. A reactivation can be effected by putting the desiccant into an oven for a period of 60 to 120 minutes at approx. 100 to 180°C. A successful reactivation is visible when the colour changes from transparent to blue.
- 2.5.4 The function of the hygrometer is based on an expansion and contraction of a metal spiral at different humidity values. Therefore a maintenance of the hygrometer in regular intervals is recommended. This can be done by exposing it to air fully saturated with water. Just put a humid – not wet – cloth around the desiccator. After approx. 30 minutes the hygrometer will show a value of about 100 % atmospheric humidity. If this value is not indicated, the hygrometer is to be re-adjusted by means of a slotted screw driver on its back. Now remove the humid cloth. After another 60 to 180 minutes the hygrometer will again indicate the correct values of the atmospheric humidity.

## 2.6 Disturbances / Repairs

The following information to clear possible disturbances shall help you to recognize or to eliminate possible reasons for the disturbances by yourselves. Should there be a disturbance you cannot handle, please do not hesitate to contact us. Of course we will try to help you in this way. If then we do not find a solution to your problem, we ask you to return the cleaned desiccator into the original packing to us. Then we will check if a cost-effective repair is possible.

### Objection

Humidity into the desiccator does not sink

### Reason and Help

hygrometer damaged  
door is not closed tightly  
the object into the desiccator is too humid  
the SILICA-GEL is saturated (can be seen by a change of colour ) – dry or change it  
too less of SILICA-GEL in use  
the air does not circulate enough

## 2.7 Waste disposal

The desiccators as well as their packing are made of recyclable materials. Please take further information about the materials used from the following specifications.

## 2.8 Specifications

### 2.8.1 Desiccator

Cat.-No. :	V 1871-07/ V 1879-07
Material of plate :	PMMA (polymethylmethacrylate)
Temperature limits :	-20° to +70°C
Material of frame :	aluminium
Material of lock :	PP black, magnetic tape
Dimensions :	W x H x D, inner : 260 x 480 x 320 mm W x H x D, outer : 304 x 520 x 375 mm
Weight :	approx. 7 kg
Max. number of shelves :	26 pieces
Max. weight of shelves :	3 kg / total 12 kg

### 2.8.2 Desiccator


Cat.-No. :	V 1873-07
Material of plate :	PMMA (polymethylmethacrylate)
Temperature limits :	-20° to +70°C
Material of frame :	aluminium

- Material of lock : PP black, magnetic tape  
Dimensions : W x H x D, inner : 260 x 480 x 320 mm  
W x H x D, outer : 304 x 520 x 375 mm  
Weight : approx. 15 kg  
Max. number of shelves : 26 pieces  
Max. weight of shelves : 3 kg / total 12 kg
- 2.8.3 Desiccator  
Cat.-No. : V 1899-07  
Material of plate : PMMA (polymethylmethacrylate)  
Temperature limits : -20° to +70°C  
Material of frame : aluminium  
Material of lock : PP black, magnetic tape  
Dimensions : W x H x D, inner : 260 x 480 x 320 mm  
W x H x D, outer : 304 x 520 x 375 mm  
Weight : approx. 7 kg  
Max. number of shelves : 26 pieces  
Max. weight of shelves : 2 kg / total 12 kg
- 2.8.4 Desiccator  
Cat. No. : V 1910-07  
Material of plate : ESDV  
Temperature limits : -20° to +70°C  
Material of frame : aluminium  
Material of lock : PP black, magnetic tape  
Dimensions : W x H x D, inner : 260 x 480 x 320 mm  
W x H x D, outer : 304 x 520 x 375 mm  
Weight : approx. 7 kg  
Max. number of shelves : 26 pieces  
Max. weight of shelves : 2 kg / total 12 kg
- 2.8.5 Shelves including a pair of rails  
Cat.No. : V 1860-67  
Material of shelves : PMMA (polymethylmethacrylate)  
Material of guides : aluminium  
Temperature limits : -20° to +70°C  
Dimensions : 320 x 258 x 4 mm (L x W x H)  
Weight : approx. 0,370 kg  
max. load limit : 3 kg
- 2.8.6 Desiccant tray  
Cat. No. : V 1861-07  
Material : MF (melamine-formaldehyde)  
Temperature limits : -40° to +110°C  
Dimensions (outer) : 283 x 223 x 30 mm (L x W x H)  
Weight : approx. 0,170 kg
- 2.8.7 Hygrometer with fastening device  
Art.-No. : V 1863-07  
Type : analogue  
Material of fastening device: Aluminium  
Indicating range : 0 - 100 %, max. deviation  $\pm 5$  %  
Weight : approx. 0,025 kg

**2.9 Spare parts / Accessories**

- 2.9.1 1 replacement shelf including a pair of rails  
Art.-No. V 1860-07
- 2.9.2 Replacement tray  
Art.-No. V 1861-07
- 2.9.3 Replacement hygrometer with fastening device  
Art.-No. V 1863-07
- 2.9.4 Replacement desiccant (SILICA-GEL)  
Art.-No. V 1895-04 500 ml bottle  
Art.-No. V 1895-08 1000 ml bottle
- 2.9.5 Replacement shelf made of stainless steel for temperatures up to 150°C  
Art.-No. V 1860-51

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