

# FastGene® apTaq HotStart Polymerase

## Technical Data Sheet

### Product Description

The FastGene® apTaq HotStart Polymerase uses an aptamer based HotStart activation technology. The aptamer based HotStart inhibits any polymerase activity at temperatures below 45 °C and has an immediate activation.

### Product Applications

The FastGene® apTaq HotStart Polymerase is ideally suited for:

- Fast PCR
- Routine PCR
- PCR using complex templates
- SNP Analysis
- Any standard PCR application for which a hot start formulation of a high-quality thermostable DNA polymerase is required.

### Shipping and Storage

FastGene® apTaq HotStart Polymerase kits are shipped on ice packs. Upon arrival, store kit components at -20 °C, in a constant-temperature freezer. When stored under these conditions and handled correctly, full activity of the kit is retained until the expiry date indicated on the kit label.

### Handling

Always ensure that the product has been fully thawed and mixed before use. Reagents may be stored at 4 °C for short-term use (up to 1 month). Return to -20 °C for long-term storage. Please make sure to centrifuge briefly each tube prior usage.

### Kit Codes and Components

LS34s LS34	FastGene® apTaq HotStart Polymerase FastGene® apTaq HotStart Polymerase	Sample 500 Units
---------------	--	---------------------

### Related Products

LS20	FastGene® Taq DNA Polymerase	100 Units
LS21	FastGene® Taq DNA Polymerase	500 Units
LS22	FastGene® Taq DNA Polymerase	2000 Units
LS23	FastGene® HotStart TAQ DNA Polymerase	100 Units
LS24	FastGene® HotStart TAQ DNA Polymerase	250 Units
LS25	FastGene® HotStart TAQ DNA Polymerase	1000 Units
LS26	FastGene® TAQ Ready Mix PCR Kit	50 x 50µl rxns
LS27	FastGene® TAQ Ready Mix PCR Kit	250 x 50µl rxns
LS28	FastGene® Optima	250 Units
LS29	FastGene® Optima HotStart ReadyMix	500 x 25 µl rxns
LS30	FastGene® Optima	50 Units
LS31	FastGene® Optima HotStart ReadyMix	20 x 25 µl rxns

### Direct PCR

LS05	DNAreleasy Advance	10 preps
LS06	DNAreleasy Advance	50 preps

### Quick Notes

- Replace Taq for Standard PCR, with no protocol change.
- Buffer is supplied at 5x concentration with separate magnesium for optimization.
- Products can be cloned in T-overhang cloning vectors.

### Optimizing MgCl<sub>2</sub>

Aptamer-based HotStart polymerase and/or reactions with increased amounts of template (e.g., genomic DNA), primer and nucleotides might require higher MgCl<sub>2</sub> concentrations. A separate 50 mM MgCl<sub>2</sub> solution is supplied with the enzyme. Please adjust the MgCl<sub>2</sub> concentration according to the table below:

Final conc. in 50 µl reaction	2 mM	2.5 mM	3 mM	3.5 mM	4 Mm
Volume of 50mM MgCl <sub>2</sub> Stock	2 µl	2.5 µl	3 µl	3.5 µl	4 µl

# FastGene® apTaq HotStart Polymerase

## Technical Data Sheet

### FastGene® apTaq Protocol

FastGene® apTaq HotStart Polymerase can be used to replace any commercial Taq DNA polymerase in an existing protocol. To allow the most seamless integration of FastGene® apTaq HotStart Polymerase into existing protocols, be sure to match reaction conditions, particularly the MgCl<sub>2</sub>, primer and enzyme concentrations, as closely as possible.

#### Step 1: Prepare the PCR master mix

- Ensure that all reagents are properly thawed and mixed.
- Prepare a PCR master mix containing the appropriate volume of all reaction components common to all or a subset of the reactions to be performed.
- Calculate the required volumes of each component based on the following table:

Component	50 µl rxn <sup>1</sup>	Final conc.
PCR-grade water	Up to 50 µl	N/A
5 U/µl of FastGene® apTaq	0.5 µl <sup>2</sup>	2.5 U <sup>2</sup>
5X FastGene® apTaq Buffer	10 µl	1X
50 mM MgCl <sub>2</sub>	2 µl <sup>3</sup>	2 mM <sup>3</sup>
dNTP Mix (10 mM each)	1 µl	0.2 mM each
Forward Primer (10 µM)	2.5 µl	0.5 µM
Reverse Primer (10 µM)	2.5 µl	0.5 µM
Template DNA <sup>4</sup>	As required	As required

<sup>1</sup> Reaction volumes of 10 – 50 µl are recommended. For volumes smaller than 50 µl, scale reagents down proportionally.

<sup>2</sup> For small to midrange products, half the concentration of enzyme can be used. For GC-rich and other difficult templates, higher enzyme concentrations (up to 5 U per 50 µl reaction) may be required.

<sup>3</sup> A final MgCl<sub>2</sub> concentration of 2 mM is sufficient for most standard applications. For assays that do not perform well with 2 mM MgCl<sub>2</sub>, the optimal MgCl<sub>2</sub> concentration for each primer–template combination should be determined empirically.

<sup>4</sup> ≤100 ng for genomic DNA; ≤1 ng for less complex DNA (e.g. plasmid, lambda).

#### Step 2: Set up individual reactions

- Transfer the appropriate volume of PCR master mix, template and primer to individual PCR tubes/wells or a PCR plate.
- Cap or seal individual reactions, mix and centrifuge briefly.

#### Step 3: Run the PCR

Step	Temperature	Duration	Cycles
Initial denaturation	95 °C	1 min <sup>1</sup>	1
Denaturation	95 °C	30 sec	20-40 <sup>3</sup>
Annealing <sup>2</sup>	T <sub>m</sub> + 2 °C	15-30 sec	
Extension	72 °C	1 min/kb	
Final extension (optional) <sup>4</sup>	72 °C	1 min/kb	1
Store	4 – 10 °C	HOLD	1

<sup>1</sup> Initial denaturation for 1 min at 95 °C is recommended for most assays. For GC-rich targets (>65% GC), 5 min at 95 °C may be used.

<sup>2</sup> An annealing temperature 2 °C higher than the calculated melting temperature (T<sub>m</sub>) of the primer set is recommended as a first approach. If low yields and/or nonspecific amplification is obtained, an annealing temperature gradient PCR is recommended to determine the optimal annealing temperature of the primer pair (gradient should be at least between -2°C and +3°C).

<sup>3</sup> 35 cycles are sufficient for most assays. A higher number of cycles may be necessary for assays requiring higher sensitivity, while lower cycle numbers can be used if the template copy number is high.

<sup>4</sup> Final extension should be included if PCR products are to be cloned into TA cloning vectors.

For information on product use limitations and licenses:

<https://www.nippongenetics.eu/en/terms-and-conditions/>

For technical support please contact:

[info@nippongenetics.eu](mailto:info@nippongenetics.eu)