AwEasy

Short Instruction Manual

- Put the samples into the disposable sample cup and put the lid on it. The probe and the sample must stay in the same room for acclimation (normally room temperature)
- Avoiding warming up the probe and sample through body heat. Do not hold too long in your hands.

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1 GENERAL DESCRIPTION

The AwEasy is a standalone measurement head for water activity and temperature. It can be operated as standalone device together with the free AwEasy App for Android and iOS.

For more information please visit www.processsensing.com or www.rotrotronic.com or scan the QR code (also on the measurement head), which gives you the direct access to the Rotronic Online Manual with further explanations.


2 MEASUREMENT ARRANGEMENT

Before charging the AwEasy put the correct socket adapter on the power adapter (PS1006). Then connect the USB-A to USB-C cable with the USB-A plug to the power adapter and with the USB-C to the wireless charging plate (AwEasy-Charger). To disconnect the device from the power supply, disconnect the power adapter from the socket. **NOTE:** The power adapter must always be easily accessible.

WARNING: The power supply can only be used up to an ambient temperature of 45°C. See also the temperature limits for the battery and charging in the **TECHNICAL DATA**.

3 STARTING THE APP

Please go to the App Store (iOS) or the Play Store (Android) and download the AwEasy App. After the installation of the AwEasy App you can then start the application. By starting up the App for the first time, setup a user and register with your email address. The email address needs to be verified over a link send by a email. **NOTE:** Please make sure before using the device, that it's fully charged. For turning on or off the AwEasy press and hold the touch button on the top (below the display) for at least 10 seconds. To add the AwEasy measurement head to your device list, go into *Measurements* and click on the top right the "+*" to add a new AwEasy device. The App will automatically search for the AwEasy measurement head

As soon as the measurement head is connected it will be shown in the *Probes* overview with a green arrow on the left side. If there is no active connection with the measurement head, it will be shown as disconnected (with no green arrow).

4 STARTING THE MEASUREMENT

To start a measurement, click on the specific connected measurement head and click on the bottom of the measurement overview the *Add new measurement* button.

Add new measurement

After that the *Measurement Settings* menu will open. To start a new measurement, it needs to be filled out at least the *Measurement Name* and the *Batch Number*. *Content* and *Notes* are optional for the measurement. In the end the measurement method has to be chosen. Rotronic recommends for the **AWQ measurement a dwell time of at least 4 minutes and a temperature stability of 0.1°C/min**. For the standard **AW measurement** Rotronic recommends an **AW stability of 0.0001 aw/min** and a **temperature stability of 0.01°C**. Optional it's also possible to set an *AW Setpoint*. After setting the *Measurement Settings* it is then possible to start the measurement by clicking the *Start* button. As soon as the measurement starts a moving green line will appear and the *Start* button will change to *Stop* (to stop the current measurement).

NOTE: Do not use the AwEasy direct after charging. The charging generates heat in the battery, which can affect the measurement. Please wait at least 1 – 2 hours.

5 STANDALONE MODE

To use the standalone mode, it must be set up first over the AwEasy App. Go from the Main Menu into *Settings*. In the *Settings* menu the *File Name* and the *Batch Number* needs to be defined. In there you can select different types like *Date*, *Serial Number* etc. for the generated *Batch Number* or *File Name*. Most important here is that in the *Batch Number* the *Increment Number* is selected as type. Next go into *Measurement* and select the desired AwEasy in the probes overview. Set up the measurement as described in **STARTING MEASUREMENT** above and click on *Save*.

After that, the AwEasy measurement head can be used without a Smartphone. The measurement data from the AwEasy will be automatically loaded to Smartphone as soon as it is reconnected. To add a new measurement and to start or stop a measurement just double tap the touch button. The following icons will appear on the AwEasy display:

Add new measurement	Start measurement	Stop measurement
+	▶	■


6 GENERAL GUIDELINES

To make accurate and reproducible measurement, follow the following general guidelines to measure water activity:

- Measure water activity only in rooms with constant temperature; Ideal are rooms with small daily and yearly temperature variations. Avoid heating of the measuring equipment by sun radiation or radiators, draught, body radiation etc.


7 SAMPLE CUP

Put the sample into the sample cup. To avoid cross-contamination do not fill more than up to the fill line.

Put the probe immediately on the sample cup after having filled in the sample to avoid humidity exchange with the ambient. The measurement time may vary depends on the probe, humidity and temperature value are ready to read when the values are stable for at least 1 – 2 minutes.

8 DIMENSIONS

9 TECHNICAL DATA

AwEasy	
Sensor	HYGROMER WA-1
Accuracy @ 23 ±5 °C	± 0.008 aw / 0.8 %rh / ±0.1 K
Operating conditions	0...70°C / 0...100%rh
Storage conditions	-10...70°C / 0...100%rh
Ambient temperature for charging	0...40°C (charging does not start when the temperature exceeds the limits)
Battery temperature for charging	0...60°C (charging stops when the temperature exceeds the limits)
Long-term stability	< 1 %rh / year (normal conditions)
Parameters shown	aw, %rh, °C, °F
Status light	Yes
AW Quick Mode	Yes
Display	1.2" LCD color display
Button	Multifunctional touch button
Trend display	Yes
Clock	Real time clock provides a time stamp for every measurement
Batterie capacity	500 mAh
Power supply	5 V DC 1.5A 7.5 W max (AC adapter)
Power supply plugs	US / EU & UK
USB cable	USB-A to USB-C
Firmware update	Direct over the App
Protocols (values)	Date, time, aw / rh / °C / °F, aw setpoint, serial no., etc.
Operating frequencies	2.4...2.5GHz (Bluetooth LE) / 140.75-174.75kHz (Charger)
Radio-frequency power transmitted	max 0.5nW (by charger)
Acoustic alarm	No
Calculations	No psychrometric calculations available
Dimensions	Ø68 x 65 mm
Housing material	Aluminum, stainless steel & PC
Weight	600 g
CE conform	CE 2014/30/EU
IP protection rating	IP65

10 FCC & IC & MIC DECLARATION

This equipment (FCC: 2AG4)-AW EYIC: 30048-AWEY) has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

Per FCC 15.19(a)(3) and (a)(4) This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference that may cause undesired operation. The Federal Communications Commission (FCC) warns the users that changes or modifications to the unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Per RSS-Gen, Section 8.4 This device (Aw Easy Charger) complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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