


DDN 9

NEUTRAL DISINFECTANT DETERGENT FOR THE
PRE-DISINFECTION OF MEDICAL DEVICES



■ AREAS OF APPLICATION

- Operating room, endoscopy, ophthalmology and contactology

■ PRODUCT ADVANTAGES

- Proven efficacy on *Pseudomonas aeruginosa* biofilm
- Instrumentation integrity guaranteed for 72 hours in immersion-bath
- Compatible with a wide range of materials
- Can be used in ultrasonic cleaners
- Effective on enveloped viruses

Fb Franklab®
notre expertise l'Ultra-Propreté



0.5%

CE
0459



INSTRUCTIONS FOR USE

Preparation of the solution:

- Wear appropriate Personal Protective Equipment (refer to the Safety Data Sheet)
- Dilute the required volume of DDN 9 to 5 litres of mains water and homogenise the solution (see dilution table)

dosage for 5L	Concentration of the solution		
	0,5%	2%	5%
	25 ml	100 ml	250 ml

Use :

- Immerse the soiled medical devices completely in the solution
 - Brush the medical devices, and swab them if necessary
 - Contact time: according to the desired antimicrobial activity (see disinfectant properties)
 - Rinse thoroughly
- Renew the bath for each use (or follow the local recommendations in force)

DISINFECTING PROPERTIES

Activity	Standard	Concentration	Contact time
Bactericidal	EN 13727	0,5%	5 minutes
	EN 14561	0,5%	15 minutes
	EN 13697	0,5%	5 minutes
	EN 1276	0,25%	5 minutes
Yeasticidal	EN 13624	0,5%	5 minutes
	EN 14562	0,5%	15 minutes
	EN 1650	0,25%	5 minutes
	EN 13697	0,5%	5 minutes
Virucidal	EN 14476 PRV Virus (hepatitis B model)	0,5%	5 minutes
	EN 14476 DBV Virus (hepatitis C model)	0,5%	5 minutes
	EN 14476 HIV Virus	0,5%	15 minutes
	EN 14476 Herpes Virus	0,5%	5 minutes
	EN 14476 Vaccinia Virus	0,5%	15 minutes
	EN 17111 Vaccinia Virus	0,5%	15 minutes
Tuberculocidal	EN 14348	2%	60 minutes
		5%	15 minutes
	EN 14563	2%	60 minutes
		5%	15 minutes

PRECAUTIONS OF USE

See (SDS 010449)

FEATURES

- Form: **Liquid**
- Aspect : **Opalescent**
- pH at 0.5%: **7**
- Expiration date: **2 years**

COMPOSITION

- Non-ionic surfactants
- Sequestrants
- Quaternary ammonium

PACKAGING

- Box of 250 doses of 25 ml.....ref. **1044947**
- Box of 12 x1 L.....ref. **1044916**
- Box of 4 x5 L.....ref. **1044920**



Hazardous product - Follow precautions of use. Always read label and product information before use. Medical device Class IIb.