

AxyPrep™ Mag DyeClean

(Sequencing Dye-Terminator Removal System)

Summary

Removing excess dye terminator is an essential step prior to Sanger sequencing. Carryover of the excess dye into the sequencing reactions may result in dye blobs resulting in inaccurate results. The AxyPrep™ Mag DyeClean kit utilizes a unique paramagnetic bead-based purification system for Sanger sequencing reaction dye terminator clean-up. The protocol is simple and comprises binding, washing and elution steps which, if desired, can be performed directly in the thermal cycling plate. This kit requires no centrifugation or filtration steps making it more amenable to automation.

The kit is suitable for use with various dye terminator chemistries including Big Dye versions 1.1 and 3.1 and DYEnamic ET.

| Product Highlights |
|--|
| Simple clean-up process completed in 25 minutes/96 samples |
| Automation compatible |
| Streamline manual throughput using IMAG™ |
| No centrifugation or filtration steps required |
| Easily scalable from tubes, through 96 and 384 well microplates |
| High signal to noise ratio resulting in long reads and high QV20+ scores |

Overview of Process:

- Bind sequencing extension products to magnetic beads, and separate on magnet plate.
- Wash beads to remove unincorporated dyes, nucleotides, salts, and other contaminants.
- Elute DNA using aqueous buffer.

AxyPrep MAG DyeClean Kits

The volume of AxyPrep Mag DyeClean required per sequencing reaction depends on plate format **ONLY**. Please refer to the charts below to determine how many clean-ups each kit can provide.

| Reaction volume | 5 mL (MAG-DYCL-5) | 50 mL (MAG-DYECL-50) | 250 mL (MAG-DYECL-250) |
|-----------------|----------------------|-------------------------|---------------------------|
| 96 Well Format | 500 | 5000 | 25000 |
| 384 Well Format | 1000 | 10000 | 50000 |

Materials Supplied in the Kit

- ✓ AxyPrep Mag DyeClean_beads only

- Store at 4°C upon arrival, for up to 12 months.
- DO NOT FREEZE.

Materials to be provided by the User:

Consumables & Hardware:

| Name | Recommended Model | Recommended Vendor and P/N |
|-----------------------------|---|---|
| 96-well PCR reaction plate | 96-well round/ flat bottom microtiter plate. Plate selection depends on the PCR reaction volume | Corning, Inc., www.corning.com # 3797, 96 well round bottom # 3591, 96 well flat bottom # 3957, 0.5 mL v bottom 96 # 3365, 360 µL round 96 # 3364, 360 µL flat 96 # 3371, 96 clear pro |
| | 96-well cycling plate | Axygen, PCR-96-FS-C, PCR-96M2-HS-C, www.axxygen.com |
| 384-well PCR reaction plate | 384 well cycling plate | Axygen, PCR-384M2-C, www.axxygen.com |
| PCR Plate Seals | Easy Peel Heat Sealing Foil | Axygen, MF-111, www.axxygen.com |
| Liquid handling robotics | Compatible with open platform robotics | Contact Axygen Biosciences Technical support for compatible AxyPrep Mag methods and accessories to your automation |
| multichannel hand pipette | AxyPet | Single, 8 and 12 Multichannel |

Reagents to be supplied by the user:

- ✓ 85% Ethanol made from non-denatured ethanol
 - 25 mL of 85% ethanol per 96 well plate is required
 - For best results it is recommended that the 85% ethanol is prepared no more than 3 days prior to use and stored in a tightly capped container.
- ✓ Elution Buffer: Reagent grade water or 0.1 mM EDTA (pH 8.0) – dependent upon sequencing instrument is available
 - 0.1 mM EDTA (pH 8.0) is used to lower the signal in cases where the signal is too strong for certain sequencing instruments
 - Reagent grade water may be used to provide maximum signal
 - Please refer to Elution step table on page 5

IMAG™ Handheld Magnetic Separation Devices Selection Guide:

(Not provided in the Kit – please contact your Corning Axygen supplier for details on IMAG™)



The IMAG™ handheld Magnetic devices have been designed and optimized for different AxyPrep Mag protocols. These Magnets address different volumes for the tubes and plate types shown below.

Tube based:

| Protocol | Manufacturer | Part number | Plate description | Plate Material | Part Number |
|-------------------------|--------------|--------------|-------------------------------------|----------------|-----------------|
| AxyPrep Mag Kits | Axygen | SCT-050-SS-C | 0.5 mL Self Standing Screw cap tube | Polypropylene | IMAG-12T |
| | Axygen | SCT-150-SS-C | 1.5 mL Self Standing Screw cap tube | Polypropylene | |
| | Axygen | SCT-200-SS-C | 2.0 mL Self Standing Screw cap tube | Polypropylene | |

Plate based:

| Protocol | Manufacturer | Part number | Plate description | Plate Material | Part Number |
|-------------------------|--------------|---------------|---------------------|----------------|-----------------|
| AxyPrep Mag Kits | Corning | 3364 | 96 flat 360 µL | Polypropylene | IMAG-96P |
| | Corning | 3591 | 96 flat bottom | Polystyrene | |
| | Corning | 3365 | 96 round 360 µL | Polypropylene | |
| | Corning | 3371 | 96 clear pro round | Polypropylene | |
| | Corning | 3797 | 96 round bottom | Polystyrene | |
| | Corning | 3957 | 96 v bottom 0.5 mL | Polypropylene | |
| | Axygen | PCR-96-FS-C | 96 PCR full skirt | Polypropylene | |
| | Axygen | PCR-96M2-HS-C | 96 PCR half skirt | | |
| | Corning | 3959 | 96 round bottom 1mL | | |
| | Corning | 3961 | 96 round bottom 2mL | | |

Procedure for 96 Well Format:

1. PREPARATION STEP:

- a. Briefly vortex the AxyPrep Mag DyeClean paramagnetic beads to fully re-suspend the beads prior to use. For optimal results, the reagent should be at **room temperature** and appear homogenous and consistent in color before use.

2. BINDING STEP:

- a. Add 10 µL of AxyPrep Mag DyeClean to each sample. This is a fixed volume and is independent of the sequencing reaction volume.
- b. Add freshly prepared, 85% ethanol, to each sample according to the table below. Pipette mix 7 times or until the solution is homogenous. It is critical that the solution is homogenous for complete binding of the sequencing products to the magnetic beads.

| AxyPrep Mag DyeClean for 96 well plates | | |
|---|----------------------------------|----------------------------|
| Sequencing Reaction Volume (µL) | AxyPrep Mag DyeClean Volume (µL) | Volume of 85% Ethanol (µL) |
| 5 | 10 | 31 |
| 10 | 10 | 41 |

Note: For a given reaction, the volume of 85% Ethanol can be determined from the following equation:

$$\text{Volume of 85\% Ethanol} = 2.077 \times (\text{AxyPrep Mag DyeClean Volume (10 } \mu\text{L)} + \text{Sequencing Reaction Volume})$$

Observation: If the sample is not well mixed the ethanol will float to the top of the sample mixture, and the AxyPrep Mag DyeClean will sink to the bottom

- c. Place the sample plate onto a 96 well magnetic plate for 3 – 5 minutes or until solution is clear. The magnetic beads will form a tight pellet on the side of the well if using the IMAG™ or a crescent on the side of the well for other magnetic plates.

3. WASHING STEP:

All subsequent washes are performed while the plate is situated on the magnet.

- a. **Aspirate the cleared solution slowly from the beads and discard.**
The solution contains excess fluorescent dye and contaminants therefore it is critical to remove as much as possible without disturbing the beads.
- b. **Add 100 µL of 85% ethanol to each sample. Incubate for at least 30 seconds to allow the beads to resettle before continuing to the next step.**

It is not necessary to mix or re-suspend the beads during this step.

- c. **Completely remove the ethanol without disturbing the beads.**

d. Repeat steps b and c for a total of two 85% ethanol washes.

e. Let the samples air-dry for 10 minutes at room temperature.

The sample plate can be situated on or off the magnet while drying.

Note: Excessive drying can lead to degradation of the fluorescent dye bound to the DNA. The beads will appear cracked if over-dried.

4. ELUTION STEP:

- a. To elute remove the plate from the magnet, add 40 μ L of elution buffer (see chart below) and incubate the plate for 5 minutes at room temperature. Elution of the sequencing products is rapid therefore it is not necessary for the beads to go back into solution for complete recovery.

Notes: The suggested elution buffers are either reagent grade water or 0.1 mM EDTA (pH 8.0). Selecting the appropriate elution buffer will depend on the multiple variables including sensitivity of the sequencing detector, the amount of BigDye used per sequencing reaction and the type of template to be sequenced. Using water as an elution buffer will typically provide the maximum signal, while EDTA may be used to lower the signal. A lower signal may be desired if the fluorescent dye yields a signal intensity beyond the capabilities of detection for the instrument. Please use the following table as a general guideline for choosing an elution buffer.

| | ABI 3100 / 3130 | ABI 3700 | ABI 3730 |
|-------------------------------------|------------------------------|---------------------|---------------------|
| >2 μ L BigDye with PCR Products | 0.1 mM EDTA | 0.1 mM EDTA | 0.1 mM EDTA |
| <2 μ L BigDye with PCR Products | 0.1 mM EDTA | Di H ₂ O | 0.1 mM EDTA |
| >2 μ L BigDye with Plasmids | Plasmids 0.1 mM | Di H ₂ O | 0.1 mM EDTA |
| <2 μ L BigDye with Plasmids | Plasmids Di H ₂ O | Di H ₂ O | Di H ₂ O |

- b. Place the reaction plate back on to the magnet to separate the beads from the eluate. Incubate at room temperature for 3 – 5 minutes or until the solution becomes clear. The sequencing product is now in the eluate.
- c. Transfer 35 μ L of the clear eluate into a new plate for loading on the detector.
- d. Leave 5 μ L – 10 μ L of the solution in the plate to prevent transfer of beads into the detector. Residual beads may interfere with the sequencing instrument function.

Storage: Seal samples and store at 4 °C, for up to 24 hours, prior to loading. If samples will not be loaded within 24 hours, store at -20 °C. Samples are stable at -20 °C for at least 1 month.

Procedure for 384 Well Format:

1. PREPARATION STEP:

- a. Briefly vortex the AxyPrep Mag DyeClean paramagnetic beads to fully re-suspend the beads prior to use. For optimal results, the reagent should be at **room temperature** and appear homogenous and consistent in color before use.

2. BINDING STEP:

- a. Add 5 µL of AxyPrep Mag DyeClean to each sample. This is a fixed volume and is independent of the sequencing reaction volume.
- b. Add freshly prepared, 85% ethanol, to each sample according the table below. Pipette mix 7 times or until the solution is homogenous. It is critical that the solution is homogenous for complete binding of the sequencing products to the magnetic beads.

| AxyPrep Mag DyeClean for 384 well plates | | |
|--|----------------------------------|----------------------------|
| Sequencing Reaction Volume (µL) | AxyPrep Mag DyeClean Volume (µL) | Volume of 85% ethanol (µL) |
| 5 | 5 | 14.3 |
| 10 | 5 | 21.4 |

Note: For a given reaction, the volume of 85% Ethanol can be determined from the following equation:

$$\text{Volume of 85\% Ethanol} = 1.428 \times (5 \mu\text{L} + \text{Sequencing Reaction Volume})$$

Observation: If the sample is not well mixed the ethanol will float to the top of the sample mixture, and the AxyPrep Mag DyeClean will sink to the bottom

- c. Place the sample plate onto a magnetic plate for 3 – 5 minutes or until solution is clear. The magnetic beads will form a crescent on the side of the well.

5. WASHING STEP:

All subsequent washes are performed while the plate is situated on the magnet.

- a. **Aspirate the cleared solution slowly from the beads and discard.**
The solution contains excess fluorescent dye and contaminants therefore it is critical to remove as much as possible without disturbing the beads.

- b. **Add 30 µL of 85% ethanol to each sample. Incubate for at least 30 seconds to allow the beads to resettle before continuing to the next step.**

It is not necessary to mix or re-suspend the beads during this step.

- c. **Completely remove the ethanol without disturbing the beads.**

d. Repeat steps b and c for a total of two 85% ethanol washes.

e. Let the samples air-dry for 10 minutes at room temperature.

The sample plate can be situated on or off the magnet while drying.

Note: Excessive drying can lead to degradation of the fluorescent dye bound to the DNA. The beads will appear cracked if over-dried.

6. ELUTION STEP:

- a. To elute remove the plate from the magnet, add 15 - 30 μ L of elution buffer (see chart below) and incubate the plate for 5 minutes at room temperature. Elution of the sequencing products is rapid therefore it is not necessary for the beads to go back into solution for complete recovery.

Notes: The suggested elution buffers are either reagent grade water or 0.1 mM EDTA (pH 8.0). Selecting the appropriate elution buffer will depend on the multiple variables including sensitivity of the sequencing detector, the amount of BigDye used per sequencing reaction and the type of template to be sequenced. Using water as an elution buffer will typically provide the maximum signal, while EDTA may be used to lower the signal. A lower signal may be desired if the fluorescent dye yields a signal intensity beyond the capabilities of detection for the instrument. Please use the following table as a general guideline for choosing an elution buffer.

| | ABI 3100 / 3130 | ABI 3700 | ABI 3730 |
|-------------------------------------|------------------------------|---------------------|---------------------|
| >2 μ L BigDye with PCR Products | 0.1 mM EDTA | 0.1 mM EDTA | 0.1 mM EDTA |
| <2 μ L BigDye with PCR Products | 0.1 mM EDTA | Di H ₂ O | 0.1 mM EDTA |
| >2 μ L BigDye with Plasmids | Plasmids 0.1 mM | Di H ₂ O | 0.1 mM EDTA |
| <2 μ L BigDye with Plasmids | Plasmids Di H ₂ O | Di H ₂ O | Di H ₂ O |

- b. Place the reaction plate back on to the magnet to separate the beads from the eluate. Incubate at room temperature for 3 – 5 minutes or until the solution becomes clear. The sequencing product is now in the eluate.
- c. Transfer the clear eluate into a new plate for loading on the detector.
- d. Leave 2 – 5 μ L of the solution in the plate to prevent transfer of beads into the detector. Residual beads may interfere with the sequencing instrument function.

Storage: Seal samples and store at 4 °C, for up to 24 hours, prior to loading. If samples will not be loaded within 24 hours, store at -20 °C. Samples are stable at -20 °C for at least 1 month.

Please contact Axygen Biosciences for sales support at: axgsales@corning.com and for technical support at: axgsupport@corning.com

* The PCR process is covered by patents owned by Roche Molecular Systems, Inc., and F. Hoffman-La Roche, Ltd.

All trademarks are property of their respective owners.