**LANGE** 

The Lange logo is a red square containing a white stylized symbol that resembles a target or a crosshair.

DOC022.98.90250

**sensION™ + MM340**

Manual  
Bedienungsanleitung  
Manuel d'utilisateur

07/2010, Edition 1

07/2010, Heft 1

07/2010, Edition 1

ENGLISH

DEUTSCH

FRANÇAIS

CONTACT

## Index

|  |    |
|--|----|
| <b>1. Technical Specifications</b> .....           | 2  |
| <b>2. General information</b> .....                | 3  |
| 2.1. Safety information .....                      | 3  |
| 2.1.1. Use of hazard information .....             | 3  |
| 2.1.2. Precautionary labels .....                  | 3  |
| 2.2. General product information .....             | 3  |
| 2.3. Instrument composition .....                  | 3  |
| <b>3. Installation</b> .....                       | 4  |
| 3.1. Assemblage .....                              | 4  |
| 3.2. Connections .....                             | 5  |
| <b>4. Operation</b> .....                          | 5  |
| 4.1. Description .....                             | 5  |
| 4.2. Start-up .....                                | 6  |
| 4.3. Factory configuration .....                   | 6  |
| 4.4. General diagram .....                         | 7  |
| 4.5. pH calibration .....                          | 8  |
| 4.6. Other options for calibration .....           | 9  |
| 4.7. pH measurements .....                         | 10 |
| 4.8. Re-programming pH measurement .....           | 11 |
| 4.9. Printing data of pH measurements .....        | 12 |
| 4.10. mV measurements .....                        | 14 |
| 4.11. Ion selective measurements .....             | 15 |
| 4.12. Data Logger .....                            | 16 |
| 4.13. System .....                                 | 17 |
| 4.14. Temperature re-adjustment .....              | 18 |
| 4.15. User buffers .....                           | 18 |
| 4.16. Isopotential, pHiso .....                    | 19 |
| 4.17. Recognized buffers .....                     | 19 |
| <b>5. Sending data via the RS 232 C</b> .....      | 20 |
| <b>6. Maintenance</b> .....                        | 20 |
| 6.1. Clean the instrument .....                    | 20 |
| 6.2. Sensor cleaning .....                         | 20 |
| <b>7. Warning messages</b> .....                   | 21 |
| <b>8. Accessories and spares</b> .....             | 21 |
| <b>9. Warranty, liability and complaints</b> ..... | 21 |
| <b>Contact Information</b>                         |    |

# 1. Technical Specifications

---

## Measuring ranges:

pH -2.00 to 19.99  
mV  $\pm 2000$   
ISE  $10^{-6}$  M to  $10^{-1}$  M  
Temp -20 to 150 °C (-4 to 302 °F)

## Resolution

pH 0.1/0.01/0.001  
mV 0.1/1  
ISE programmable  
Temp 0.1 °C (0.1 °F)

## Measuring error ( $\pm 1$ digit)

pH  $\leq 0.002$   
mV  $\leq 0.2$   
Temp  $\leq 0.2$  °C ( $\leq 0.4$  °F)

## Reproducibility ( $\pm 1$ digit)

pH  $\pm 0.001$   
mV  $\pm 0.1$   
Temp  $\pm 0.1$  °C ( $\pm 0.1$  °F)

## Automatic temperature compensation

Manual entrance.  
With Pt 1000 temperature probe (A.T.C.).  
With NTC 10 K $\Omega$  probe.  
Isopotential pH programmable, standard value 7.00.

## pH calibration

With 1, 2 or 3 buffers selectable inside the range.  
Automatic recognition of technical and DIN 19266 buffers and up to 3 user's buffers.  
Calibration at any pH value (indirect calibration).  
Theoretical calibration.  
Manual introduction of the calibration parameters.  
Programmable calibration validity between 0 h and 7 days.  
Automatic recalibration warning.  
Automatic rejection of electrodes in poor conditions.

## pH calibration criteria

Asymmetry potential  $\pm 58$  mV.  
Slope 41 to 62 mV/pH, at 25 °C (77 °F) (sensitivity 70 to 105 %).

## mV calibration

Automatic recognition of redox standard 220 mV at 25 °C (77 °F).  
Calibration at any mV value.  
Manual introduction of the calibration parameters.

## ISE calibration

2 to 5 programmable standards.  
Selectable calibration units mol/l, M, g/l, mg/l, ppm, mmol/l, mM,  $\mu$ mol/l or %.  
Up to 4 calibration curves for each channel.

## Temperature readjustment

Correction of the temp. probe deviation (A.T.C.) at 25 °C (77 °F) and 85 °C (185 °F).

## Measuring mode

By stability (selectable stability criteria).  
Continuous.  
By time.

## Data Logger

Storage capacity up to 340 readings, by channel.  
Storage the last 9 pH calibrations.

## Languages

English, German, Spanish, French, Italian, Portuguese.

## Display

Graphic, backlit liquid crystal, 128 x 64 dots.

## Inputs and outputs

2 Combined or indicator electrodes, BNC connector (Imp.  $>10^{12}\Omega$ ).  
2 Reference electrodes, banana connector.  
1 A.T.C. type Pt 1000 (or NTC 10 K $\Omega$ ), banana or telephonic connector.  
2 Stirrers, RCA connector.  
RS232C for printer or PC, telephonic connector.  
External PC key board, mini DIN connector.

## Ambient conditions

Working temperature 5 to 40 °C (41 to 104 °F).  
Storage temperature -15 to 65 °C (5 to 149 °F).  
Relative humidity  $< 80\%$  (not condensed).

## Power supply

External plug-in power supply 100-240 VAC 0.4 A 47-63 Hz.

## Materials

Enclosure, ABS and PC. Keypad, PET with protective treatment.

## Physical parameters

Weight: 1100 g.  
Size: 350 x 200 x 110 mm.

## 2. General information

As a result of constant improvements to our products sometimes differences may exist between this manual and the instructions supplied with the instrument.

### 2.1. Safety information

Please, read carefully this information before installing and using the instrument !

Pay attention to all danger and caution statements.

#### 2.1.1. Use of hazard information

##### **DANGER**

Indicates a potentially or imminently hazardous situation that, if not avoided, will result in death or serious injury.

##### **WARNING**

Indicates a potentially or imminently hazardous situation that, if not avoided, may result in death or serious injury.

##### **CAUTION**

Indicates a potentially hazardous situation that, if not avoided, may result in minor or moderate injury.

**Important note:** Indicates a situation that, if not avoided, could lead to damage to the instrument. Important information that requires special emphasis.

**Note:** Information that supplements points in the main text.

#### 2.1.2. Precautionary labels

Read carefully all labels and tags attached to the instrument.

Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of life equipment to the Producer for disposal at no charge to the user.

**Note:** For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.

### 2.2. Description

The **sensION™** *MM340* is an instrument with 2 measurement channels to measure pH, mV or concentration with a selective electrode.

### 2.3. Instrument composition

| Version         | Electrode       | Accessories | Manual |
|-----------------|-----------------|-------------|--------|
| LPV2210.98.0002 | —               | ✓           | ✓      |
| LPV2210T.98.002 | LZW5010T.97.002 | ✓           | ✓      |
| LPV2214T.98.002 | LZW5014T.97.002 | ✓           | ✓      |
| LPV2211T.98.002 | LZW5011T.97.002 | ✓           | ✓      |
| LPV2221T.98.002 | LZW5021T.97.002 | ✓           | ✓      |

##### **Accessories:**

- Buffer solutions, 135 ml flasks.
- Magnetic stirrer.
- Calibration flasks, with magnetic bar inside.
- Electrolyte for the electrode, KCl 3M, 50 ml flask.
- Plug-in power supply.
- Electrode's stick and support.

## 3. Installation

### DANGER

Only qualified personnel, following the local security regulations, should conduct the tasks described in this section of the manual.

### 3.1. Assemblage

1. Unpack the instrument.

2. Verify that the shift is in the same position than in the picture.

3. Drive the stick into the hole.

4. Turn the shift in a counterclockwise until the stick is well fixed.


Place the o-ring at 135-140 mm from the base in order to avoid the shock of the sensor against the bottom of the glass.

5. Place the electrode's support on the stick by softly pressing the button.

Holder for 3 sensors.  
PN= LZW9321.99

6. Place the electrode on the electrode's support. Insert the electrode's cable through the groove.

#### Note: Stirrer replacement

Before replacing the stirrer, disconnect it from the rear panel and release the cable from the bottom of the instrument.

1. Turn the stirrer to release or fix from base.

2. Pull up the stirrer and replace it. Place the cable of the new stirrer through the bottom groove of the instrument.


## 3. Installation

### 3.2. Connections

#### DANGER

Make sure the power adapter matches your local power supply.

1. ST.1, Magnetic stirrer, Channel 1.
2. Combined pH electrode (or indicator), Channel 1.
3. Reference electrode (separated electrodes), Channel 1.
4. Temperature probe for the electrodes version T, Channel 1.
5. Separated temperature probe (telephone connector), Channel 1.
6. ST.2, Magnetic stirrer, Channel 2.
7. Combined pH electrode (or indicator), Channel 2.
8. Reference electrode (separated electrodes), Channel 2
9. PC keyboard, mini DIN connector.
10. RS-232 for printer or PC, telephonic connector.
11. Plug-in power supply (12 V).

## 4. Operation

### 4.1. Description

#### Keys

Switch on / Switch off.

Confirms the option displayed. When there is more than one option, accepts the one marked in negative.

Move backwards in the menus.

Selection between the different options displayed.  
Selection of numeric values.

Change between Channel 1 and Channel 2.  
Move forwards / backwards a digit when entering a numeric value.

Show the next display in those menus that have more than one.

## 4. Operation


### 4.2. Start-up

Connect to the rear panel of the instrument:


- the electrode to the corresponding channel.
- the temperature probe (A.T.C.) integrated in the pH electrode or as a separated sensor
- the magnetic stirrer supplied with the instrument to ST 1.
- the plug-in power supply

Follow the next steps:

For the first time...

Successive times...

**Note:**

If you have a 2nd. stirrer it must be connected to Channel 2 (ST2). Afterwards it must be configured, see System page 17.

### 4.3. Factory configuration

The **sensION™ + MM340** is an instrument with 2 measurement channels to measure pH, mV or concentration with a selective electrode. The factory configuration is as follows:

#### Channel 1

Units: pH  
Resolution 0.01 pH.  
Measurements: By stability with standard criterion.  
In continuous pressing twice the key .  
Calibration in 1, 2 or 3 points, with technical buffers.  
Calibration frequency, every 24 hours.

#### Channel 2

Units: ISE  
Measurements: By stability with standard criterion.  
In continuous pressing twice the key .  
Curve number 1.  
Five points calibration, units mol/l.  
Standards: 0.00001, 0.0001, 0.001, 0.01, 0.1.

#### Common to Channel 1 and 2

Automatic sample coding.  
Data sending, report type, etc. in accordance with the start up configuration.

#### Re-programming

To perform any modification in the default programming, for example measuring units of the channels, different measuring modes, calibration with other type of buffers, etc. see corresponding chapters.


#### Stirrers

In one standard instrument, Channels 1 and 2 will work with stirrer number 1 (ST1).  
If you have a 2nd stirrer it must be connected to Channel 2 (ST2).  
Afterwards it must be configured, see System page 17.  
The second stirrer will be assigned to Channel 2.

## 4. Operation


### 4.4. General diagram

This is a general view over the way to access to all possibilities of the instrument:

#### Example of report

|                                |                        |
|--------------------------------|------------------------|
| INSTRUMENT CONDITIONS          |                        |
| Monday 30 July 2010            | SN:70602 V1.0<br>15:50 |
| -----                          |                        |
| CHANNEL 1                      |                        |
| Units .....                    | pH                     |
| Resolution .....               | 0.01                   |
| Measurement .....              | By stability           |
| Criteria .....                 | standard               |
| Type of calibration .....      | Technical buffers      |
| Frequency of calibration ..... | every 24 h.            |
| CHANNEL 2                      |                        |
| Units .....                    | ISE                    |
| Measurement .....              | By stability           |
| Criteria .....                 | standard               |
| Curve number .....             | 1, F-                  |
| Units measurement .....        | mol/l                  |
| Standard 1 .....               | 0.00001                |
| Standard 2 .....               | 0.0001                 |
| Standard 3 .....               | 0.001                  |
| Standard 4 .....               | 0.01                   |
| Standard 5 .....               | 0.1                    |
| SYSTEM                         |                        |
| Stirrer number 2 .....         | non                    |
| Codification .....             | automatic              |
| Data send .....                | printer                |
| Type of printer .....          | thermal                |
| Type of report .....           | reduced                |
| Samples by report .....        | several                |
| Readjust temperature .....     | NO                     |

**Note:**  
The > < sign on the display indicates flashing.

#### Notes:

1. Channel 1 and Channel 2 can measure pH, mV or ISE with no difference. Measuring units can be selected in System, see page 17.
2. The system configuration, sample codification, data output, type of report, etc, are common parameters for the two channels.
3. Data Logger. The measurements will be stored in the Data Logger of the corresponding channel.

## 4. Operation

### 4.5. pH calibration

To measure pH correctly, the electrode and the instrument set need to be calibrated regularly with buffer solutions. This compensates deviations of asymmetry potential and slope of the electrode over time.

The calibration with buffer solutions corrects variations in the electrode response.

The instrument allows one-, two- and three- point calibration. The calibration parameters remain in memory until a new calibration is performed.

The standard calibration is performed with technical buffers and standard stability criterion. In addition, the instrument offers other options for calibration.

It is highly recommended to apply the same stability criteria to measure and calibrate.

#### One-point calibration

This type of calibration is acceptable when measuring pH values similar to those of the buffer solution.

It corrects the asymmetry potential of the electrode.

#### Two-point calibration

This is the most common calibration.

It is recommended to start with pH 7 as first buffer. pH 4 or pH 9 buffers can be used as a second buffer, depending on the measured zone (acid or alkaline).

The two-point calibration compensates the asymmetry potential and the electrode's sensitivity (slope).

#### Three-point calibration

This type of calibration is recommended when measurements are made across the whole pH scale.

As a first point pH 7 is recommended.

For the second and third points, two of the remaining values must be chosen.

The three-point calibration pH 4.00, 7.01 and 9.21 (at 25°C, 77°F) compensates the asymmetry potential and electrode's sensitivity in the acid and alkaline zone.

#### Calibration with technical buffers

Buffers pH 4.00, 7.01 and 9.21 (at 25°C, 77°C) are supplied together with the instrument.


#### Preparation

Fill the flasks with the corresponding buffer solution (40 ml).

These buffers can be used to perform more than 10 calibrations.

Every calibration flask is supplied with a magnetic bar.

Diagram for calibration

#### Important note:

The pH value of the buffer appearing on the display depends on the temperature of calibration. See buffers pH values as a function of the temperature on page 19.

## 4. Operation

### 4.5. Other options for calibration

#### Automatic buffers recognition

Technical buffers

pH 2.00, 4.01, 7.00, 9.21 and 10.00 at 25°C (77°F).

Buffers DIN 19266

pH 1.679, 4.006, 6.865, 7.000, 9.180, 10.012 and 12.454.

User's buffers

Suitable when the two above mentioned options will not be used. The **sensION™+ MM340** allows introducing the table pH/Temp. for the desired buffers. See System, page 17.

#### Calibration to X value

Manual re-adjustment of the measured pH at any value of the scale. The pH readjustment is equivalent to calibration with only one buffer.

#### Calibration data introduction

Manual introduction of calibration data for an electrode calibrated previously, for example in calibration laboratory.

#### Theoretical calibration

Selecting "theoretical calibration" the calibration data of the electrode in use are replaced by the theoretical values at 25°C (77°F).

This step is prior to the obtention of a calibration certificate.

#### Calibration frequency

Select between 0 h and 7 days. The instrument by default is programmed for daily calibration and indicates the remaining time to the new calibration.

Selecting 0 h the automatic re-calibration warning is deactivated.


#### Stability criteria.

Rapid criterion (variation < 0.02 pH in 6 s).

Standard criterion (variation < 0.01 pH in 6 s).

Strict criterion (variation < 0.002 pH in 6 s).

It is highly recommended to apply the same stability criteria to measure and calibrate.

**Important note:**  
Once the calibration is re-programmed, it will be performed with the new programming.

See Stability Criteria and View mV on page 10.

The instrument goes to MEASURE menu. The introduced calibration data will be used.

The instrument goes to MEASURE menu. The theoretical calibration will be used.

The instrument goes to the option "CALIBRATE pH" and it must be recalibrated.

## 4. Operation

### 4.7. pH measurements

#### Measurement by stability

This is the default measuring mode for the **sensION™ + MM340**. Meanwhile the reading is in evolution, it appears on the display momentarily to orient the user. When the electrode signal remains invariable during a certain period of time, the reading is fixed on the display, the reading is stable.

Criterion of "rapid" stability: 0.02 pH (1 mV) in 6 seconds.

Criterion of "standard" stability: 0.01 pH (0.5 mV) in 6 s.

Criterion of "strict" stability: 0.002 pH (0.1 mV) in 6 seconds.

If in a certain time the reading is not stable, the instrument automatically passes to continuous measurement (the measured value is displayed all the time).

To select another stability criterion or to change measuring parameters see "Measure re-programming", page 11.

**Magnetic stirrer:** the instrument controls automatically the start/stop of the stirrer.

#### Continuous measurement

The instrument displays the value obtained by the electrode all the time. These data can be stored or printed with programmable interval of time. To select this type of measurement see "Re-programming pH measurement".


From the standard measuring mode (by stability), pressing twice the key  the measured value is shown on the display at any moment but without option for data acquisition.

#### Measurement by time


The measurement is shown on the display after programmed time. To select this measuring mode see "Re-programming pH measurement".

#### Manual samples coding


With **sensION™ + MM340** it is possible to introduce a specific code (maximum 15 characters) per each sample using an external PC key board or bar code. See "System", page 17.

#### Manual introduction of the temperature

#### With connected temperature probe.


## 4. Operation

### 4.8. Re-programming pH measurement

|                        | Default programming (standard)  | Other options  |
|------------------------|---|--|
| <b>Resolution</b>      | <ul style="list-style-type: none"> <li>0.01 pH</li> </ul>   | <ul style="list-style-type: none"> <li>0.1 pH</li> <li>0.01 pH</li> <li>0.001 pH</li> </ul>  |
| <b>Measuring mode</b>  | <ul style="list-style-type: none"> <li>By stability.                             <ul style="list-style-type: none"> <li>Standard criterion (variation &lt; 0.01 pH in 6 s).</li> </ul> </li> <li>In continuous.                             <ul style="list-style-type: none"> <li>The result is stored when the measurement has finished.</li> <li>The result is printed when the measurement has finished.</li> </ul> </li> </ul> | <ul style="list-style-type: none"> <li>By stability.                             <ul style="list-style-type: none"> <li>Rapid criterion (variation &lt; 0.02 pH in 6 s).</li> <li>Standard criterion (variation &lt; 0.01 pH in 6 s).</li> <li>Strict criterion (variation &lt; 0.002 pH in 6 s).</li> </ul> </li> <li>In continuous. Selection of the time interval for:                             <ul style="list-style-type: none"> <li>Data storage.</li> <li>Data printing.</li> </ul> </li> <li>By time. Selection of the time interval for:                             <ul style="list-style-type: none"> <li>Data storage.</li> <li>Data printing.</li> </ul> </li> </ul> |
| <b>View mV</b>         | <ul style="list-style-type: none"> <li>No</li> </ul>  | <ul style="list-style-type: none"> <li>Simultaneous display of pH and mV concerned.</li> </ul>   |
| <b>Limits</b>          | <ul style="list-style-type: none"> <li>No</li> </ul>  | <ul style="list-style-type: none"> <li>Higher and lower limits. Acoustic warning when the measurement is out of limits, in the printed ticket appears "A" beside the measured value.</li> </ul>  |
| <b>pH isopotential</b> | <ul style="list-style-type: none"> <li>pHiso 7.00</li> </ul>  | pHiso 7.00 / Manual introduction / Value calculation.  |

#### "Re-programming" flow chart

#### Important note:

Once the measurement has been re-programmed, it will be performed in accordance with the new parameters.

It is highly recommended to apply the same stability criteria to measure and calibrate.


Data storage in the Data Logger. It coincides with the printing and PC sending. Selecting 0 store/print at the end of the measurement.

## 4. Operation

### 4.9. pH measurement with a connected printer

#### 4.9.1. Printer configuration

The printer can be configured at the start-up of the instrument or afterwards according to:


#### 4.9.2. Calibration

Once the calibration has finished the calibration data will be automatically printed.

The report can be "reduced", "standard" or "GLP" in accordance with the selected option in "System", see page 17.

##### Reduced report


Example of calibration reports:

##### Reduced report.

| REPORT of CALIBRATION        |             | CANAL 1  |         |
|------------------------------|-------------|----------|---------|
|                              |             | SN:70802 | V1.10   |
| Thursday, 14 February 2010   |             | 08:53    |         |
| Technical Buffers            | Slope mV/pH | Sens. %  | A.P. mV |
| 4.01, 7.00                   | -59.2       | 100.3    | -0.1    |
|                              |             | T. 'C    | t. s    |
|                              |             | 25.3     | 6       |
| Calibration frequency, 24 h. |             |          |         |

##### Standard and GLP report

Only if work with users was selected, see "System" page. 17.

##### Standard and GLP report.

| Laboratory Test                |             | REPORT of CALIBRATION      |         | CANAL 1 |      |
|--------------------------------|-------------|----------------------------|---------|---------|------|
|                                |             | SN:70802                   | V1.10   |         |      |
|                                |             | Thursday, 14 February 2010 |         | 08:53   |      |
| pH ELECTRODE                   |             |                            |         |         |      |
| Code:                          | 50 14T      |                            |         |         |      |
| Serial n.:                     | 7890        |                            |         |         |      |
| pHiso.(pH 7):                  | pH 7.00     |                            |         |         |      |
| CALIBRATED 14-02-2010 08:53:00 |             |                            |         |         |      |
| Technical Buffers              | Slope mV/pH | Sens. %                    | A.P. mV | T. 'C   | t. s |
| 4.01, 7.00                     | -59.2       | 100.3                      | -0.1    | 25.3    | 6    |
| Calibration frequency, 24 h.   |             |                            |         |         |      |
| User: Pedro Martinez           |             |                            |         |         |      |


## 4. Operation


### 4.9.3. Measurement

Once the measurement has finished, the corresponding report will be printed automatically.  
The report can be "reduced", "standard" or "GLP" in accordance with the selected option in "System", see page 17.

#### Reduced report (factory programming)

#### Standard and GLP report

#### Sample report (see System)

**Several:** The results from the different samples will be printed one by one until the option MEASURE is quitted.  
**One:** A complete report per sample.

#### Voluntary printing and report repetition (copies)

During measurement in continuous, pressing the key the measured value at that moment is printed. If the measurement is by stability pressing the key a copy of the printed report can be obtained.  
(If a PC keyboard is used, it must be pressed F1 key).

Example of reports:

#### Reduced report.

```

REPORT of MEASUREMENTS
                                SN:70802 V1.10
Thursday, 14 February 2010      08:53
-----
CONDITIONS                      CHANNEL 1
STABILITY MEASUREMENTS STANDARD CRITERIA
-----
Sample      pH      Ts'C    Time St%
-----
000001      4.01  23.7   00:08  25
000002      4.01  23.8   00:11  25
-----
    
```

#### Standard report.

```

Laboratory
Test
=====
REPORT of MEASUREMENTS
                                SN:70802 V1.10
Thursday, 14 February 2010      08:53
-----
pH ELECTRODE
Code:                               50 14T
Serial number:                       7890
pHiso.(pH 7):                         pH 7.00
CALIBRATED 14-02-2010 08:50:37
-----
CONDITIONS                      CHANNEL 1
STABILITY MEASUREMENTS STANDARD CRITERIA
-----
Sample      pH      Ts'C    Time St%
-----
000001      4.01  23.7   00:08  25
000002      4.01  23.8   00:11  25
-----
User: Pedro Martinez
    
```

#### GLP report.

```

Laboratory
Test
=====
REPORT of MEASUREMENTS
                                SN:70802 V1.10
Thursday, 14 February 2010      08:53
-----
pH ELECTRODE
Code:                               50 14T
Serial number:                       7890
pHiso.(pH 7):                         pH 7.00
CALIBRATED 14-02-2010 08:50:37
Technical      Sens. A.P. T. t. St.
Buffers        %      mV 'C s %
-----
4.01, 7.00    100.3 -0.1 25 6 25
-----
CONDITIONS                      CHANNEL 1
STABILITY MEASUREMENTS STANDARD CRITERIA
-----
Sample      pH      Ts'C    Time St%
-----
000001      4.01  23.7   00:08  25
000002      4.01  23.8   00:11  25
-----
User: Pedro Martinez
    
```

## 4. Operation

### 4.10. mV measurement


The instrument factory configuration is:

Channel 1, pH.

Channel 2, ISE.

Anyone of the 2 channels can measure pH, mV or ISE.

The measurement units are selected in System.


Such amended measures, it is necessary to save the changes.

#### Measurement conditions

**Resolution:** 1 mV / 0.1 mV

**Measuring mode:** By stability, In continuous or By time.

For more details, see "pH measurement", page 10.


#### mV measurement with printer

See pH measurement with printer, page 12.

#### Important note:

Measuring mV the temperature is not shown. The temperature is displayed only if a temperature probe is connected.

### mV calibration

## 4. Operation

### 4.11. Measurement with ion selective electrode

#### Curve number 1

**sensION™ + MM340** can store in memory up to 4 calibration curves for each channel.

#### Measurement by stability

It is the standard measurement mode of **sensION™ + MM340**. See details in “pH measurement”, page 10.

Criterion of “rapid” stability: variation < 0.5 mV in 6 s.

Criterion of “standard” stability: variation < 0.1 mV in 10 s.

Criterion of “strict” stability: variation < 0.05 mV in 10 s.

It is highly recommended to apply the same criteria to measure and calibrate.

#### Measurement in continuous or by time


See details in “pH measurement”, page 10.

**Note:** Anyone of the 2 channels can measure pH, mV or ISE. The measurement units are selected in System.


#### ISE measurement with printer

See pH measurement with printer in page 12.

#### ISE Calibration

#### Programming ISE calibration.

#### Reprogramming ISE measurement conditions.

#### Curve number

1 to 4: The **sensION™ + MM340** can store in memory up to 4 calibration curves for each channel.

**Stability criteria:** see above “Measurement by stability” It is highly recommended to apply the same criteria to measure and calibrate.

#### Units

mol/l, M, g/l, mg/l, ppm, mmol/l, mM, µmol/l, %: Concentration units of the calibration standards. The measurement will be performed in the same units..

#### Name of the ion

5 characters to relate the curve with the ion to measure .

#### Number of standards: 2 to 5


#### Standards concentration

It is recommended to introduce the standard concentration in increasing order, starting with the lower concentration. The instrument will dispose them automatically.

## 4. Operation

### 4.12. Data Logger

The **sensION™ + MM340** automatically memorizes the last 340 measurements by each Channel, the current calibrations, the sensor's history and the two Channels measuring conditions.  
All these data can be viewed at any moment on display, printed or sent to a PC.

### Report examples

| MEASUREMENT DATA    |       |      |       |     | Channel 1     |  |
|---------------------|-------|------|-------|-----|---------------|--|
|                     |       |      |       |     | SN:70602 V1.0 |  |
| Monday 30 July 2010 |       |      |       |     | 15:50         |  |
| Sample              | pH    | °C   | Time  | St% |               |  |
| 15-03-07 14:36      |       |      |       |     |               |  |
| 38510               | 10.15 | 23.5 | 00:15 | 30  |               |  |
| 15-03-07 14:40      |       |      |       |     |               |  |
| 38511               | 9.85  | 23.2 | 00:22 | 30  |               |  |
| 15-03-07 14:46      |       |      |       |     |               |  |
| 38512               | 10.23 | 22.8 | 00:18 | 30  |               |  |
| 16-03-07 08:36      |       |      |       |     |               |  |
| Pozo1               | 7.48  | 23.1 | 00:25 | 30  |               |  |

| pH CALIBRATION DATA                   |       |      |        |      |              |  | Channel 1     |  |
|---------------------------------------|-------|------|--------|------|--------------|--|---------------|--|
|                                       |       |      |        |      |              |  | SN:70602 V1.0 |  |
| Monday 30 July 2010                   |       |      |        |      |              |  | 15:50         |  |
| CALIBRATED 29-07-2010 15:42:51 22.5°C |       |      |        |      |              |  |               |  |
| Buffers                               | Slope | Sens | PAsym. | Time |              |  |               |  |
| Technical                             | mV/pH | %    | mV     | s    |              |  |               |  |
| 4.01,7.00                             | 58.16 | 99.0 | -1.5   | 18   |              |  |               |  |
| 7.00,9.21                             | 57.91 | 98.5 | -1.5   |      |              |  |               |  |
|                                       |       |      |        |      | Stirring 30% |  |               |  |
| Calibration frequency, every 24 h     |       |      |        |      |              |  |               |  |


| INSTRUMENT CONDITIONS    |       | SN:70602 V1.0     |  |
|--------------------------|-------|-------------------|--|
| Monday 30 July 2010      |       | 15:50             |  |
| CHANNEL 1                |       |                   |  |
| Units                    | ..... | pH                |  |
| Resolution               | ..... | 0.01              |  |
| Measurement              | ..... | By stability      |  |
| Criteria                 | ..... | standard          |  |
| Type of calibration      | ..... | Technical buffers |  |
| Frequency of calibration | ..... | every 24 h.       |  |
| CHANNEL 2                |       |                   |  |
| Units                    | ..... | ISE               |  |
| Measurement              | ..... | By stability      |  |
| Criteria                 | ..... | standard          |  |
| Curve number             | ..... | 1, F              |  |
| Units measurement        | ..... | mol/l             |  |
| Standard 1               | ..... | 0.00001           |  |
| Standard 2               | ..... | 0.0001            |  |
| Standard 3               | ..... | 0.001             |  |
| Standard 4               | ..... | 0.01              |  |
| Standard 5               | ..... | 0.1               |  |
| SYSTEM                   |       |                   |  |
| Stirrer number 2         | ..... | non               |  |
| Codification             | ..... | automatic         |  |
| Data send                | ..... | printer           |  |
| Type of printer          | ..... | thermal           |  |
| Type of report           | ..... | reduced           |  |
| Samples by report        | ..... | several           |  |
| Readjust temperature     | ..... | NO                |  |

| ELECTRODE HISTORY       |               |               | Channel 1         |  |
|-------------------------|---------------|---------------|-------------------|--|
|                         |               |               | SN:70602 V1.0     |  |
| Monday 30 July 2010     |               |               | 15:50             |  |
| Installed:              |               |               | 01-03-07 15:39:32 |  |
| In service:             |               |               | 06:30             |  |
| Number of measurements: |               |               | 216               |  |
|                         | Maximum value | Minimum value |                   |  |
| pH                      | 13.75         | 1.15          |                   |  |
| mV                      | --            | --            |                   |  |
| T (°C)                  | 38.2          | 14.3          |                   |  |

Standard and GLP report: If this type of report is selected, there will be data about the electrode, report header and user's name in the reports with calibration data and electrode's history.

## 4. Operation

### 4.13. System

#### Measurement units. pH, mV or ISE.

Measurement units of channels.  
If the measurement of Channel 2 is pH, it can be used the temperature of Channel one in the temperature compensation.

#### Sample ID:

**Automatic:** A consecutive number per each sample.  
**Manual:** The user introduces by an external PC keyboard a specific code per each sample before starting the measurement.

**Data output:** Select "deactivated" if there is neither a printer nor a connected PC. If a "printer" is selected the following options will appear:

#### Printer data output, options:

#### Only with Standard and GLP report

**Users:** The name will appear on the printed report (17 characters), up to 10 users.

**Report header.** Two lines of 40 characters for the company name, etc, that will appear in the printed report.

**Sensor identification:** Model and serial number of the sensor used in the 2 channels. If the measurement unit is ISE it is possible to introduce up to 4 different sensors, associated to the calibration curves. These data will appear on the printed report.

#### Stirrer 2.

When the second stirrer is connected automatically is assigned to Channel 2.

#### Temperature readjustment.


See page 18.

#### Stirling. Stirrer activation.

## 4. Operation

### 4.14. Temperature readjustment

The **sensION™ + MM340** allows the correction of the deviation presented by a temperature probe (integrated in the electrode or separated sensor) at 25° and 85°C (77 °F and 185 °F). In this way the pH-meter can be used as a precise thermometer. The temperature probe is connected to Channel 1.

### 4.15. User buffers

#### Application


When the buffers used in the calibration differ from the technical and the DIN 19266,

#### **sensION™ + MM340**

The instrument offers the possibility to program the table pH/Temp. for up to 3 buffers.

#### Table programming:

- Select "User Buffers" following the diagram that appears in "Other pH calibration options", page 9.

#### Important note:

- The temperature related to the pH values of buffer 2 and 3 will be assigned automatically in accordance with the values introduced for buffer 1.

## 4. Operation

### 4.16. Isopotential (pHiso)

#### Definition

Potential (mV) of one pH electrode that does not alter with the temperature. This is the intersection point of the different calibration lines at different calibration temperature. Normally, this value corresponds to pH 7. Nevertheless, in reality the value is slightly different.

#### Application

It is recommended to determine the pH of the real isopotential ( $\neq 7$ ) when after calibration at room temperature, the measurements will be performed at different temperature, and in addition a very high precision is required.

#### sensION™ + MM340

The instrument offers three possibilities for the application of pHiso:

**pH 7:** A value generally accepted by all the manufacturers.


**Introducing the real value:** If the pHiso is known, it can be introduced by the user.

**Value calculation,** by the instrument.

#### Procedure:

- Calibrate the instrument in 2 or 3 points with the buffers at room temperature, following the common procedure.
- Heat the buffers. Their temperature must differ with 10 or more degree centigrade from the calibration temperature.
- Select "Calculate" from the menu "Measuring conditions" (see pH measurement re-programming on page 11) and follow the instructions of the pH-meter.
- The instrument will calculate the pHiso of the electrode in use and will memorize and apply it in the corresponding calculation for the temperature compensation in the measurements.

#### Example for electrode behavior

- a) Electrode calibration line at room temperature
- b) A calibration line with pHiso=7 to apply when measure at 70°C (158 °F)
- c) Real calibration line of the electrode at 70°C (158 °F)
- P1) Classic isopotential, pHiso 7
- P2) Real isopotential of the electrode, pHiso 6.

**Error:** Applying the classic isopotential (pHiso7) the pH measurement in buffer 5 at 70°C (122 to 158 °F) will be 5.15.

### 4.17. Recognized buffers

#### Buffer solutions

The **sensION™ + MM340** has in memory the following data (pH values as a function of the temperature).

#### Technical buffers (DIN 19267)

| °C        | °F        | pH          |             |             |             |              | mV         |
|-----------|-----------|-------------|-------------|-------------|-------------|--------------|------------|
| 0         | 32        | 2.01        | 4.01        | 7.12        | 9.52        | 10.30        | --         |
| 10        | 50        | 2.01        | 4.00        | 7.06        | 9.38        | 10.17        | 245        |
| 20        | 68        | 2.00        | 4.00        | 7.02        | 9.26        | 10.06        | 228        |
| <b>25</b> | <b>77</b> | <b>2.00</b> | <b>4.01</b> | <b>7.00</b> | <b>9.21</b> | <b>10.01</b> | <b>220</b> |
| 30        | 86        | 2.00        | 4.01        | 6.99        | 9.16        | 9.96         | 212        |
| 40        | 104       | 2.00        | 4.03        | 6.97        | 9.06        | 9.88         | 195        |
| 50        | 122       | 2.00        | 4.06        | 6.97        | 8.99        | 9.82         | 178        |
| 60        | 140       | 2.00        | 4.10        | 6.98        | 8.93        | 9.76         | 160        |
| 70        | 158       | 2.01        | 4.16        | 7.00        | 8.88        | --           | 142        |
| 80        | 176       | 2.01        | 4.22        | 7.04        | 8.83        | --           | --         |
| 90        | 194       | 2.01        | 4.30        | 7.09        | 8.79        | --           | --         |

#### Buffer solutions DIN 19266

| °C        | °F        | pH           |              |              |              |              |               |               |  |
|-----------|-----------|--------------|--------------|--------------|--------------|--------------|---------------|---------------|--|
| 5         | 32        | 1.668        | 4.004        | 6.951        | 7.087        | 9.395        | 10.245        | 13.207        |  |
| 10        | 50        | 1.670        | 4.000        | 6.923        | 7.059        | 9.332        | 10.179        | 13.003        |  |
| 20        | 68        | 1.675        | 4.001        | 6.881        | 7.016        | 9.225        | 10.062        | 12.627        |  |
| <b>25</b> | <b>77</b> | <b>1.679</b> | <b>4.006</b> | <b>6.865</b> | <b>7.000</b> | <b>9.180</b> | <b>10.012</b> | <b>12.454</b> |  |
| 30        | 86        | 1.683        | 4.012        | 6.853        | 6.987        | 9.139        | 9.966         | 12.289        |  |
| 40        | 104       | 1.694        | 4.031        | 6.838        | 6.970        | 9.068        | 9.889         | 11.984        |  |
| 50        | 122       | 1.707        | 4.057        | 6.833        | 6.964        | 9.011        | 9.828         | 11.705        |  |
| 60        | 140       | 1.723        | 4.085        | 6.836        | 6.968        | 8.962        | --            | 11.449        |  |
| 70        | 158       | 1.743        | 4.126        | 6.845        | 6.982        | 8.921        | --            | --            |  |
| 80        | 176       | 1.766        | 4.164        | 6.859        | 7.004        | 8.885        | --            | --            |  |
| 90        | 194       | 1.792        | 4.205        | 6.877        | 7.034        | 8.850        | --            | --            |  |

## 5. Sending data via the RS 232 C

### Specifications

Baud Rate: 9600 bps

Word length: 8 bits

Stop bits: 2 bits

Parity: None

### How to activate communication

At the start-up of the instrument or from "System", see page 17.

### Data sending to a printer (thermic or dot printer)

40 columns printers.

Thermic: Compatible EPSON (ESC/POS emulation)


Dot-Matrix: Compatible CITIZEN

### Data sending to PC, Hyperterminal.

Data sending through the RS232C in printer format.

With the standard application Hyperteminal from Windows, the data can be captured in the same format as the printer reports, see pages 12 and 13. They can be exported to an Excel program, as well.

See PC connection cables at page 21 (Replacement parts).

### Data sending to PC, ComLabo software

**ComLabo** (PN LZW8999.99). ComLabo Software allows controlling from a PC several modules, pH-meters, conductivity meters, automatic burettes, Sampler, etc.

**ComLabo Easy** (PN LZW8997.99). ComLabo Easy Software allows obtaining data of pH-meters and conductivity meters from a PC.

### Main features shared

Windows based software.

Data Logger. All calibration and measuring data are stored in the PC.

Sample code assignment from PC

Data export to EXCEL.

## 6. Maintenance

### 6.1. Clean the instrument

**Important note:** Never use cleaning such as turpentine, acetone or similar products to clean the instrument including the display and accessories.

Only clean the housing and accessories using a soft, damp cloth. Mild soap solution may also be used. Dry the cleaned parts carefully with a soft cotton cloth.

### CAUTION

When using chemicals or solvents, comply with the instructions of the producer and all local safety regulations.


### 6.2. Sensor cleaning

See the manual of the sensor in use.



## 7. Warning messages

### During calibration

### During measurement

Measuring by stability.  
Rapid criterion.

Measuring by stability.  
Standard criterion.

Measuring by stability.  
Strict criterion.

ISE: Measured value superior to the highest standard used in calibration.

ISE: Measured value inferior to the lowest standard used in calibration.

## 8. Accessories and spares

| Code          | Description   |
|---------------|---|
| LZW9110.98    | Set of 3 x 50 ml printed flasks for benchtop pH-meters calibration  |
| LZW9322.99    | Magnetic stirrer  |
| LZW9400.99    | Standard redox solution of 220 mV, 250 ml flask.  |
| LZW9463.99    | pH buffer solution 4.01, 250 ml flask.  |
| LZW9464.98    | pH buffer solution 7.00, 250 ml flask.  |
| LZW9465.99    | pH buffer solution 9.21, 250 ml flask.  |
| LZW9500.99    | Electrolytic solution, KCl 3M, 250 ml flask.  |
| LZW9502.99    | Electrolytic solution, KCl + glycerine, 250 ml flask.   |
| LZW8201.99    | Thermal printer.  |
| LZW8200.99    | Dot-impact printer.   |
| LZW9135.99    | RS 232 C cable, DB-9.   |
| LZW9135USB.99 | RS 232 C cable, USB.  |
| LZW8999.99    | Software "ComLabo" for intercommunication benchtop instruments with a PC. Kit with CD, USB cable and adapter.   |
| LZW8997.99    | Software "ComLabo Easy" for data acquisition instruments from a PC desktop. Kit with CD, USB cable and adapter. |
| LZW9319.99    | 2nd Magnetic stirrer, with sensors holder   |
| LZW9321.99    | Three-sensor holder.  |

For electrodes please ask for specific brochure.

## 8. Warranty, liability and complaints

---

The manufacturer warrants that the product supplied is free of material and manufacturing defects and undertakes the obligation to repair or replace any defective parts at zero cost. The warranty period for instruments is 24 months.

With the exclusion of the further claims, the supplier is liable for defects including the lack of assured properties as follows: all those parts that can be demonstrated to have become unusable or that can only be used with significant limitations due to a situation present prior to the transfer of risk, in particular due to incorrect design, poor materials or inadequate finish will be improved or replaced, at the supplier's discretion. The identification of such defects must be notified to the supplier in writing without delay, however at the latest 7 days after the identification of the fault. If the customer fails to notify the supplier, the product is considered approved despite the defect. Further liability for any direct or indirect damages is not accepted.

If instrument-specific maintenance and servicing work defined by the supplier is to be performed within the warranty period by the customer (maintenance) or by the supplier (servicing) and these requirements are not met, claims for damages due to the failure to comply with the requirements are rendered void.

Any further claims, in particular claims for consequential damages cannot be made.

Consumables and damage caused by improper handling, poor installation or incorrect use are excluded from this clause.

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NEDERLAND B.V.**