


How to Collect Blood Using the


PAXgene® Blood ccfDNA Tube

For Circulating Cell-Free DNA (ccfDNA)




Required Items:

- 1a. Ensure that the PAXgene® Blood ccfDNA Tube (A) is at room temperature (15–25 °C) prior to use and is properly labeled with specimen identification.
- 1b. The PAXgene® Blood ccfDNA Tube should be the last tube drawn in the phlebotomy procedure.




Venipuncture:

2. Using a blood collection set such as the BD Vacutainer® UltraTouch™ Push Button (B) and the One-Use Holder (C), perform venipuncture using your institution's recommended procedure for standard venipuncture



Blood Collection:

- 3a. Hold the PAXgene® Blood ccfDNA tube vertically below the blood donor's arm during blood collection.
- 3b. Allow at least 10 seconds for a complete blood draw to take place. Ensure that the blood has stopped flowing into the tube before removing the tube from the holder.



After Blood Collection:

- 4a. Gently invert the PAXgene® Blood ccfDNA Tube 8 times.
- 4b. Store the PAXgene® Blood ccfDNA Tube upright at (2–30 °C) for up to seven days or 37 °C for up to one day until centrifugation and plasma processing.†

www.PreAnalytiX.com

Before Blood Collection



After Blood Collection



Ref # Product Description

Ref #	Product Description
768165	PAXgene® Blood ccfDNA Tube
367364	BD Vacutainer® UltraTouch™ Push Button Blood Collection Set
364815	BD Vacutainer® One-Use Holder

PAXgene® Blood ccfDNA Tubes and BD blood collection products are distributed by BD and its distributors. Visit www.preanalytix.com and www.bdbiosciences.com/eu/cbp for more information.

†Refer to PAXgene Blood ccfDNA Tube product circular.