
# SteriMax

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*Annealing without Flame*

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Double Tube System  
Ceramic / Quartz glass

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
**Holding device  
for WLD-TEC  
inoculation  
loop holders only  
suitable for all angles**

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**Universal  
support for all  
inoculation  
loop holders  
for angular pos. 1 - 3**

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
**Tray**  
for left side, made  
of stainless steel  
Art.-Nr./Art.-No.  
6.000.341

## Annealing without Flame!

The electrical sterilization system SteriMax is most suitable for the sterilization of inoculation loops, needles and instruments. It is ideal for laboratories and safety cabinets where the use of gases and open flames is not permitted.

### Robust


The unique double tube system, which is made of robust quartz glass and wear-resistant ceramic, reaches its optimal temperature of 900°C (1650°F) after only a few minutes. To sterilize the inoculation loop, simply insert it into the ceramic annealing tube. After 5 - 7 seconds, the inoculation loop is sterilized<sup>1</sup>.

A back side spray protector is integrated for more safety when handling with pathogenic material. Even clogging is not an issue. The ceramic annealing tube can be removed easily for in-depth cleaning. A 100% stainless steel design makes the SteriMax resistant to even the most extreme laboratory conditions.


### Flexible and safe

The sterilization tube can lock itself into place, using the specially designed angular adjusting device. The low and stable housing facilitates ergonomic operation; the unique design protects the working surface from contamination. Exceptional passive safety feature: After prolonged use, the residual heat display signals a hot surface of the sterilization tube in order to protect the user from burns.

- **Stainless steel design**
- **Long - lasting double tube system**
- **Sterilization period 5-7 seconds<sup>1</sup>**
- **Self - locking angular adjustment**
- **Two sizes available**

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**Quartz glass heating element  
wear-persistent**

**The range:**

with 7 angular adjustments  
Residual heat display  
1 inoculation loop holder and inoculation loop  
Tray for 3 instruments, screwdriver  
2-year warranty

**Double Tube System:**

**Diameter Ø 15 mm, length 118 mm:**

for 230 V mains voltage **Art.-Nr. 5.001.000**  
for 120 V mains voltage **Art.-Nr. 5.002.000**  
for 100 V mains voltage **Art.-Nr. 5.003.000**

**Double Tube System:**

**Diameter Ø 20 mm, length 161 mm:**

for 230 V mains voltage **Art.-Nr. 5.101.000**  
for 120 V mains voltage **Art.-Nr. 5.102.000**  
for 100 V mains voltage **Art.-Nr. 5.103.000**

**Technical Data:**

Double Tube System  
Ø 15 mm: length 118 mm  
Ceramic annealing tube Ø 15 mm  
Quartz glass heating element Ø 21 mm  
Ø 20 mm: length 161 mm  
Ceramic annealing tube Ø 20 mm  
Quartz glass heating element Ø 26.6 mm  
7  
Angular adjustments:  
Sterilization temperature: 900 °C - 950 °C  
(1650 °F - 1742 °F)  
Power supply Ø 15 mm: 160 W  
Power supply Ø 20 mm: 225 W  
Heating up time: 10 - 16 minutes  
Sterilization time: 5 - 7 seconds<sup>1</sup>  
Heat transfer: Average warming of a stainless  
steel surface, located above in  
motionless air and at an ambient  
temperature of 20 °C (131 °F)  
Distance: 20 cm: 55 °C (131 °F)  
30 cm: 38 °C (100 °F)  
40 cm: 32 °C (90 °F)  
Safety features: residual heat display  
Casing: stainless steel, UV and  
solvent resistant  
Weight Ø 15 mm: 960 g  
Measurements Ø 15 mm: 175 x 129 x 139 mm (w x h x d)  
Weight Ø 20 mm: 1060 g  
Measurements Ø 20 mm: 175 x 129 x 169 mm (w x h x d)  
Licenses: CE: EN 61326-1, EN 61010-1,  
EN 61010-2-010  
EU guidelines: 2014/30/EU,  
2014/35/EU, 2011/65/EU

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