Sicherheits-Laborgasbrenner

Safety Laboratory Gas Burner

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with the enclosed description fulfills the requirements of the Low Voltage Directive 2006/95/EC (valid to 19. April 2016) the Low Voltage Directive 2014/35/EU (valid from 20. April 2016), the regulations of the EMC Directive 2014/30/EC and the eco design Directive 2009/125/EC.

Hereby, we certify that our products, regardless of the production location, RoHS compliant and fulfill the directive 2011/65/EC (revised version: directive 2002/95/EC).

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*A. Wiegner*

Firmenstempel / Company stamp

Armin Wiegner  
 Vice President Research & Development

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2.2 EN 61010-2:010-2014  
 Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 1: General requirements  
 Part 2: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 3: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 4: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 5: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 6: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 7: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 8: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 9: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 10: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 11: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 12: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 13: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 14: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 15: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 16: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 17: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 18: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 19: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 20: Safety requirements for electrical equipment for measurement, control, and laboratory use  
 Part 20:10: Particular requirements for laboratory equipment for the heating of materials

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


Altengottern, 20.04.2016


- 1 - Function knob
- 2 - Gas / Air adjustment
- 3 - LED Standard
- 4 - LED Start-Stop
- 5 - LED Button
- 6 - LED Burner head hot / BHC
- 7 - Burner head
- 7a - Flame orifice

- 8 - Connector for foot pedal
- 9 - Power connector
- 10 - Gas inlet  
R 1/4" L (left hand thread)
- 11 - Holding device for inoculation loop holder

- 12 - Burner head screw
- 13 - Monitor electrode
- 14 - Ignition electrode


- 15 - Active nozzle
- 16 - Nozzle holder for alternative gas
- 17 - Cover of the burner shaft
- 18 - Retaining screw for cover of the burner shaft

- 15a - Position screw for the cover
- 19 - Guide slots for tilt adjustment
- 20 - [Fuego basic RF] Function LED and programming key

- 21 - [Fuego basic RF] Battery compartment (2 x AA, Mignon)
- 21a - [Fuego basic RF] Transport insulation for the batteries

## Table of contents

Use.....	22
Safety Precautions.....	22
The range.....	24
<b>1. Set up</b> .....	25
1.1 Choice of nozzle.....	25
1.2 Gas connection.....	25
1.3 Mains power supply.....	26
1.4 Foot pedal connection (tethered).....	26
<b>2. Operation</b> .....	26
2.1 On-Off switch, operating mode.....	26
2.2 Flame regulation.....	26
2.3 [Fuego basic RF] Radio foot pedal activation and function.....	27
2.4 [Fuego basic RF] Function LED on the radio foot pedal.....	27
2.5 [Fuego basic RF] Foot pedal programming.....	28
2.6 Switch-off the unit.....	28
<b>3. Application programs</b> .....	28
3.1 Button Start-Stop.....	28
3.2 Pedal Standard.....	29
3.3 Pedal Start-Stop.....	29
<b>4. Safety symbols and safety functions</b> .....	29
4.1 Residual heat display.....	29
4.2 Burner Head Control (BHC).....	29
4.3 Automatic unit switch off.....	29
<b>5. Error displays</b> .....	29
5.1 Ignition failure.....	30
5.2 Flame failure.....	30
5.3 Overtemperature.....	30
5.4 Burner head assembly monitor.....	30
5.5 Burner Head Control (BHC).....	30
5.6 Foot pedal connection control.....	30
<b>6. Cleaning and sterilizing</b> .....	31
6.1 Burner head disassembly and cleaning.....	31
6.2 Burner shaft cleaning.....	31
<b>7. Turbo flame</b> .....	32
<b>8. Tilt adjustment</b> .....	32
<b>9. [Fuego basic RF] Changing the batteries</b> .....	32
<b>10. Warranty</b> .....	32
Technical data.....	33
Troubleshooting guide.....	34
Declaration of Conformity.....	36

Read these instructions carefully to familiarize yourself with the product. Please retain these operating instruction for future reference.

## Use

**Safety laboratory gas burners for heating and flame sterilizing.  
Ideal for use in cleanroom workbenches and laboratory.**

All users who have been assigned to use this device must have read and understood these operating instructions or have been instructed by an expert user so that this device can be used safely without causing danger.

**ATTENTION: OPERATE THIS LABORATORY GAS BURNER UNDER  
CONSTANT SUPERVISION ONLY!**

## Safety Precautions

- On unpacking the unit, check for possible transportation damages. Do not operate the unit if damages are visible.
- After use or for any longer period of time without attendance, turn the main gas supply off and turn off the gas burner at the function knob (1).
- All gas connections must be adequately tightened (left-hand thread). Ensure gas proofness with suitable test equipment. DO NOT seal up the thread of the gas connection (10) of the laboratory gas burner with Teflon tape, etc.
- BEFORE using the device carefully check the gas feed tube for leaks. Check this even if the device has been installed by your distributor. To do this, carry out all the procedures mentioned in these operating instructions (see paragraph 1.).
- In the event that gas can be smelled: immediately turn off the gas supply to the device. Extinguish any open flames. Pull out the mains plug. Check all gas connections for gas proofness. If the smell of gas persists, the appropriate authorities must be notified (janitor, gas utility company, Fire Brigade).  
**LEAKING GAS CAN CAUSE A FIRE OR AN EXPLOSION. THIS MAY RESULT IN SEVERE INJURIES, FATAL ACCIDENTS AND DAMAGE TO PROPERTY.**
- The device can be dangerous if operated or used in an incorrect manner by untrained staff.
- An incorrect gas connection may create a hazard. Observe the installation instructions in the manual.

- Do not store spare or unconnected gas cartridges / gas bottles in the vicinity of this device.
- Even in an apparently empty gas cartridge / gas bottle, some gas may still remain. Gas cartridges / gas bottles should be transported and stored accordingly. Empty gas cartridges should be properly disposed.
- Do not use the device if there is a smell of gas or if there is a leak.
- NEVER try to loosen or unfasten gas connections while the gas supply is turned on and the device is in operation.
- Pay attention to your relevant rules for using liquid gas.
- Only use DVGW safety tubings with thread or tubing connectors. Check the condition of the tube / hose frequently. Depending upon type of tube / hose, hose clamps are required.
- NEVER use an open flame to look for leaks.
- Do not smoke if you are searching for leaks.
- Keep hands or other parts of the body away from the burner orifice (7a).
- Do not operate the unit near flammable liquids or hazardous materials.
- Unattended operation of the unit is not permissible.
- Always work in a well-ventilated area.
- Note that the burner orifice (7, 7a) remains hot after the flame has been extinguished. Do not touch. Can cause burns.
- Allow sufficient time for flame orifice (7a) to cool down prior to cleaning, disinfecting, servicing or transport. Ensure that the unit and the gas supply are turned off.
- Because of the connectors at the back of the unit the backside should not be sterilized with a flame.
- Allow sufficient time for burner head (7) to cool down prior to disassembling.
- Operate the unit with assembled burner head (7) only.
- After cleaning the burner head (7) allow sufficient time to dry before assembling again.
- Do not allow anything to fall into the flame orifice (7a).
- Before mounting a nozzle check the O-Ring (22). Replace the O-Ring if damaged or worn (see paragraph 1.1).
- When working with this device, always wear protective glasses.

## The range

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Art.-No. 8.201.000

### Fuego SCS basic RF

Art.-No. 8.203.000-RF (only available in Europe)

- [Fuego basic] Tethered foot pedal (R7)
- [Fuego basic RF] Radio frequency foot pedal (battery operated) (R8)
- 3 standard-programs for button (function knob) and foot pedal
- SCS (Safety Control System)
- BHC (Burner Head Control)
- Removable and decomposable burner head
- Tilt mechanism, right / left (R1)
- Holding device for 3 inoculation loop holders (11)
- Nozzles for natural gas, propane/butane gas
- Turbo flame
- Wrench 17 mm (R2) for gas connection
- Screwdriver (R3) for burner head and cover of the burner shaft
- Tubing connector with swivel nut for 10 mm pipe spout tubing (R4)
- Power connection with 4 adapter (R5), global (level 6)
- [Fuego basic RF] 2 Mignon batteries (AA), (R6), already inserted in the battery compartment of the foot pedal (R8)
- Instruction manual
- 2-year warranty

Fig. G

## 1. Set up

Initial operation is to be carefully carried out as described in the following paragraphs.

Failure to observe the instruction manual **may give rise to hazards from leaks and / or bursts of flame.**

### 1.1 Choice of nozzle

The unit is shipped with the nozzle for natural gas (N, Ø 0.8 mm) installed. The nozzle for propane / butane (P, Ø 0.5 mm) can be found in the nozzle holder (16). The nozzle on the device must be fitted or the gas type in use, in which case proceed as follows: On the bottom of the device, remove the nozzle P for propane / butane by unscrewing it from the nozzle holder (16) with a coin (1 Cent coin) or the screwdriver (R2). Then remove the active nozzle N for natural gas (15) (preinstalled) and exchange the nozzles. Tighten both nozzles again after changing over. The optional nozzle B (Ø 0.45 mm) is only for usage with the adapter CV 360.

Each time before mounting a nozzle, check the O-ring (22) for damage. **Replace O-ring if damaged, Art. No.: 8.000.010**

### 1.2 Gas connection

Connect the device via the gas inlet (10) to the gas supply, if the operating pressure for natural gas is 18 - 25 mbar, for propane / butane 47.5 - 57.5 mbar. Only use DVGW or other gas approved safety tubing with thread or tubing connectors. Check the condition of the tube / hose regularly.

Depending on the type of hose, hose clips and/or the hose clamps and the swivel nut included (R4) are to be used.

A DVGW tested or other gas approved pressure regulator (50mbar) must be used for liquid gas. Pay attention to your relevant rules for using liquid gas.

If using gas cartridges, use only an original WLD-TEC gas cartridge adapter! Observe the operating instructions in question for the gas cartridge adapter.

**DO NOT seal up the thread of the gas inlet (10) of the laboratory gas burner with Teflon tape, etc.**

All gas connections must be adequately tightened with a wrench SW17 (R2) (included in delivery). Note: left-hand thread. Check gas proofness with testing equipment or soapy water. To do this with the gas supply turned on, wet the gas


Fig. 1.1/1

Fig. 1.1/2

Fig. 1.2/1

connections with test solution or soapy water and look for bubbles. If bubbles form or if one bubble becomes larger, there is a leak. Immediately turn off the gas supply. Check all gas supply connections or gas joints (e.g.: at the gas cartridge / gas bottle, at the laboratory gas burner or at an additional gas leak protection). In the event of leaks, tighten up the gas joints with a wrench and test again for gas proofness.

**Never use an open flame to look for leaks.**

### 1.3 Mains power supply

Plug the power cable into the socket **(9)** at the back of the unit or into the socket of the tethered foot pedal. The mains power supply unit included is designed for a voltage from 100 - 240 V, 50/60 Hz. The mains power supply unit may only be connected if these values comply with the intended electricity supply. Before using it, push the corresponding socket adapter onto the mains power supply plug.

### 1.4 Foot pedal connection (tethered)

Insert the connection cable of the foot pedal into the socket **(8)** at the back of the unit.

If no pedal is connected, the LED's "Standard" **(3)** and "Start-Stop" **(4)** will flash alternately.

## 2. Operation

This section describes the operation of the unit with the function knob and the gas / air adjustment, as well as the functioning of the [\[Fuego basic RF\] radio foot pedal](#).

### 2.1 On-Off switch, operating mode

Switch the unit on by a short push on the function knob **(1)**. It can be turned off by a long push (over 2 seconds) on the function knob **(1)**. By turning the function knob the foot pedal operating mode "Standard" or "StartStop" or the operating mode "Button" can be chosen. The corresponding LED's will light up.

### 2.2 Flame regulation

The flame can be varied in size and intensity by turning the gas / air knob **(2)**. Regulating the size of the flame is done by turning the gas /air adjustment. Regulating the air intake for the flame intensity is carried out by **pulling outwards and at the same time rotating** the gas / air adjustment.

When using for the first time or changing the gas type, turn the gas / air adjustment knob two revolutions to the left. Then pull the gas / air adjustment outwards and simultaneously turn it one revolution to the left to open the air intake.

## 2.3 [Fuego basic RF] Radio foot pedal activation and function

**Switch the burner off prior to operation to prevent inadvertent activation of the burner.**

Remove the transport insulation from the batteries **(21a)** prior to first use. The two mignon batteries (AA) have already been inserted into the battery compartment of, the wireless foot pedal **(21)**.

The LED's Standard **(3)** and Start-Stop **(4)** will flash alternately until a wireless foot pedal has been detected within range (2 - 5 metres).

Once the transport insulation has been removed, switch burner on again. The wireless foot pedal does not have a mains switch, but is automatically switched on when operated. The foot pedal also deactivates itself after a certain time once the laboratory gas burner has been switched off.

The radio system can be used without changing any further settings as the laboratory gas burner and the wireless foot pedal have already been programmed to each other **(gas burner and foot pedal programming, see paragraph 2.5)** by the manufacturers.

**The laboratory gas burner is started as soon as the foot pedal is operated and the LED's (3) and (4) stop blinking. Attention: Flame!**

The corresponding LED for the foot pedal program **(3)** or **(4)** that has been selected will also light up continuously in green and signals that the wireless foot pedal is connected to the appliance. (In cases the alternate flashing continues after operation and the foot pedal is not being detected, **(see paragraph 2.5.)**)

The LED's "Standard" **(3)** and "Start-Stop" **(4)** of the laboratory burner will flash alternately as soon as the wireless foot pedal is no longer within range (2 – 5 metres) or the radio signal is interrupted. In this case the flame will be switched off automatically for safety reasons.

As an alternative, the burner can also be operated using a tethered foot pedal (optional): In this case, plug the foot pedal cable into the foot pedal connector **(8)** on the back of the burner. The burner will detect the foot pedal and the alternate flashing of the LEDs Standard **(3)** and Start-Stop **(4)** will cease as soon as the foot pedal has been plugged in.

## 2.4 [Fuego basic RF] Function LED on the radio foot pedal

The function LED **(20)** flashes at different intervals to indicate different operational status.

**No LED lights up:**

The foot pedal is deactivated / switched off.

**Brief, faint flashes:**

Fully functional, unimpeded radio transmission.

**Brief, intense flashes:**

Renewed radio transmission at weak radio signal strength or signal overlap.

**Brief, intense flashes at rapid set intervals:**

The RF stick is out of range or has been switched off. If no burner is detected after 40 flashes, the foot pedal will be switched off automatically.

**Intermittent lighting up of the LED, duration approx. 1 sec:**

Batteries are low (see paragraph 9).

**Slight glow:**

Programming in progress (see paragraph 2.5).

## 2.5 [Fuego basic RF] Foot pedal programming

The laboratory burner and foot pedal can be newly programmed to each other if the signal is interrupted or in cases of parallel use of several radio systems.

For programming of the radio channel, place the laboratory gas burner close to the wireless foot pedal (max. 30 centimetres). Programming has a restricted range to avoid inadvertent re-programming of other burners. (In cases where undesired re-programming of a neighbouring burner occurs, switch it off during the programming procedure.)

### Programming procedure:

**Adhere to the sequence for the programming procedure. Close the gas regulator (2) by turning it clockwise to avoid inadvertent activation of the flame.**

1. Switch off the laboratory gas burner.
2. Switch on foot pedal by operating briefly or press again briefly.
3. Switch on the laboratory gas burner.
4. Shortly afterwards, press the programming key **(20)** on the foot pedal for approx. 1 second to start the programming procedure. The function LED **(20)** will now start to glow. This procedure can last for between 0.5 and 30 seconds. As soon as the laboratory burner and the foot pedal have been successfully programmed, the glowing is switched off and the function LED **(20)** switches to brief, faint flashes.

The laboratory burner and the foot pedal have already been programmed to each other by the manufacturers.

## 2.6 Switch-off the unit

The unit can be turned off by pushing the function knob **(1)** for more than 2 seconds. [Fuego basic RF] The radio foot pedal is deactivated automatically.

**To depressurize gas hose, shut off the gas supply and activate the burner again until the flame extinguishes in order to burn the residual gas. Afterwards turn off the gas burner at the function knob (1).**

## 3. Application programs

This chapter explains the various application programs ("Button Start-Stop", "Pedal Standard" und "Pedal Start-Stop").

### 3.1 Button Start-Stop

The flame is ignited by operation of the function knob **(1)**. The flame is extinguished after renewed actuation of the function knob **(1)**. Additionally the flame is automatically extinguished when the burning timer has expired after 60 min.

### 3.2 Pedal Standard

The flame is ignited by operation of the foot pedal. The foot pedal remains depressed for the duration of use. The flame is extinguished upon release of the pedal.

### 3.3 Pedal Start-Stop

The flame is ignited by operation of the foot pedal. The flame is extinguished after renewed actuation of the foot pedal. Additionally the flame is automatically extinguished when the burning timer has expired after 60 min. Alternatively the flame can be extinguished by a short push on the function knob.

While no foot pedal is connected, only Button "Start-Stop" is available.

## 4. Safety symbols and safety functions

The safety LED's appear during usage to warn the user of potential hazards.

### 4.1 Residual heat display

The LED „Burner Head Hot / BHC“  (6) indicates a hot burner head. As long as the LED shines, DO NOT TOUCH the burner head. This can cause serious burns. Even after switching-off the unit, the residual heat LED remains "on" till the burner head is cooled down.

**Disconnecting the power supply or removing the power cord will clear the residual heat display even if the burner head is still hot.**

### 4.2 Burner Head Control (BHC)

If the burner head is clogged, the amber LED "Burner Head HOT / BHC" (6) will flash you have to clean the burner head (**see paragraph 6.1**), otherwise the maximum burning time in the operating mode "Button" and Pedal "Start-Stop" is limited to 30 seconds (**see paragraph 5.5**). If burning times longer than 30 seconds are required in case of a clogged burner head, the operating mode "Pedal Standard" can be used without time limit.

### 4.3 Automatic unit switch off

The unit switches itself off automatically after 4 hours if the flame has not been lit in this period. All indicated malfunctions are automatically switched off after 4 hours, too (**see paragraph 5**). For further operation, switch the unit on again.

## 5. Error displays

If there is an error during operation, eventually the gas valve of the burner automatically closed and different flashing LEDs indicate the error.

All error displays can be reset by a long push (> 2 seconds) on the function knob (1). In case of overtemperature the unit needs to be cooled down and in case of burner head assembly monitor, the burner head needs to be reinstalled before a reset is possible.

## 5.1 Ignition failure

### **Green LED "Button" (5), "Standard" (3) or "Start-Stop" (4) blinks 2x**

This signal appears and indicates a malfunction if the flame fails to ignite after 7 seconds. In case of ignition failure check the burner head (7) for possible clogging, check the correct input pressure of the gas supply and verify that the correct nozzle is installed. In case of this malfunction the gas valve of the burner will be shut off automatically.

Nozzle N: natural gas, 18 - 25 mbar

Nozzle P, B: propane / butane gas, 47.5 - 57.5 mbar

## 5.2 Flame failure

### **Green LED "Button" (5), "Standard" (3) or "Start-Stop" (4) blinks 3x**

This signal indicates a malfunction if the flame is extinguished by external factors and fails to reignite within 5 sec. In case of flame failure check the burner head (7) for possible clogging and verify the correct input pressure of the gas supply. In case of this malfunction the gas valve of the burner will be shut off automatically.

## 5.3 Overtemperature

### **Green LED "Button" (5), "Standard" (3) or "Start-Stop" (4) blinks 4x**

This signal indicates a malfunction if the interior temperature has exceeded 70 °C. At a normal room temperature with normal air circulation the unit is suited for continuous operation. In case of overtemperature increase the air ventilation or change the operation site. In case of this malfunction the gas valve of the burner will be shut off automatically.

## 5.4 Burner head assembly monitor

### **Green LED "Button" (5), "Standard" (3) or "Start-Stop" (4) blinks 5x**

This message indicates that the burner head is removed. Further operation is possible after the burner head is reinstalled.

## 5.5 Burner Head Control (BHC)

### **Orange LED "Burner Head HOT / BHC" (6) flashes**

This signal indicates that the time limit (30 seconds) is turned on in operating mode "Start-Stop" and "Button" due to a clogged burner head. For cleaning the burner head see paragraph 6.1.

## 5.6 Foot pedal connection control

### **LED "Standard" and "Start-Stop" are blinking alternately**

Tethered foot pedal is not connected with the burner. Plug in the foot pedal.

[Fuego basic RF] Wireless foot pedal out of range or signal interrupted.

Switch on wireless foot pedal or move within range (see paragraph 2.3).

If necessary, re-programme the wireless foot pedal (see paragraph 2.5).

## 6. Cleaning and sterilizing

Allow sufficient time for burner orifice (**7, 7a**) to cool down before disassembling or cleaning the burner head. Check the unit is disconnected and that the gas supply is turned off at the mains. The burner can be cleaned with customary commercial disinfectants. Additionally, it is possible to remove the burner head and to clean it separately. The stainless steel and glass construction allow 100% UV-radiation sterilization and short time surface flame sterilization.

Because of the connectors at the back of the unit the backside should not be sterilized with a flame!

### 6.1 Burner head disassembly and cleaning

Allow sufficient time for burner orifice (**7,7a**) to cool down before disassembling or cleaning the burner head. Check the unit is turned off, that the gas supply is turned off at the mains. Clean the burner head with customary commercial disinfectants, sterilize it in an autoclave or wash it in a dishwasher. To remove the burner head proceed as follows:

Unscrew the burner head screw (**12**) completely with the included screwdriver. Turn approx. 8 revolutions to the left.

Now remove the burner head from the device by pulling it upwards. Reinstallation is performed in the reverse sequence. The dismantled burner head can be even dismantled into the individual components for in-depth cleaning: Unscrew both screws (**12a**) with the included screwdriver (**R3**) and take off the base plate (**12b**) of the burner head which was fixed by the two screws (**12a**). After the base plate is removed both electrodes can be pulled out for separate cleaning. Reinstallation is performed in the reverse sequence.

**i** When dismantling the burner head completely the sealing ring placed around the burner head screw (**12**) could dropout. Ensure that the sealing ring is placed around the burner head screw (**12**) when reassembling.


Fig. 6.1/1

### 6.2 Burner shaft cleaning

Unscrew the screw (**18**) completely with the included screwdriver (**R3**). Take off the cover (**17**) of the burner shaft.

Now the burner shaft can be cleaned or solid substances which have fallen into the unit can be removed. Reinstallation is performed in the reverse sequence.

Take care that the notch of the cover fits to the screw (**18a**).

Fig. 6.2/1

## 7. Turbo flame

If the cover of the burner shaft (17) is removed the flame is extremely firm and consistent. To take off the cover of the burner shaft unscrew the screw (18) completely with the included screwdriver (R3). With an open burner shaft the intensity of the flame cannot be adjusted by the air knob (2) any longer. During the use of the turbo flame most of the needed air is taken inside through the open burner shaft. Remounting the cover of burner shaft (see paragraph 6.2).

## 8. Tilt adjustment

Insert the tilt adjustment (R1) into the slots (19) at the bottom of the unit. The tilt adjustment can be used to the left or right side to protect the burning chamber from contamination when working with liquids.


Fig. 8./1

## 9. [Fuego basic RF] Changing the batteries

Change the batteries in the wireless foot pedal as soon as the function LED (20) signals "batteries low". (intermittent lighting up of the LED, duration approx. 1 sec.)

**Switch off the laboratory gas burner while changing the batteries, to avoid inadvertent activation of the laboratory burner through unintended pressing of the foot pedal while inserting the batteries.**

Change both mignon batteries (AA) in the battery compartment of the wireless foot pedal (21). Pay attention to correct insertion in accordance with the batteries' polarity (see label in the battery holder of the foot pedal).

## 10. Warranty

All Fuegos are covered under our two-year manufacturer warranty against any manufacture defects in material and workmanship. The WLD-TEC warranty guarantees all Fuegos under normal usage conditions and does not cover any damages as a direct result of user misuse or/and abuse. The warranty is void upon any unauthorized servicing, disassembly or modifications.

[Fuego basic RF] In respect of parts to wear and tear (e.g. batteries) this warranty shall be valid for six months from the date of purchase.

## Technical data

Technology:	Microprocessor
<b>Programs</b>	
Foot pedal:	Standard (flame during pressed foot pedal) Start-Stop with timer, 60 min
Button:	Start-Stop with timer, 60 min

EN

### Safety features

Safety Control System (SCS) with gas safety cut off:	ignition and flame control, temperature monitor, burner head clogging and assembly monitor (BHC)
Automatic unit switch off:	4 h
Residual heat display:	indicates a hot burner head

### Gas supply and consumption

Gas supply:	1/4" left + filter
Gas types:	natural gas E/LL, 18 - 25 mbar, liquid gas 47.5 - 57.5 mbar, II <sub>2ELL3B/P</sub>
Connected load:	95 l/h natural gas (LL), 70 g/h liquid gas
Continuous cartridge operation: (approx.)	<i>Campingaz</i> : CV 360 - 40 min, C 206 - 170 min, CP 250 - 210 min, CV 470 - 370 min, CG 1750 - 150 min, <i>Guilbert</i> : Express 444 / 445 - 50 min <i>Coleman</i> : C 250 - 160 min, C 500 - 320 min

### Temperatures

Flame temperature:	1350 °C on liquid gas / 1300 °C on natural gas (E)
Temperature threshold level:	1 kW liquid gas, 1 kW natural gas

### Electrical

Power consumption:	2 VA (stand by max. 0.1 VA)
Power connection:	100 - 240V / 50/60Hz / max. 0.3 A; 9V DC / 1.3A Level 6 (stand by power consumption max. 0.1W)

[Fuego basic RF] Power supply  
foot pedal:

2 x AA (Mignon batterie)

### Mechanical

Casing and operating controls:	stainless steel / glass, UV and solvent resistant
Burner head:	removable and decomposable, stainless steel
Cover of the burner shaft:	Ø 23 mm, with drains
Measurements (w x h x d):	103 x 49 x 130 mm
Weight:	700 g

### [Fuego basic RF]

#### Radio Frequency System

Range:	2 - 5 m
Frequency:	2400 - 2483.5 MHz
Output power:	< 10 mW
Channels:	253
Batteries foot pedal:	2 x AA (Mignon batterie)
Life time batteries:	approx. 2000 h

### Licences

DIN-DVGW Reg.-No.:	NG-2211AS0167
CE:	EN 61326-1, EN61010-1, EN61010-2-010 [Fuego basic RF] EN 301489, EN 300440-1, EN 300440-2 2014/30/EU, 2014/35/EU, 2011/65/EU [Fuego basic RF] 2014/53/EU (RED)
EU guidelines:	

## Troubleshooting guide

- **The green LED "Standard" or "Start-Stop" does not light up**

Check for correct connection and specification of the power adapter.  
Ensure that the original power adapter is used.

Specifications: 9 V / DC, 1.3 A

Polarity: 

---

- **The tethered foot pedal does not work**

Check for correct connection of the foot pedal. Ensure that the foot pedal socket and plug is not twisted or broken.

---

- **No Flame**


In case of ignition or flame failure check if the burner head is clogged.  
Verify the input pressure of the used gas. Ensure that the correct nozzle is installed in the unit.

Nozzle N: natural gas, 18 - 25 mbar

Nozzle P, B: propane / butane gas, 47.5 - 57.5 mbar

---

- **LED "Burner head HOT / BHC" flashes / Inspection of the burner head (clogging)**

Take care that there are no liquids or other substances at the Bypass (area between the inner and the outer ring)

Especially remove substances in the marked areas at the electrodes. If there are contaminants in this area, the flame cannot encircle the electrodes correctly.

Clean this areas with a brush. The burner head can be cleaned with customary commercial disinfectants, or it can be sterilized in an autoclave or washed in a dishwasher.

---

- **In operating mode "Button" or Pedal "Start-Stop" flame burns 30 seconds, only**

BHC time limit is active, LED "Burner head hot / BHC" is flashing.

The burner head is clogged and must be cleaned.

**(see paragraph 4.2, 5.5 & 6.1).**

---

- **Flame too small / large / soft**

Check the position of the air and gas adjustment.

Check if the correct nozzle is installed.

Nozzle N: natural gas, 18 - 25 mbar

Nozzle P, B: propane / butane gas, 47.5 - 57.5 mbar

Check if the drilling of the active nozzle is blocked (**see paragraph 1.1**). Unscrew the active nozzle. If the drilling is blocked clean with a brush or compressed air.

- **No ignition spark / LED Burner head "HOT / BHC" flashes but the burner head is clean**

Check if the ceramic electrodes are in good condition. In some cases the electrodes may break. To check move the ends of the electrodes. If they are not moving they should be okay. If they are moving more than 0.5 mm the electrodes are broken. The electrodes can be dismantled and changed by the user (**see paragraph 6.1**).

- **The burner shuts-off due to overtemperature frequently**

In case of overtemperature increase the air ventilation or change the operation site.

- **Green LED "Button" (5), "Standard" (3) or "Start-Stop" (4) blinks 2x**  
Ignition failure (**see paragraph 5.1**).

- **Green LED "Button" (5), "Standard" (3) or "Start-Stop" (4) blinks 3x**  
Flame failure (**see paragraph 5.2**).

- **Green LED "Button" (5), "Standard" (3) or "Start-Stop" (4) blinks 4x**  
Overtemperature (**see paragraph 5.3**).

- **Green LED "Button" (5), "Standard" (3) or "Start-Stop" (4) blinks 5x**  
This message indicates that the burner head is removed or not mounted correctly. Reinstall burner head and reset the unit by a long push on the function knob (**see paragraph 5.4**).

- **Amber LED "Burner head HOT / BHC" (6) is on permanently**  
Residual heat display is active. Attention: DO NOT TOUCH the burner head (**see paragraph 4.1**).

- **LEDs "Standard" (3) and "Start-Stop" (4) are blinking alternately**  
Tethered foot pedal not connected to the burner (**see paragraph 1.4**).  
[Fuego basic RF] Wireless foot pedal out of range or signal interrupted (**see paragraph 2.3**). In cases of frequent occurrence of this error, check the function LED on the wireless foot pedal (**20**) and programme a new channel or replace weak batteries, as appropriate (**see paragraph 2.4, 2.5 or 9**).

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Wir, the manufacturer, hereby confirm, *that the product,*

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*with the enclosed description fulfills the requirements of the Low Voltage Directive 2006/95/EC (valid to 19. April 2016) the Low Voltage Directive 2014/35/EU (valid from 20. April 2016) the regulations of the EMC Directive 2014/30/EC and the eco design Directive 2009/125/EC.*

*Hereby, we certify that our products, regardless of the production location, RoHS compliant and fulfill the directive 2011/65/EC (revised version: directive 2002/95/EC).*

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*A. Wegener*

Firmenstempel / Company stamp

Armin Wegener  
 Vice President Research & Development

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