ME&DI

INSTRUCTIONS FOR USE

ANSELL SURGICAL Cat IIa GLOVES

Available in 28 languages

v.3 2025-06

Please click on the box below to select your preferred language.

North America

US Legal Manufacturer/Fabricant:
Ansell Healthcare Products LLC
2301 Robb Drive,
Reno, NV 89523, USA
+1 732 345 5400
+1 800 800 0445

Latin America and Caribbean:

Ansell Commercial México S.A. de C.V.
Blvd. Bernardo Quintana No. 7001-C,
Q7001 Torre II, Suites 1304,
1305 y 1306. Col. Centro Sur, C.P. 76079
Querétaro, Qro. México
+52 442 296 2050

UK Importer / Responsible Person:

Nitritex Ltd
Ground Floor, 15 Kings Court,
Willie Snaith Road, Newmarket,
Suffolk, CB8 7SG, United Kingdom

Europe, Middle East and Africa Region

Ansell Healthcare Europe NV
Boulevard International 55
Brussels, B-1070
Belgium

Brazil

Ansell Brazil Ltda.
Rua das Figueiras 474 - 4º Andar,
Santo André, SP - Brazil
CEP 09080-300
CNPJ: 03.496.778/0001-21

Australia

Authorised Sponsor ANZ: Ansell Limited
Level 3, 678 Victoria Street,
Richmond, Vic, 3121
+61 1800 337 041
+61 1800 803 578

Middle East and Africa (EMEA)

Ansell Healthcare Europe NV
Boulevard International 55
Brussels, B-1070
Belgium

Asia Pacific

Ansell Global Trading Center (M) Sdn Bhd
Prima 6, Prima Avenue,
Block 3512, Jalan Teknokrat 6
63000 Cyberjaya, Selangor, Malaysia
+60 3 8310 6688
+60 3 8318 6699

Switzerland

Authorised Representative for Switzerland:
Cosanum AG
Brandstrasse 28
8952 Schlieren
Switzerland

www.ansell.com

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. © 2025 Ansell Limited. All Rights Reserved.

CONTENTS

Please click on the content below to navigate to your selected page.

DOMINIQUE DUTSCHER

SYMBOLS & PICTOGRAMS

EN ISO 21420:2020

(01)

EN 388:2016
+A1:2018

A B C D E

(02)

EN ISO 374-5: 2016

VIRUS

(03)

EN ISO 374-1:2016+A1:2018
TYPE A, B OR C

ABCDEFGHIJKLMNOPS

(04)

EN 421:2010

(05)

(06)

(07)

(08)

(09)

(10)

(11)

(12)

(13)

(14)

(15)

(16)

(17)

(18)

(19)

(20)

(21)

EN ANSELL SURGICAL CAT IIa GLOVES

A. USE: This Instruction for Use is to be used in combination with the specific information that appears on the first packaging. Gloves are designed as a two way barrier protection against cross contamination at the surgical site and also to protect the hands mainly against chemical splash risks and comply with the standards shown by the relevant pictograms. If a serious incident occurs with the glove, please report it to Ansell Healthcare Europe NV and the Competent Authority where it was being used.

EXPLANATION OF SYMBOLS & PICTOGRAMS (01) EN ISO 21420:2020 / ISO 15223-1:2021 – Please read the Instructions for Use, prior to using the products, or contact Ansell for more information. If a level X is mentioned under any of the pictograms, this means this test is not applicable and glove is not designed and therefore not to be used for this specific hazard. **(02) EN 388:2016+A1:2018 – A B C D** E = Protection from mechanical risks: A: Abrasion resistance (performance levels 0 to 4) B: Blade cut resistance (performance levels 0 to 5) If level X is indicated for this property the TDM as per E is the reference performance result for cut resistance C: Tear resistance (performance levels 0 to 4) D: Puncture resistance (performance levels 0 to 4) E: TDM ISO EN 13997 cut resistance (performance levels A to F) Warning! The performances (A to E) claimed for the gloves are based on tests performed on the palm area of the gloves only. For gloves with two or more layers, these overall performance levels may not necessarily reflect the performance of the glove's outermost layer. For gloves where the palm and back and cuff are different, mechanical protection is only applicable to the palm of the glove. **(03) EN ISO 374-5:2016** – Protection against bacteria, fungi and virus. **(04) EN ISO 374-1:2016+ A1:2018 TYPE A, B OR C** – Protection against chemical hazards: Type A = chemical breakthrough time > 30 minutes against at least 6 chemicals as per list defined in EN ISO 374-1:2016 Type B = chemical breakthrough time > 30 minutes against at least three chemicals as per list defined in EN ISO 374-1:2016 Type C = chemical breakthrough time > 10 minutes against at least one test chemical as per list defined in EN ISO 374-1:2016 (no code underneath the pictogram) A = methanol – B = acetone – C = acetonitrile – D = dichloromethane – E = carbon disulfide – F = toluene – G = diethylamine – H = tetrahydrofuran – I = ethyl acetate – J = n-heptane – K = sodium hydroxide, 40% – L = sulphuric acid, 96 % – M = nitric acid, 65% – N = acetic acid, 99% – O = ammonia, 25% – P = hydrogen peroxide, 30% – S = hydrofluoric acid, 40% – T = formaldehyde, 37%.

B. WARNING! Chemical resistance data has been assessed under laboratory conditions and relates only to the chemical tested. It can be different if used in a mixture. For gloves equal or longer than 400 mm, the chemical resistance data is based on samples taken, 80 mm from the end of the cuff. The data may not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. Check that the gloves are suitable for the intended use. Conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the chemical due to changes in physical properties. Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Chemical permeation data, tested per EN16523-1:2015+A1:2018, and degradation data, tested per EN ISO 374-4:2019, are available upon request and/or via ansell.com, through the Ansell product page/downloads criteria/CE-certified Chemical Permeation test data. **(05) EN 421:2010** – Protection against radioactive contamination. **(06) CE** = Product is compliant and certified to the requirements of the European Regulations on Personal Protective Equipment 2016/425 and Medical Devices 2017/745. For compliance as a PPE: The CE mark is followed by a four digit code which refers to the identification number of the Notified Body that is in charge of the PPE category III conformity assessment (Module C2 or D), for products to protect against serious risks. To obtain the EU-or UK Conformity Declaration to PPE Regulation, please go to: www.ansell.com/regulatory **(07) MEDICAL DEVICE** = For compliance as a Medical Device product, the CE and UKCA mark is followed by a four digit code which refers to the identification number of the Notified / Approved Body used For EU: I.D. 2797 – BSI Netherlands. For UK: I.D. 0086 – BSI UK. To obtain the EU-or UK Conformity Declaration to Medical Devices Regulation, please visit the product page on www.ansell.com **(08) UKCA** = Product is compliant and, where applicable, certified to the requirements of the Regulation 2016/425 on personal protective equipment, as amended to apply in GB and UK Medical Devices Regulation 2002 (SI 2002 No. 618 as amended) PPE Type examination certificate (Module B) and Supervised product checks (Module C2) or Conformity to type based on quality assurance of the production process (Module D) by Centexbel Belgium (I.D. 0493), Technologiepark 70, B-9052 Zwijnaarde, for EU: Type-examination certificate (Module B) and conformity to type certificates based on Supervised product checks (Module C2) or Conformity to type based on quality assurance of the production process (Module D) for CE marking are used as the basis for applying a UKCA. In some cases, Type-examination certificate (Module B) by Satra Technology Centre, Wyndham Way, Telford Way, Kettering, Northamptonshire, NN16 8SD, UK. For more details, please consult UK Conformity Declaration. **(09) SINGLE USE ONLY (10) A SINGLE STERILE BARRIER SYSTEM, STERILISED USING IRRADIATION (11) UA.TR.116** = Medical device Class IIa under Technical Regulation on medical devices, approved by Decree of Cabinet of Ministers of Ukraine No. 753 on October 2nd., 2013 under conformity assessment procedure by the Notified Body "UNI-CERT" LLC, with reference (UA.TR.116) **(12) MADE WITH NATURAL RUBBER LATEX (13) NOT MADE WITH NATURAL RUBBER LATEX (14) EN 455-3:2022 – CONTAINS POTENTIAL TYPE IV CHEMICAL ALLERGENS (15) DO NOT USE IF PACKAGE IS DAMAGED (16) KEEP DRY (17) KEEP AWAY FROM SUNLIGHT (18) AVOID OZONE (19) LOT NUMBER (20) MANUFACTURING DATE (21) EXPIRY DATE** For more detailed information on the product's performance, please consult Ansell.

C. PRECAUTIONS FOR USE 1. Before usage, inspect the gloves for any defects or imperfections such as holes, pinholes and tears. If the gloves are ripped or punctured during use, dispose of them immediately. If in doubt, do not use the gloves, get a new pair. 2. Do not reverse the gloves. 3. It is essential to keep all chemicals from contact with the skin, even if they are thought to be harmless. Ensure the chemicals cannot enter via the cuff. 4. Contaminated gloves should be cleaned or washed or wiped dry before removal. Avoid touching contaminated surfaces with bare hands. 5. Gloves should not come in contact with a naked flame nor for protection against heat. 6. Gloves shall not be used for protection against ionising radiation nor for use in containment enclosures. 7. Contents of the pouch are sterile unless the pouch is opened or damaged. Do not use if the pouch is opened or damaged. 8. For Medical Use - Lifetime of the Glove – For normal use, following routine inspection, Ansell recommends Surgical gloves be changed every 60-90 minutes. For gloves used with chemotherapy, we recommend they are changed every 30 minutes, or immediately if damaged. **PROPER DONNING & DOFFING:** How to don gloves: 1 - Remove one glove from the package and inspect it to be sure no pinholes or tears are present. 2 - If gloves are

ambidextrous, they can be worn on either hand. If not, align the glove's fingers and thumb with the proper hand before donning. 3 - Insert five fingers into the cuff and pull the cuff over the wrist. 4 - Check for a secure fit around the fingers and palm. The cuff should fit snugly around the wrist. How to doff gloves: 1 - Grasp the outside edge of the glove near the wrist. 2 - Peel the glove away from the hand, turning it inside out. Hold it in the opposite gloved hand. 3 - Slide an ungloved finger under the wrist of the remaining glove, being careful not to touch the outside of the glove. 4 - Peel the remaining glove off from the inside, creating a "bag" containing both gloves. Discard.

D. INGREDIENTS / HAZARDOUS INGREDIENTS: Some gloves might contain ingredients which are known to be a possible cause of allergies in sensitised persons, who may develop irritant and/or allergic contact reactions. If allergic reactions should occur, obtain medical advice immediately. Warning: If gloves contain natural latex, this would be mentioned on the packaging. In that case, **THIS PRODUCT MAY CAUSE ALLERGIC REACTIONS, INCLUDING ANAPHYLACTIC RESPONSES** to sensitised people. For more information, please contact Ansell.

E. CARE INSTRUCTIONS STORAGE: Keep away from direct sunlight; store in a dry place at 25°C to a maximum of 40°C and keep in the original packaging. Keep away from ozone sources. If gloves are properly stored, as indicated above, they won't lose their performances and won't change the glove characteristics significantly. If gloves could be affected by ageing or storage, the expiry date is mentioned on the packaging materials.

CLEANING: The gloves are single-use only and not designed to be laundered nor to be reused. Re-use of the glove could cause post contamination and postoperative complication. Cleaning and re-sterilization has not been validated for this product by the manufacturer. Product integrity may be compromised by any cleaning or re-sterilization process used.

F. DISPOSAL: After use, the product should be incinerated or disposed as per the clinical waste disposal protocol. Used products which have been in contact with chemicals or contaminated with infectious or other hazardous materials should be disposed after each working shift and not reused. They should also be disposed once they show any signs of degradation during usage, such as tearing, holes, discoloration and weakening of the gloves.

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