ME&SU

INSTRUCTIONS FOR USE

ANSELL SINGLE USE GLOVES

Available in 32 languages

v.1 2024-08

Please click on the box below to select your preferred language.

DOMINIQUE DUTSCHER

Europe, Middle East and Africa Region

Ansell Healthcare Europe NV
Boulevard International 55,
Brussels, B-1070, Belgium
T +32 2528 74 00

Ansell Brazil Ltda.

Rua das Figueiras, 474 – Quarto Andar
Bairro Jardim – Santo André,
CEP 09080-300, SP, Brasil
T +55 11 3356 3100

UK Importer / Responsible Person:

Nitritex Ltd
Ground Floor, 15 Kings Court,
Willie Snaith Road, Newmarket,
Suffolk, CB8 7SG, United Kingdom

www.ansell.com

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. © 2024 Ansell Limited. All Rights Reserved.

CONTENTS

Please click on the content below to navigate to your selected page.

DOMINIQUE DUTSCHER

SYMBOLS & PICTOGRAMS

EN ISO 21420:2020

(01)

EN 388:2016
+A1:2018

A B C D E

(02)

EN ISO 374-5:2016

VIRUS

(03)

EN ISO 374-1 : 2016 + A1: 2018
Type A, B or C

ABCDEFGHIJKLMNPOST

(04)

ISO 18889:2019

X

(05)

EN 421:2010

(06)

(07)

(08)

(09)

(10)

(11)

(12)

(13)

(14)

(15)

(16)

(17)

(18)

(19)

(20)

(21)

(22)

(23)

(24)

(25)

(26)

CA XX.XXX

(27)

STERILE R

(28)

LOT

(29)

BATCH

(30)

(31)

(32)

(33)

USE: This Instruction for Use is to be used in combination with the specific information that appears on the first packaging. Gloves are designed as a single use two way barrier protection against cross contamination and also protect the hands mainly against chemical splash risks and comply with the standards shown by the relevant pictograms.

EXPLANATION OF SYMBOLS & PICTOGRAMS THAT MAY APPEAR ON GLOVES/PACKAGING: (01) EN ISO 21420:2020 / ISO 15523-1:2016

– Please read the Instructions for Use, prior to using the products, or contact Ansell for more information. If a level X is mentioned under any of the pictograms, this means this test is not applicable and glove is not designed and therefore not to be used for this specific hazard.

(02) EN 388: 2016 + A1: 2018 – A B C D E – Protection from mechanical risks – A: Abrasion resistance (performance levels 0 to 4) – B: Blade cut resistance (performance levels 0 to 5) If level X is indicated for this property the TDM as per E is the reference performance result for cut resistance – C: Tear resistance (performance levels 0 to 4) – D: Puncture resistance (performance levels 0 to 4) – E: TDM ISO EN 13997 cut resistance (performance levels A to F) **Warning!** The performances (A to E) claimed for the gloves are based on tests performed on the palm area of the gloves only. For gloves with two or more layers, these overall performance levels may not necessarily reflect the performance of the glove's outermost layer. For gloves where the palm and back and cuff are different, mechanical protection is only applicable to the palm of the glove.

(03) EN ISO 374-5: 2016 – Protection against bacteria, fungi and virus. **(04) EN ISO 374-1: 2016 + A1: 2018**

TYPE A, B OR C – Protection against chemical hazards: Type A = chemical breakthrough time > 30 minutes against at least 6 chemicals as per list defined in EN ISO 374-1:2016. Type B = chemical breakthrough time > 30 minutes against at least three chemicals as per list defined in EN ISO 374-1:2016. Type C = chemical breakthrough time > 10 minutes against at least one test chemical as per list defined in EN ISO 374-1:2016 (no code underneath the pictogram) A = methanol, B = acetone, C = acetonitrile, D = dichloromethane, E = carbon disulfide, F = toluene, G = diethylamine, H = tetrahydrofuran, I = ethyl acetate, J = n-heptane, K = sodium hydroxide, 40%, L = sulphuric acid, 96 %, M = nitric acid, 65%, N = acetic acid, 99%, O = ammonia, 25%, P = hydrogen peroxide, 30%, S = hydrofluoric acid, 40%, T = formaldehyde, 37%. **(05) ISO 18889: 2019** – Protection against pesticides – Glove suitable when the potential risk is relatively low. These gloves are not suitable for use with concentrated pesticide formulations and/or for scenarios where mechanical risks exist. The pesticide shall not have the possibility to penetrate between the garment sleeve and the glove. If the overlap is less than approximately 50 mm between the glove and the sleeve, a glove with a longer length should be used. For pesticide resistance data, the duration of the test is not based on actual use time since the permeation test is an accelerated test in which the surface of the specimen is in constant contact with the test chemical. Although the duration of the exposure may be for a longer period during field application with a dilute formulation, the entire surface is not in constant contact with the test chemical. **Warning!** Chemical resistance data has been assessed under laboratory conditions and relates only to the chemical tested. It can be different if used in a mixture. For gloves equal or longer than 400 mm, the chemical resistance data is based from samples taken 80 mm from the end of the cuff. The data may not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. Check that the gloves are suitable for the intended use. Conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the chemical due to changes in physical properties. Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Chemical permeation data, tested per EN 16523-1:2015+A1:2018, and degradation data, tested per EN ISO 374-4:2019, are available upon request and/or via ansell.com, through the Ansell product page/downloads criteria/CE-certified Chemical Permeation test data.

(06) EN 421:2010 – Protection against radioactive contamination. **(07) CE** – Product is compliant and certified to the requirements of the European Regulation on Personal Protective Equipment 2016/425. When MD is mentioned with the CE mark, this is also a Class I product under the Medical Device Regulation 2017/745. **(08) UKCA** – Product is compliant, where applicable, certified to Regulation 2016/425 on Personal Protective Equipment, as amended to apply in GB. When **(09) MD** is mentioned the product is compliant and certified to the requirements of the UK Medical Devices Regulation 2002 (SI 2002 No. 618 as amended) The CE mark is followed by a four digit code which refers to the identification number of the Notified Body that is in charge of the PPE category III conformity assessment (Module C2 or D), for products to protect against serious risks. PPE Type examination certificate (Module B) and Supervised product checks (Module C2) or Conformity to type based on quality assurance of the production process (Module D) by: For EU: Centexbel Belgium (I.D. 0493), Technologiepark 70, B-9052 Zwijnaarde. For Great Britain: Type-examination certificate (Module B) and conformity to type certificates based on Supervised product checks (Module C2) or Conformity to type based on quality assurance of the production process (Module D) for CE marking are used as the basis for applying a UKCA. In some cases, Type-examination certificate (Module B) by Satra Technology Centre, Wyndham Way, Telford Way, Kettering, Northamptonshire, NN16 8SD, UK. For more details, please consult UK Conformity Declaration. For Medical products if the CE and UKCA mark is followed by a four digit code which refers to the identification number of the Notified / Approved Body this means the product is a medical device class I sterile under the European Regulation 2017/745. For EU : I.D. 2797 – BSI Netherland For UK: I.D. 0086 – BSI UK To obtain the EU- or UK Conformity Declaration, please go to: www.ansell.com/regulatory **(10) FOOD PICTOGRAM** – Suitable for contact with foodstuffs. Products carrying this pictogram are in conformity with the European Regulations 1935/2004 and 2023/2006 as well as with all applicable National Regulations for Food-contact materials. **(11) KOHSA** – Product is compliant and certified to the requirements of the Korean Occupational Health & Safety Act legislation for PPE. **(12) UKRSEPRO** – Certified in Ukraine **(13) NON-STERILE (14) NOT MADE WITH NATURAL RUBBER LATEX (15) MADE WITH NATURAL RUBBER LATEX (16) SINGLE USE ONLY (17) EN 455-3:2023** – Contains potential Type IV chemical allergens. **(18) KEEP AWAY FROM SUNLIGHT (19) KEEP DRY (20) AVOID OZONE (21) STORAGE TEMPERATURE BETWEEN 10°C AND 40°C (22) KEEP AWAY FROM HEAT AND RADIOACTIVE SOURCES (23) DO NOT USE IF PACKAGE IS DAMAGED (24) EXPIRY DATE (25) MANUFACTURED DATE (26) MANUFACTURER (27) CA XX.XXX** – Certificate of Approval, as certified to the requirements of the Brazilian Regulation (whereas xx.xxxx refers to the certificate number). **(28) STERILIZED USING IRRADIATION (29) LOT NUMBER (30) BATCH NUMBER (31) ANSELL INTERNAL NUMBER (32) CATALOG NUMBER**

For more detailed information on the product's performance, please consult Ansell. **(33) SIRIM CERTIFICATION MARK** – Product is compliant and certified to the requirements of the Malaysia SIRIM QAS International Sdn. Bhd.

PRECAUTIONS FOR USE: 1. Before usage, inspect the gloves for any defects or imperfections such as holes, pinholes and tears. If the gloves are ripped or punctured during use, dispose of them immediately. If in doubt, do not use the gloves, get a new pair. 2. It is essential to keep all chemicals from contact with the skin, even if they are thought to be harmless. Ensure the chemicals cannot enter via the cuff. Remove the glove immediately if contaminated by a concentrated spill of pesticides. 3. Contaminated gloves should be cleaned or washed or wiped dry before removal. Avoid touching contaminated surfaces with bare hands. 4. Gloves should not come in contact with a naked flame nor used for protection against heat. 5. Gloves shall not be used for protection against ionising radiation nor for use in containment enclosures. 6. Gloves suitable for contact with foodstuffs may show some migration against specific foodstuffs. Please obtain advice from Ansell or consult the Ansell Food Conformity declaration to know if specific restrictions apply and for which specific foodstuffs the gloves can be used. 7. If gloves are being used in explosive environments, please ensure they meet the EN 16350 requirements. Persons wearing these gloves should be properly earthed, e.g. by wearing adequate footwear & clothing. **Warning!** The gloves shall not be unpacked, opened, adjusted or removed whilst in flammable or explosive atmospheres. The electrostatic properties of the gloves might be adversely affected by ageing, wear, contamination and damage and might not be sufficient for oxygen enriched flammable atmospheres where additional assessments are necessary 8. For Sterile Products – Contents of the pouch are sterile unless the pouch is opened or damaged. Do not use if the pouch is opened or damaged. 9. For Medical Use - Lifetime of the Glove – For normal use, following routine inspection, Ansell recommends Examination gloves to be changed every 15 minutes or during patient care if moving from a contaminated body site to another body site and after caring for a patient.

PROPER DRESSING & DOFFING; How to don gloves: 1. Remove one glove from the package and inspect it to be sure no pinholes or tears are present. 2. If gloves are ambidextrous, they can be worn on either hand. If not, align the glove's fingers and thumb with the proper hand before donning. 3. Insert five fingers into the cuff and pull the cuff over the wrist. 4. Check for a secure fit around the fingers and palm. The cuff should fit snugly around the wrist. **How to doff gloves:** 1. Grasp the outside edge of the glove near the wrist. 2. Peel the glove away from the hand, turning it inside out. Hold it in the opposite gloved hand. 3. Slide an ungloved finger under the wrist of the remaining glove, being careful not to touch the outside of the glove. 4. Peel the remaining glove off from the inside, creating a "bag" containing both gloves. Discard.

INGREDIENTS / HAZARDOUS INGREDIENTS: Some gloves might contain ingredients which are known to be a possible cause of allergies in sensitised persons, who may develop irritant and/or allergic contact reactions. If allergic reactions should occur, obtain medical advice immediately. **Warning!** If gloves contain natural latex, this would be mentioned on the packaging. In that case, **THIS PRODUCT MAY CAUSE ALLERGIC REACTIONS, INCLUDING ANAPHYLACTIC RESPONSES** to sensitised people. For more information, please contact Ansell.

CARE INSTRUCTIONS: STORAGE: Keep away from direct sunlight; store in a cool dry place and keep in the original packaging. Keep away from ozone sources. If gloves are properly stored, as indicated above, they won't lose their performances and won't change the glove characteristics significantly. If gloves could be affected by ageing or storage, the expiry date is mentioned on the packaging materials. **CLEANING:** The gloves are single-use only and not designed to be laundered nor to be reused. Re-use of the glove could cause post contamination and postoperative complication. Cleaning and re-sterilization has not been validated for this product by the manufacturer. Product integrity may be compromised by any cleaning or re-sterilization process used. **DISPOSAL:** Used products which have been in contact with chemicals or contaminated with infectious or other hazardous materials such as residual pesticides should be disposed after each working shift and not reused. They should also be disposed once they show any signs of degradation during usage, such as tearing, holes, discoloration and weakening of the gloves. Dispose of according to Local Authority Regulations. Landfill or incinerate under controlled conditions. If product has been used in a clinical setting after use, the product should be incinerated or disposed as per the clinical waste disposal protocol.

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