



# Viral RNA / DNA isolation

## User manual

NucleoMag<sup>®</sup> VET

July 2019 / Rev. 07

**MACHERY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)



DOMINIQUE DUTSCHER SAS

## Contact MN

### Germany and international

#### **MACHEREY-NAGEL GmbH & Co. KG**

Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

Tel.: +49 24 21 969-0

Toll-free: 0800 26 16 000 (Germany only)

Fax: +49 24 21 969-199

E-mail: [info@mn-net.com](mailto:info@mn-net.com)

### Technical Support Bioanalysis

Tel.: +49 24 21 969-270

E-mail: [tech-bio@mn-net.com](mailto:tech-bio@mn-net.com)

## USA

#### **MACHEREY-NAGEL Inc.**

2850 Emrick Blvd. · Bethlehem, PA 18020 · USA

Tel.: +1 484 821 0984

Toll-free: 888 321 6224 (MACH)

Fax: +1 484 821 1272

E-mail: [sales-us@mn-net.com](mailto:sales-us@mn-net.com)

## France

#### **MACHEREY-NAGEL SARL à associé unique**

1, rue Gutenberg · 67722 Hoerd · France

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)

## Switzerland

#### **MACHEREY-NAGEL AG**

Hirsackerstr. 7 · 4702 Oensingen · Switzerland

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)

[www.mn-net.com](http://www.mn-net.com)

DOMINIQUE DUTSCHER SAS

## Table of contents

1 Components	4
1.1 Kit contents	4
1.2 Material to be supplied by user	5
2 Product description	6
2.1 The basic principle	6
2.2 Kit specifications	6
2.3 Magnetic separation systems	6
2.4 Adjusting the shaker settings	7
2.5 Handling of beads	8
2.6 Elution procedures	8
3 Storage conditions and preparation of working solutions	9
4 Safety instructions	10
5 Protocols	13
5.1 Preparation of sample materials	13
5.2 Isolation of viral RNA/DNA and bacterial DNA from blood, tissue homogenates, serum, plasma, other body fluids and washes	14
5.3 Detailed protocol for KingFisher® Flex	18
5.4 Isolation of viral RNA (PRRS) from porcine sperm samples	19
6 Appendix	20
6.1 Troubleshooting	20
6.2 Ordering information	21
6.3 Product use restriction/warranty	22

# 1 Components

## 1.1 Kit contents

NucleoMag® VET		
REF	1 x 96 preps 744200.1	4 x 96 preps 744200.4
NucleoMag® B-Beads	2 x 1.25 mL	10 mL
Lysis Buffer VL1	30 mL	100 mL
Binding Buffer VEB	110 mL	3 x 110 mL
Wash Buffer VEW1	75 mL	300 mL
Wash Buffer VEW2	75 mL	300 mL
Elution Buffer VEL	30 mL	125 mL
Carrier RNA*	400 µg	4 x 400 µg
Carrier RNA Buffer	500 µL	4 x 500 µL
Proteinase K (lyophilized)*	75 mg	3 x 75 mg
Proteinase Buffer PB	8 mL	15 mL
User manual	1	1

\* For preparation of working solutions and storage conditions see section 3.

## 1.2 Material to be supplied by user

Product	REF	Pack of
<b>Separation plate for magnetic beads separation,</b> e.g., Square-well Block (96-well block with 2.1 mL square-wells)	740481 740481.24	4 24
<b>Lysis tubes for incubation of samples and lysis,</b> e.g., Rack of Tubes Strips (1 set consists of 1 Rack, 12 Strips with 8 tubes (1.2 mL wells) each, and 12 Cap Strips)	740477 740477.24	4 sets 24 sets
<b>Elution plate for collecting purified nucleic acids,</b> e.g., Elution Plate U-bottom (96-well 0.3 mL microtiterplate with 300 µL u-bottom wells)	740486.24	24
e.g., Elution Plate Flat-bottom (96-well 0.3 mL microtiterplate with 300 µL flat-bottom wells)	740673	20
<b>For use of kit on KingFisher® Flex instrument:</b> KingFisher® Accessory Kit A (Square-well Blocks, Deep-well Tip Combs, Elution Plates for 4 x 96 NucleoMag® VET preps using KingFisher® Flex platform)	744950	1 set
<b>For use of kit on KingFisher® Duo/Duo Prime instrument:</b> KingFisher® Duo Accessory Kit (KingFisher® Deep-well Blocks, KingFisher® Duo 12 Tip Combs, KingFisher® Duo Elution Strips for 8 x 12 NucleoMag® VET preps using KingFisher® Duo/Duo Prime platform)	744952	1 set

### Reagents:

- 80 % ethanol

## 2 Product description

### 2.1 The basic principle

The **NucleoMag® VET** kit is designed for the isolation of viral DNA or RNA from cell-free body fluids such as serum or plasma, blood or homogenized tissue sample suspensions. This kit provides reagents and magnetic beads for isolation of 96 samples from 100–200 µL. The procedure is based on the reversible adsorption of nucleic acids to paramagnetic beads under appropriate buffer conditions. Sample lysis is achieved by incubation with a Lysis Buffer VL1 containing chaotropic ions supported by Proteinase K digestion. For binding of nucleic acids to the paramagnetic beads, Binding Buffer VEB and the NucleoMag® B-Beads are added to the lysate. After magnetic separation, the paramagnetic beads are washed to remove contaminants and salts using Wash Buffers VEW1 and VEW2 and 80 % ethanol. Residual ethanol from previous wash steps is removed by airdrying. Finally, highly pure viral RNA/DNA is eluted with low salt Elution Buffer VEL or water. Purified viral RNA/DNA can directly be used for downstream applications. The **NucleoMag® VET** kit can be used either manually or automated on standard liquid handling instruments or automated magnetic separators.

### 2.2 Kit specifications

**NucleoMag® VET** is designed for rapid manual and automated small-scale preparation of viral RNA/DNA from cell-free body fluids such as serum or plasma samples, blood samples or homogenized tissue suspensions. The kit is designed for use with NucleoMag® SEP magnetic separator plate (see ordering information) or other magnetic separation systems (see section 2.3). Manual time for the preparation of 96 samples is about 120 minutes. The purified RNA/DNA can be used directly as template for RT-PCR, PCR, or any kind of enzymatic reactions.

**NucleoMag® VET** allows easy automation on common liquid handling instruments or automated magnetic separators. The actual processing time depends on the configuration of the instrument and the magnetic separation system used. Typically, 96 samples can be purified in less than 120 minutes using the NucleoMag® SEP on the automation platform.

The NucleoMag® VET kit is intended for use by professional users such as technicians and physicians experienced and trained in molecular biological techniques including experience with swabs and other potentially infectious, veterinary sample materials.

### 2.3 Magnetic separation systems

For use of **NucleoMag® VET**, the use of the magnetic separator NucleoMag® SEP is recommended. Separation is carried out in a Square-well Block (see ordering information, section 6.2). The kit can also be used with other common separators.

Magnetic separator	Separation plate or tube
NucleoMag® SEP (MN REF 744900)	Square-well Block (MN REF 740481)
Tecan Te-MagS™	1.5 mL tubes without lid (Sarstedt)

### Static magnetic pins

Separators with static magnetic pins, for example, NucleoMag® SEP (for manual use and for use on liquid handling workstations): This type of separator is recommended in combination with a suitable microplate shaker for optimal resuspension of the beads during the washing and elution steps. Alternatively, beads can be resuspended in the buffer by pipetting up and down several times. For fully-automated use on liquid handling workstations, a gripper tool is required, the plate is transferred to the magnetic separator for separation of the beads and transferred to the shaker module for resuspension of the beads.

### Movable magnetic systems

Separators with moving magnetic pins: Magnetic pins/rods are moved from one side of the well to the other and vice versa. Beads follow this movement and are thus pulled through the buffer during the wash and elution steps. Separation takes place when the system stops.

### Automated separators

Separators with moving magnets: Magnetic beads are transferred into suitable plates or tubes. Beads are resuspended from the rod-covered magnets. Following binding, washing or elution beads are collected again with the rod-covered magnets and transferred to the next plate or tube.

## 2.4 Adjusting the shaker settings

When using a plate shaker for the washing and elution steps, the speed settings have to be adjusted carefully for each specific separation plate and shaker to prevent cross-contamination from well to well. Proceed as follows:

### Adjusting shaker speed for binding and wash steps:

- Load 600 µL dyed water to the wells of the separation plate. Place the plate on the shaker and start shaking with a moderate speed setting for 30 seconds. Turn off the shaker and check the plate surface for small droplets of dyed water.
- Increase speed setting, shake for an additional 30 seconds, and check the plate surface for droplets again.
- Continue increasing the speed setting until you observe droplets on top of the separation plate. Reduce speed setting, check again, and use this setting for the washing step.

### Adjusting shaker speed for the elution step:

- Load 100 µL dyed water to the wells of the collection plate and proceed as described above.

## 2.5 Handling of beads

### Distribution of beads

A homogeneous distribution of the magnetic beads to the individual wells of the separation plate is essential for a high well-to-well consistency. Therefore, before distributing the beads, make sure that the beads are completely resuspended. Shake the storage bottle well or place it on a vortexer shortly. Premixing magnetic beads with the binding buffer allows easier homogenous distribution of the beads to the individual wells of the separation plate. During automation, a premix step before aspirating the beads/binding buffer mixture from the reservoir is recommended to keep the beads resuspended.

### Magnetic separation time

Attraction of the magnetic beads to the magnetic pins depends on the magnetic strength of the magnetic pins, the selected separation plate, distance of the separation plate from the magnetic pins, and the volume to be processed. The individual times for complete attraction of the beads to the magnetic pins should be checked and adjusted on each system. It is recommended using the separation plates or tubes specified by the supplier of the magnetic separator.

### Washing the beads

Washing the beads can be achieved by shaking or mixing. In contrast to mixing by pipetting up and down, mixing by shaker or magnetic mixing allows simultaneous mixing of all samples. This reduces the time and number of tips needed for the preparation. Resuspension by pipetting up and down, however, is more efficient than mixing by a shaker or magnetic mix.

Method	Resuspension efficiency	Speed	Number of tips needed
Magnetic mix	+	++	Low
Shaker	++	++	Low
Pipetting	+++	+*	High

## 2.6 Elution procedures

Purified viral RNA/DNA can be eluted directly with the supplied Elution Buffer VEL. Elution can be carried out in a volume of  $\geq 50 \mu\text{L}$ . It is essential to cover the NucleoMag® Beads completely with elution buffer during the elution step. The volume of dispensed elution buffer depends on the magnetic separation system (e.g., the position of the pellet inside the separation plate). For efficient elution, the magnetic bead pellet should be resuspended completely in the elution buffer. For some separators, high elution volumes might be necessary to cover the whole pellet.

\* 8-channel pipetting device



### 3 Storage conditions and preparation of working solutions

*Attention:*

VL1, VEB, VEW1, VEW2 and the Carrier RNA Buffer contain chaotropic salt (e.g. guanidine hydrochloride and/or sodium perchlorate) which can form highly reactive compounds when combined with bleach (sodium hypochlorite)! DO NOT add bleach or acidic solutions directly to the sample preparation waste. Wear suitable protective clothing, gloves and safety goggles!

- All components of the **NucleoMag® VET** kit should be stored at room temperature (18–25 °C) and are stable for up to one year.
- All buffers are delivered ready to use.

Before starting any **NucleoMag® VET** protocol, prepare the following:

- **Proteinase K:** Before first use of the kit, add 3.35 mL Proteinase Buffer PB to each vial of the **lyophilized Proteinase K**. Dissolved Proteinase K solution should be stored in aliquots at -20 °C.
- **Carrier RNA:** Before first use of the kit, add 500 µL Carrier RNA Buffer to each vial **lyophilized Carrier RNA**. Store dissolved Carrier RNA solution in aliquots at -20 °C.






NucleoMag® VET		
REF	1 x 96 preps 744200.1	4 x 96 preps 744200.4
Proteinase K (lyophilized)	1 vial (75 mg) Add 3.35 mL Proteinase Buffer	3 vials (75 mg/vial) Add 3.35 mL Proteinase Buffer to each vial
Carrier RNA (lyophilized)	1 vial (400 µg) Add 500 µL Carrier RNA Buffer	4 vials (400 µg/vial) Add 500 µL Carrier RNA Buffer to each vial

## 4 Safety instructions

The following components of the **NucleoMag® VET** kits contain hazardous contents.

*Wear gloves and goggles and follow the safety instructions given in this section.*

Only harmful features do not need to be labeled with H and P phrases up to 125 mL or 125 g. *Mindergefährliche Eigenschaften müssen bis 125 mL oder 125 g nicht mit H- und P-Sätzen gekennzeichnet werden.*

Component	Hazard contents	GHS symbol	Hazard phrases	Precaution phrases
Inhalt	Gefahrstoff	GHS-Symbol	H-Sätze	P-Sätze
VL1	Guanidine hydrochloride 36–50 % <i>Guanidinhydrochlorid 36–50 %</i> CAS 9003-98-9	 WARNING ACHTUNG	302, 319	264W, 280sh, 301+312, 330
VEB	Sodium perchlorate 15–40 % + ethanol 35–55 % <i>Natriumperchlorat 15–40 % + Ethanol 35–55 %</i> CAS 7601-89-0, 64-17-5	 DANGER GEFAHR	226, 302	210, 264W, 301+312, 330
VEW1, VEW2	Sodium perchlorate 15–40 % + ethanol 20–35 % <i>Natriumperchlorat 15–40 % + Ethanol 20–35 %</i> CAS 7601-89-0, 64-17-5	 WARNING ACHTUNG	226, 302	210, 264W, 301+312, 330
Carrier RNA Buffer	Guanidinium thiocyanate 30–45 % <i>Guanidiniumthiocyanat 30–45 %</i> CAS 593-84-0	 WARNING ACHTUNG	302, 412	264W, 273, 301+312, 330
Proteinase K	Proteinase K, lyophilized 90–100 % <i>Proteinase K, lyophilisiert 90–100 %</i> CAS 39450-01-6	 DANGER GEFAHR	315, 319, 334	261sh, 280sh, 342+311

**Hazard phrases**

- H 226 Flammable liquid and vapour.  
*Flüssigkeit und Dampf entzündbar.*
- H 302 Harmful if swallowed.  
*Gesundheitsschädlich bei Verschlucken.*
- H 315 Causes skin irritation.  
*Verursacht Hautreizungen.*
- H 319 Causes serious eye irritation.  
*Verursacht schwere Augenreizung.*
- H 334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.  
*Kann bei Einatmen Allergie, asthmaartige Symptome oder Atembeschwerden verursachen.*
- H 412 Harmful in contact with skin.  
*Schädlich für Wasserorganismen, mit langfristiger Wirkung.*

**Precaution phrases**

- P 210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
*Von Hitze, heißen Oberflächen, Funken, offenen Flammen sowie anderen Zündquellenarten fernhalten. Nicht rauchen.*
- P 261sh Avoid breathing dust/vapors.  
*Einatmen von Staub/Dampf vermeiden.*
- P 264W Wash with water thoroughly after handling.  
*Nach Gebrauch mit Wasser gründlich waschen.*
- P 273 Avoid release to the environment.  
*Freisetzung in die Umwelt vermeiden.*
- P 280sh Wear protective gloves/eye protection.  
*Schutzhandschuhe/Augenschutz tragen.*
- P 301+312 IF SWALLOWED: Call a POISON CENTER / doctor if you feel unwell.  
*BEI VERSCHLUCKEN: Bei Unwohlsein GIFTINFORMATIONSZENTRUM/Arzt anrufen.*
- P 330 Rinse mouth.  
*Mund ausspülen.*
- P 342+311 If experiencing respiratory symptoms: Call a POISON CENTER/doctor.  
*Bei Symptomen der Atemwege: GIFTINFORMATIONSZENTRUM/Arzt anrufen.*



The symbol shown on labels refers to further safety information in this section.  
*Das auf Etiketten dargestellte Symbol weist auf weitere Sicherheitsinformationen dieses Kapitels hin.*

When working with the NucleoMag® VET kit wear suitable protective clothing (e.g., lab coat, disposable gloves, and protective goggles). For more information consult the appropriate Material Safety Data Sheets (MSDS available online at <http://www.mn-net.com/msds>).

Caution: Guanidin hydrochloride in Lysis Buffer VL1, sodium perchlorate in buffer VEB, VEW1, VEW2 and guanidine thiocyanate in Carrier RNA Buffer can form highly reactive compounds when combined with bleach! Thus, do not add bleach or acidic solutions directly to the sample preparation waste.

The waste generated with the NucleoMag® VET kit has not been tested for residual infectious material. A contamination of the liquid waste with residual infectious material is highly unlikely due to strong denaturing lysis buffer and Proteinase K treatment but it cannot be excluded completely. Therefore, liquid waste must be considered infectious and should be handled and discarded according local safety regulations.

**Disposal**

Dispose hazardous, infectious or biologically contaminated materials in a safe and acceptable manner and in accordance with all local and regulatory requirements.

## 5 Protocols

### 5.1 Preparation of sample materials

#### a) Blood and serum/plasma samples

A sample volume of 100–200 µL blood can be used. Do not use higher volumes. When using less than 200 µL samples, adjust with PBS buffer to 200 µL.

#### b) Tissue samples

Homogenize tissue samples. Typically 5–10 mg sample material can be homogenized in 400 µL PBS buffer using a bead based homogenizer. If necessary, higher amounts of sample material can be used (up to 25 mg). It should be considered that the copurified total nucleic acid may cause inhibition in the subsequent PCR assays. After homogenization of the tissue, centrifuge and use up to 200 µL clear supernatant for the purification protocol. If using less than 200 µL, adjust with PBS buffer to a final volume of 200 µL.

For isolation of viral RNA:

Tissue samples can be also disrupted in a buffer containing chaotropic salt (e.g., Buffer RA1, see ordering information 6.2) and beta-mercaptoethanol or TCEP reducing agent (see ordering information 6.2).

#### c) Swab samples

Incubate the swabs with PBS, sodium chloride, or cell culture medium for 30 min with shaking. Remove and squeeze out the swab. Proceed with 200 µL of the particle-free buffer or medium for purification protocol.

#### d) Feces

Mix 1 volume of feces (e.g., 500 µL) with an equal volume of PBS buffer. Mix vigorously by vortexing for 1 min. Allow the particles to settle down or centrifuge with low speed (e.g., at 500 x g). For difficult-to-lyse bacteria, mechanical disruption (e.g., treatment using suitable glass beads) may be required. Take the supernatant and use 200 µL for the purification protocol.

#### e) TRIzol® lysis

For sample materials such as semen, a TRIzol® lysis may be required. Homogenize 10–30 mg tissue or up to 250 µL blood with 1 mL TRIzol® reagent to manufacturer's instructions. After phase separation by centrifugation, remove aqueous, colorless (upper) phase (approximately 400 µL). For further processing, start with step 2 of the purification protocol by mixing 400 µL of the aqueous phase with 600 µL Buffer VEB and 20 µL NucleoMag® B-Beads.

#### f) Milk samples

Usually a sample volume higher than 200 µL is used. Typically, 1 mL of a normal milk sample is centrifuged (e.g., 11.000 x g for 3 min). Discard the supernatant and resuspend the pellet in 400 µL PBS. Proceed with 200 µL sample input for the purification protocol.

Sour milk samples need an additional pretreatment step. Therefore, incubate the sour milk particles/chunk in an appropriate amount of lysis buffer for 1–3 h at 56 °C (ideally shaking). Pellet residual particles and proceed with 400 µL of the lysate with the step 2 of the purification protocol.

## 5.2 Isolation of viral RNA/DNA and bacterial DNA from blood, tissue homogenates, serum, plasma, other body fluids and washes

### Preparation of sample material





The standard protocol is related to a volume of 200 µL (homogenized) sample. For the preparation of different sample materials (e.g., tissue, swabs, feces), please see the indications at section 5.1.











### Protocol at a glance

- For hardware requirements refer to section 2.3.
- For detailed information on each step see page 16.

#### Before starting the preparation:

- Check that Proteinase K and Carrier RNA were prepared according to section 3.

1	Lyse sample	<b>200 µL (homogenized) sample</b> <b>20 µL Proteinase K</b> <b>4 µL Carrier RNA</b> <b>180 µL VL1</b>  <b>Mix</b>  <b>RT, 15 min</b>	
2	Bind nucleic acid to NucleoMag® B-Beads	<b>20 µL B-Beads</b> <b>600 µL VEB</b>  <b>Mix by shaking for 5–10 min at RT</b> <i>(Optional: Mix by pipetting up and down)</i>  <b>Remove supernatant after 2 min separation</b>	    

3	Wash with VEW1	Remove Square-well Block from NucleoMag® SEP	
		600 µL VEW1	
		Resuspend: Shake 1 min at RT	
4	Wash with VEW2	Remove Square-well Block from NucleoMag® SEP	
		600 µL VEW2	
		Resuspend: Shake 1 min at RT	
5	Wash with 80 % ethanol	Remove Square-well Block from NucleoMag® SEP	
		600 µL 80 % ethanol	
		Resuspend: Shake 1 min at RT	
6	Air-dry magnetic beads	Air-dry 10 min at RT	
7	Elute RNA / DNA	Remove Square-well Block from NucleoMag® SEP	
		50-100 µL VEL	
		Shake 5 min at RT <i>(Optional: Mix by pipetting up and down)</i>	
		Separate 2 min and transfer RNA / DNA into elution plate / tubes	

## Detailed protocol

This protocol is designed for magnetic separators with static pins (e.g., NucleoMag® SEP) and suitable plate shakers. It is recommended using a Square-well Block for separation (see ordering information, section 6.2). Alternatively, isolation of RNA/DNA can be performed in reaction tubes with suitable magnetic separators. This protocol is for manual use and serves as a guideline for adapting the kit to robotic instruments.

### 1 Lyse sample

Pre-dispense **20 µL Proteinase K** and **200 µL of sample** to a suitable reaction tube. Add **180 µL Buffer VL1** to the reaction tube. Optional: add **4 µL of the Carrier RNA** stock solution to the reaction tube. Mix well by repeated pipetting up and down and incubate at **room temperature** for **15 min** with shaking. Alternatively, lysis step can be performed in Tube Strips (see ordering information).

Following the lysis incubation, spin down to collect any sample from the lysis tube lids and transfer each lysate to the wells of a Square-well Block.

### 2 Bind nucleic acid to magnetic beads

Add **20 µL resuspended B-Beads** and **600 µL Buffer VEB** to the lysed sample.

Mix by pipetting up and down 6 times and **shake** for **5 min** at **room temperature**. Alternatively, when processing the kit without a shaker, pipette up and down 10 times and incubate for 5 min at room temperature.

NucleoMag® B-Beads and Buffer VEB can be premixed.

*Be sure to resuspend the NucleoMag® B-Beads before removing them from the storage bottle. Vortex storage bottle briefly until a homogenous suspension has been formed.*

Separate the magnetic beads against the side of the wells by placing the Square-well Block on the NucleoMag® SEP a magnetic separator. Wait at least **2 min** until all the beads have been attracted to the magnets. Remove and discard supernatant by pipetting.

*Do not disturb the attracted beads while aspirating the supernatant.*

### 3 Wash with VEW1

Remove the Square-well Block from the NucleoMag® SEP magnetic separator. Add **600 µL Buffer VEW1** and resuspend the beads by shaking until the beads are resuspended completely (1–3 min). Alternatively, resuspend beads completely by repeated pipetting up and down.

Separate the magnetic beads by placing the Square-well Block on the NucleoMag® SEP magnetic separator. Wait at least **2 min** until all the beads have been attracted to the magnet. Remove and discard supernatant by pipetting.



**4 Wash with VEW2**

Remove the Square-well Block from the NucleoMag® SEP magnetic separator. Add **600 µL Buffer VEW2** and resuspend the beads by shaking until the beads are resuspended completely (1–3 min). Alternatively, resuspend beads completely by repeated pipetting up and down.

Separate the magnetic beads by placing the Square-well Block on the NucleoMag® SEP magnetic separator. Wait at least **2 min** until all the beads have been attracted to the magnet. Remove and discard supernatant by pipetting.

---

**5 Wash with 80 % ethanol**

Remove the Square-well Block from the NucleoMag® SEP magnetic separator. Add **600 µL 80 % ethanol** and resuspend the beads by shaking until the beads are resuspended completely (1–3 min). Alternatively, resuspend beads completely by repeated pipetting up and down.

Separate the magnetic beads by placing the Square-well Block on the NucleoMag® SEP magnetic separator. Wait at least **2 min** until all the beads have been attracted to the magnet. Remove and discard supernatant by pipetting.

---

**6 Air-dry magnetic beads**

Air-dry the magnetic bead pellet for **10 min at room temperature**.

---

**7 Elute RNA / DNA**

Remove the Square-well Block from the NucleoMag® SEP magnetic separator. Add desired volume of **Buffer VEL (50–100 µL)** to each well of the Square-well Block and resuspend the beads by shaking **5 min at room temperature**. Alternatively, resuspend beads completely by repeated pipetting up and down and incubate for **5 min at 56 °C**.

Separate the magnetic beads by placing the Square-well Block on the NucleoMag® SEP magnetic separator. Wait at least **2 min** until all the beads have been attracted to the magnets. Transfer the supernatant containing the purified nucleic acids to either elution plates or tube strips (see ordering information).

---

### 5.3 Detailed protocol for KingFisher® Flex

*Note:* The required method file 'NucleoMag® VET' for the instrument is available at Technical Support Bioanalysis (tech-bio@mn-net.com).

**Important:** Always prepare deep-well block with samples first and add reagents exactly in the order as given below.

#### Before starting the preparation:

- Check that Proteinase K and Carrier RNA were prepared according to section 3.
- KingFisher® Accessory Kit A (see ordering information)

#### 1 Prepare sample /lysis plate (part I)

Dispense **20 µL Proteinase K solution** to each well of the 96-well deep-well block. Add **200 µL blood sample /homogenized tissue sample** to each well of the 96-well deep-well block, mix by pipetting up and down. Add **180 µL Buffer VL1** and mix by pipetting up and down 3 times.

*Optional:* Shake at 1,000 rpm for 15 min at room temperature.

Continue with the preparation of the wash and elution plates before adding magnetic beads and binding buffer to the sample plate.

#### 2 Prepare wash and elution plates

Wash plates:

Fill **600 µL Buffer VEW1** to each well of an empty Thermo 96-well deep well plate.

Fill **600 µL Buffer VEW2** to each well of an empty Thermo 96-well deep well plate.

Fill **600 µL 80 % ethanol** to each well of an empty Thermo 96-well deep well plate.

Elution plate:

Fill **100 µL Buffer VEL** to each well of an empty Thermo 200 µL 96-well plate.

#### 3 Prepare sample /lysis plate (part II)

Add **20 µL B-Beads** and **575 µL buffer VEB** to each well of the sample /lysis plate.

**4 Run purification protocol on instrument**

Start the isolation of nucleic acids on the KingFisher® Flex instrument.

Start the method file 'NucleoMag® VET'.

Insert plates as indicated on the KingFisher® instrument display.

Method starts with a mixing step (combined lysis and binding step) after setting up the last plate to the instrument.

---

**5 Remove eluted nucleic acids**

The instrument stops after the final elution step. Follow the instructions on the instrument's display and unload the plates from the instrument.

Purified RNA/DNA can be used for further PCR based analysis.

---

**5.4 Isolation of viral RNA (PRRS) from porcine sperm samples**

**Before starting the preparation:**

- The additional Lysis Buffer RA1 is required (see ordering information 6.2).
- 

**1 Precipitation**

Centrifuge 1 mL sperm sample for 4 min at 12.000 x *g*.

Discard the supernatant after centrifugation.

---

**2 Lyse sample**

Add 400 Lysis Buffer RA1 and mix by pipetting.

Incubate for 10 min at 70° C.

---

**3 Clear lysate**

Centrifuge the lysed sample for 1 min. at 15.000 x *g*

Use 400 µL of the cleared lysate and proceed with step 2 of standard protocol, see section 5.2.

---

## 6 Appendix

### 6.1 Troubleshooting

Problem	Possible cause and suggestions
Poor yield / low sensitivity	<i>Insufficient elution buffer volume</i>
	<ul style="list-style-type: none"> <li>• Beads pellet must be covered completely with elution buffer.</li> </ul>
	<i>Insufficient performance of elution buffer during elution step</i>
	<ul style="list-style-type: none"> <li>• Remove residual buffers during the separation steps completely. Remaining buffers decrease the efficiency of following wash and elution steps.</li> </ul>
Poor yield / low sensitivity	<i>Beads dried out</i>
	<ul style="list-style-type: none"> <li>• Do not let the beads dry as this might result in lower elution efficiencies.</li> </ul>
	<i>Aspiration of attracted bead pellet</i>
Poor yield / low sensitivity	<ul style="list-style-type: none"> <li>• Do not disturb the attracted beads while aspirating the supernatant, especially when the magnetic bead pellet is not visible in the lysate.</li> </ul>
	<i>Aspiration and loss of beads</i>
	<ul style="list-style-type: none"> <li>• Time for magnetic separation too short or aspiration speed too high.</li> </ul>
Low purity / low sensitivity	<i>Insufficient washing procedure</i>
	<ul style="list-style-type: none"> <li>• Use only the appropriate combinations of separator and plate, for example, Square-well Block in combination with NucleoMag® SEP.</li> <li>• Make sure that beads are resuspended completely during the washing procedure. If shaking is not sufficient to resuspend the beads completely mix by repeated pipetting up and down.</li> </ul>
Poor performance of nucleic acids in downstream applications	<i>Carry-over of ethanol from wash buffers</i>
	<ul style="list-style-type: none"> <li>• Be sure to remove all of the 80 % ethanolic wash solution from the final wash, as residual ethanol interferes with downstream applications.</li> </ul>
Poor performance of nucleic acids in downstream applications	<i>Ethanol evaporation from wash buffers</i>
	<ul style="list-style-type: none"> <li>• Close buffer bottles tightly, avoid ethanol evaporation from buffer bottles as well as from buffer filled in reservoirs. Do not reuse buffers from buffer reservoirs.</li> </ul>

*Time for magnetic separation too short*

Carry-over of beads

- Increase separation time to allow the beads to be completely attracted to the magnetic pins before aspirating any liquid from the well.

*Aspiration speed too high (elution step)*

- High aspiration speed during the elution step may cause bead carry-over. Reduce aspiration speed for elution step.

**6.2 Ordering information**

Product	REF	Pack of
NucleoMag® VET	744200.1 744200.4	1 x 96 preps 4 x 96 preps
NucleoMag® SEP	744900	1
Square-well Blocks	740481 740481.24	4 24
Self adhering PE Foil	740676	50 sheets
Rack of Tube Strips (set consists of 1 Rack, 12 Tube Strips with 8 tubes each, and 12 Cap Strips)	740477 740477.24	4 sets 24 sets
Elution Plate U-bottom	740486.24	24
KingFisher® Accessory Kit A (set consists of Square-well Blocks, Deep-well Tip Combs, set consists of Elution Plates for 4 x 96 NucleoMag® VET preps using KingFisher® Flex platform)	744950	1 set
KingFisher® Duo Accessory Kit (set consists of KingFisher® Deep-well Blocks, KingFisher® Duo 12 Tip Combs, KingFisher® Duo Elution Strips for 8 x 12 NucleoMag® VET preps using KingFisher® Duo/Duo Prime platform)	744952	1 set
Buffer RA1 (60 mL)	740961	60 mL
Reducing Agent TCEP	740395.107	107 mg

Visit [www.mn-net.com](http://www.mn-net.com) for more detailed product information.

### 6.3 Product use restriction/warranty

**NucleoMag® VET** kit components are intended, developed, designed, and sold FOR RESEARCH PURPOSES ONLY, except, however, any other function of the product being expressly described in original MACHEREY-NAGEL product leaflets.

MACHEREY-NAGEL products are intended for GENERAL LABORATORY USE ONLY! MACHEREY-NAGEL products are suited for QUALIFIED PERSONNEL ONLY! MACHEREY-NAGEL products shall in any event only be used wearing adequate PROTECTIVE CLOTHING. For detailed information please refer to the respective Material Safety Data Sheet of the product! MACHEREY-NAGEL products shall exclusively be used in an ADEQUATE TEST ENVIRONMENT. MACHEREY-NAGEL does not assume any responsibility for damages due to improper application of our products in other fields of application. Application on the human body is STRICTLY FORBIDDEN. The respective user is liable for any and all damages resulting from such application.

DNA/RNA/PROTEIN purification products of MACHEREY-NAGEL are suitable for IN VITRO-USES ONLY!

ONLY MACHEREY-NAGEL products specially labeled as IVD are also suitable for IN VITRO-diagnostic use. Please pay attention to the package of the product. IN VITRO-diagnostic products are expressly marked as IVD on the packaging.

IF THERE IS NO IVD SIGN, THE PRODUCT SHALL NOT BE SUITABLE FOR IN VITRO-DIAGNOSTIC USE!

ALL OTHER PRODUCTS NOT LABELED AS IVD ARE NOT SUITED FOR ANY CLINICAL USE (INCLUDING, BUT NOT LIMITED TO DIAGNOSTIC, THERAPEUTIC AND/OR PROGNOSTIC USE).

No claim or representations is intended for its use to identify any specific organism or for clinical use (included, but not limited to diagnostic, prognostic, therapeutic, or blood banking). It is rather in the responsibility of the user or - in any case of resale of the products - in the responsibility of the reseller to inspect and assure the use of the DNA/RNA/protein purification products of MACHEREY-NAGEL for a well-defined and specific application.

MACHEREY-NAGEL shall only be responsible for the product specifications and the performance range of MN products according to the specifications of in-house quality control, product documentation and marketing material.

This MACHEREY-NAGEL product is shipped with documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. MACHEREY-NAGEL's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Supplementary reference is made to the general business terms and conditions of MACHEREY-NAGEL, which are printed on the price list. Please contact us if you wish to get an extra copy.

There is no warranty for and MACHEREY-NAGEL is not liable for damages or defects arising in shipping and handling (transport insurance for customers excluded), or out of accident or improper or abnormal use of this product; defects in products or

components not manufactured by MACHEREY-NAGEL, or damages resulting from such non-MACHEREY-NAGEL components or products.

MACHEREY-NAGEL makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, REPRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO MACHEREY-NAGEL PRODUCTS.

In no event shall MACHEREY-NAGEL be liable for claims for any other damages, whether direct, indirect, incidental, compensatory, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of MACHEREY-NAGEL products to perform in accordance with the stated specifications. This warranty is exclusive and MACHEREY-NAGEL makes no other warranty expressed or implied.

The warranty provided herein and the data, specifications and descriptions of this MACHEREY-NAGEL product appearing in MACHEREY-NAGEL published catalogues and product literature are MACHEREY-NAGEL's sole representations concerning the product and warranty. No other statements or representations, written or oral, by MACHEREY-NAGEL's employees, agent or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

Product claims are subject to change. Therefore please contact our Technical Service Team for the most up-to-date information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Applications mentioned in MACHEREY-NAGEL literature are provided for informational purposes only. MACHEREY-NAGEL does not warrant that all applications have been tested in MACHEREY-NAGEL laboratories using MACHEREY-NAGEL products. MACHEREY-NAGEL does not warrant the correctness of any of those applications.

Last updated: 07/2010, Rev. 03

Please contact:

MACHEREY-NAGEL GmbH & Co. KG

Tel.: +49 24 21 969-270

tech-bio@mn-net.com

---

Trademarks:


KingFisher is a registered trademark of Thermo Fisher Scientific

NucleoMag is a registered trademark of MACHEREY-NAGEL GmbH & Co KG

Te-MagS is a trademark of Tecan Group Ltd., Switzerland

TRIzol is a registered trademark of Molecular Research Center, Inc.

All used names and denotations can be brands, trademarks, or registered labels of their respective owner – also if they are not special denotation. To mention products and brands is only a kind of information (i.e., it does not offend against trademarks and brands and can not be seen as a kind of recommendation or assessment). Regarding these products or services we can not grant any guarantees regarding selection, efficiency, or operation.



**Plasmid DNA**

**Clean up**

**RNA**

**DNA**

**Viral RNA and DNA**

**Protein**

**High throughput**

**Accessories**

**Auxiliary tools**

DOMINIQUE DUTSCHER SAS





[www.mn-net.com](http://www.mn-net.com)

**MACHEREY-NAGEL**



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

DE/International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: [info@mn-net.com](mailto:info@mn-net.com)

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)

US:

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: [sales-us@mn-net.com](mailto:sales-us@mn-net.com)

A044359/0801