


## Instructions for Use

CELLlevator™ - Stacking Device for CELLdisc™ (Item-No. 878071)



## Intended Use

The CELLlevator™ facilitates easy and secure stacking of two CELLdisc™ devices.

## Product Details

The CELLlevator™ consists of three identical segments with indentations and centering ribs which assembled form a triangle.

The corrugated side of the stripe with the Geiner Bio-One logo must always be on the outside of each part and the final stacking device.

→ Assembly guide on the reverse side →

Feature	Properties
Material	Polypropylene
Max. Loading Capacity	8 kg (refers to the assembled CELLlevator™)
Autoclavability	Yes (120 °C, 2 bar), max. 3 times*
TP-range	4°C-37°C
Sterility	No

\* Autoclaving the CELLlevator™ more than 3 times is not recommended. After each auto-clave cycle the integrity of the segments has to be analysed by the user.

For further information please visit our website [www.gbo.com](http://www.gbo.com) or contact us.


### Greiner Bio-One:

<b>Austria</b> office@at.gbo.com	<b>Belgium</b> info@be.gbo.com	<b>Brazil</b> office@br.gbo.com	<b>China</b> office@cn.gbo.com	<b>France</b> accueil@gbo.com	<b>Germany</b> info@de.gbo.com	<b>Hungary</b> office@hu.gbo.com
<b>Italy</b> office@it.gbo.com	<b>Japan</b> info@jp.gbo.com	<b>Netherlands</b> info@nl.gbo.com	<b>Portugal</b> info@vacuette.pt	<b>Spain</b> info@es.gbo.com	<b>UK</b> info@uk.gbo.com	<b>USA</b> office@us.gbo.com




greiner bio-one


## CELLlevator™ - Assembling Guide



**1** Join segment (1) with segment (2) using the indentations. Centering ribs at the side of the indentation guide the segments together.



**2** Join segment (3) accordingly with the indentation of segment (1) from above.



**3** Lift segment (2) and insert into the remaining indentation from segment (3). Slight bending of the segments during the insertion is acceptable.



**4** The final assembled CELLlevator™ should form a uniform triangle.



**5** The CELLlevator™ must be positioned within the circumferential ring of the top plate of the CELldisc™.



**6** The other CELldisc™ can now be positioned on top of the CELLlevator™. Ventilation of the lower CELldisc™ is assured by the notches on the segments whereas planarity of the upper CELldisc™ is assured by the elevations on the segments.



A video animation of the assembling process is available on our website:  
<https://shop.gbo.com/en/germany/products/bioscience/cell-culture-products/celldisc/878071.html>