1.10044.0001

MQuant® Cyanide Test CN⁻

1. Method

Cyanide ions react with a chlorinating agent to form cyanogen chloride, which in turn reacts with 1,3-dimethylbarbituric acid in the presence of pyridine to form a violet dye (König reaction). The cyanide concentration is measured **semi-quantitatively** by visual comparison of the reaction zone of the test strip with the fields of a color scale.

2. Measuring range and number of determinations

Measuring range / color-scale graduation	Number of determinations
1 - 3 - 10 - 30 mg/l CN ⁻	100

3. Applications

This test measures only cyanide ions (free cyanide).

Sample material:

Wastewater, especially from the electroplating industry

4. Influence of foreign substances

This was checked individually in solutions with 3 and 0 mg/l CN⁻. The determination is not yet interfered with up to the concentrations of foreign substances given in the table. Cumulative effects were not checked; such effects can, however, not be excluded.

Concentrations of foreign substances in mg/l			
Ag⁺	1	Fe²⁺	1000
Al³⁺	1000	Fe³⁺	1000
Ba²⁺	1000	Hg⁺	1
Br⁻	5	Hg²⁺	1
Ca²⁺	1000	I⁻	5
Cd²⁺	1000	K⁺	1000
Cl⁻	1000	Mg²⁺	1000
Co²⁺	1000	MnO₄⁻	50
CrO₄²⁻	50	Na⁺	1000
Cu²⁺	1	NH₄⁺	1000
		Ni²⁺	1000
		NO₂⁻	50
		NO₃⁻	1000
		Pb²⁺	1000
		PO₄³⁻	1000
		S²⁻	100
		SCN⁻	1
		SO₄²⁻	1000
		Zn²⁺	1000

5. Reagents and auxiliaries

Please note the warnings on the packaging materials!

The test strips and test reagents are stable up to the date stated on the pack when stored closed at +15 to +25 °C.

Package contents:

Tube containing 100 test strips
1 bottle of reagent CN-1
1 bottle of reagent CN-2
1 red dosing spoon
1 test vessel

Other reagents:

MQuant® Universal indicator strips pH 0 - 14, Cat. No. 1.09535
Sodium hydroxide solution 1 mol/l Titripur®, Cat. No. 1.09137
Sulfuric acid 0.5 mol/l Titripur®, Cat. No. 1.09072
Cyanide standard solution CRM Certipur®, 1000 mg/l CN⁻, Cat. No. 1.04695

6. Preparation

- Samples containing more than 30 mg/l CN⁻ must be diluted with distilled water.
- **The pH must be within the range 6 - 7.** Adjust, if necessary, with sodium hydroxide solution or sulfuric acid.

7. Procedure

Rinse the test vessel several times with the pre-treated sample.

Pretreated sample (15 - 30 °C)	5 ml	Fill the test vessel to the 5-ml mark.
Reagent CN-1	1 level red dosing spoon	Add and dissolve by swirling.
Reagent CN-2	5 drops ¹⁾	Add and swirl.

Immediately immerse the reaction zone of the test strip in the measurement sample **for 30 sec.**

Remove the strip, allow excess liquid to run off via the long edge of the strip onto an absorbent paper towel, and, **within 10 sec.**, determine with which color field on the label the color of the reaction zone coincides most exactly.

Read off the corresponding result in mg/l CN⁻.

¹⁾ **Hold the bottle vertically while adding the reagent!**

Notes on the measurement:

- The color of the reaction zone may continue to change after the specified reaction time has elapsed. This must not be considered in the measurement.
- If the color of the reaction zone is equal to or more intense than the darkest color on the scale, repeat the measurement using **fresh**, diluted samples until a value of less than 30 mg/l CN⁻ is obtained.

Concerning the result of the analysis, the dilution (see also section 6) must be taken into account:

$$\text{Result of analysis} = \text{measurement value} \times \text{dilution factor}$$

8. Method control

To check test strips, test reagents, and handling: Dilute the cyanide standard solution with distilled water to 10 mg/l CN⁻ and analyze as described in section 7.

Additional notes see under

www.sigmaaldrich.com/qa-test-kits

9. Notes

- **Reclose** the reagent bottles and **the tube containing the test strips immediately after use.**
- **Dispose of chemical waste in accordance with the local regulations.**
- Rinse the test vessel **with distilled water only.**

The life science business of Merck operates as MilliporeSigma in the U.S. and Canada.

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. Merck, Supelco, Sigma-Aldrich, and MQuant are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

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