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# 1. Safety

## 1.1 General

- These instructions facilitate the safe and efficient handling of the equipment.
- Before operating the device for the first time the operator must have read and understood the operating manual.
- Adherence to all the specified safety instructions and instructions for actions in this operating manual is a fundamental requirement for working safely.
- In addition, all local accident prevention guidelines and general safety regulations must be observed.
- Store these operating instructions in a well accessible location.

## 1.2 Safety Symbols

 <b>DANGER</b>	DANGER - Indicates an imminently hazardous situation, which, if not avoided, will result in death.
 <b>WARNING</b>	WARNING - Indicates a potentially hazardous situation, which, if not avoided, can result in death or serious injury.
 <b>CAUTION</b>	CAUTION - Indicates a potentially hazardous situation, which, if not avoided, can result in minor or moderate injury.
<b>NOTICE</b>	NOTICE – Is used to address practices not related to physical injury.
	Safety alert symbol is a symbol that indicates a general hazard. Failure to follow these instructions could lead to personal injury or damage to equipment.
	This symbol indicates a potential danger caused by hot surface.
	This symbol refers to jeopardizing persons by electric power.

### 1.3 General cause of risks

#### Product- and media safety

 <b>DANGER</b>	<p><b>1. Explosion hazard!</b> <b>Danger of death due to explosion!</b></p> <ul style="list-style-type: none"><li>• Do not use the device in an unsafe environment and especially not in explosive atmospheres. It is not EX-protected.</li><li>• With substances capable of forming an explosive mixture, appropriate safety measures must be applied, e.g. working under a fume hood.</li></ul>
 <b>WARNING</b>	<p><b>1. Fire Hazard!</b> <b>Burns caused by fire!</b></p> <ul style="list-style-type: none"><li>• The max. working temperature may not be over or just under the flash point of the media. The hotplate surface can act as an ignition source.</li><li>• Beware of hazards due to:<ul style="list-style-type: none"><li>- chemical reactions</li><li>- flammable materials</li><li>- combustible media with a low boiling temperature.</li></ul></li></ul> <p><b>2. Risk of accident and health hazard!</b> <b>Risk of injury caused by hazardous substances!</b></p> <ul style="list-style-type: none"><li>• All relevant national regulations regarding work safety and accident prevention must be observed when working with hazardous, toxic or pathogenic substances.</li></ul>
 <b>CAUTION</b>	<p><b>1. Hot surfaces!</b> <b>Contact with the skin causes burns!</b></p> <ul style="list-style-type: none"><li>• Do not touch the heating surface or upper housing parts during operation.</li><li>• Pay attention to the residual heat on the heating plate after switching off.</li><li>• Do not directly touch hot media, substances or objects.</li><li>• Avoid unsafe standing position of vessels or incorrect container size.</li></ul> <p><b>2. Compromise against to the device safety!</b> <b>Danger due to functional impairment!</b></p> <ul style="list-style-type: none"><li>• Avoid vapors and gases that might attack electrical/ mechanical components of the device or use a properly working suction system.</li></ul>

## Electrical Safety

**Dangerous electrical voltages!**  
**Contact with dangerous voltages can result in life threatening electric shock.**

- Do not use the device if:
  - it is damaged
  - cable is damaged.
- Use the device only with the protective grounding conductor connected.
- Never insert objects through slots and vent holes.
- Liquids must not be poured directly onto heating surface.
- Arrange the cables so that the electrical wires do not come in to contact with high-temperature parts.
- Only a qualified electrician may open the device.

### **1.4 Intended use**

The hotplate can be used for heating up **non-explosive** materials, substances and objects.

When using flammable substances and mixtures safety precautions need to be taken into account, e.g working in a fume cupboard. Liquids must not be poured directly onto heating surface.

Use the device only as described in these Operating Instructions. Any other usage is considered to be improper and may lead to material damages and even personal injuries. Any selfconversion and alteration undertaken on the device is not permitted. The manufacturer does not bear any liability for damages caused due to improper use.

Use the device in indoor environments like that of a laboratory of research, teaching, industry or trade.

## **2. Operating the device**

### **2.1 Set up**

The devices have been designed for continuous operation, even at the maximum temperature and during continuous operation, the housing will heat up only to a slight extent. This allows them to be set up also on heat-sensitive surfaces. The area needs to be even, solid, nonslip and clean. Stable, horizontal placement can be achieved by rotating the appliance feet.

**Caution: Electric shock hazard!** Unwind the main cable completely. Make sure that the line is safe and do not touch any hot parts of the device.

Compare specified voltage on hotplate label with available mains voltage. The 230V units are ready for connection up to 3300W. Only earthed connectors with protective contact are permissible.

The electrical connection of devices with more than 3300W power input must only be carried out by a certified electrician. The 400V-devices have a 5-pole cable, whereby wire colours are yellow/green basically the protective conductor, and blue basically the wire for the neutral conductor. The neutral conductor is needed for the control voltage of the electronic (230V). The remaining three cores can be connected arbitrarily to the 400V network.

### **2.2 Switching on**

The devices have a single bottom operation. As soon as the bottom has left the zero marking, the yellow control lamp illuminates.

### **2.3 Select heating power**

An energy regulator switches the heating element on and off in controlled periodical intervals. The numbers on the scale indicate the power in %-Watt and should help to find established values once more.

**Notice:** Odor and smoke formation is typical for new devices.

- Please use device under a fume hood or next to an open window and heat up at maximum for about one hour for first time use.

### **2.4 Setting of desired temperature (except HD0-HD3, SD7-SD9)**

The required temperature is adjusted with the °C knob. As well as the thermostat, the devices are equipped with an infinitely adjustable energy regulator. This enables the heating power to be adapted to the heat requirement. Reduce the output to the minimum required for the application. Excessive heating power results in overheating - despite a switched-off thermostat (automatic procedure) - furthermore it causes a large control bandwidth.

### **2.5 Automatic safety cutoff**

These heating devices are equipped with an additional safety temperature limiter; this switches the device permanently OFF at all poles if the safety

temperature is exceeded (control range 30-110°C at approx. 130-140°C, control range of 50-300°C + 100-370°C at approx. 390°C). The protection can only be reactivated after cooling down and troubleshooting. This is done by pushing the red button on the rear of the device (next to the power cable).

## 2.6 Switching off

**Caution: Residual heat!** Please note, that hotplate will remain hot after turning it off.

The unit is switched off by turning the power controller all the way to the zero position. The yellow control lamp goes out.


## 3. Temperature behaviour of sandbaths

The data shown below are guidelines and depend on the quality and graining of the sand.

Since sand is a poor heat conductor, warming towards the surface of the sand happens very slowly. This measurement was made after 90 minutes of heating up.

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As a rule of thumb: 10 % fall in temperature per 5 mm sand.

## 4. Failures

In case of failures or damages, disconnect the device from the mains immediately and allow it to cool.

### Error messages

Error	Cause	Correction
The device heats up, but the set temperature is not reached.	Heating power is set too low.	Increase heating power.
The device does not heat.	Heating element or electronics are defective.	Have the unit repaired by the manufacturer or replace unit by a qualified electrician.
The safety temperature limiter is released.	Contact or thermostat is defective.	Replace unit by a qualified electrician or have the unit repaired by the manufacturer.

### Repair advice to the qualified electrician

**Caution: Danger of fire and of etching!** The capillary tubes must not be damaged, the min. bending radius is 5 mm.

If you experience an error that is not listed in this table, please contact your supplier. Please do not try to eliminate it by yourself. Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

## 4. Maintenance and cleaning

**Caution: Risk of electric shock and risk of burns!** Disconnect device from main power and let it cool down before starting to do maintenance work or cleaning.


Apart from occasional cleaning there is no maintenance necessary. Do not allow moisture to enter the unit during cleaning.


## 5. Disposal



**The meaning of the symbol of the crossed-out bin:** For a sustainable and professional disposal of electronic devices the user is obligated to dispose legacy equipment separately from general waste. The device marked with the above symbol must be given to an appropriate facility for recycling.

In case of doubt, obtain information on environmentally sound disposal from the local municipal authority or special disposal companies.

Ersatzteile	Best.-Nr./ Order No.	Spare parts
Energierегler (alle Typen) Schaltschütz (außer HB 110 + 300) Kontrolleuchte grün Kontrolleuchte gelb	EGO 50.57 AB 2-3 KL 9 gr KL 9 ge	Energy regulator (all types) Relay (except HB 110 + 300) Control lamp, green Control lamp, yellow
Thermostate Temperaturbereich 30...110 °C Temperaturbereich 50...300 °C Temperaturbereich 100..370 °C Begrenzer HD..., HE 1, SD...	T 55.10/110 T 55.10/300 T 55.04/370 T 55.13/370	Thermostats Temp. range 30...110°C Temp. range 50...300°C Temp. range 100..370°C Limiter HD .., HE 1, SD..
Heizung incl. Isolierung für HD 0 - HT 03 HD 1 - HT 13, SD 7 - ST 73 HE 1- TH 13 HD 2 - HT 23, SD 8 - ST 83 HD 3 - HT 33, SD 9 - ST 93 HB 110 HB 300	HZ 3030 HZ 3535 HZ 4429 HZ 5035 HZ 5843 HBZ 1.0 HBZ 2.0	Heating incl. Insulation for HD 0 - HT 03 HD 1 - HT 13, SD 7 - ST 73 HE 1 - TH 13 HD 2 - HT 23, SD 8 - ST 83 HD 3 - HT 33, SD 9 - ST 93 HB 110 HB 300
<b>Bei Ersatzteilbestellung bitte Typ und Serien-Nr. des Gerätes an- geben.</b>		<b>Please state serial number and type of the device in any spare order.</b>

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**Harry Gestigkeit GmbH**