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# Safety Information for Vacuum Equipment

Hazards described in this document are graded by use of the following formatting:

- **Danger!** Hazardous situation which, if not avoided, will result in death or serious injury as well as in damage to equipment and environment.

- ▶ **Warning!** Hazardous situation which, if not avoided, could result in death or serious injury as well as in damage to equipment and environment.

- ▷ **Caution!** Hazardous situation which, if not avoided, could result in minor or moderate injury as well as in damage to equipment.

Notice. Not related to personal injury. Disregarding could result in damage to equipment.

## 1. General information:

- ▶ Do not use any damaged or malfunctioning equipment.

- ▶ Do not modify the equipment without authorization. Please read the manual and observe its instructions. This safety information is part of the manual. Keep the manual complete and accessible to personnel at all times!

- ▷ Use the mounted handle or the recessed grips when moving the equipment.

## 2. Intended use / Ambient conditions:

- The equipment is not approved for operation in potentially explosive atmospheres and not approved for the pumping of potentially explosive atmospheres; exception (ATEX): see rating plate.

- The equipment is not suitable for pumping unstable substances and substances which react explosively under impact (mechanical stress) and/or when being exposed to elevated temperatures without air.

- The equipment is not suitable for pumping self-inflammable substances, substances which are inflammable without air, and explosive substances.

- The equipment is not approved for operation below ground.

- ▶ Do not use the equipment or any system parts on humans or animals.

- ▶ Prevent any part of the human body from coming into contact with vacuum.

- ▶ The individual components are to be con-

nected, combined and operated only according to their design or with original manufacturer's spare parts and accessories.

- ▶ Admissible ambient and gas inlet temperatures during operation: +10°C to +40°C (rotary vane pumps: +12°C to +40°C). Monitor the maximum temperatures if installing the equipment in a cabinet or a housing. Make sure ventilation is adequate. Install an external automatic ventilation system if necessary.

- ▶ Do not aspirate particles and dust.

- ▶ Do not pump liquids (except BVC systems).

- ▶ Do not use the equipment to generate pressure (except ME 4R (NT)).

- ▶ If pumping hot process gases, make sure that the maximum permitted gas inlet temperature is not exceeded.

- ▷ Do not pump substances which may form deposits inside the equipment. Check the pump chamber regularly and clean if necessary.

- ▷ Ensure that the materials of the wetted parts are compatible with the pumped substances. Consider interactions and chemical reactions of the pumped media.

- ▷ Adopt suitable measures when using the equipment outdoors, at altitudes of more than 2000 m above mean sea level, or amid conductive pollution or external condensation.

- ▷ Do not operate this product near flames.

Use the equipment and all system parts only as intended, that is, for generation, measurement, and control of vacuum in vessels designed for that purpose. The BVC systems are intended for aspiration of liquids as well.

Use the condenser (exhaust waste vapor condenser, immission condenser, dry ice condenser, EK Peltronic) only for condensation of vapors at the pump outlet or inlet.

## 3. Installation and operation of the equipment:

- Equipment with power plug must be plugged into an outlet that is properly grounded. Power cords must be undamaged and in compliance with regulations. Failure to connect to ground may result in deadly electrical shock. Equipment without fuse: provide a slow blow fuse according to their current draw.

- Never operate this equipment if it has a damaged cord or plug.
- Equipment without power plug has to be wired by an electrician.

- ▶ Keep the electrical power cord away from heated surfaces.
- ▷ Check the power source for voltage, phase, and frequency (see rating plate).
- ▷ Check that the voltage selection switch, if present, is set correctly. Unplug the equipment before setting the voltage selection switch! Note: If the equipment is switched on with wrong voltage selection, the equipment may be damaged!

Single-phase motors are protected by a self-hold thermal cutout in the winding. Three-phase motors are protected by a circuit-breaker isolating all poles in case of overload. Separate from power to reset.

- ▷ If the supply voltage is below nominal voltage, the lock of the thermal cutout may be impaired, and the pump may restart on its own after cooling down sufficiently. Take appropriate precautions if this may lead to a dangerous situation.

A temperature sensor at the circuit board protects the motor of pumps with 24V DC motor (also unit PC 3001). In case of excess temperature, the pump is switched off.

- Prevent the release of dangerous, toxic, explosive, corrosive, noxious or polluting fluids, vapors and gases. If needed, install an appropriate collecting and disposal system and take protective action for personnel, equipment, and environment.
- Attention: Dangerous and/or explosive mixtures may form if outside air and certain pumped substances react in the equipment or at the equipment's outlet.
- You must take suitable precautions to prevent any formation of explosive mixtures in the housing. In case, e.g., of a diaphragm failure, mechanically generated sparks, hot surfaces or static electricity may ignite these mixtures. Use inert gas for gas ballast or venting, if necessary.
- Drain appropriately or otherwise remove any potentially explosive mixtures at the outlet of the equipment (use antistatic lines, follow explosion-protection guidelines), or dilute them with inert gas to non-explosive concentrations.

- ▶ Comply with applicable regulations when disposing of chemicals. Take into consideration that chemicals may be contaminated. Take adequate precautions to protect people from the effects of dangerous substances. Use appropriate protective clothing and safety goggles.

- ▶ The residual leak rate of the equipment might render possible an exchange of gas, albeit extremely slight, between the environment and the vacuum system. In case of a leak, pumped substances might be released into the environment or into the pump housing or motor.

- ▶ Use only original manufacturer's spare parts and accessories. Otherwise the safety and performance of the equipment, as well as the electromagnetic compatibility of the equipment might be reduced.

- ▶ The CE mark or the cTÜVus mark (see rating plate) may be voided if not using original manufacturer's spare parts.

- ▷ Pay attention to the safety symbol "Hot surfaces" on the equipment. Adopt suitable measures to prevent any danger arising from hot surfaces. Provide a suitable contact guard if necessary.

- ▷ Provide a firm, level platform for the equipment. Check that the system which you are going to evacuate is mechanically stable. Check that all fittings are secure. Ensure a stable position for the pump without any mechanical contact other than the pump feet. Ensure that hose connections are gas-tight.
- ▷ Note: Flexible elements will shrink when evacuated.

- ▶ Failure of the equipment (e.g., due to power interruption), failure of connected components or of parts of the supply (e.g., coolant), or a change of process parameters as well as the switching of valves or pumps must not be allowed to lead to a dangerous situation under any circumstances.

Ensure that the application will always turn in a safe status. Provide appropriate protective measures to allow for the possibility of failure and malfunction. Comply with all applicable and relevant safety requirements (regulations and guidelines). Implement the required actions and adopt suitable safety measures.

- ▷ Ensure that no foreign objects, no parts of your clothing, hair or fingers can be caught or drawn in at the inlet of the equipment. Never insert fingers or drop any other object into the inlet or outlet.

Use only hoses at the equipment's inlet and outlet with an inner diameter at least as large as the diameter of the equipment's tubing. Do not transfer forces (via hoses, tubing) onto the equipment's connections.

Allow the equipment to equilibrate to ambient temperature if you bring it from cold environment into a room prior to operation. Notice if there is water condensation on cold surfaces. Keep a minimum distance of 5 cm between the cooling fan and surrounding items.

Check fan regularly for dust/dirt and clean if necessary. Do not place the equipment on soft surfaces (e.g., rubber foam) during operation. This may reduce or block the fans' air supply. Do not cover the equipment.

- ▶ Due to the high compression ratio, the equipment may generate overpressure at the outlet. Check pressure compatibility with system components at the outlet.
- ▶ Do not permit any uncontrolled pressurizing. For example, a blocked exhaust line, closed exhaust isolation valve, or clogged silencer creates a risk of bursting! Always provide a free and pressureless exhaust outlet.
- ▶ Check compatibility with the maximally permitted pressure at inlet and outlet ports and at gas connections as well as with pressure differences (see manual). Do not operate the equipment with overpressure at the inlet.
- ▶ Avoid overpressure of more than 1.2 bar absolute in the event that inert gas is connected to the pump, to the gas ballast or to a venting valve. The inert gas has to be dry and clean.

Do not start the pump if the pressure difference between inlet and outlet exceeds 1.1 bar at maximum.

Prevent the backflow of condensates.

### 3.1 Systems with separator / exhaust waste vapor condenser:

- Cold surfaces can absorb a great amount of gas, which may expand abruptly in case of warming. This may lead to unmanageable overpressure in the system, creating a risk of bursting!

- ▶ Check the overpressure safety relief device at the exhaust waste vapor condenser at appropriate intervals.

- ▶ Ensure that the coolant outlet tubing is always unobstructed. Install an optional

coolant valve only in the supply line of the exhaust waste vapor condenser.

- ▶ There must be no damage on the glass surface. Do not use damaged components.
- ▶ Ensure sufficient air admittance to the EK Peltronic condenser.

Secure coolant hoses to prevent their accidental slipping.

Check liquid level in catchpots regularly and drain condensate in time to avoid overfilling.

### 3.2 Systems with cold trap:

See 3.1

Suitable coolants: e.g., Liquid nitrogen (LN<sub>2</sub>) or acetone-dry ice (CO<sub>2</sub>) mixtures.

- Prevent the formation of explosive or flammable mixtures. Oxygen may condense at the temperature of liquid nitrogen. Use inert gas for venting, if necessary.

- ▶ Check the coolant level of the cold trap at appropriate intervals. Coolant may overflow unexpectedly, e.g., in case of large amounts of gas.

- ▶ Comply with all applicable safety measures and requirements when using cryogenic coolants. Use only transport receptacles intended for coolants. Ensure pressure compensation between coolant receptacle and atmosphere at any time. Wear safety goggles and protective gloves.

### 3.3 Systems with vacuum measurement / vacuum control:

- ▶ At pressures above approximately 1060 mbar (795 Torr) the pressure reading becomes incorrect. Release pressure immediately to avoid a risk of bursting!

- ▶ Comply with maximum permissible pressures at the pressure transducer/vacuum gauge.

- ▶ You must take suitable precautions to prevent any dangerous situation from arising if the controller starts a pump, switches a coolant valve or opens a venting valve. Attention: If the controller is set to "Autostart: On", the process will start immediately after a power failure without pressing any further key. It is the user's responsibility to ensure that no dangerous status of the system due to the automatic start-up can occur and to provide appropriate safety measures.

Position pressure transducers, vacuum gauges and vacuum lines in such a way that

condensate cannot flow towards the pressure transducer or vacuum gauge.

#### 3.4 Systems with gas ballast:

- ▶ Make sure that air/gas intake through the gas ballast valve can never lead to hazardous, explosive or otherwise dangerous mixtures. If in doubt, use inert gas or close the gas ballast. Ensure that gas ballast line connections are gas tight.
- ▶ Pumping at high inlet pressure may lead to overpressure at the gas ballast valve. Pumped gases or condensate might be expelled if the valve is open.

#### 3.5 Systems with oil-sealed pumps:

- ▶ Check the pump's oil level and the oil's condition before starting the pump and subsequently at regular intervals, top up as necessary. Use only oil of the recommended type. Use special oils if pumping acid vapors, oxygen or strong oxidizing agents. Specifications are only valid if using standard oil. Adopt precautionary measures when handling pump fluids, lubricants, and solvents (e.g., safety-clothing, protective goggles).

#### 3.6 Aspiration systems BVC:

- ▶ Never use a damaged collecting bottle or expose it to vacuum.
- ▶ Use the BVC system only in combination with a working hydrophobic filter.
- ▶ During the aspiration process the collecting bottle of the BVC system must not contain any disinfecting solutions which corrode the materials of the pumping unit. In particular, no disinfectants may be used which release chlorine or oxygen radicals. Maximum liquid level in collecting bottle: approximately 80 %.

#### 3.7 Compressor pump ME 4R (NT):

- ▶ Do not generate overpressure in vessels other than in those designed for that purpose. Check the mechanical stability of the pressure vessel to avoid a risk of bursting! Install an overpressure safety relief device in the pressure system, if necessary. Note: Comply with maximum permissible pressure at the outlet (see manual). Flexible elements will expand when pressurized.

## 4. Maintenance and repair:

Ensure that maintenance is done only by suitably trained and supervised technicians.

- Before starting maintenance, disconnect the electrical power cord.
- Wait two minutes after isolating the equipment from power to allow the capacitors to discharge.
- Never operate the equipment if covers or other parts are disassembled. Ensure that the equipment cannot be operated accidentally.
- Note: The equipment may be contaminated with noxious or otherwise dangerous process chemicals, which have been pumped during operation. Ensure that the equipment is completely decontaminated or cleaned before maintenance commences.

- ▶ Vent the equipment before starting maintenance. Isolate the equipment and other components from the vacuum system. Allow sufficient cooling of the equipment. Drain condensate, if applicable. Motor capacitors are wear parts. Wear parts have to be replaced regularly. Motor capacitors have to be replaced by an electrician.

Return of equipment for repair only with fully completed Health and safety clearance form.

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# Declaration of Conformity – China RoHS 2

## DECLARATION OF CONFORMITY – China RoHS 2

VACUUBRAND GMBH + CO KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(+VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

### Environmentally Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmentally Friendly Use Period for VACUUBRAND products is 40 years.

此表格是按照SJ/T 11364-2014中规定所制定的。  
This table is created according to SJ/T 11364-2014.

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Declaration of Conformity – China RoHS 2

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
**Note:** Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。  
 O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。  
 X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

除上表所示信息外, 还需声明的是, 这些部件并非是有意图用铅 (Pb)、汞 (Hg)、铬 (Cd)、六价铬 (Cr(+VI))、多溴联苯 (PBB) 或多溴二苯醚 (PBDE) 来制造的。  
 Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (Cr+VI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by VACUUBRAND may enter into further devices (e.g., rotary evaporator) or can be used together with other appliances (e.g., usage as booster pumps).  
 With these products and appliances in particular, please note the EFUP labeled on these products. VACUUBRAND will not take responsibility for the EFUP of those products and appliances.

Place, date: Wertheim, 06/04/2020

  
 (Dr. F. Gitmans)  
 Managing Director

i.A.   
 (Dr. A. Wollschläger)  
 Regulatory Affairs Manager

**VACUUBRAND GMBH + CO KG**  
 Alfred-Zippe-Str. 4  
 97877 Wertheim  
 Germany

Tel.: +49 9342 808-0  
 Fax: +49 9342 808-5555  
 E-Mail: [info@vacuubrand.com](mailto:info@vacuubrand.com)  
 Web: [www.vacuubrand.com](http://www.vacuubrand.com)

**DOMINIQUE DUTSCHER SAS**

**DOMINIQUE DUTSCHER SAS**

**DOMINIQUE DUTSCHER SAS**

**VACUUBRAND®**

**Manufacturer:**

VACUUBRAND GMBH + CO KG  
Alfred-Zippe-Str. 4  
97877 Wertheim  
GERMANY

[www.vacuubrand.com](http://www.vacuubrand.com)

**Phone:**

- Head office: +49 9342 808-0
- Sales: +49 9342 808-5550
- Service: +49 9342 808-5660

Fax: +49 9342 808-5555

Email: [info@vacuubrand.com](mailto:info@vacuubrand.com)

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Copyright 2020