## Gebrauchsanleitung TPP - TubeSpin® Bioreaktor / Instruction for Use TPP TubeSpin® Bioreactor

### Anleitung

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### Instruction

The TubeSpin® Bioreactors are intended for performing media optimization and process development experiments of suspension cells in high numbers, cultivation in shakers at 37 °C temperature and centrifugation.

- *Introduction of medium: put tube into appropriate rack, fill medium as with the classic 50 ml centrifuge tube design using a pipette and allowing the medium to flow down the sides of the tube, please note the recommended filling volume.*
- *Closure: screw down filter screw cap tightly. Holes of partly different size above a PTFE membrane serve as sterile barrier enable optimal gas exchange with the surrounding atmosphere; and minimize liquid loss by evaporation.*
- *Cultivation: mass transfer is achieved by shaking in an appropriate incubation shaker (temperature, CO<sub>2</sub> and humidity controlled; regulation of shaking frequency).*
- *Shaker: Please refer to the relevant shaker information, for ex. [www.kuhner.ch](http://www.kuhner.ch)*
- *Centrifugation: use the applicable centrifuge rotor or adapter.*
- *Maximum centrifugal acceleration force (gx): can be influenced by following criteria: use of organic solutions, temperature below 0 °C, centrifuging without centrifuge adapter, etc.*

Information to the shanking culture can be found on [www.shakingtechnology.com](http://www.shakingtechnology.com).

### Caution

Consult the maximum centrifugal acceleration force (gx or RCF) as well as the necessity to use the applicable centrifuge adapter. Alcohol < 90 % plus rubbing can dissolve and/or blur the white color print on the tube. The TubeSpin® Bioreactors should only be inserted/withdrawn into/from the rack when the shaker is at a completely halt.

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**Instruction**

Racks are intended for stable hold f. ex. during manipulation and storage of tubes.

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**Rack # 99017, Rack „L“ # 99019**

- Working temperature: +121 to -190 °C, autoclavable.
- Identification: alpha numerical inscription
- Storage: can be disassembled and stored saving space.
- Tube size: fits similarly designed tubes 15 ml and 50 ml.

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**Rack # 99013**

- Working temperature: +121 to -190 °C, autoclavable.
- Identification: alpha-numerical inscription
- Tube size: fits 30 x TPP TubeSpin® Bioreactor 600.

Abmessungen	Measurements	99013	99017	99019
Länge mm	Length mm	366	209	209
Breite mm	Width mm	146	98	172
Höhe mm	Height mm	70	62	62
Röhrchen pro Gestell	Tubes per rack	3	18	30
Material	Materials		PP	
Einsatz-Temperatur	Working Temperature		+121 to -190 °C	