

Arcis DNA Sample Prep Kit

(Bulk Kit)

UFL002B Arcis DNA Sample Prep 50 rxn



Instructions for use


1. General Information (Bulk Kit)

Arcis DNA Sample Prep is a ready to use kit comprising two pre-filled tubes enabling pre-analytical processing of a variety of sample types including whole blood, saliva, bacteria, plant or animal cells or tissues.

In 3 minutes, with no prior sample preparation, the Arcis DNA Sample Prep kit allows you to go from cells to downstream nucleic acid investigations without the need for isolation or purification.

2. Materials Provided

| Material Provided | Quantity | Number of Preps |
|----------------------|----------|-----------------|
| Tube 1: Lysis Buffer | 1 Tube | 50 |
| Tube 2: Wash buffer | 1 Tube | |



3. Storage Conditions

Tubes are shipped and stored at room temperature. The bulk kit format has been designed for high-throughput genomic applications. We advise that the buffers are not exposed to the environment for extended periods and/or continual short exposures, if this occurs then unused buffer not used within 48 hours after opening should be disposed of. If required the buffers can be pre-aliquoted and then stored for later use.

4. Applications

The product can be used on virtually any DNA containing sample. The nucleic acid released has been successfully applied to molecular biology techniques including PCR, sequencing and cloning.

The Arcis DNA Sample Prep Kit accelerates pre-analytical processing in the following targets:

- Human Samples: including whole blood (fresh/frozen, with and without heparin or EDTA), plasma, urine and buccal swabs
- Microbiological samples: including bacteria (E.coli, S.aureus, P.aeruginosa, K. pneumoniae), viruses (HBV/HCV) and parasites (plasmodium)

5. Applications

The Arcis DNA Sample Prep Kit is particularly suited to the following areas:

- Laboratory-free situations including point of care and field-based diagnostics, near-patient testing, and integrated sample-processing workflows.
- Forensic Analysis
- Veterinary research
- Genomics and biomarker analysis

Instructions for Use continued

6. Typical Protocol- Rapid processing of whole blood

Ensure samples have thawed completely before starting this procedure.

- 6.1. Add 30µl of whole blood to 150µl of Reagent 1 (or scale up for larger sample volume).
- 6.2. Incubate for one minute at room temperature.
- 6.3. Take 10µl of the above mixture and combine with 40µl of Reagent 2 (or scale up for larger sample volume).
- 6.4. Add appropriate volume into PCR master mix (e.g. 5µl per 25µl reaction).

For applications where enhanced sensitivity is required, the following protocol modifications can be used:

Samples can be mixed with Reagent 1 at a 1:4 or 1:3 ratio to reduce sample dilution (See Table 1).

Samples that have been processed in step 6.1 can be added to Reagent 2 at 1:3, 1:2 or 1:1 ratio to reduce sample dilution (See Table 2).

Table 1: Processing samples in Reagent 1

| Sample Volume (µl) | Reagent 1 Volume (µl) | Ratio |
|--------------------|-----------------------|-------|
| 30 | 120 | 1:4 |
| 60 | 240 | 1:4 |
| 90 | 360 | 1:4 |
| 30 | 90 | 1:3 |
| 60 | 180 | 1:3 |
| 90 | 270 | 1:3 |

Table 2: Washing samples in Reagent 2

| Extract from Tube 1 (µl) | Reagent 2 Volume (µl) | Ratio |
|--------------------------|-----------------------|-------|
| 10 | 30 | 1:3 |
| 10 | 20 | 1:2 |
| 10 | 10 | 1:1 |

**This product is for laboratory research use only.
CAUTION: Not for diagnostic use.**

7. Supplier Contact Details

Cole-Parmer Ltd

Beacon Road, Stone,
Staffordshire,
ST15 0SA,
United Kingdom

Tel: +44 (0)1785 812121
Fax: +44 (0)1785 810405
Email: cpsales@coleparmer.com

