Instructions for use
Instructions d'utilisation
Instrucciones de uso
Istruzioni per l'uso
Instrukcja obsługi

NADAL[®] Giardia-Crypto Combo Test (test cassette)

REF 521007

Version A 1.0

2012-09-07

CE

IVD

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Carl-Peschken-Straße 9
47441 Moers
Germany

Friedenstraße 32
D-93053 Regensburg
Germany

Tel.: 0049 941 29010-0
Fax: 0049 941 29010-50

info@nal-vonminden.com
www.nal-vonminden.com

Directors:
Sandra von Minden
Roland Meißner
Thomas Zander

Commercial reg.: HRB 5679
VAT-ID-Nr.: DE 189 016 086
Tax-No: 244/133/00130

IBAN: DE98750500000840017073
BIC: BYLADEM1RBG
Bank: Sparkasse Regensburg

Point of Care Tests

Lab Assay Products

Lab Services

About us

EN	Instructions for use.....	3
FR	Instructions d'utilisation.....	6
E	Instrucciones de uso.....	9
I	Istruzioni per l'uso.....	12
PL	Instrukcja obsługi.....	15
	Symbols.....	19

1. Intended Use

The NADAL® Giardia-Crypto Combo test cassette by nal von minden is a qualitative chromatographic immunoassay rapid test for the detection of *Giardia* and *Cryptosporidium* antigens in stool samples. The test is intended for professionals in in-vitro diagnostic use only.

2. Clinical Significances

Cryptosporidium is a gastrointestinal protozoan parasite that infects humans and animals, causing a disease known as intestinal, tracheal and pulmonary cryptosporidiosis.

Giardia intestinalis, also called *Giardia lamblia*, is a gastrointestinal protozoan parasite causing a clinical picture that ranges from the symptom-free carrier to cases of acute or chronic diarrheic disease, especially in children. Both parasites are found throughout the world, particularly in areas with deficient sanitary conditions.

As does *Cryptosporidium*, *Giardia lamblia* infects through ingestion of its cysts, which can remain viable even in chlorinated water. When cysts reach the duodenum, each one releases up to 4 flagellated trophozoites (invasive form) that colonize the cells in the duodenum and proximal jejunum. The life cycle of *Cryptosporidium* is very short; each generation can develop and mature in 12-24 hours. This, together with its self-infectious character, leads to the intestinal tract being infected by a great number of parasites in only a few days. Diarrhoea, weight loss and abdominal cramps are the clinical symptoms of this disease, which can produce an electrolyte imbalance in immunosuppressed individuals. Pulmonary and tracheal cryptosporidiosis is associated with cough and fever, often accompanied by intestinal alterations. These symptoms generally develop 4-6 days after infection, although they can appear at any moment at 2-10 days.

In the case of *Giardia lamblia*, symptoms appear 1-3 weeks after contagion, following an incubation period of approximately 8 days. The most frequent presentation is characterized by diarrhoea, weight loss, colic-type abdominal pain and arrested growth and development. Symptom onset may be sudden or gradual; the disease can be self-limited or produce severe diarrhoea with poor intestinal absorption. In children, digestive function alterations are common. There are persistent symptoms, such as poor sugar absorption (particularly disaccharides like sucrose and lactose), fats and liposoluble vitamins. Lactose intolerance is common following acute infection and relapses or re-infection can appear. In rare cases, the trophozoites invade the intestinal wall, causing inflammatory diarrhoea and fever. Whether they produce toxic substances has not been documented.

3. Principle of the Test

The NADAL® Giardia-Crypto Combo test cassette is based on the immunological capture of colored microparticles as they pass through a membrane on which the monoclonal antibody has been immobilized.

NADAL® Giardia-Crypto Combo test cassette use monoclonal *Cryptosporidium parvum*- and *Giardia lamblia* specific antibodies that detect all forms of the parasites' life cycle.

The tests are based on the use of three microspheres: blue microspheres which are covalently linked to a monoclonal anti-*Cryptosporidium parvum* antibody, red microspheres which are covalently linked to a monoclonal anti-*Giardia lamblia* antibody, and green microspheres which are used as test control.

The parasites present in stool samples react with the latex particles which are coated with monoclonal antibodies specific against the antigen. This *latex particles/antibodies/parasite* complex migrates through a chromatographic process towards the reaction area. In this area, there are anti-*Cryptosporidium* and anti-*Giardia* antibodies, which react with the *latex particles/antibodies/parasite* complex. This reaction leads to the appearance of a blue line in the case of Crypto and a red line in the case of Giardia. These lines are used to interpret the result, at ten minutes' room-temperature incubation.

4. Reagents and Materials Supplied

The kit has the following components:

- 10 NADAL® Giardia-Crypto Combo test cassettes
- 10 Plastic tubes containing buffer solution
- 1 instruction sheet

Required Materials (not included in the kit):

- Stop watch or timer
- Pipettes that can be graduated (for liquid stool samples)
- Pipette tips (for liquid stool samples)

5. Storage & Stability

The NADAL® Giardia-Crypto Combo test kit can be stored at temperatures between 2°C and 30°C (it can be stored in a refrigerator). The expiration date is printed on the foil pouch and the outer package label.

6. Warnings and Precautions

- Patient samples (stool) may contain infectious agents and should be treated and discarded as potentially dangerous biological materials.
- The buffer contains sodium azide as an antimicrobial agent. Avoid direct contact with skin and mucosa. Discard appropriately. Do not use the buffer if there are indications of contamination or precipitation.

- Do not store or prepare food, eat, drink or smoke in the area where the reagents and samples are handled.
- Wear disposable gloves when handling the samples. Wash your hands thoroughly when you have finished working.
- Do not exchange components from kits with different LOT numbers.
- Use all the reagents only *in vitro*.
- Before using, let all kit components and samples reach room temperature, because cold reagents and/or samples can reduce test functionality. About 20-30 minutes are usually sufficient for reaching room temperature.
- Do not use kit components after the expiration dates.
- If the package is broken, the product can still be used as long as none of the components have been damaged.
- It is important to add the correct sample amount (4-5 drops of suspension). If less than indicated is used, the chromatography may not occur because the sample may not reach the reaction area; if there is too much, brown lines may appear instead of red and blue ones.
- The used product should be discarded as indicated by current legislation.
- Do not use the test if a colored line appears in the result area before you start to use it.
- It is important to take the appropriate amount of stool sample: 30-50 mg of solid stool or approx. 100 ul of liquid stool. An excessive sample amount prevents correct chromatography.


7. Specimen Collection and Preparation

Stool samples should be collected in an appropriate sterile container as soon as possible after symptom onset. The sample should be representative, if possible. The samples can be stored in the refrigerator (at approximately 4°C) for 1-2 days before being analyzed. To conserve samples for a longer period, store them in a freezer at -20°C, with no further handling. In this case, the sample must be completely defrosted, taken to room temperature and homogenized before analysis.


8. Procedure of the Test

Before testing: please let reach the tests and specimen room temperature (15 – 30°C). After opening the foil pouch the tests should be performed quickly, as they are sensitive to humidity.

1. Unscrew the vial cap carefully, so as not to spill the extraction buffer.


2. Collect a sufficient amount of fecal matter (30-50 mg) with the end of the applicator. Take sample from three different sample sites at least, in order to get a sample as much representative as possible.

If the stool is liquid, collect 100 microliters with the help of a pipette and transfer to the vial.

3. Introduce the applicator with the specimen into the vial. Screw the cap on tightly and **shake** vigorously to ensure a homogenous mixture.

4. Take the reaction device out of the aluminum pouch. Discard the small bag of desiccant; it is not used during the test, as its only purpose is to protect the test against moisture.

Break the upper end of the vial.

5. Add **4 or 5 drops** to the sample area of the reaction device (round window marked with an arrow). Do not add solid particles along with the liquid.

6. Wait exactly **10 minutes** and then read and interpret the results.

Any line that may appear after 10 minutes due to the nature of the sample has no diagnostic value.

9. Interpretation of the Results

To interpret the results please look at the line pattern that has formed in the result window marked with C und T1 and T2.

Positive Result:

Cryptosporidium:

A **GREEN** control line C and a separate **BLUE** test line T1 appear in the result window. The intensity can vary according to the antigen concentration present.

Giardia:

A **GREEN** control line C and a separate **RED/PINK** test line T2 appear in the result window. The intensity can vary according to the antigen concentration present.

Negative Result:

There is only a single **GREEN** control line C in the result window. This line should always appear.

Invalid Result:

The test is **INVALID** if no green control line appears because the procedure was not followed correctly, the reagents have deteriorated or because an incorrect amount of sample was added. Repeat the test with a new test cassette.

Note:

Any line that may appear after 10 minutes due to the nature of the sample has no diagnostic value.

Final diagnosis should not be based on only the result from one test. The diagnosis should be established by the correlation of test results with other appropriate data and with clinical symptoms.

10. Quality control

The test is invalid if no green line appears, whether because the test was not performed properly or the reagents have deteriorated. If this happens, repeat the analysis, following the work protocol indicated in these instruction sheets closely.

WARNING: Including a control with an established result is recommended to ensure that the data obtained are correct.

11. Limitations

- Final diagnosis should not be based on only the result from one test. The diagnosis should be established by the correlation of test results with other appropriate data and with clinical symptoms.
- The NADAL® Giardia-Crypto Combo test should be used only for detecting *Cryptosporidium parvum* and *Giardia lamblia* antigens in stool samples.
- This is a qualitative test; no quantitative interpretation of the result should be made with respect to the intensity of the positive line.
- A large number of samples were evaluated to ensure proper test functioning. The correlation of results with other techniques (ELISA) was excellent. However, interferences in the test performance cannot be excluded.
- If there is too much sample, brown lines may appear instead of red and blue ones, or there may be no test and control lines because chromatography is not produced correctly. The brown lines have no diagnostic value. If this happens, repeat the test with a smaller amount of stool sample or dilute the previously-prepared extract.
- No cross reaction with other substances has been observed during test evaluation. A negative result does not totally exclude the possibility of *Cryptosporidium parvum* and *Giardia lamblia* infection. The importance of the results should be evaluated in relation to the patient's clinical symptoms.
- The analysis of some samples can produce lines of imprecise color, caused by negative samples in the majority of the cases. The test should be repeated if lines of indeterminate color appear. If the same result is obtained again, the analysis should be performed with another analytic method.

12. Performance Characteristics

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

Assay 1 - A CE-marked ELISA product is taken as the product of reference. NADAL® Giardia-Crypto Combo test was evaluated against ELISA for detecting *Cryptosporidium* antigens.

The results were as follows:

		Ridascreen Crypto ELISA		Total
		+	--	
NADAL® Crypto	+	19	0	19
test cassette	--	8	53	61
Total		27	53	80

Concordance percent: 100% [(19+53)/80] = 90,0 %.

Estimated sensitivity concordance: 100% [19/(19+8)] = 70,4 % *

Estimated specificity concordance: 100% [53/(0+53)] > 99,9 %

*The 8 divergent samples were analyzed by PCR to determine the presence of *Cryptosporidium* DNA, with a negative result.

Assay 2 - Parasite detection via microscope was taken as the technique of reference. NADAL® Giardia-Crypto Combo test was evaluated against detecting *Cryptosporidium* antigens under the microscope.

The results were as follows:

		Microscopy		Total
		+	--	
NADAL® Crypto	+	15	0	15
test cassette	--	0	54	54
Total		15	54	69

Concordance percent: 100% [(15+54)/69] > 99,9 %.

Estimated sensitivity concordance: 100% [15/(15+0)] > 99,9 %

Estimated specificity concordance: 100% [54/(0+54)] > 99,9 %

Assay 3 - A CE-marked ELISA product is taken as the product of reference. NADAL® Giardia-Crypto Combo test was evaluated against ELISA for detecting *Giardia* antigens.

The results were as follows:

		Ridascreen Giardia ELISA		Total
		+	--	
NADAL® Giardia	+	26	0	26
test cassette	--	1	53	54
Total		27	53	80

Concordance percent: 100% [(26+53)/80] = 98,8 %.

Estimated sensitivity concordance: 100% [26/(26+1)] = 96,3 % *

Estimated specificity concordance: 100% [53/(0+53)] > 99,9 %

Assay 4 - Parasite detection via microscope was taken as the technique of reference. NADAL® Giardia-Crypto Combo test was evaluated against detecting *Giardia* antigens under the microscope.

The results were as follows:

		Microscopy		Total
		+	--	
NADAL® Giardia	+	15	1	16
test cassette	--	0	53	53
Total		15	54	69

Concordance percent: 100% [(15+53)/69] > 98,5 %.

Estimated sensitivity concordance: 100% [15/(15+0)] > 99,9 %

Estimated specificity concordance: 100% [53/(1+53)] > 98,1 %

REPEATABILITY- INTRA-ASSAY PRECISION

Five replicas of each concentration of the sensitivity curve were assayed with one lot, obtaining the same results.

REPRODUCIBILITY

INTER-DAY PRECISION

Using 1 product lot, 4 replicas of the sensitivity curve were performed over 10 consecutive days. Only a difference less than a 1:2 dilution was found, which is acceptable and tolerable for immunochromatography technique.

INTER-LABORATORY PRECISION

Three different laboratory technicians assayed those same samples, maintaining elevated precisions and concordances. Only a difference of a 1:2 dilution was found, which is acceptable and tolerable for immunochromatography assay.

INTER-LOT PRECISION

A sensitivity curve was performed in duplicate with 3 product lots. The analyses were performed by the same person on the same day. Only a difference less than a 1:2 dilution was found, which is acceptable and tolerable for immunochromatography technique.

HOOK EFFECT

No known concentration standard was available. However, the analyses of the sensitivity curve based on the maximum concentration possible of the internal standard did not show the existence of any Hook effect.

INTERFERING SUBSTANCES

The substances indicated in the table, at the concentrations specified, did not interfere with the result. One lot was used to carry out the study.

Metronidazole	1,2 mg/ml
Racecadotril	0,72 mg/ml
Loperamide	24 µg/ml
Braun Atropine	1,68 µg/ml
Cimetidine	1,92 mg/ml
Omeprazole	0,144 mg/ml
Neomycin	6,62 mg/ml
Ampicillin	7,2 mg/ml

Ibuprofen	5,76 mg/ml
Acetylsalicylic Acid	4,8 mg/ml
Sucrose	8 mg/ml
Blood	10% (v/v)
Metronidazole	1,2 mg/ml

INTERFERING MICROORGANISMS

The microorganisms indicated below did not cause any interference in the result.

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13. Références

1. Clark D.P. "New insights into human Cryptosporidiosis". *Clinical Microbiology Reviews* (1999). Vol 12, n° 14, p. 554-563.
2. Garcia L.S. y Shimizu R.Y. "Detection of *Giardia lamblia* and *Cryptosporidium parvum* antigens in human fecal specimens using the ColorPAC combination rapid solidphase qualitative immunochromatographic assay". *Journal of Clinical Microbiology* (2000). Vol 38, n° 3, p. 1267-1268.
3. Sharp S.E. et al. "Evaluation of the Triage Micro Parasite Panel for detection of *Giardia lamblia*, *Entamoeba histolytica/Entamoeba dispar*, and *Cryptosporidium parvum* in patient stool specimens". *Journal of Clinical Microbiology* (2001). Vol 39, n° 1, p. 332-334.
4. Katanik M.T. et al. "Evaluation of ColorPac Giardia/Cryptosporidium rapid assay and ProSpecT Giardia/Cryptosporidium microplate assay for detection of *Giardia* and *Cryptosporidium* in fecal specimens". *Journal of Clinical Microbiology* (2001). Vol 39, n° 12, p. 4523-4525.

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13. Bibliografia

1. Clark D.P. "New insights into human Cryptosporidiosis". *Clinical Microbiology Reviews* (1999). Vol 12, n° 14, p. 554-563.
2. Garcia L.S. y Shimizu R.Y. "Detection of *Giardia lamblia* and *Cryptosporidium parvum* antigens in human fecal specimens using the ColorPAC combination rapid solidphase qualitative immunochromatographic assay". *Journal of Clinical Microbiology* (2000). Vol 38, n° 3, p. 1267-1268.
3. Sharp S.E. et al. "Evaluation of the Triage Micro Parasite Panel for detection of *Giardia lamblia*, *Entamoeba histolytica/Entamoeba dispar*, and *Cryptosporidium parvum* in patient stool specimens". *Journal of Clinical Microbiology* (2001). Vol 39, n° 1, p. 332-334.
4. Katanik M.T. et al. "Evaluation of ColorPac Giardia/Cryptosporidium rapid assay and ProSpecT Giardia/Cryptosporidium microplate assay for detection of Giardia and Cryptosporidium in fecal specimens". *Journal of Clinical Microbiology* (2001). Vol 39, n° 12, p. 4523-4525.

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nal von minden GmbH