



Amersham Cy3B mono-reactive dye

Reagents for the labeling of biological compounds with Cy3B monofunctional dyes

Product Booklet

Table of Contents

1	Introduction	3
2	Components	5
3	Background	5
4	Example protocols	13
5	References	21

1 Introduction

Product codes

PA63100

PA63101

PA63106

PA63120

PA63121

PA63126

PA63130

PA63131

PA63136

Important

Read these instructions carefully before using the products.

Intended use

The products are intended for research use only, and shall not be used in any clinical or *in vitro* procedures for diagnostic purposes.

Safety

For use and handling of the products in a safe way, refer to the Safety Data Sheets.



CAUTION

These dyes are intensely colored and very reactive. Care should be exercised when handling the dye vials to avoid staining clothing, skin, and other items.

Note: *This article contains example protocols. It is the responsibility of the user to design and optimise protocols that are appropriate for the target compound that is being labeled. These example protocols utilise chemicals that may be hazardous, and should only be performed by appropriately qualified and trained persons.*

Storage

Store refrigerated at 2–8°C in the dark. Do not use if desiccant capsule in foil pack is either pink or green. Aqueous solutions of Cy3B NHS esters and maleimides are readily hydrolyzed back to the free acid. Therefore, do not store aliquots of aqueous Cy3B esters or maleimide solutions. Use immediately and discard residues. Aliquots of Cy3B NHS esters, hydrazides and maleimide in anhydrous DMSO are more stable and may be stored at -20°C, but for no longer than 2 weeks. Aliquots of Cy3B hydrazide in aqueous solutions appear to be fairly stable at neutral pH, therefore, it is recommended that aliquots of aqueous Cy3B hydrazide solutions can be stored for up to 2 weeks at -20°C.

2 Components

- Foil packs, each containing 1, 5, or 25 mg dried dye.
- Product specification sheet with instructions for using the dye.
- Reconstitute the material to 1 mg/ml in anhydrous DMF/ DMSO. Reconstituted material may be stored for up to 2 weeks at -20°C in aliquots to avoid repeat freeze thaw cycles.

3 Background

Cyanine reagents have been shown to be useful as fluorescent labels for biological compounds (1,2). These dyes are intensely fluorescent and highly water soluble, providing significant advantages over other existing fluorophores (3).

The CyDye™ are fluorescent cyanine compounds that produce an intense signal easily detected using appropriate detection equipment. The CyDye supplied here are monofunctional NHS-esters, maleimides or hydrazides and are provided in a dried, pre-measured form ready for the labeling of compounds containing free amino groups, thiol groups or carbonyl groups. The Cy3B dye is orange fluorescing cyanine that produces an intense signal easily detected using most rhodamine filter sets in the appropriate instrumentation.




Fig 1. Cy3B NHS ester monofunctional dye




Fig 2. Cy3B maleimide monofunctional dye




Fig 3. Cy3B hydrazide monofunctional dye

Labeling with Cy3B NHS esters

The most convenient and widely used functional group for the labeling of peptides and proteins is the primary amino group provided by the ϵ -amino group of lysine or the N-terminal amino group. Lysine is a relatively common amino acid and most proteins will have at least one. In many cases one or more lysine residues will be accessible to labeling reagents (5). The most useful reaction for labeling at amino groups is acylation. For maximum convenience, stable active esters that may be stored as solid materials, particularly NHS esters, have been extensively used over many years for the acylation of amino groups. The labeling of proteins is generally performed in an aqueous buffer; hydrolysis of the NHS ester is the major competing reaction of the acylation reaction. The rate of hydrolysis can be increased by raising the pH and by using dilute solutions of proteins. The pH affects the balance of the NHS ester hydrolysis rate versus the rate of reaction with primary amines. However, pH values between 7 and 9 are commonly used for most protein labeling reactions, together with phosphate, bicarbonate/carbonate and borate buffers. Others may be employed, but they should not contain a source of primary or secondary amines, e.g. Tris. These general principles apply to labeling reactions when using Cy3B NHS esters.

Labeling with CyDye hydrazides

A convenient and widely used functional group for the labeling of free carbonyl groups of glycoproteins and carbohydrates is the hydrazide. For maximum convenience, stable active hydrazide may be stored as solid materials. The labeling of glycoproteins and carbohydrates is generally performed in an aqueous buffer; hydrolysis of the hydrazide group is the major competing reaction.

Labeling with CyDye maleimides

A convenient and widely used functional group for the labeling of free thiol groups of peptides proteins, and oligonucleotides is the maleimide. Cysteine is the most common amino acid residue containing a thiol reactive group. For maximum convenience, stable active maleimides may be stored as solid materials.

Altering the protein concentration and reaction pH will change the labeling efficiency of the reaction. The optimal pH for the reaction of maleimides is near 7.0. In the pH range 7.0–7.5 the protein thiol groups are sufficiently nucleophilic so that they almost exclusively react with the maleimide dye in the presence of the more numerous protein amines, which are protonated and relatively unreactive.

Labeling of antibodies and biologically active proteins

When labeling antibodies or other proteins with NHS esters, maleimides or hydrazides the optimum conditions have to be established experimentally. The extent of labeling to give maximum fluorescence between the different CyDye vary and should be taken into account, as shown by Gruber et al (2); it will be necessary to optimize the ratio of Cy3B NHS ester to protein and pH to give the final dye to protein (D/P) ratio that is required.

Waggoner and co-workers (4) investigated the labeling of antibodies using CyDye NHS esters. They found that the brightest antibodies had D/P ratios between 4 and 12; at higher D/P ratios self-quenching was observed. In a more recent study by the University of Linz, it has been found that a D/P ratio of 2:1 gives the brightest signal. The general comments given above about labeling with NHS esters, maleimides and hydrazides should be taken into consideration when the labeling protocol is being designed. Biologically active proteins will vary greatly in terms of their properties (size, morphology, solubility etc.) and these should be taken into account. These properties may affect the choice of separation method of the labeled protein from free dye. Each case has to be considered on its own merits. Methods of separating excess free dye from labeled antibody other than gel filtration (e.g. dialysis) may be used.

At Cytiva, an anti-glutathione-S-transferase (GST) polyclonal antibody has been labeled using different CyDye NHS ester to antibody ratios. Ratios of 1:1, 5:1, 10:1 and 20:1 gave final D/P ratios of 0.28:1, 1.16:1, 2.3:1 and 4.6:1 respectively; these results are in general agreement with data previously reported by Waggoner and co-workers (4).

The scale of the labeling reaction is another factor that will affect the degree of labeling obtained; with small-scale labelings (100 µg or less) poor recoveries obtained during the purification processing can be a significant problem. The dye should be accurately aliquoted in anhydrous DMSO solution for use in small scale labelings. It is important that the biological properties of the labeled protein are maintained, and there should be some way of determining this. Information in the literature may provide guidance on the particular protein being used. Generally, higher degrees of labeling are more likely to have an effect on the biological properties of the protein. In some cases, a lysine residue accessible to the labeling reagent may be critical for the biological properties of the protein.

Protocol 1 has been designed for the preparation of Cy™ - labeled IgG antibodies. It is designed to utilise CyDyes in the quantities provided in the bulk pack sizes, 1–25 mg, and label protein to a final molar dye/protein (D/P) ratio between 4 and 12. This assumes an average protein molecular weight of 155 K Daltons. Other D/P ratios can be obtained by using different amounts of protein and/or dye.

Note: *The following materials and procedures have been optimized for IgG antibodies. Other proteins may also be readily labeled, however, choice of buffers, separation media, and technique may need to be varied in order to produce optimal results.*

Altering the protein concentration and reaction pH will change the labeling efficiency of the reaction. Optimal labeling with NHS esters generally occurs at pH 9.3. Proteins have been successfully labeled with this dye at a pH as low as 7.3, however, labeling times must be significantly longer at lower pH. Higher protein concentrations usually increase labeling efficiency. Solutions of up to 10 mg/ml protein have produced good conjugation reactions.

Conjugation of dye to antibody

Empirically we have determined the reaction molar stoichiometry required to produce a dye protein ratio in the range 6–12. This is based on mW values of,

Antibody (Ab) = 155 K Daltons and, Cy3B NHS ester = 766
Cy3B Hydrazide = 688

These figures are used to determine an adjusted labeling factor, relating the required stoichiometry by weight, which can then be used to determine the quantity of Ab (in mg) required for labeling with the amount of dye reagent (also in mg).

Adjusted dye labeling factor.

It is also necessary to calculate a dye purity factor which takes into account the % chromophore and NHS ester content of a specific batch of reagent. This information can be found on the *Batch Analysis sheet* supplied with the product.

For Example

% Chromophore =

% NHS ester by HPLC =

$$\text{Dye Purity Factor} = \frac{(\% \text{ Chromophore})}{100} \frac{(\% \text{ NHS ester HPLC})}{100}$$

Example

% Chromophore = 95%

% NHS ester by HPLC = 82.9%

$$\text{Dye Purity Factor} = (95/100) (82.9/100) = 0.787$$

Now calculate the amount of antibody to use for the preparation based on the amount of dye in milligrams, the adjusted dye labeling factor (this figure is product specific please see table at the beginning of the dye calculation section) and the dye purity NHS ester to be used.

$$\text{Quantity of Antibody} = \frac{\text{mgs of dye} \times \text{dye purity factor}}{\text{Dye labeling factor}}$$

Example

$$\text{Quantity of Antibody} = \frac{8.62 \times 0.787}{0.071} = 95.5 \text{ mgs}$$

4 Example protocols

Labeling of 1 mg Goat IgG Antibody with Cy3B NHS ester

Antibody Labeling

Step	Action
1	Prepare a solution of Goat IgG antibody at 1 mg/ml in 0.1 M sodium borate buffer at pH 8.3.
2	Prepare a 10 mg/ml solution of Cy3B monofunctional NHS ester (MW 761) in Dimethylsulfoxide (DMSO)/ Dimethylformamide (DMF) [1:1] by adding 300 μ l DMSO/DMF solution to 3 mg Cy3B mono NHS ester. This will give the desired ratio of CyDye NHS ester to antibody of 3:1.
3	Add this solution gradually to the antibody solution while stirring. Stir the solution for a further 60 minutes at room temperature in the dark.
4	The reaction is stopped by the addition of 1.3 ml of Sodium Dihydrogen orthophosphate buffer pH << 7
5	Purify the labeled antibody by gel chromatography.

Purification of Labeled Antibody

Step	Action
1	Prewash a PD-10 column with 25 ml buffer pH 7.5.
2	Load sample onto column.

Step Action

- 3 Add a further 2 × 1 mls buffer.
- 4 Elute the labeled IgG from the column with water (2 ml).
- 5 The unbound dye is the eluted from the column with a further 10 ml water.

- IgG concentration is measured at 280 nm
$$[\text{IgG}] = \{A_{280} - k_{\text{corr dye}} * A_{\text{dye max}}\} / 210\,000$$
- Dye concentration is measured by $A_{\text{dye max}}$
$$[\text{dye}] = (A_{\text{dye max}} / e_{\text{dye}})$$

$$D/P_{\text{final}} = [\text{Cy3B}] / [\text{antibody}]$$

Note:

- a. Any Cy3B NHS ester solution should be used immediately. Do not aliquot Cy3B NHS esters in aqueous solutions for storage, always use immediately and discard any residues.*
 - b. This protocol is designed to give one label per antibody. As stated above, two dye labels per antibody gives the best signal. To achieve this it is suggested that the dye to antibody in the labeling reaction should be increased to 6:1.*
-

Labeling of an Antibody with Cy3B Hydrazide

Material supplied

- Foil packs, each containing dried dye to label 1 mg of antibody

- Product specification sheet with instructions for using the dye

Material required but not supplied

- Conjugation buffer: phosphate buffer pH 7.0

Note: *All solvents and buffers should be degassed. This protocol has been designed for the preparation of Cy3B-labeled IgG antibodies. It is designed to label 1 mg. This assumes an average protein molecular weight of 155 K Daltons.*

Oxidation of antibody

Step	Action
------	--------

- | | |
|---|--|
| 1 | Dissolve the antibody in coupling buffer (10 mg/ml). Place 1 mg of antibody to be labeled (0.1 ml) in a 0.5 ml microfuge tube and flush with nitrogen. |
| 2 | Add 150 μ l Sodium Periodate solution (0.1 M) to the antibody, seal tightly and mix thorough for 30–45 minutes in the absence of light at ambient temperature. |
| 3 | While the oxidation reaction is incubating, equilibrate the gel filtration column in coupling buffer.
Place the periodate reaction on the column: collect and combine the fractions containing antibody |

Coupling of the antibody with a Cy3B hydrazide

Step	Action
------	--------

- | | |
|---|---|
| 1 | Dissolve the Cy3B hydrazide in coupling buffer to afford a 0.02 M solution. |
| 2 | Combine equal volumes of antibody to Cy3B hydrazide solution and gently mix for 3 hours at ambient temperature or overnight at 4°C. |
| 3 | Purify the product using a 10 ml gel filtration column that has been pre-treated with PBS. Elute the product with PBS. |
| 4 | Isolate the first colored fractions and scan sample from 250–800 nm against a PBS blank. |
| 5 | Record the wavelength at the maximum absorbance 548 nm and the 280 nm absorbance. |
| 6 | Determine the dye/protein ratio using the calculation below. |

$$\text{Dye/protein} = \frac{(A_{\text{max}}) (\text{Antibody Extinction Coefficient})}{((A_{280}) - (\text{correction factor})(A_{\text{max}})) (\text{CyDye Extinction Coefficient})}$$

Correction factor for Cy3B = 0.08

The dye/protein ratio may vary depending on the antibody and equivalence of CyDye used. Increased reaction time could also lead to saturated labeling. Incubation times and dye/protein ratios must be determined for each biomolecule used.

Step Action

Note:

Unconjugated dye can also be separated from the labeled antibody by dialysis. Dialysis does not give as efficient and rapid separation as gel filtration. It is recommended that gel filtration be used whenever possible.

Labeling of an oligonucleotide with Cy3B Maleimide

Material supplied

- Foil packs, each containing dried dye to label 1 mg of antibody
- Product specification sheet with instructions for using the dye

Material required but not supplied

- 0.17 M phosphate buffer, pH 8.0
- 0.04 M Dithiothreitol (DTT) (USB 15397)
- NAP™-10 column (Sephadex™ G-25 DNA Grade contains 0.15% Kathron CG as preservative)
- 0.2 M Sodium phosphate, pH 7.2
- Anhydrous DMSO
- 3 M Sodium acetate
- 0.2 M Sodium phosphate, pH 7.2, containing 0.04 M DDT

Note: 1. *pH will change the labeling efficiency of the reaction. Optimal labeling occurs at pH 7.2 in degassed buffer*

2. *To minimize oxidation of the thiol, carry out in an oxygen free environment, i.e under nitrogen, using degassed solvents/buffers*

This protocol was used to label ~120 nmol of thiol modified oligonucleotide with Cy3B maleimide reagent.

Reduction of oligonucleotide disulphide linkage

Step	Action
------	--------

- | | |
|---|--|
| 1 | Incubate oligonucleotide in 0.17 M phosphate buffer, pH 8.0, containing 0.04 M DTT for 16 hours at ambient temperature to reduce the disulphide linkage. |
| 2 | Remove DTT and thiol by-products using a NAP-10 column and following manufacturer's instructions. During this procedure, sterile de-ionized water should be used in place of the recommended equilibration buffer. |
| 3 | Monitor each eluent fraction for the presence of DNA at 260 nm. |
| 4 | Isolate the appropriate fractions and concentrate in vacuo. |

Conjugation of a Cy3B to an oligonucleotide

Step	Action
------	--------

- | | |
|---|--|
| 1 | To the flask add 375 μ l sterile, de-ionized water and 125 μ l of labeling buffer (0.2 M sodium phosphate, pH 7.2, containing 0.04 M DDT) mix gently and transfer to screw top vial. |
|---|--|

Step Action

- 2** To 1 mg Cy3B maleimide reagent add 120 μ l anhydrous DMSO, then transfer 100 μ l to reaction vial and incubate, with gentle mixing, at ambient temperature for 120 minutes. The dye/oligonucleotide ratio during labeling is ~20:1.
 - 3** Separation of the product from the free Cy3B. Precipitate oligonucleotide by adding 60 ml 3 M sodium acetate buffer, pH 5.2, and 1.65 ml ice cold ethanol and incubate on dry ice or at -80°C for 60 minutes.
 - 4** Centrifuge in a microfuge for 15 minutes at 14 000 rpm.
 - 5** Remove supernatant and wash pellet with 70% (v/v) ice-cold ethanol.
 - 6** Dissolve pellet in 500 μ l sterile, de-ionized water and purify labeled fraction by HPLC. Monitor at the absorbance wavelength of Cy3B, 559 nm.
-

Monofunctional dye characteristics table

	Cy3BNHSEster	Cy3BHydrazide	Cy3B Maleimide
Formula weight	771	688	796
Absorbance max	559 nm	559 nm	559 nm
Extinction coefficient in max Methanol	130000 M ⁻¹ cm ⁻¹	130000 M ⁻¹ cm ⁻¹	130000 M ⁻¹ cm ⁻¹
Emission max	570 nm	570 nm	570 nm
Quantum yield in Methanol	>0.7	>0.7	>0.7




Fig 4. Cy3B dye absorption and fluorescence spectra

5 References

1. Mujumdar, R.B., Ernst, L.A., Mujumdar, S.R., Lewis, C.J. and Waggoner, A.S. Cyanine dye labeling reagents: sulfoindocyanine succinimidyl esters. *Bioconj. Chem.*, 4(2), 105-111 (1993).
2. Gruber, H.J., Hahn, C., D., Kada, G., Riener, C., K., Harms, G., S., Ahrer, W., Dax, T., G., Knaus, H-G. Anomalous fluorescence enhancement of Cy3 and Cy3.5 versus anomalous fluorescence loss of Cy5 and Cy7 upon covalent linking to IgG and noncovalent binding to Avidin. *Bioconjugate Chem.* 11, 696-704 (2000).
3. Wessendorf, M.W. and Brelje, T.C. *Histochemistry*, 98(2), 81-85, (1992).
4. Southwick, P.L., Ernst, L.A., Tauriello, E.W., Parker, S.R., Mujumdar, R.B., Mujumdar, S.R., Clever, H.A., and Waggoner, A.S. Cyanine dye labeling reagents - carboxymethylindocyanine succinimidyl esters. *Cytometry*, 11(3), 418-430, (1990).
5. "Introduction to Molecular Biology" G.H.Haggis, D.Michie, A.R.Muir, K.B.Roberts & PMB Walker, Longmans (Bristol) Green & Co. LTD, 1965.
6. At the University of Linz, an independent labeling study on a range of dye reagents has recently been concluded. The complete results of this study will be available from Cytiva at a later date.



cytiva.com

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

Amersham™, Cy, CyDye, NAP and Sephadex are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

The purchase of CyDye products includes a limited license to use the CyDye products for internal research and development but not for any commercial purposes. Cy and CyDye are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva. A license to use the Cy and CyDye trademarks for commercial purposes is subject to a separate license agreement with Cytiva. Commercial use shall include:

1. Sale, lease, license or other transfer of the material or any material derived or produced from it.
2. Sale, lease, license or other grant of rights to use this material or any material derived or produced from it.
3. Use of this material to perform services for a fee for third parties, including contract research and drug screening.

If you require a commercial license to use the Cy and CyDye trademarks, please contact LSlicensing@cytiva.com.

All other third-party trademarks are the property of their respective owners.

© 2020–2021 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

PA631001PL AE V:5 03/2021