

Datasheet

# Cryopan I

## Serum-free Freezing Medium

| Product   | Description   | Catalogue-No. | Size   |
|-----------|---|---------------|--------|
| Cryopan I | Serum-free freezing medium without human and animal components, includes DMSO | P07-92010     | 10 ml  |
|           |   | P07-92050     | 50 ml  |
|           |   | P07-92100     | 100 ml |
|           |   | P07-92500     | 500 ml |

### Product description

Cryopan I is a serum-free, ready-to-use freezing medium for animal and human cells (adherent and suspension cells). It contains no human or animal components and is especially suitable for cells from serum-free culture.

### Storage conditions

Storage conditions: - 20°C  
 Stability: 2 years from date of production  
 Filling: 10 ml, 50 ml, 100 ml, 500 ml, other sizes on request

### Composition


Cryopan I consists of a chemically defined and optimised mixture of salts, sugar, DMSO and additional antifreeze-substances. It contains no animal and human components.

### Suitability

Cryopan I is for the cryoconservation of human and animal cells.

### Special Advantages

Due to its serum-free formulation, Cryopan is especially suitable for the conservation of serum-free cultured cells. The optimised composition guarantees a high cell viability after the thawing process. Best results are attained for all cell types including primary human cells and human cell lines. As Cryopan I contains no animal or human components, a risk of contamination with viruses, BSE or disturbing proteins is excluded.



### Instructions for Use

### 1. Freezing cells with Cryopan I

For optimal results only vital cells in the log-growth phase should be used.

- Thaw Cryopan I and store till using at 2-8°C.
- Trypsinate adherent cells, transfer the cells into the culture medium, stop the trypsin activity with trypsin inhibitor and centrifuge.
- Discard the supernatant and wash the cell pellet in PBS (without Ca<sup>2+</sup>/Mg<sup>2+</sup>).
- After an additional centrifugation step (100 – 200 g, 5 - 10 minutes) transfer the cells into PBS and determine cell number and cell vitality by trypan blue cell viability staining.
- Centrifuge the cells again and discard the supernatant.
- Transfer the cells into the cold Cryopan I (5x10<sup>5</sup> - 2x10<sup>6</sup> cells/ml Cryopan I).
- Suspend the cells carefully mixing the suspension by repeated pipetting until there are no more cell clumps.
- Refill the cell suspension into labelled cryotubes (0,5 - 1,5ml/tube).
- Freeze the cells in an automatically or manually controlled cryo freezer. The optimal freezing rate is approximately 1°C/minute.
- Alternatively, put the tubes for 15 minutes into a refrigerator, so that the freezing medium can penetrate into the cells. After this step freeze the tubes at -20°C for 2 hours box and put them into the vapour phases of liquid nitrogen over night.
- Store cryotubes in a cryotank with liquid nitrogen.

### 2. Thawing cells

- Remove the cryotubes from the cryotank and thaw them as soon as possible in warm water (< 1 minute).
- Disinfect the exterior of the cryo tubes with alcohol and convict the cells under sterile conditions into a centrifuge tube with 10 ml growth medium and mix it carefully.
- Centrifuge the cells (150 -200 g, 5 - 10 minutes).
- Discard the supernatant and recover the cells into the designated culture medium. Determine the cell viability by an appropriate method, like FACS or trypan blue cell viability staining.

### **Technical support**

For technical support, questions or remarks please contact your local PAN-Biotech partner or the technical department of PAN-Biotech via email ([info@pan-biotech.com](mailto:info@pan-biotech.com)) or phone +49-8543-601630.

---

FOR RESEARCH USE ONLY! Not approved for human or animal diagnostic or therapeutic procedures.