

# Amersham AMDEX streptavidin-HRP

## Product Specification Sheet

### Introduction

#### Product code

RPN4401

#### Important

Read these instructions carefully before using the products.

#### Intended use

The products are intended for research use only, and shall not be used in any clinical or *in vitro* procedures for diagnostic purposes.

The product may not be used to perform or offer commercial services of any kind and may not be incorporated into a diagnostic kit of any type intended for resale.

#### Safety

For use and handling of the products in a safe way, refer to the Safety Data Sheets.

#### Storage

Store, in the form supplied, at 2-8°C.

#### Expiry

Please see label.

#### Specifications

<b>Description</b>	Streptavidin conjugated to multiple HRP molecules on a dextran backbone.
<b>Buffer</b>	Stabilized in 50% TRIS Buffer and 50% HRP Conjugate Stabilizer
<b>Conjugate concentration</b>	Approximately 1 mg HRP/mL.
<b>Batch number</b>	Please see label.
<b>Volume</b>	500 µL.
<b>Dilution</b>	For ELISA, it is recommended that the conjugate should be diluted in the range 1:4000 to 1:8000.

#### Use

The conjugate is supplied for use in ELISA systems.

It can be used in direct substitution for a standard streptavidin HRP conjugate, without the need to alter existing protocols.

It is important to optimise the dilution by setting up a range test. Dilution recommendations are given overleaf, but are issued for guidance only.

**Note:** *It is essential that the conjugate is diluted in a suitable buffer.*

For optimum results the following buffer is suggested:

Starting buffer: PBS pH 7.2, 0.1% BSA, 0.1% Tween™.

Optimise depending on application and signal: noise ratio to PBS pH 7.2, up to 3% BSA, up to 1% Tween.

### Troubleshooting

**Table 1. Symptom:** No signal

Diagnosis	Remedy
One or more reagents omitted.	Re-run the assay, taking care over additions.
One or more reagents used too dilute.	Re-run the assay using higher concentrations.

**Table 2. Symptom:** High non-specific binding (NSB)

Diagnosis	Remedy
Improperly washed plate.	Re-run the assay, taking special care over the wash steps (Thorough rinsing with a wash-bottle is preferable to using an automatic plate-washer. It is important that plates are banged firmly to remove all liquid from wells after each step).
One or more reagents used too concentrated.	Re-run the assay using lower concentrations.
Insufficient blocking agents present.	Re-run the assay in buffers with higher concentrations of Tween and BSA.

**Table 3. Symptom:** Optical density (OD) readings too high.

Diagnosis	Remedy
If standard readings within expected range, samples are too concentrated.	Re-run the assay after diluting the samples in a suitable buffer.
If all readings too high, substrate incubation too long.	Re-run the assay, allowing less time between additions of substrate and stop solutions.
Improperly washed plate.	Re-run the assay (see <i>Table 2</i> above).
Amdex conjugate too strong.	Dilute the conjugate further.

**Table 4. Symptom:** OD readings too low.

Diagnosis	Remedy
Samples too dilute to read on standard curve.	Samples outside assay range. Re-assay after concentrating down.
All readings (including standards) too low; substrate incubation too short.	Re-run the assay, allowing more time between additions of substrate and stop solutions.
Amdex conjugate too dilute.	Re-run the assay with the conjugate at a high concentration.

**Table 5. Symptom:** Poor standard curve shape.

Diagnosis	Remedy
Errors in diluting standards.	If the problem hampers reading of the assay, re-run with fresh standards.
Improperly washed plate.	Re-run the assay (see <a href="#">Table 2, on page 1</a> ).

cytiva.com

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

Amersham™ is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Tween is a trademark of Croda Group of Companies.

All other third-party trademarks are the property of their respective owners.

© 2020-2021 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](http://cytiva.com/contact)

29257029 AE V:4 01/2021

