



**illustra**

Nucleon Phytopure

Genomic DNA Extraction Kits

Product Booklet

# Table of Contents

1	Introduction .....	3
2	Components of the system .....	4
3	Description .....	4
4	Principle steps in PhytoPure extraction .....	7
5	Critical parameters .....	7
6	Additional equipment and solutions required .....	8
7	Protocols .....	8
8	Troubleshooting guide .....	17
9	Additional information .....	20
10	DNA extraction products .....	21

# 1 Introduction

## Product codes

RPN8510

RPN8511

## Important

Read these instructions carefully before using the products.

## Intended use

The products are intended for research use only, and shall not be used in any clinical or *in vitro* procedures for diagnostic purposes.

## Safety

All chemicals should be considered as potentially hazardous. For use and handling of the products in a safe way, refer to the Safety Data Sheets.



### CAUTION

The protocol requires the use of:

- Chloroform: carcinogen.
- Cat 3: harmful, irritant.
- Mercaptoethanol: toxic.
- Isopropanol: flammable.
- Ethanol: highly flammable.

## Storage

Store at room temperature.

## Stability

Nucleon™ kit components are stable for up to 18 months, 3 months once opened, when stored under the recommended conditions. Performance is consistent when stored under the recommended conditions using the recommended procedures.

## Expiry

For expiry date please refer to outer packaging label.

## 2 Components of the system

*illustra™ Nucleon PhytoPure, plant and fungal DNA extraction kits (RPN8510/8511)* consist of:

Component	Product Code	
	0.1 g samples RPN8510	1.0 g samples RPN8511
Reagent 1	31 mL	245 mL
Reagent 2	11 mL	85 mL
PhytoPure Resin	6 mL	12 mL

## 3 Description

Nucleon PhytoPure systems from Cytiva for extracting DNA from plant and fungal samples are capable of producing high yields of high quality DNA in a fraction of the time taken by conventional methods. Not only is it more efficient but is also considerably simpler.

While most plant DNA extraction techniques are effective in removing proteins, they are much less successful with polysaccharides. Polysaccharides are very common contaminants in plant DNA extracts, and often result in difficult-to-handle, 'slimy' DNA pellets. This problem is compounded as polysaccharides, particularly those of an anionic nature, can be inhibitory to the further enzymatic analysis of the DNA.

Nucleon PhytoPure DNA extraction systems have been developed specifically to solve these problems. The system protocol is rapid and eliminates the requirement for phenol or cetyltrimethylammonium bromide (CTAB). After breaking of the cell wall, the cells are lysed in a reagent containing potassium SDS which is known to form complexes with proteins and polysaccharides. Chloroform is then added along with the specially modified Nucleon PhytoPure proprietary resin. This resin covalently binds polysaccharides, resulting in a high quality DNA preparation at the end of the protocol.

## **Role of Nucleon PhytoPure resin in plant DNA extraction.**

The specially modified Nucleon PhytoPure resin suspension performs two vital functions during the extraction procedure.

Nucleon PhytoPure resin particles contain free boric acid ( $-B(OH)_2$ ) groups. This resin is reactive towards polysaccharides through a well established mechanism whereby the boric acid reacts with 1,2 dihydroxy compounds to yield a cyclic boric acid ester.


The polysaccharides are bound by Nucleon PhytoPure resin particles and therefore removed from the sample.

Nucleon PhytoPure resin forms a semi-solid stratum between a lower organic and upper DNA-containing aqueous phase during the extraction process. This facilitates DNA recovery without repeated partitions, ensuring high yields of high quality DNA suitable for further analysis by restriction enzyme digestions, RAPD and AFLP.

**Plant and Filamentous Fungi species from which DNA has been successfully extracted using PhytoPure.**

<i>Arabidopsis</i>	<i>Helianthus annuus</i>	<i>Pisum sativum</i>
<i>Araucaria araucaria</i>	<i>Helianthus tuberosus</i>	<i>Rhododendron spp</i>
<i>Beta vulgaris</i>	<i>Hevea brasiliensis</i>	<i>Salix spp</i>
<i>Brassica oleracea</i>	<i>Humulus lupulus</i>	<i>Solanum tuberosum</i>
<i>Brassica napus</i>	<i>Irvingia gabonensis</i>	<i>Sorghum seeds</i>
<i>Capsicum annuum</i>	<i>Lolium</i>	<i>Sphagnum spp</i>
<i>Capsicum frutescens</i>	<i>Lotus japonicus</i>	<i>Spinacea oleracea</i>
<i>Cereals (barley, maize, rye, wheat)</i>	<i>Lupins albus</i>	<i>Swietenia macrophyla</i>
<i>Cocos nucifera</i>	<i>Lycopersicon esculentum</i>	
<i>Cryptomeria</i>	<i>Malus spp</i>	
<i>Eucalyptus globulis/ grandis</i>	<i>Musa spp</i>	<i>Aspergillus niger</i>
<i>Fragaria x ananassa</i>	<i>Nicotiana</i>	<i>Mortierella alpina</i>
<i>Fucus</i>	<i>Phaseolus vulgaris</i>	<i>Colletotrichum gloeosporioides</i>
<i>Gossypium hirsutum</i>	<i>Pinus sylvestris</i>	<i>Septoria nodorum</i>

## 4 Principle steps in PhytoPure extraction



## 5 Critical parameters

- Carry out steps 1–3 as quickly as possible.
- Use chloroform stored at  $-20^{\circ}\text{C}$  as this has been found to increase the efficiency of removal of complexed proteins/polysaccharides.
- Spin speed of 1300 g allows the resin to form a barrier at the organic/aqueous interface, facilitating the removal of the upper phase without danger of contamination.

Using a higher spin speed may cause the resin to spin to the bottom of the tube, and more care must be taken in removing the upper phase.

## 6 Additional equipment and solutions required

### Equipment

- Microcentrifuge
- Bench top centrifuge
- Assorted range of high precision pipettes
- Pasteur pipettes
- Polypropylene centrifuge tubes
- Water bath
- Rotary mixer

### Solutions

- Chloroform, AnalaR grade or similar (stored at -20°C)
- Isopropanol, AnalaR grade or similar
- Ethanol, AnalaR grade or similar
- RNase (if required)
- Mercaptoethanol (if required)
- Liquid nitrogen / dry ice

## 7 Protocols

### **illustra Nucleon PhytoPure for small samples, 0.1 g (RPN8510)**

The protocol requires the use of cold isopropanol, 70% ethanol and chloroform stored at -20°C.



## Breaking of the cell wall

Step	Action
------	--------

- |   |  |
|---|--|
| 1 | Add three volumes of dry ice (or liquid nitrogen if preferred) to 0.1 g (fresh weight) of plant tissue which has been frozen at -20°C. |
|---|--|

**Note:**

*The volumes presented in the protocol are for fresh weight of tissue. If dried tissue is to be used, please reduce the weight of tissue to be extracted by approximately a factor of 5. If the user requires to extract larger amounts of tissue, the volumes in the protocol can be scaled-up proportionately to the increased weight of tissue. However, PhytoPure resin should be scaled up according to tube-size ie:*

Tube size (mL)	Volume of resin (µL)
1.5	100
5	200
10	300

- |   |  |
|---|--|
| 2 | Grind the tissue in the dry ice (or liquid nitrogen) to yield a free flowing powder.       |
| 3 | Transfer the powder, using a chilled spatula, to a suitable polypropylene centrifuge tube. |

## Cell Lysis

Step	Action
------	--------

- |   |   |
|---|---|
| 1 | Add 600 $\mu\text{L}$ of Reagent 1, ensuring that all the reagent ingredients are fully dissolved. (Optional incubations with mercaptoethanol and RNase should be performed at this stage, see note below). |
|---|---|

**Note:**

*This protocol has been found not to require mercaptoethanol as a constituent. If the user prefers to include this reagent, add it at this stage, to a concentration in Reagent 1 of 10 mM. The extracted DNA may contain small amounts of RNA. Should RNA-free DNA be required, perform an RNase digestion at this step by adding RNase to a concentration of 20  $\mu\text{g}/\text{mL}$  after the addition of Reagent 1. Incubate at 37°C for 30 minutes.*

- |   |   |
|---|---|
| 2 | Mix thoroughly with a spatula.  |
| 3 | Add 200 $\mu\text{L}$ of Reagent 2.   |
| 4 | Invert several times until a homogeneous mixture is obtained.   |
| 5 | Incubate at 65°C in a shaking water bath for 10 minutes. Alternatively, regular manual agitation during the incubation is sufficient. |
| 6 | Place sample on ice for 20 minutes.   |

## DNA extraction

Step	Action
------	--------

- |   |   |
|---|---|
| 1 | Remove sample from ice and add 500 $\mu$ L of chloroform which has been stored at $-20^{\circ}\text{C}$ . |
|---|---|

**Note:**

*Chloroform at  $-20^{\circ}\text{C}$  is most effective in aiding the removal of complexed proteins/polysaccharides.*

- |   |   |
|---|---|
| 2 | Add 100 $\mu$ L of Nucleon PhytoPure DNA extraction resin suspension. |
|---|---|

**Note:**

*Ensure the resin is fully suspended by vigorous shaking immediately before use. It is important to ensure that the resin bottle contains equal proportions of resin to buffer. Allow the resin to settle and add sterile TE buffer if necessary.*

- |   |  |
|---|--|
| 3 | Shake on a tilt shaker for 10 minutes, at room temperature. Alternatively, regular manual agitation during the incubation is sufficient. |
|---|--|

- |   |                                      |
|---|--------------------------------------|
| 4 | Centrifuge at 1300 g for 10 minutes. |
|---|--------------------------------------|

- |   |   |
|---|---|
| 5 | Without disturbing the Nucleon resin suspension layer, transfer (using a pipette) the upper DNA containing phase, above the brown resin layer, into a fresh tube. |
|---|---|

Step	Action
------	--------

**Note:**

*The upper phase may appear green and cloudy but does not affect the quality of the DNA. An additional spin at speeds greater than 1300 g may be used to clarify the transferred aqueous phase.*

## DNA precipitation

Step	Action
------	--------

1 Add an equal volume of cold isopropanol.

2 Gently invert the tube until DNA precipitates.

**Note:**

*Precipitated DNA may be hooked out at this stage, using a heat-sealed pipette. This DNA does not require a 70% ethanol wash and should be placed directly into TE buffer or sterile water.*

3 Centrifuge at a minimum of 4000 g for 5 minutes to pellet the DNA.

4 Wash the DNA pellet with cold 70% ethanol.

5 Centrifuge at a minimum of 4000 g for 5 minutes to pellet the DNA.

6 Discard the supernatant.

7 Air-dry the DNA pellet for 10 minutes. Remove any remaining ethanol droplets from the tube.

8 Resuspend the DNA in TE buffer or water as required.

---

**Step Action**

---

**Note:**

*An RNase digestion can be performed at this stage if required (see note in [Cell Lysis on page 10](#)).*

---

## illustra Nucleon PhytoPure for large samples, 1.0 g (RPN8511)

The protocol requires the use of cold isopropanol, 70% ethanol and chloroform stored at -20°C.

### Breaking of the cell wall

---

**Step Action**

---

- 1 Add three volumes of dry ice (or liquid nitrogen if preferred) to 1.0 g (fresh weight) of plant tissue which has been frozen at -20°C.

**Note:**

*The volumes presented in the protocol are for fresh weight of tissue. If dried tissue is to be used, please reduce the weight of tissue to be extracted by approximately a factor of 5. If the user requires to extract larger amounts of tissue, the volumes in the protocol can be scaled-up proportionately to the increased weight of tissue. However, PhytoPure resin should be scaled up according to tube-size ie:*

Tube size (mL)	Volume of resin (µL)
1.5	100
5	200

---

## Step Action

---

Tube size (mL)	Volume of resin ( $\mu$ L)
10	300

---

- 2 Grind the tissue in the dry ice (or liquid nitrogen) to yield a free flowing powder.
  - 3 Transfer the powder, using a chilled spatula, to a suitable polypropylene centrifuge tube.
- 

## Cell Lysis

### Step Action

---

- 1 Add 4.6 mL of Reagent 1, ensuring that all the reagent ingredients are fully dissolved. (Optional incubations with mercaptoethanol and RNase should be performed at this stage, see note below)

**Note:**

*This protocol has been found not to require mercaptoethanol as a constituent. If the user prefers to include this reagent, add it at this stage, to a concentration in Reagent 1 of 10 mM. The extracted DNA may contain small amounts of RNA. Should RNA-free DNA be required, perform an RNase digestion at this step by adding RNase to a concentration of 20  $\mu$ g/ml after the addition of Reagent 1. Incubate at 37°C for 30 minutes.*

- 2 Mix thoroughly with a spatula.

<b>Step</b>	<b>Action</b>
-------------	---------------

---

- |          |   |
|----------|---|
| <b>3</b> | Add 1.5 mL of Reagent 2.  |
| <b>4</b> | Invert several times until a homogeneous mixture is obtained.   |
| <b>5</b> | Incubate at 65°C in a shaking water bath for 10 minutes. Alternatively, regular manual agitation during the incubation is sufficient. |
| <b>6</b> | Place sample on ice for 20 minutes.   |
- 

## **DNA extraction**

<b>Step</b>	<b>Action</b>
-------------	---------------

---

- |          |   |
|----------|---|
| <b>1</b> | Remove sample from ice and add 2 mL of chloroform which has been stored at -20°C.<br><br><b>Note:</b><br><i>Chloroform at -20°C is most effective to aid the removal of complexed proteins/polysaccharides.</i>   |
| <b>2</b> | Add 200 µL of Nucleon PhytoPure DNA extraction resin suspension.<br><br><b>Note:</b><br><i>Ensure the resin is fully suspended by vigorous shaking immediately before use. It is important to ensure that the resin bottle contains equal proportions of resin to buffer. Allow the resin to settle and add sterile TE buffer if necessary.</i> |

<b>Step</b>	<b>Action</b>
-------------	---------------

---

- |          |   |
|----------|---|
| <b>3</b> | Shake on a tilt shaker for 10 minutes, at room temperature. Alternatively, regular manual agitation during the incubation is sufficient.                          |
| <b>4</b> | Centrifuge at 1300 g for 10 minutes.  |
| <b>5</b> | Without disturbing the Nucleon resin suspension layer, transfer (using a pipette) the upper DNA containing phase, above the brown resin layer, into a fresh tube. |

**Note:**

*The upper phase may appear green and cloudy but does not affect the quality of the DNA. An additional spin at speeds greater than 1300 g may be used to clarify the transferred aqueous phase.*

---

## **DNA precipitation**

<b>Step</b>	<b>Action</b>
-------------	---------------

---

- |          |  |
|----------|--|
| <b>1</b> | Add an equal volume of cold isopropanol.       |
| <b>2</b> | Gently invert the tube until DNA precipitates. |

**Note:**

*Precipitated DNA may be hooked out at this stage, using a heat-sealed pipette. This DNA does not require a 70% ethanol wash and should be placed directly into TE buffer or sterile water.*

- |          |  |
|----------|--|
| <b>3</b> | Centrifuge at a minimum of 4000 g for 5 minutes to pellet the DNA. |
|----------|--|



Step	Action
4	Wash the DNA pellet with cold 70% ethanol.
5	Centrifuge at a minimum of 4000 g for 5 minutes to pellet the DNA.
6	Discard the supernatant.
7	Air-dry the DNA pellet for 10 minutes. Remove any remaining ethanol droplets from the tube.
8	Resuspend the DNA in TE buffer or water as required.

**Note:**

*An RNase digestion can be performed at this stage if required (see note in [Cell Lysis on page 14](#)).*

## 8 Troubleshooting guide

Problem	Possible cause	Remedy
1. Extracted DNA fails to amplify by PCR or digest with restriction enzymes.	1. The DNA is impure requiring modifications of the procedure to reduce the contamination.	<p>1.1. Ensure that the plant tissue is fully ground at the start of the protocol.</p> <p>1.2. Ensure full resuspension of the ground material in Reagent 1 and that Reagent 2 is thoroughly mixed with the lysate. FULLY dissolved material is essential for maximum yield and purity.</p>

Problem	Possible cause	Remedy
2. On isopropanol precipitation, a large pellet forms which is obviously not comprised of DNA alone.	2. The DNA is impure requiring modifications of the procedure to reduce the contamination.	2. Check the amount of plant material being extracted does not exceed that recommended for each protocol. The ratio of reagent volumes to tissue is important. Do not exceed the stated amounts of tissue. Either reduce the amount of starting material or increase reagent volumes proportionately.
3. On resuspension of the DNA pellet in water or TE buffer the resulting solution is viscous and perhaps cloudy.	3. The DNA is impure requiring modification of the procedure to reduce the contamination.	3. A post extraction clean-up on the resuspended DNA solution at the end of the protocol may be performed. Add PhytoPure resin to a final concentration of 1%, mix for a few minutes and spin out the PhytoPure resin.
4. Large amounts of cellular debris present.	4. The system may be overloaded or complete lysis has not occurred.	4.1. See the solutions above relating to reagent mixing and amount of tissue being extracted.

Problem	Possible cause	Remedy
		<p>4.2. Perform an additional chloroform extraction after cooling on ice for 20 minutes. Centrifuge at 1300 g for 10 minutes and retain the aqueous upper phase. Re-extract with chloroform and PhytoPure resin as per the protocol. This will remove large amounts of cell debris prior to the addition of PhytoPure resin which would otherwise become entrapped in the debris preventing polysaccharide binding.</p> <p>4.3. Centrifuge at greater than 1300 g after PhytoPure resin addition and subsequent shaking as per protocol. This will cause PhytoPure resin to spin to the bottom of the tube rather than form a barrier at the interface. Care must therefore be taken when removing the aqueous upper phase.</p>

Problem	Possible cause	Remedy
5. DNA pellet appears brown.	5. Phenol / tannin mediated oxidation occurred.	5. In the wide variety of plant materials tested mercaptoethanol has not been found to be necessary. However, in plant tissues with very high levels of phenol / tannins the addition of 2-mercaptoethanol to the lysis buffers (Reagent 1 & 2 combined) to a final concentration of 10 mM may alleviate this problem.

---

## 9 Additional information

### Calculation of centrifugal force

To ensure that Nucleon protocols are universally applicable to all centrifuges, centrifugal force is expressed as g rather than rpm values. To convert rpm to g please refer to the rotor manufacturers manual. If this is not available use the formula illustrated below.

$$g = 1.12r (\text{rpm}/1000)^2$$

$$\text{rpm} = 1000\sqrt{g/1.12r}$$

r = maximum radius of the rotor in mm

## 10 DNA extraction products

### **illustra Nucleon HT,**

for hard tissue, and paraffin sections, 50 preparations of up to 25 mg per prep.

**RPN8509**

### **illustra Nucleon PhytoPure,**

for plant and fungal DNA extraction kit, 50 preparations of 0.1 g.

**RPN8510**

### **illustra Nucleon PhytoPure,**

for plant and fungal DNA extraction kit, 50 preparations of 1.0 g.

**RPN8511**

### **illustra Nucleon BACC1,**

for 50 preparations of 1 mL whole blood or cultured cells ( $1$  to  $3 \times 10^6$ ).

**RPN8501**

### **illustra Nucleon BACC2,**

for 50 preparations of 10 mL of whole blood or cultured cells ( $3 \times 10^6$  to  $1 \times 10^7$ ).

**RPN8502**

### **illustra Nucleon BACC3,**

for 50 preparations of 10 mL of whole blood.

**RPN8512**



cytiva.com

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

illustra is a trademark of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

Nucleon is a trademark of Tepnel Life Sciences PLC.

All other third-party trademarks are the property of their respective owners.

© 2020–2021 Cytiva

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

29645877 AA V:4 02/2021