

# illustra ExoProStar

## Enzymatic PCR and Sequence reaction clean

### Instruction

#### Introduction

#### Product code

29010216

#### Important

Read these instructions carefully before using the products.

#### Intended use

The products are intended for research use only, and shall not be used in any clinical or *in vitro* procedures for diagnostic purposes.

#### Safety


For use and handling of the products in a safe way, refer to the Safety Data Sheets.

#### Recommended Protocol

Step	Action
1	Remove the illustra™ Alkaline Phosphatase and Exonuclease I from the freezer and keep on ice whilst preparing the reaction.
2	Take a 5 µL Aliquot of the completed PCR reaction mix.
3	Add 1 µL of illustra Alkaline Phosphatase to the reaction mix.
4	Add 1 µL of illustra Exonuclease I to the reaction mix.
5	Incubate at 37°C for 15 minutes.
6	Incubate at 80°C for 15 minutes to inactivate the enzymes.

The PCR product is now ready for use in sequencing or other downstream applications as required.

This reaction is scalable for larger quantities of PCR reaction mix. Simply add Alkaline Phosphatase and Exonuclease I in proportion to the volume of PCR reaction mix required.



## cytiva.com

Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.

ExoProStar™ and illustra are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

All other third-party trademarks are the property of their respective owners.

© 2020–2021 Cytiva

Any use of Name of the software is subject to Cytiva Standard Software End-User License Agreement for Life Sciences Software Products. A copy of this Standard Software End-User License Agreement is available on request.

All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.

For local office contact information, visit [cytiva.com/contact](https://cytiva.com/contact)

29010216 AC V:7 02/2021

